

Extreme X8



Rev4.:10/2023


ORDER FORM

Errors and immissions excepted..

Billing Address Or Customer-ID		Shipping Address		☐ Identical To Billing Address
C.-ID:		C.-ID:		
Name:		Name:		
Street:		Street:		
Town/PC:		Town/PC:		
E-Mail:		E-Mail:		
Phone:		Phone:		
Order Information		User Weight		
☐ Order ☐ Calculation ☐ Built-4-Me®		Art. No.:		
Purchase Order No.:		User weight: kg		
Reimbursement Code:				
Order date: Day Month Year		Base Price £12,495		

Client Measurements

a) Top Of Head To Seat		mm
b) Top Of Shoulder To Seat		mm
c) Underarm To Seat		mm
d) Elbow To Seat		mm
e) Top Of Pelvis To Seat		mm
f) Shoulder Width		mm
g) Trunk Width		mm
i) Hip Width		mm
j) Actual Sitting Depth		mm
m) Foot Length		mm
n) Lower Leg Length		mm
o) Actual Seat Width		mm



Technical Data

User Weight Capacity	155 kg	Ground Clearance	100 mm (depends on the configuration)
Base Width	710 mm	Drive Wheels & Castors	4 x 355 mm (14") Drive Wheels
Chassis length	1030 mm	Charger	220 V/10 A
Drive Control	R-Net (120A Power Module)	Battery Size	80 Ah (108 Ah optional)
Drive Motors	10 km/h, 4 x 700 W, 4-Pole Motors	Turning radius	1415 - 1470 mm
Standard configuration	Including MPS back system, MPS seat system, maintenance-free 80 Ah batteries and 10 A charger		
Tyres	Maximum tyre pressure depending on user weight from 0.25 to 0.55 bar (see table in instruction manual)		

Drive Unit**Power Wheelchair**

☒ MAG010006	EXTREME X8	£ 12,495
---	-------------------	----------

Seating System

☐ MAGMPSSEATVINYL	MPS Seat in VINYL	£ 0
☐ MAGMPSSEATFABRIC	MPS Seat in FABRIC	£ 0

Section A – MPS Seat Size**Seat width**

☐ MAGSW405MM	405 mm	£ 0
☐ MAGSW455MM	455 mm	£ 0
☐ MAGSW510MM	510 mm	£ 0
☐ MAGSW560MM	560 mm	£ 0
☐ MAGSW610MM	610 mm	£ 0

Seat depth

☐ MAGSD405MM	405 mm	£ 0
☐ MAGSD455MM	455 mm	£ 0
☐ MAGSD510MM	510 mm	£ 0

Seat Tilt Modules

☐ MAG2852-KIT	Fixed Frame, 3° tilt (no power functions)	£ 0
☐ MAG3502-KIT	Fixed Frame, 3° tilt (only lighting)	£ 0

Seat Height

☐ MAG6552	19" (470mm)	£ 0
☐ MAG3323-PI-2	20" (495mm)	£ 0
☐ MAG3323-PI-4	21" (520mm)	£ 0

BACKREST**Recline Option**

Recline option for MPS seat and MPS back		£
☐ MAG3900-S-3901-F	Manual Recline (90° - 170°)	£ 0

MPS Back

☐ MAGMPSBACK	MPS back	£ 0
-------------------------------------	----------	-----

MPS Back UPH

☐ MAGFABRICBR	Fabric	£ 0
☐ MAGVINYLBR	Vinyl	£ 0

MPS Backrest – Sizes

€

Backrest width Backrest Height Minus Backrest Width Must Exceed 50 mm

☐ MAGBW405MM	405 mm	XXX,-
☐ MAGBW455MM	455 mm	XXX,-
☐ MAGBW510MM	510 mm	XXX,-
☐ MAGBW560MM	560 mm	
☐ MAGBW610MM	610 mm	

Backrest height Backrest Height Minus Backrest Width Must Exceed 50 mm

☐ MAGBH535MM	535 mm	XXX,-
☐ MAGBH585MM	585 mm	XXX,-

Push Handles

£

☐ MAG0057-PI	MPS Push Rail	135
☐ MAGPRHA	MPS Height adjustable handles	
☐ MAG0056-WOKOR-2	MPS Peg handles	

Backrest Bag Hooks

£

☐ MAG5724	MPS Bag hooks (pair) on backrest	35
----------------------------------	----------------------------------	----

HEADREST**Headrest**

£

☐ MAG2770VI	MPS Headrest, MPS Mount	0
------------------------------------	-------------------------	---

Headrest UPH

£

☐ MAGFABRICH	Fabric	0
☐ MAGVINYLHR	Vinyl	0

ARMREST**Armrest – Types**

£

☐ MAG3572	Armrest, Single Post	0
☐ MAG3994	Armrest, Double Post	160
☐ MAGFLEXIARM	Armrest, Flip Back	575
☐ MAG6503-6502	Armrest, Flip Back, with Locked Position	550
☐ MAG3331-6500	Armrest, Flip Back, with Locking Post	600
☐ MAG3727-6501	Armrest, Flip Back, with Locking Post and Drop Forward	700

Armrest – Heights

☐ MAGAH220MM	220 mm	☐ MAGAH260MM	260 mm	☐ MAGAH300MM	300 mm	☐ MAGAH340MM	340 mm	☐ MAGAH380MM	380 mm
☐ MAGAH230MM	230 mm	☐ MAGAH270MM	270 mm	☐ MAGAH310MM	310 mm	☐ MAGAH350MM	350 mm	☐ MAGAH390MM	390 mm
☐ MAGAH240MM	240 mm	☐ MAGAH280MM	280 mm	☐ MAGAH320MM	320 mm	☐ MAGAH360MM	360 mm		
☐ MAGAH250MM	250 mm	☐ MAGAH290MM	290 mm	☐ MAGAH330MM	330 mm	☐ MAGAH370MM	370 mm		

Armpad – Style

£

☐ MAGAP75-WOKOR-2	75 mm armpads	0
☐ MAG3971-3972	MPS armpads	80
☐ MAG3712	MPS flip-out armpads	320
☐ MAGAP50-2	50 mm armpads	55
☐ MAG6711	Waterfall armpads	160
☐ MAG3775	Gutter armpads	160
☐ MAG6750	Gutter armpads with elbow support pad	160
☐ MAG6751	Gutter armpads with hand pad	160
☐ MAG6752	Gutter armpads with hand pad and elbow support pad	160

REV1 (01/01/2026)

Armpad – UPH			£
☐	MAGAPVINYL	Vinyl	0
☐	MAGAPFABRIC	Fabric	0

Armrest Bag Hooks			£
☐	MAGHOOKA2RH	Right Side Armrest bag hooks	35
☐	MAGHOOKA2LH	Left Side Armrest bag hooks	35

POSITIONING

Lap Belt							£
☐	MAG4973-MPS-4974-R	Standard belt	500 mm	600 mm	700 mm	800 mm	0
☐	MAG4973-4974-PBPAD	Standard belt with 2 posture pads	500 mm	600 mm	700 mm	800 mm	100
☐	MAG6540	4 Point Belt	500 mm	600 mm	700 mm	800 mm	180
☐	MAG6540-2XPBPAD-4	4 Point Belt with 2 posture pads	500 mm	600 mm	700 mm	800 mm	185
☐	MAG4730	Retractable Belt					110
	☐ MAGLBLH	Left sided buckle					0
	☐ MAGLBRH	Right sided buckle					0

Harness			£
☐	MAGHARP	Harness standard	160
☐	MAGHARQ	Harness pediatric	160

Thoracic Support (pad and mount)			£
☐	MAG1506	MPS Adjustable Laterals	530

Hip Pads - In Vinyl			£
☐	MAGHPKIT	Adjustable Hip Pads	270

LEGREST

Legrest – Choice			£
☐	MAGCPLEGREST	Fixed Centre Mount Legrest	0
☐	MAG3528	Standard Swing away Legrest	0

Footplate			£
☐	MAG4055FIX	Centre Post: 1 piece fixed footplate	0
☐	MAG4055FIX-2	Centre Post: 2 piece fixed footplate	0
☐	MAG4055	Centre Post: 1 piece flip up footplate	0
☐	MAG4055SFP	Centre Post: 1 piece flip up, SHORT	0
☐	MAG4055-2	Centre Post: 2 piece flip up footplate	0
☐	MAG4055FP-2	Centre Post: 2 piece flip up, SHORT	0
☐	MAGFLIPUP	Swing-Away: 2 piece flip up	0

Footplate Options		£
☐ MAG4312	Swing Away: Footplate Lock	55
☐ MAG5323	Centre Post: Metal edges	55
☐ MAG4260	Abduction Block	55

Footplate Height							
☐ MAGLL300MM	300 mm	☐ MAGLL340MM	340 mm	☐ MAGLL380MM(2)	380 mm	☐ MAGLL420MM (2)	420 mm
☐ MAGLL310MM	310 mm	☐ MAGLL350MM	350 mm	☐ MAGLL390MM(2)	390 mm	☐ MAGLL430MM (2)	430 mm
☐ MAGLL320MM	320 mm	☐ MAGLL360MM	360 mm	☐ MAGLL400MM(2)	400 mm		
☐ MAGLL330MM	330 mm	☐ MAGLL370MM	370 mm	☐ MAGLL410MM(2)	410 mm		

Calf Support		£
☐ MAGCALFSUPPORT	Calf Support Not with MAGCPOWER	135

Knee Pads		£
☐ MAG4247	Knee Pads for Fixed Centre Mount	125
☐ MAG4660	Knee Pads for Standard Swing-Away Legrests	125
☐ MAG4104	Knee Pads for Elevating Swing-Away Legrests	125

CONTROLLER

Controller – R-Net		£
☐ MAGPG-JSMLEDL	Remote Joystick, R-Net, LED	0
☐ MAGPG-JSM-PI	Remote Joystick, R-Net, Colour LCD	160
☐ MAGPG-JSM2IRL	Remote Joystick, R-Net, Colour LCD Advanced, Infrared	315
☐ MAGPG-JSMIRBT	Remote Joystick, R-Net, Colour LCD Advanced, Infrared & Bluetooth	530

Joystick Position		£
☐ MAGJSRIGHT	Joystick Right	0
☐ MAGJSLEFT	Joystick Left	0
☐ MAGPGCYPGC2.03725	Interchangeable to Left and Right	215
☐ MAG3330JSM-RH	Joystick Mounted as Attendant Control - Right	315
☐ MAG3330JSM-LH	Joystick Mounted as Attendant Control - Left	315

Swingback Options		£
☐ MAG3083-KIT-OS	Swingback Outside (Magnets)	0
☐ MAG3083-KIT-IS	Swingback Inside (Magnets)	0
☐ MAG3082-KIT-OS	Swingback Outside (Locking Latch)	125
☐ MAG3082-KIT-IS	Swingback Inside (Locking Latch)	125

Joystick Mounting Position		£
☐ MAG3713	Height Adjustable Joystick Stem	105
☐ MAG3729SB	Quick Release Stem	0

Dual Control		£
☐ MAG3330AM-KIT	Dual Control	290
Dual Control Mounting Location		
☐ MAGACRH	Mounting - RHS	
☐ MAGACLH	Mounting - LHS	

Joystick Knob Option		£
☐ MAGKNOBSTD	Standard joystick knob	0
☐ MAGPG-KS80	Stick S80	80
☐ MAGPG-KCHIN	Chin C38	80
☐ MAGPG-KB32	Ball B32	80
☐ MAGPG-KB38	Ball B38	80
☐ MAGPG-KM48	Mushroom M48	80
☐ MAGPG-KTBAR	T-Bar T80	80
☐ MAGPG-KSB70	Soft Ball SB70	80

5 Switch Box		£
☐ MAG158709	5 Switch Box (Button)	320
☐ MAG158708	5 Switch Box (Toggle)	320

5 Switch Box Location		£
☐ MAG5430	Same side as joystick	0
☐ MAG5431	Opposite side to joystick	0
☐ MAG5433	Opposite side to joystick + socket	0
☐ MAG5432	As Attendant	300

Omni Control		£
☐ MAG5436	Omni 2 Input Control Module	1170
☐ MAGPG-OMNISCANNER	Omni 2 Input Control Module (Scanning Mode)	1170

Omni Control Mounting Location		£
☐ MAGOMNILH	Left	0
☐ MAGOMNIRH	Right	0

BASE OPTIONS

Tyre		£
☐ MAG7050	Pneumatic	0
☐ MAG7058	Neverflat for Extreme X8	225
☐ MAG7056	Extreme X8 Puncture Res Kit	600
X8 Drive wheel rim		
☐ MAG6671	Silver rims	0
☐ MAG6670	Black rims	0

LED Lights		£
☐ MAG5227	Omit lights	0
☐ MAG2857	Lighting mounted on seat	0
☐ MAG6530	Lighting mounted on base	0

Colour Selection		£
☐ MAGRED	Flame Red	0
☐ MAGCOWLBL	Ocean Blue	0
☐ MAGBLACK	Stealth Black	0
☐ MAGORANGE	Desert Tangerine	0
☐ MAGYELLOW	Sunrise Yellow	0
☐ MAGMGREEN	Forest Green	0
☐ MAGCOWLFC	Forest Camo	580
☐ MAGCOWLAC	Army Camo	580
☐ MAGCOWLFL	Flames Fury	580

Batteries		£
☐ MAG6668	Battery, 70 AH	0
☐ MAG6669	Battery, 90 AH	275

Chargers		£
☐ MAG160104	Charger, 10ah	0

Extreme X8 Base Options		£
☐ MAG2899	Crash tested tie downs	0
☐ MAG5420	X8 lifting rings	100
☐ MAG2944	Steering lock	165
☐ MAG4384	Extreme X8 stability roller	165

Accessories		£
☐ MAG4471	RAM® Cup holder	55
☐ MAG4627	Magic Mobility Steel cup holder (1)	105
☐ MAG4628	Magic Mobility Steel cup holder (Right sided) (1)	105
☐ MAG4705	Accessory charger (1x12V + 2xUSB) - side rail mounted, same side as joystick	160
☐ MAG6629	Accessory charger (2xUSB) - side rail mounted, opposite to joystick	160
☐ MAG4801	Accessory charger (2xUSB) - armrest mounted, same side as joystick	160
☐ MAG6594	Accessory charger (2xUSB) - armrest mounted, opposite to joystick	160
☐ MAG5360	Loud horn	55
☐ MAG4834	A4 Accessory bag (mounted on the backrest)	105
☐ MAG4834RH	A4 Armrest bag (mounted on the right armrest)	105
☐ MAG4834LH	A4 Armrest bag (mounted on the left armrest)	105
☐ MAG4834-2	A4 Armrest bags (mounted on the left and right armrest)	215
☐ MAG5320RHI	A5 Armrest bag (mounted internally, right armrest)	80
☐ MAG5320LHI	A5 Armrest bag (mounted internally, left armrest)	80
☐ MAG5320RHE	A5 Armrest bag (mounted externally, right armrest)	80
☐ MAG5320LHE	A5 Armrest bag (mounted externally, left armrest)	80
☐ MAG48345320I	A4 & A5 Armrest bags (A4 joystick side, A5 internal facing)	160
☐ MAG48345320	A4 & A5 Armrest bags (A4 joystick side, A5 external facing)	160
☐ MAGCOVER	Wheelchair cover	0

Tray Tables		£
☐ MAG2677	Tray Table, Swing Away, Small	335
☐ MAG7700	Tray Table, Swingover with Magnet	335

Comments For Adaption



Extreme X8

£

Base Price

Total Of Additional Equipment

Packaging And Transport

Exclusive VAT

VAT

Total Price

Date

Signature

Day

Month

Year



Sunrise Medical GmbH

Kahlbachring 2-4

69254 Malsch/HD

Tel: +49 (0) 7253-9800

Fax: +49 (0) 7253-980111

E-Mail: kundenservice@sunrisemedical.de

www.SunriseMedical.de

Dimension Matrix MPS-SEAT

◀ BACK



Armrest Bag

◀ BACK



Mounts under armrest
(unless otherwise requested)

If swingback joystick is selected it is generally more suitable to fit accessory bag to opposite side

© Magic 1901

Colors

◀ BACK



Camo Forest

MAG5527



Desert Tangerine

MAG5523



Flame Red

MAG5522



Flames Fury

MAG5526



Forest Green

MAG5521



Ocean Blue

MAG5528



Stealth Black

MAG5525



Sunrise Yellow

MAG5524

Drive Wheels Rims & Special Tyre Options

◀ BACK

MAG5924 Silver rims
MAG5927 Black rims

MAG5927




MAG5927 Extreme X8 Puncture Res Kit




MAG4384 Extreme X8 stability roller

Seat Tilt Modules

[◀ BACK](#)

Headrest

[◀ BACK](#)

MPS HEADREST

Armrests

← BACK

ARMREST TYPE (1)



STANDARD

Note: 25mm of extra adjustment is added to bottom



DOUBLE STEM

For bariatric users



FLEXI

Flip-up, paralleling


Note: Depending on requested armrest height the flexi arms will be assembled using slightly different parts. Therefore the range for each armrest will not be the full 220-440mm

Requested Height	Actual Range (Flexi on MPS back)	Requested Height	Actual Range (Flexi on Rehab back)
220-330mm	220-370mm	220-330mm	180-360mm
340-370mm	250-400mm	340-440mm	275-440mm
380-440mm	330-440mm		


ARMPAD ADJUSTMENTS (ALL ARMREST TYPES)



FLEXI ARMREST ADJUSTMENTS



ARMPADS




MPS




MPS Flip-out


Allows the arm pads to be in tight to the user while still allowing flip up (Flexi arms only)




75mm



50mm



Waterfall



Gutter


If MPS, 50mm, 75mm, gutter or waterfalls armrests are selected, the length of armpad fitted will generally be based on the seat depth. See below table for details. If specific armpad length is required please state this on order form

SEAT DEPTH ▶	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"	
50mm Arm Pad	255	255	305	305	305	305	305	305	380	380	380	380	380	← Armpad length
75mm Arm Pad	255	255	305	305	305	305	305	305	380	380	380	380	380	
Gutter Arm Pad	255	255	305	305	305	305	305	305	380	380	380	380	380	
Waterfall Arm Pad	255	255	305	305	305	305	305	305	380	380	380	380	380	
MPS Arm Pad	-	-	-	-	355	355	405	405	460	460	460	460	460	← Armpad length

magic
mobility

Joystick Mounting

◀ BACK



Swingback Joystick Release Options

Magnet	Latch
- Requires more force but no dexterity to release - Can potentially release if wheelchair is being driven erratically or over rough terrain - Releases more easily if impacted	- Requires minimal force but some dexterity to release - Will stay locked at all times unless latch is pushed - Less likely to release if impacted, potentially causing more damage - Can be locked in fully swing-back position (optional) - Custom levers can be screwed to latch to make releasing it easier for clients with less hand movement



JOYSTICK STEM OPTIONS



HEIGHT ADJUSTABLE JOYSTICK STEM



QUICK RELEASE MOUNT

Joystick Knobs

◀ BACK



Footplates

◀ BACK

FOOTPLATE SELECTION



1 PIECE FLIP UP



2 PIECE FLIP UP




1 PIECE FIXED

Fastening bolts can be removed to allow footplate to flip up if required in the future



2 PIECE FIXED

Can adjust each footplate to a different angle




SHORT FOOTPLATES




© Magis Mobility 2016

Legrest MK2 Footplate Lock




Knee Pads


◀ BACK



**Centre mount legrest
(fixed and power elevating)**
3 axes of adjustment




Power or gas swingaway legrest
Fixed kneepad attached to legrest pivot



Standard swingaway legrest
3 axes of adjustment

Rear Boot Plate

◀ BACK



Cup Holder

◀ BACK

