

FRONTIER V6 REHAB SEATING




ORDER FORM

Billing Address Or Customer-ID		Shipping Address		☐ Identical To Billing Address
C.-ID:		C.-ID:		
Name:		Name:		
Street:		Street:		
Town / PC:		Town / PC:		
E-Mail:		E-Mail:		
Phone:		Phone:		
Order Information		User Weight		
☐ Order ☐ Calculation ☐ Built-4-Me®		Max. Load: 182 kg (without seat lift and seat inclination) Max. Load: 155 kg (with seat lift and seat tilt)		
Purchase Order No.:		User Weight: kg		
Art.-Nr.: 24100000				

Client Measurements

a) Top Of Head To Seat		mm
b) Top Of Shoulder To Seat		mm
c) Underarm To Seat		mm
d) Elbow To Seat		mm
e) Top Of Pelvis To Seat		mm
f) Shoulder Width		mm
g) Trunk Width		mm
i) Hip Width		mm
j) Actual Sitting Depth		mm
m) Foot Length		mm
n) Lower Leg Length		mm
o) Actual Seat Width		mm



Technical Specifications

Max. User Weight	182 kg (without seat lift and seat inclination) 155 kg (with seat lift and seat tilt)	Ground Clearance	100 mm (depends on the configuration)
Base Width	710 mm	Drive Wheels	14"
Base Length	1150 mm	Castor Wheels	9"
Electronics	R-Net	Charger	220 V/10 A
Drive Motors	10 km/h, 2 x 800 W, 4-Pole Motors	Batteries	74 Ah
Seat Height	435 mm	Turning Radius	2000 mm
Static Stability	10°	Dynamic Stability	10°
Tyres	Maximum tyre pressure depending on user weight from 0.25 to 0.55 bar (see table of operating instructions)		

DRIVE BASE

Power Wheelchair

☒ MAG010002	FRONTIER V6 - MID-WHEEL DRIVE - ALL TERRAIN	£ 13,500
---	--	----------

Colour Selection

☐ MAGRED	Flame Red (Red)	£ 0
☐ MAGSBLUE	Space Blue (Dark Blue)	£ 0
☐ MAGBLACK	Black Onyx (Black)	£ 0
☐ MAGNGREEN	Lucent Lime (Neon Green)	£ 160
☐ MAGMGREEN	Forest Green (Green)	£ 160
☐ MAGAGREEN	Army Green (Olive Green)	£ 160
☐ MAGYELLOW	Sunrise Yellow (Yellow)	£ 160
☐ MAGORANGE	Burnt Orange (Orange)	£ 160
☐ MAGHPINK	Hot Pink (Pink)	£ 160
☐ MAGSILVER	Silver Streak (Silver)	£ 160
☐ MAGCOWLBL	Blaze Blue (Light Blue)	£ 160
☐ MAGDVIOLET	Ultra Violet (Violet)	£ 160

Base Options

☒ MAGWC19X8TDK	Transport Eyelets	£ 0
--	-------------------	-----

Drive Wheels - Size

☒ MAGWH14BSP-2	Drive Wheel, 14"	£ 0
--	------------------	-----

Drive Wheels - Rims

☐ MAGWH145BPS-PI-X	Rims, Silver	£ 0
☐ MAGWH145BPB-PI-X	Rims, Black	£ 0

Drive Wheels - Tyres

☐ MAGDTPNEUMATIC	Drive Wheel Tyres, Air, Lug Profile, Black	£ 0
☐ MAGNFLAT-2	Neverflat For V6 (Puncture-Proof Tyre)	£ 235
☐ MAG4707-4	Puncture Resistance Kit (Puncture-Proof Tyres, Insert Between Tyre And Tube)	£ 600

Castor Wheels

☐ MAGWH2XXX-2	Castor Wheel Tyres, Air, 9 "	£ 0
☐ MAGWH2XXX-4	Castor Wheel Tyres, Puncture Proof, 9 "	£ 0

LED Lights & Indicators

☒ MAG2855	LED light and indicator according to StVZO, mounted on the frame	£ 0
---	--	-----

Splash Guard

☒ MAG3199	Splash Guard, Drive Wheel For Frontier V6	£ 0
---	---	-----

SEATING SYSTEM

Seat Modules

		£
☐ MAGFIXEDTILT	Fixed Tilt, Preset to 3°	0
☐ MAG2851-KIT	Powered Tilt 0° - 50°	1,950
☐ MAG2850-KIT	Powered Lift & Tilt (30cm lift & 0° - 45° Tilt)	2,850
☐ MAG4711-KIT	Powered Tilt -13° - 45° Tilt	3,350
☐ MAG4710-KIT	Powered Lift & Tilt (30cm lift & -13° - 45° Tilt)	4,450

Seat Height

		£
☐ MAG2622-4	435 mm	0
☐ MAG3323-PI-2	460 mm (1)	0
☐ MAG3323-PI-4	485 mm (1)	0

Seating System

		£
☐ MAGREHABSEATPAN	Rehab seat	0

Seat Sizes

Seat Width

☐ MAGSW305MM	305 mm	☐ MAGSW430MM	430 mm	☐ MAGSW560MM	560 mm		0
☐ MAGSW330MM	330 mm	☐ MAGSW455MM	455 mm	☐ MAGSW585MM	585 mm		0
☐ MAGSW355MM	355 mm	☐ MAGSW480MM	480 mm	☐ MAGSW610MM	610 mm		0
☐ MAGSW380MM	380 mm	☐ MAGSW510MM	510 mm				0
☐ MAGSW405MM	405 mm	☐ MAGSW535MM	535 mm				0

Seat Depth

☐ MAGSD12MM	305 mm (2)	☐ MAGSD430MM	430 mm	☐ MAGSD560MM	560 mm		0
☐ MAGSD13MM	330 mm (2)	☐ MAGSD18MM	455 mm	☐ MAGSD585MM	585 mm		0
☐ MAGSD14MM	355 mm (2)	☐ MAGSD480MM	480 mm	☐ MAGSD610MM	610 mm		0
☐ MAGSD15MM	380 mm (2)	☐ MAGSD510MM	510 mm				0
☐ MAGSD405MM	405 mm	☐ MAGSD535MM	535 mm				0

SEAT

Seat Cushions

		£
☐ MAGWOSEATCUSH	Omit Cushion	0
☐ MAGCMRC	MPS Seat Cushion, Removable (3)	0
☐ MAGPADDEDSC	Seat Cushion Flat, Upholstered	235

Seat Cushion Covers

		£
☐ MAGSUFABRIC	Seat cushion, cover, fabric	0
☐ MAGSUVINYL	Seat cushion, cover, vinyl	0

(1) Note: Seat height kits are required for lower leg lengths of 390 - 430 mm

(2) Not available with SW 430mm to SW 610mm

(3) Available in sizes from SW 430mm

BACKREST

Backrest for Rehab Seat

		£
☐ MAGREHABCANES	Rehab back	0
☐ MAGREHABFOLD	Rehab back, folds down to the front (1)	375

Backrest Frame Type

		£
Combination rehab seat & rehab back		
☐ MAG3920S-3921-F	Tool Angle Adjustable Backrest, 0 ° - 15 ° (0 °, 7.5 °, 15 °) (1)	0
☐ MAG3924-S-3925-F	Manual Recline Backrest, 0 ° - 170 ° (stepless)	395
☐ MAG3928-S-3929-F	Powered Recline Backrest, 0° - 160°	1,175
☐ MAG3932-S-3933-F	Antishear Powered Recline Backrest, 0° - 155°	1,875
Combination rehab seat & rehab back, foldable to the front		
☐ MAG3922-S-3923-F	Tool Angle Adjustable Backrest, 0 ° - 15 ° (0 °, 7.5 °, 15 °) (1)	0
☐ MAG3926-S-3927-F	Manual Recline Backrest, 0 ° - 170 ° (stepless)	395
☐ MAG3930-S-3931-F	Powered Recline Backrest, 0° - 160°	1,175
☐ MAG3934-S-3935-F	Antishear Powered Recline Backrest, 0° - 155°	1,875

Back Cushion

		£
☐ MAGNOINSERT	No Backrest Cushion Required	0
☐ MAGMPSINSERT	Backrest Cushion, MPS (see Section A)	0
☐ MAGBACKRESTINSERT	Backrest Cushion, Rehab (see Section B)	535

SECTION A - MPS Back Sizes

		£
Backrest Width		
☐ MAGBW405MM	405 mm (For Seat Width 405, 430 and 455 mm)	0
☐ MAGBW430MM	430 mm (For Seat Width 405, 430, 455 and 480 mm)	0
☐ MAGBW455MM	455 mm (For Seat Width 405, 430, 455, 480 and 510 mm)	0
☐ MAGBW480MM	480 mm (For Seat Width 430, 455, 480, 510 and 535 mm)	0
☐ MAGBW510MM	510 mm (For Seat Width 455, 480, 510, 535 and 560 mm)	0
☐ MAGBW535MM	535 mm (For Seat Width 480, 510, 535, 560 and 585 mm)	0
☐ MAGBW560MM	560 mm (For Seat Width 510, 535, 560, 585 and 610 mm)	0
☐ MAGBW585MM	585 mm (For Seat Width 535, 560, 585 and 610 mm)	0
☐ MAGBW610MM	610 mm (For Seat Width 560, 585 and 610 mm)	0
Backrest Height		
☐ MAGBH430MM	430 mm	0
☐ MAGBH480MM	480 mm	0
☐ MAGBH535MM	535 mm	0
☐ MAGBH585MM	585 mm	0
☐ MAGBH635MM	635 mm	0

SECTION A - MPS Backrest Cover

		£
☐ MAGFABRICBR	MPS back cushion cover, Fabric	0
☐ MAGVINYLBR	MPS back cushion cover, Vinyl	0

(1) In combination with the armrest, flexible, the back can be adjusted from 0 ° - 17 ° (in 4 steps)

SECTION B - Rehab Backrest Cushion Sizes		£
☐ MAGBRMRLS	Large, 430 mm Height	0
☐ MAGBRMRL	Large, 480 mm Height	0
☐ MAGBRMRLXT	Large, 530 mm Height	0
☐ MAGBRMRMS	Medium, 430 mm Height	0
☐ MAGBRMRM	Medium, 480 mm Height	0
☐ MAGBRMRMXT	Medium, 530 mm Height	0
☐ MAGBRMRJ	Small, 350 mm Height	0
☐ MAGBRMRJXS	Small, 410 mm Height	0
SECTION B - Rehab Backrest Cover		£
☐ MAGFABRICBR	Rehab back cushion cover, Fabric	0
☐ MAGVINYLBR	Rehab back cushion cover, Vinyl	0
SECTION B - Side Cushions / Side Bolsters		£
Rehab Backrest, Large & Medium Only		
☐ MAGBOLSTERS	Side Cushion, flat contour, 127 mm deep	55
☐ MAGBOLSTERSTD	Side Cushion, standard contour, 177 mm deep	80
☐ MAGBOLSTERD	Side Cushion, deep contour, 237 mm deep	105
SECTION B - Lumbar Support		£
Rehab Backrest Cushion Only		
☐ MAGLUMBARS	Shallow, Padding 15 mm	0
☐ MAGLUMBARSTD	Standard, Padding 25 mm	0
☐ MAGLUMBARD	Deep, Padding 35 mm	0
☐ MAGNOLUMBAR	Omit	0
Bracket Options		£
☐ MAG3492-PI	Back cushions, removable, with quick release	0
Push Handles		£
☐ MAG0056-WORKOR-2	Push handles, standard	0
☐ MAG4028	Stabilising bar (can also be used as a push bar)	135
☐	Omit Push Handles	0
Backrest Accessories		£
☐ MAG2040-LANDR	Backpack Hook Mounted on the Back	35
☐ MAGHOOKA2-KIT	Backpack Hook, Mounted on the Stabiliser Bar	35

HEADREST

Headrest		£
Headrest		
☐ MAG1718-V	Headrest, MPS With MPS Mount and Fabric To Match Backrest (1)	0
☐ MAGREHABSTD	Headrest Pad Rehab, Standard	160
☐ MAGREHABWIDE	Headrest Pad Rehab, Wide	220
Headrest Bracket (1)		
☐ MAG4405	Headrest bracket adjustable (multi-axis)	215
☐ MAG4406	Headrest bracket 2 axis (single stem, 2 pivot points)	165
☐ MAG4407	Headrest bracket 3 axis (single stem, 3 pivot points)	190

Headrest – Covers		£
☐ MAGFABRICLHR	Fabric Cover	0
☐ MAGVINYLHR	Vinyl Cover	0
☐ MAGDATEXHR	Dartex	60
☐ MAGCOVH	Vinyl Headrest with Canvas Cover	80
☐ MAGLEATHERHR	Leather	115

ARMRESTS

Armrest– Type		£
☐ MAG3572	Armrest, Single Post T-Bar, Standard	0
☐ MAGFLEXIARM	Armrest, Flip Back (2)	575
☐ MAG3800ARLW	Armrest, Flip Back with Locking Post	600

Armrest – Height									
☐ MAGAH220MM	220 mm	☐ MAGAH260MM	260 mm	☐ MAGAH300MM	300 mm	☐ MAGAH340MM	340 mm	☐ MAGAH380MM	380 mm
☐ MAGAH230MM	230 mm	☐ MAGAH270MM	270 mm	☐ MAGAH310MM	310 mm	☐ MAGAH350MM	350 mm	☐ MAGAH390MM	390 mm
☐ MAGAH240MM	240 mm	☐ MAGAH280MM	280 mm	☐ MAGAH320MM	320 mm	☐ MAGAH360MM	360 mm		
☐ MAGAH250MM	250 mm	☐ MAGAH290MM	290 mm	☐ MAGAH330MM	330 mm	☐ MAGAH370MM	370 mm		

Armrest – Pad		£
Right Armpad		
☐ MAG3972RH	MPS Armrest Pad	80
☐ MAGAP50R	Armrest Pad, Standard, 50 mm Wide, Right	55
☐ MAGAP75R	Armrest Pad, Standard, 75 mm Wide, Right	0
☐ MAG3772R	Amputee Armpad, Right	160
☐ MAGAPGER	Amputee Armpad with Hand Pad, Right	160
☐ MAGAPGHER	Amputee Armpad with Hand Pad and Elbow Support Pad, Right	160
☐ MAGAPGHR	Amputee Armpad with Elbow Support Pad, Right	160
☐ MAGAPWR	Waterfall Armpad, Right	160
Left Armpad		
☐ MAG3972LH	MPS Armrest Pad	80
☐ MAGAP50L	Armrest Pad, Standard, 50 mm Wide, Left	55
☐ MAG75L	Armrest Pad, Standard, 75 mm Wide, Left	0
☐ MAG3772L	Amputee Armpad, Left	160
☐ MAGAPGEL	Amputee Armpad with Hand Pad, Left	160
☐ MAGAPGHEL	Amputee Armpad with Hand Pad and Elbow Support Pad, Left	160
☐ MAGAPGHL	Amputee Armpad with Elbow Support Pad, Left	160
☐ MAGAPWL	Waterfall Armpad, Left	160

(1) Only Available with MAGMPSINSERT Backrest Cushion MPS

(2) In combination with MAGREHABFOLD Rehab back, folds down to the front, the MAGFLEXIARM option, flexibly adjustable, can be folded up, is recommended

Armrest – Covers		£
☐ MAGAPFABRIC	Armrest Cover, Fabric	0
☐ MAGAPVINYL	Armrest Cover, Vinyl	0
☐ MAGCOVA	Armrest Cover, Vinyl with Canvas Cover	110
☐ MAGDARTEX	Armrest Cover, Dartex	115
☐ MAGAPLEATHER	Armrest Cover, Leather	230

Mounting Options		£
☐ MAGARWI	25 mm Inset (1)	60
☐ MAGARWO	25 mm Offset (1)	60

Armrest Accessories		£
☐ MAGHOOKA2LH	Armrest Baghook, Right Side	35
☐ MAGHOOKA2RH	Armrest Baghook, Left Side	35

POSITIONING & ACCESSORIES

Lap Belt		£
☒ MAG4973-MPS-4974-R	Standard Posture Belt	0
☐ MAG4973-4974-PBPAD	Standard Posture Belt with 2 Pads	100
☐ MAG4730	Retractable Lap Belt	100
☐ MAGPB4P	4-Point Seat Belt	160
☐ MAGPB4P-2xPBPAD	4-Point Seat Belt with 2 Pads	185
☐ MAG4975-MPS-4976-R	Seat Belt and Mounting Brackets (Crash Tested)	290

Lap Belt Size		£
☐ MAGPB500	500 mm	0
☐ MAGPB600	600 mm	0
☐ MAGPB700	700 mm	0
☐ MAGPB800	800 mm	0

Thoracic Support		£
☐ MAGTHBKT	Standard Thoracic Support Assembly - select mount required	
Thoracic Mount - Required		
☐ MAG1506	Adjustable Brackets	530
☐ MAGTHR82	Swing-Away Brackets	770

Hip Pads - Vinyl		£
☐ MAGHPKIT	Hip Pads on Adjustable Brackets	270

Harness		£
☐ MAGHARQ	Harness Standard	160
☐ MAGHARP	Harness Paediatric	160

Chest Strap - Vinyl		£
☐ MAGSHOULDERSTRAP	Shoulder Strap	60
☐ MAGCHESTSTRAPVINYL	Chest Strap	115

(1) Only Available with MAGFLEXIARM Armrest, Flip Back and MAG3800ARLW Armrest, Flip Back with Locking Post

CONTROLS & ELECTRONICS

R-Net		£
☐ MAGPG-JSMLEDL	R-Net, LED	0
☐ MAGPG-JSM-PI-HD	R-Net, LCD Standard	160
☐ MAGPG-JSM2IRL	R-Net, LCD Advanced, Integrated Infrared	315
☐ MAGPG-JSMIRBT-HD	R-Net, LCD Advanced, Integrated Infrared & Bluetooth	530

Joystick Mounting		£
☐ MAGJSRIGHT	Right Hand	0
☐ MAGJSLEFT	Left Hand	0
☐ MAG3329JSM-RH	Joystick Mounted as Attendant Control (Right Side)	285
☐ MAG3329JSM-LH	Joystick Mounted as Attendant Control (Left Side)	285

Control Mounting		£
☐ MAG3083-KIT-IS	Swing Away Mounting, Inside (with Magnets)	0
☐ MAG3083-KIT-OS	Swing Away Mounting, Outside (with Magnets)	0
☐ MAG3082-KIT-IS	Swing Away Mounting, Inside (with Locking Latch)	125
☐ MAG3082-KIT-OS	Swing Away Mounting, Outside (with Locking Latch)	125

Joystick Stem Options		£
☐ MAGSTEMSTD	Standard	0
☐ MAG3717-KIT	Height Adjustable Joystick Stem	85
☐ MAG3729STD	Quick Release Stem	105
☐ MAG3729SB	Quick Release Stem (with swingback)	185
☐ MAG4726-KIT	Height Adjustable Drop Down Joystick (Left Hand Side)	105
☐ MAG4737-KIT	Height Adjustable Drop Down Joystick (Right Hand Side)	105

Attendant Control		£
☐ MAG3329AM-KIT	Attendant Control for Rehab Back	290

Attendant Control Mount		£
☐ MAGACRH	Attendant Control Mounted Right Hand	315
☐ MAGAC LH	Attendant Control Mounted Left Hand	315

Joystick Knob		£
☐ MAGKNOBSTD	Standard Joystick Knob	0
☐ MAGPG-KB32	Ball B32	80
☐ MAGPG-KB38	Ball B38	80
☐ MAGPG-KCHIN	Chin C38	80
☐ MAGPG-KM48	Mushroom M48	80
☐ MAGPG-KS30	Stick S30	80
☐ MAGPG-KS80	Stick S80	80
☐ MAGPG-KSB70	Soft Ball SB70	80
☐ MAGPG-KT BAR	T-Bar T80	80

FOOTRESTS

Lower Leg Length

☐ MAGLL300MM	300 mm	☐ MAGLL340MM	340 mm	☐ MAGLL380MM	380 mm	☐ MAGLL420MM	420 mm(2)(3)	
☐ MAGLL310MM	310 mm	☐ MAGLL350MM	350 mm	☐ MAGLL390MM	390 mm (1)	☐ MAGLL430MM	430 mm(2)(3)	
☐ MAGLL320MM	320 mm	☐ MAGLL360MM	360 mm	☐ MAGLL400MM	400 mm (1)			
☐ MAGLL330MM	330 mm	☐ MAGLL370MM	370 mm	☐ MAGLL410MM	410 mm (2)			

Footrest – Type

£

Centre Mount

☐ MAG4053	Centre Mount, Inc. heel loops and calf pads	0
☐ MAG2260-KIT	Powered Centre Mount, Inc. calf pads	1,600

Swing Away Legrests

☐ MAG3528-4168xx	Swing Away Legrests, Inc. calf pads	0
☐ MAG3527-XX	Manual Elevating Legrests	595
☐ MAG3526	Powered Elevating Legrests, Inc. calf pads	1,750

Footplates

£

Centre Mount Legrests

☐ MAGFLIPUPC	One-piece footplate flip-up	0
☐ MAG44172FFU	Two-piece footplate flip-up	0

Swing Away Legrests

☐ MAGFOOTPLATESA	Two-piece footplate flip-up	0
---	-----------------------------	---

Metal Edges On Footplate

£

☐ MAGMETAEDGES	Metal Edges	55
---------------------------------------	-------------	----

Foot Straps

£

☐ MAGFOOTSBU	Foot Strap with Buckle	70
☐ MAGFOOTSVEL	Foot Strap with Velcro	70

Calf Support

£

☐ MAGCALFSUPPORT	Calf Support	135
☐ MAGCPADSTD	Standard	0

Calf Support Size - Required

☐ MAGSTDCALFS	Calf Support, Standard	0
☐ MAG4167	Calf Support, Extra Deep	115

Heel Ankle Supports

£

☐ MAGHL	Heel Loops	110
☐ MAGAC	Ankle Cuffs, Standard	175

Knee Pads

£

☐ MAG4104	Knee Pads for Elevating Swingaway Legrests	125
☐ MAG4247	Knee Pads for Swingaway Legrests	125

(1) For LLL 390 - 400 mm, the seat height kit +25 mm is required (MAG3323-PI-2, surcharge see seat heights)

(2) For LLL 410 - 430 mm, the seat height kit +50 mm is required (MAG3323-PI-4, for surcharge see seat heights)

(3) Not available with MAGFIXEDTILT Fixed Tilt, Preset to 3°

ACCESSORIES

Spare Drive Wheel		£
☐ MAGWH14GS-KOR	Grey 14 Inch Hybrid Solid Spare Wheel	115
☐ MAGWH145BPB-PI	Black Rims Spare Wheel	230
☐ MAGWH145BPS-PI	Silver Rims Spare Wheel	230
☐ MAGWH145BPBNF	Neverflat, Black Rims Spare Wheel	230
☐ MAGWH145BPSNF	Neverflat, Silver Rims Spare Wheel	230
☐ MAGPR1B	Puncture Resistant, Black Rims Spare Wheel	580
☐ MAGPR1S	Puncture Resistant, Silver Rims Spare Wheel	580
Spare Castor Wheel		£
☐ MAGCASTPNE	Spare Castor Wheel, Pneumatic	75
☐ MAGCASTSLD	Spare Castor Wheel, Solid	75
Wheelchair Jack		£
☐ MAG1595	Powerchair Jack	85
Cup Holder Mount		£
☐ MAG4471	RAM Cup Holder	55
Mobile Phone Holder Mount		£
☐ MAG4472	Mobile Phone Holder	205
Tablet Holder Mount		£
☐ MAG4473	Tablet Holder	345
☐ MAG4475	Swingaway Tablet Holder	445
Ram Camera Holder		£
☐ MAG4474	Action Camera Holder	175
☐ MAG4476	Swingaway Camera Holder	345
Loud Horn		£
☐ MAG4323	Loud Horn	55
Bag		£
☐ MAG5320LH	Armrest Bag, Mounted Externally To Left Armrest	80
☐ MAG5320RH	Armrest Bag, Mounted Externally To Right Armrest	80
☐ MAG5320ACB	A4/A5 Armrest Bags, Mounted Externally Joystick Side	160
☐ MAG5321ACB	A4/A5 Armrest Bags, Mounted Internally Joystick Side	160
☐ MAG5321LH	Armrest Bag, Mounted Internally To Left Armrest	80
☐ MAG5321RH	Armrest Bag, Mounted Internally To Right Armrest	80
☐ MAGACBBR	Accessory Bag, Mounted On The Backrest	105
☐ MAGACBLH	Armrest Bag, Mounted On The Left Armrest	105
☐ MAGACBRH	Armrest Bag, Mounted On The Right Armrest	105
☐ MAGACB-2	Armrest Bags, Mounted On The Left And Right Armrests	215
Key		£
☐ MAGPG-SK	Key	60
IV Pole		£
☐ MAG4137R-1831MPSLH	IV Pole, Mounted Left	230
☐ MAG4137R-1831MPSRH	IV Pole, Mounted Right	230

Ventilator Carrier		£
☐ MAGACVC	Ventilator Carrier	460
Oxygen Carrier		£
☐ MAGACOC	O2 Carrier	345
Tray Table		£
☐ MAG2677	Swingaway Table, Small (320 x 240mm)	335
Table Material - required		
☐ MAGCLEARTAB	Tray Table, Clear	0
Table Mounting		£
☐ MAGTABLELH	Table Mount, Left	0
☐ MAGTABLERH	Table Mount, Right	0
Luggage Rack		£
☐ MAG3198	Luggage Rack (max payload 10kg)	230
Charging Socket		£
☐ MAG4705	12 V Charging Socket (1x12 V + 2x USB)	160

Comments for B4M (please ask for a quote; delivery times may vary)

Magic Mobility Frontier V6

£

Base Price

Total Of Additional Equipment

Packaging And Transport

Exclusive VAT

VAT

Total Price**Date****Signature**

Day

Month

Year

**Sunrise Medical**

Thorns Road, Brierley Hill, West Midland

DY5 2LD

Phone: +44 (0) 845 605 66 88

Fax: +44 (0) 845 605 66 89

E-Mail: enquiries@sunmed.co.ukwww.SunriseMedical.co.uk