



Helio C2

ORDER FORM - INTERNATIONAL

1- Purchasing information

P.O. number:* _____ Date: _____

Company Name:* _____

Purchasing contact:* _____

Phone:* _____

Email:* _____

Contact preference:* Phone ☐ Email ☐

Marked for: _____

2- ATP/Clinician information

ATP/Clinician:* _____

Phone:* _____

Email:* _____

Contact preference:* Phone ☐ Email ☐ * Required

3- Shipping information

Ship to:* _____

Name: _____

Address: _____

City: _____ State: _____

Zip Code: _____

Phone: _____

MODEL														
WC05	☐	Ultralight Folding Carbon Fiber Helio C2 Wheelchair (120 kg - 265 lb weight capacity)												£3,595
WC09	☐	HD Heavy Duty package (159 kg - 350 lb weight capacity)												£395
WC10	☐	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6° camber												£35
Frame														
SEAT WIDTH ³ - cm (inches)														
SW	☐	35.56 (14")	☐	38.10 (15")	☐	40.64 (16")	☐	43.18 (17")	☐	45.72 (18")	☐	48.26 (19")	☐	NC
	☐	50.80 (20") ¹										53.34 (21") ²	55.88 (22") ²	£360
SEAT DEPTH														
SD	☐	35.56 (14")	☐	38.10 (15")	☐	40.64 (16")	☐	43.18 (17")	☐	45.72 (18")	☐	48.26 (19")	☐	NC
	☐	50.80 (20")											£360	
TRANSIT TIE DOWN														
AC04	☐	Transit Tie Down (Unoccupied)												NC
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴ Short belt 50-104 cm (20" - 41")												£285
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴ Long belt 50-116 cm (20" - 46")												£285
AC35	☐	WC-19 Transit Tie Down without belt ⁴												£279

1- Available with HD option; 2- Only available with HD option

[3. Width calculator](#)

4- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option



PACKAGES							
PK01	☐	Newton Ultralight Package (include: Carbon fiber angle adjustable footplates + Newton Gravity Ultralight rear wheels)					£437
VP01	☐	Vision Package: Orange accent added to point of contact accessories including: footrest latches, "T" and flip back armrest button caps, 6" Push to lock extension and rubber inserts on Newton composite angle adjustable footplates ¹					£79

1- Please select footplates, footrest, armrest, push to lock wheel lock and the extensions in the appropriate sections; To know the compatible parts, look for this logo *VP01; Recommended with black chair; See picture of the product on page 6.

C2 FRONT Seat-to-floor Height - Please select your choice

		7.62cm (3") caster			10.16cm (4") caster			12.7cm (5") caster			15.24cm (6") caster			17.78cm (7") caster			20.32cm (8") caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
Forks / Position		STD	+ 2.54 cm 1"	+ 5.08cm 2"	STD	+ 2.54 cm 1"	+ 5.08cm 2"	STD	+ 2.54 cm 1"	+ 5.08cm 2"	STD	+ 2.54 cm 1"	+ 5.08cm 2"	STD	+ 2.54 cm 1"	+ 5.08cm 2"	STD	+ 2.54 cm 1"	+ 5.08cm 2"
4" 10.2 cm	P1	33.02cm 13" ☐	35.56cm 14" ☐	38.10cm 15" ☐	34.29cm 13.5" ☐	36.83cm 14.5" ☐	39.37cm 15.5" ☐												
	P2	34.29cm 13.5" ☐	36.83cm 14.5" ☐	39.37cm 15.5" ☐	35.56cm 14" ☐	38.10cm 15" ☐	40.64cm 16" ☐	36.83cm 14.5" ☐	39.37cm 15.5" ☐	41.91cm 16.5" ☐									
5" 12.7 cm	P1																		
	P2				35.56cm 14" ☐	38.10cm 15" ☐	40.64cm 16" ☐	36.83cm 14.5" ☐	39.37cm 15.5" ☐	41.91cm 16.5" ☐									
	P3	35.56cm 14" ☐	38.10cm 15" ☐	40.64cm 16" ☐	36.83cm 14.5" ☐	39.37cm 15.5" ☐	41.91cm 16.5" ☐	38.10cm 15" ☐	40.64cm 16" ☐	43.18cm 17" ☐	39.37cm 15.5" ☐	41.91cm 16.5" ☐	44.45cm 17.5" ☐						
	P4	36.83cm 14.5" ☐	39.37cm 15.5" ☐	41.91cm 16.5" ☐	38.10cm 15" ☐	40.64cm 16" ☐	43.18cm 17" ☐	39.37cm 15.5" ☐	41.91cm 16.5" ☐	44.45cm 17.5" ☐	40.64cm 16" ☐	43.18cm 17" ☐	45.72cm 18" ☐	41.91cm 16.5" ☐	44.45cm 17.5" ☐	46.99cm 18.5" ☐			
7" 17.8 cm	P1										41.91cm 16.5" ☐	44.45cm 17.5" ☐	46.99cm 18.5" ☐	43.18cm 17" ☐	45.72cm 18" ☐	48.26cm 19" ☐			
	P2							41.91cm 16.5" ☐	44.45cm 17.5" ☐	46.99cm 18.5" ☐	43.18cm 17" ☐	45.72cm 18" ☐	48.26cm 19" ☐	44.45cm 17.5" ☐	46.99cm 18.5" ☐	49.53cm 19.5" ☐			
	P3							43.18cm 17" ☐	45.72cm 18" ☐	48.26cm 19" ☐	44.45cm 17.5" ☐	46.99cm 18.5" ☐	49.53cm 19.5" ☐	45.72cm 18" ☐	48.26cm 19" ☐	50.80cm 20" ☐	46.99cm 18.5" ☐	49.53cm 19.5" ☐	52.07cm 20.5" ☐
	P4	41.91cm 16.5" ☐	44.45cm 17.5" ☐	46.99cm 18.5" ☐	43.18cm 17" ☐	45.72cm 18" ☐	48.26cm 19" ☐	44.45cm 17.5" ☐	46.99cm 18.5" ☐	49.53cm 19.5" ☐	45.72cm 18" ☐	48.26cm 19" ☐	50.80cm 20" ☐	46.99cm 18.5" ☐	49.53cm 19.5" ☐	52.07cm 20.5" ☐	50.80cm 20" ☐	53.3cm 21" ☐	
	P5	43.18cm 17" ☐	45.72cm 18" ☐	48.26cm 19" ☐	44.45cm 17.5" ☐	46.99cm 18.5" ☐	49.53cm 19.5" ☐	45.72cm 18" ☐	48.26cm 19" ☐	50.80cm 20" ☐	46.99cm 18.5" ☐	49.53cm 19.5" ☐	52.07cm 20.5" ☐	48.26cm 19" ☐	50.80cm 20" ☐	53.3cm 21" ☐	49.53cm 19.5" ☐	52.07cm 20.5" ☐	54.6cm 21.5" ☐
Froglegs	P1				36.83cm 14.5" ☐	39.37cm 15.5" ☐	41.91cm 16.5" ☐												
	P2	36.83cm 14.5" ☐	39.37cm 15.5" ☐	41.91cm 16.5" ☐	38.10cm 15" ☐	40.64cm 16" ☐	43.18cm 17" ☐	39.37cm 15.5" ☐	41.91cm 16.5" ☐	44.45cm 17.5" ☐									
	P3	38.10cm 15" ☐	40.64cm 16" ☐	43.18cm 17" ☐	39.37cm 15.5" ☐	41.91cm 16.5" ☐	44.45cm 17.5" ☐	40.64cm 16" ☐	43.18cm 17" ☐	45.72cm 18" ☐	41.91cm 16.5" ☐	44.45cm 17.5" ☐	46.99cm 18.5" ☐						

Height available + 0,65cm (1/4") (see stem bolt and fork combination at the end of the order form)

CASTERS				Available			
CA01-03 to CA01-08	☐	Composite wheel w/ PU tire (1")		3", 4", 5", 6", 7", 8"			NC
CA10 to CA11	☐	Pneumatic casters (1 1/4")		6", 8"			£140
CA14 to CA16	☐	Newton UltraCaster Composite (1 1/2")		4", 5", 6"			£81
CA17 to CA19	☐	Newton UltraCaster Aluminum (1 1/2")		4", 5", 6"			£140
CASTER OPTIONS							
CA12	☐	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY) ¹					£250
CA13	☐	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)					£395

1- Not available with WC-19

C2 - REAR Seat-to-floor Height

WHEEL SIZES		Height available		Specify desired height: _____
☐	20" (451) wheel	31.75cm / 12.5"	to 43.81cm / 17.25"	
☐	22" (501) wheel	34.29cm / 13.5"	to 46.35cm / 18.25"	
☐	24" (540) wheel	36.83cm / 14.5"	to 48.89cm / 19.25"	
☐	25" (559) wheel	38.10cm / 15"	to 50.16cm / 19.75"	
☐	26" (590) wheel	39.37cm / 15.5"	to 51.34cm / 20.25"	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1cm" of selected seat to floor height"

Footrest						
SWING-IN/SWING-OUT FOOTREST						
FR03-01	☐	60° Swing-in/swing-out footrest - 34.29 - 43.18cm (13.5" - 17")		Set at*:		NC
FR03-02	☐	60° Swing-in/swing-out footrest - 36.83 - 49.53cm (14.5" - 19.5")		Set at*:		NC
FR03-03	☐	60° Swing-in/swing-out footrest - 39.29 - 52.07cm (15.5" - 20.5")		Set at*:		NC
FR04-01	☐	70° Swing-in/swing-out footrest - 34.29 - 43.18cm (13.5" - 17")		Set at*:		NC
FR04-02	☐	70° Swing-in/swing-out footrest - 36.83 - 49.53cm (14.5" - 19.5")		Set at*:		NC
FR04-03	☐	70° Swing-in/swing-out footrest - 39.37 - 52.07cm (15.5" - 20.5")		Set at*:		NC
FR15-01	☐	80° Swing-in/swing-out footrest - 34.29 - 40.64cm (13.5" - 16") ⁵		Set at*:		NC
FR15-02	☐	80° Swing-in/swing-out footrest - 36.83 - 49.53cm (14.5" - 19.5") ⁵		Set at*:		NC
FR15-03	☐	80° Swing-in/swing-out footrest - 39.37 - 52.07cm (15.5" - 20.5") ⁵		Set at*:		NC
FR06	☐	90° swing-in/swing-out with front mount - Standard (no extension) ¹	19.05 - 27.94cm (7.5" - 11")	Set at*:		NC
FR07	☐	90° swing-in/swing-out with front mount - Short extension 2" ^{1, 2}	24.13 - 33.02cm (9.5" - 13")	Set at*:		NC
FR08	☐	90° swing-in/swing-out with front mount - Medium extension 4" ^{1, 2}	24.13 - 36.83cm (9.5" - 14.5")	Set at*:		NC
FOOTREST OPTIONS ³						
FR09	☐	High mount (footplate attached to the hanger with clamps) ³	12.70 - 34.29cm (5" - 13.5")	Set at*:		£70
FR10	☐	Elevating legrest with calf pad, length adjustable ⁴ Left ☐ Right ☐	35.56 - 55.88 cm (14"-22")	Set at:	INDIVIDUAL PRICE	£206
FR99	☐	Omit Footrest				NC

1- Available only with two piece adjust. footplates; 2- May interfere with front casters; 3- Not available with front mount and the elevating legrest;

4- The elevating legrest replace the footrest side selected; Not available with one piece footplate; 5- Not available with 8" caster

FOOTPLATES						
FP01	☐	Standard - Composite ¹				NC
FP02	☐	Newton adjustable angle - Composite ^{**VP01}				£180
FP03	☐	Newton adjustable angle - Aluminum				£212
FP04	☐	Newton adjustable angle - Carbon fiber with rubber inserts				£365
FP23	☐	Newton adjustable angle extra deep - Aluminum				£212
FP06	☐	One piece adjustable angle flip-up - Aluminum				£170
FP07	☐	One piece adjustable angle flip-up - Carbon fiber with rubber inserts				£380
FOOTPLATE ACCESSORIES						
FP08	☐	Length adjustable heel loop ²				£35
FP09	☐	Calf strap				£33
FP30	☐	Bodypoint Aeromesh Padded calf strap				£95

1- Not available with front mount; 2- not available with one piece footplate

SWING AWAY RESIDUAL LIMB SUPPORT				INDIVIDUAL PRICE £355		PAIR
FR27	☐	Residual limb support - flat 15.24x20.32cm (6"x8")	Left ☐ Right ☐			£710
FR28	☐	Residual limb support - flat 20.32x24.4cm (8"x10")	Left ☐ Right ☐			£710
FR29	☐	Residual limb support - flat 25.4x25.4cm (10"x10")	Left ☐ Right ☐			£710
FR30	☐	Residual limb support - flat 25.4x35.56cm (10"x14")	Left ☐ Right ☐			£710
FR32	☐	Residual limb support - contour 15.24x20.32cm (6"x8")	Left ☐ Right ☐			£710
FR33	☐	Residual limb support - contour 20.32x24.4cm (8"x10")	Left ☐ Right ☐			£710

Backrest						
FIXED ANGLE, HEIGHT ADJUSTABLE BACK CANES WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)						
BC06-912	☐	Straight Carbon fiber back cane	22.86 - 30.48cm (9" - 12")	increment of 1/2"	Set at ²	NC
BC06-1215	☐	Straight Carbon fiber back cane	30.48 - 38.10cm (12" - 15")	increment of 1/2"	Set at ²	NC
BC06-1518	☐	Straight Carbon fiber back cane	38.10 - 45.72cm (15" - 18")	increment of 1/2"	Set at ²	NC
BC06-1821	☐	Straight Carbon fiber back cane	45.72 - 53.34cm (18" - 21")	increment of 1/2"	Set at ²	NC
BC08-1315	☐	8° Bend Carbon fiber back cane	34.29 - 38.1cm (13.5" - 15")	increment of 1/2"	Set at ²	£50
BC08-1518	☐	8° Bend Carbon fiber back cane	38.1 - 45.72cm (15" - 18")	increment of 1/2"	Set at ²	£50
BC08-1821	☐	8° Bend Carbon fiber back cane	45.72 - 53.34cm (18" - 21")	increment of 1/2"	Set at ²	£50
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)						
BC09-1215	☐	Straight Carbon fiber back cane	30.48 - 38.10cm (12" - 15")	increment of 1/2"	Set at ²	£183
BC09-1518	☐	Straight Carbon fiber back cane	38.10 - 45.72cm (15" - 18")	increment of 1/2"	Set at ²	£183
BC09-1821	☐	Straight Carbon fiber back cane	45.72 - 53.34cm (18" - 21")	increment of 1/2"	Set at ²	£183
BC10-1315	☐	8° bend Carbon fiber back cane	34.29 - 38.1cm (13.5" - 15")	increment of 1/2"	Set at ²	£229
BC10-1518	☐	8° bend Carbon fiber back cane	38.1 - 45.72cm (15" - 18")	increment of 1/2"	Set at ²	£229
BC10-1821	☐	8° bend Carbon fiber back cane	45.72 - 53.34cm (18" - 21")	increment of 1/2"	Set at ²	£229
		Select back cane angle:	85° ¹ ☐ 90° ☐ 95° STD ☐ 100° ☐ 105° ☐ 110° ☐			
PUSH HANDLES OPTIONS (ONE option at a time)						
PH01	☐	Push handles				STD
PH05	☐	Fold down push handles ^{3, 5}				£140
PH06-920	☐	clamp on height adjustable push handles 22.86 -50.8cm (9" - 20") ^{3,4, 5}				£260
PH06-923	☐	clamp on height adjustable push handles 22.86 - 58.4cm (9" - 23") ^{3,4, 5}				£260
PH02	☐	Omit push handles				NC
STROLLER HANDLE AND STABILIZER BAR (only available with push handle)						
PH04	☐	Stroller handle (not available with the HD option)				£140
PH03	☐	Newton folding stabilizer bar (not available with stroller handle)				£229

1- Not available with 8 degree bend cane; 2- If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; 3- 120kg (265 lb) weight limit

4- Not compatible with swing away armrests, 8° bend back canes and amputee axle plates; 5- Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

BACK UPHOLSTERY						
BH07	☐	Slip-on Active-Mesh back upholstery				NC
BH08	☐	Adjustable tension Active-Mesh back upholstery				£150
BH99	☐	Omit Back upholstery				STD
SEAT CUSHION						
SC01	☐	5.08 cm (2") Cushion				£150
SC02	☐	7.62 cm (3") Cushion				£166
SC99	☐	Omit seat cushion				STD

Rear wheels						
ONE ARM DRIVE ¹		Only available with Newton One rear wheels, permanent axles, 0° camber				
HR22	☐	One arm drive - aluminum inner handrim ²	Left ☐	Right ☐		£1200
HR60	☐	One arm drive - plastic coated inner handrim ²	Left ☐	Right ☐		£1200
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)				NC
REAR WHEELS					Available	
RW01	☐	Mag ³			20", 22" 24"	£47
RW02	☐	Newton One - spoke wheel ⁴			20", 22" 24"	NC
RW03	☐	Newton Gravity - ultralight wheel			20", 22", 24", 25", 26"	£295
RW04	☐	Spinergy Spox	Blk STD ☐	Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐	22", 24", 25", 26"	£625
RW05	☐	Spinergy LX	Blk STD ☐	Pur ⁵ ☐ Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐ Pink ⁵ ☐ Grn ⁵ ☐ Org ⁵ ☐	22", 24", 25", 26"	£825
RW07	☐	Spinergy CLX ⁶	Blk STD ☐	Pur ⁵ ☐ Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐ Pink ⁵ ☐ Grn ⁵ ☐ Org ⁵ ☐	24", 25"	£1200
TIRES					Available	
AIRLESS						
TI01	☐	Soft Urethane 3.49cm (1 3/8") - Low Tread, black ⁹			20", 22" 24"	STD
TI02	☐	Soft Urethane 3.49cm (1 3/8") - Medium Tread, black ^{7,9}			22", 24", 25", 26"	£85
TI04	☐	Soft Urethane 2.54cm (1") - No tread, black - <i>Shox</i>			22", 24", 25", 26"	£99
TI08	☐	Pneumatic Airless insert 3.49cm (1 3/8") - Medium Tread ⁸			22", 24", 26"	£75
PNEUMATIC						
TI05	☐	Pneumatic 3.49cm (1 3/8") Medium Tread ⁹			20", 22", 24", 25", 26"	NC
TI06	☐	Pneumatic HP 2.54cm (1") - Low Tread, Puncture resistant, black - <i>SpeedLite</i>			20", 22", 24", 25", 26"	£85
TI07	☐	Pneumatic HP 2.54cm (1") - Medium Tread, Puncture resistant, black - <i>TrailBlazer</i>			24", 25", 26"	£99
TI03	☐	Pneumatic 5.08cm (2") - Large Tread, black ¹⁰			22", 24", 25"	£95

Handrims					
HANDRIM			Available		
HR01	☐	Aluminum anodized	20", 22", 24", 25", 26"		STD
HR61	☐	Aluminum black hard anodized ²	20", 22", 24", 25"		£55
HR02	☐	Plastic coated ¹	20", 22", 24", 25", 26"		£85
HR03	☐	High friction coated ^{1,2,3}	22", 24", 25", 26"		£295
HR07-HR18 and HR65-HR71	☐	Newton Air Grip ² <i>See on page 6 for color choice and price</i>	20", 22", 24", 25", 26"		
HR28	☐	Surge ^{1,2,3} <i>(1 1/2" Oval handrim w/ Gripton high friction strip)</i>	22", 24", 25", 26"		£354
HR05	☐	Surge LT ^{1,2,3} <i>(1 1/8" Oval handrim w/ Gripton high friction strip)</i>	22", 24", 25", 26"		£354
HR04	☐	Natural Fit ^{1,2,3}	22", 24", 25", 26"		£279
HR62	☐	Simi H ^{1,2,3} <i>(1 3/8" Oval handrim with high friction strip)</i>	22", 24", 25", 26"		£329
HR63	☐	Nova H ^{1,2,3} <i>(1 1/4" Oval handrim with high friction strip)</i>	22", 24", 25", 26"		£354
HR64	☐	Optimum H ^{1,2,3} <i>(1 1/8" heart-shaped handrim)</i>	22", 24", 25", 26"		£270
HR99	☐	Omit handrim			NC

Handrim option						
ASSEMBLY POSITION						
HR29	☐	Standard Handrim position	Full tab			STD
HR030	☐	Narrow Handrim position no cut	Full tab			NC
HR30	☐	Narrow Handrim position with cut ³	Short tab 2 positions			£22
HR31	☐	Super Narrow Handrim ³	Short tab 2 positions STD ☐ Short tab 1 position ☐			£22
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS						
HR23	☐	Standard grip				£79
HR27	☐	Super grip				£125
HR099	☐	Omit thumb piece				NC

AXLE ¹					
AX01	☐	Quick release axle ¹			STD
AX02	☐	Permanent axle			NC
AX03	☐	Quad release axle ¹			£112

AX04	☐	Amputee adjustable axle plate ²			£175
------	--------------------------	--	--	--	------


POSITIONING BELT					
SB01-S	☐	Velcro adjustable - Small 88.9cm (35") ¹			£40
SB01-L	☐	Velcro adjustable - Long 152.40cm (60")			£40
SB02-S	☐	Auto buckle - Small 88.9cm (35") ¹			£79
SB02-M	☐	Auto buckle - Medium 121.92cm (48")			£79
SB02-L	☐	Auto buckle - Long 152.40cm (60")			£79
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - Medium 142.24cm (56")			£79
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - Long 154.94cm (61")			£79
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - Small 129.54cm (51")			£154
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - Medium 142.24cm (56")			£154
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - Long 154.94cm (61")			£154
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - Small 129.54cm (51")			£165
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - Medium 142.24cm (56")			£165
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - Long 154.94cm (61")			£165
SB03	☐	Hardware attachment for 2 extra points (1 pair)			£15
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 73.66 (18")			£258
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 93.98cm (24")			£258
SB99	☐	Omit positioning belt			STD

1- Only available with push handles - Not available with a stroller handle or a folding stabilizer bar

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 34,3cm (13 1/2") and 50,2cm (19 3/4"); 2- Cannot be combined with

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;

Fork and Stem bolt combinations



Vision Package Accents



Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.5 to 3.75
	16	1.0 to 2.75
	17	1.0 to 3.25
22 in	18+	1.0 to 3.75
	16	1.0 to 2.0
	17	1.0 to 2.75
24 in	18+	1.0 to 3.75

Notes :