



4/A, A. de Curtis St - 80040 Cercola (NA) - Italy
Tel. +39 81 5551946 - Fax +39 81 5552507
www.neatech.it - info@neatech.it

NEATECH™




Dealer _____
User _____
Order number _____ Date _____
Shipping address _____

Prices of each row are unitary and market retail prices without VAT.
If not indicated during the filling of the order form, the product will be configured in the standard version.
Some combinations of options may result incompatible. For any request of customization contact the manufacturer.

TECHNICAL INFORMATION

Maximum user weight: **125 kg**
Wheelchair weight: **chassis 7 kg - wheel 3 kg**
Aluminum chassis
Ø400 mm pneumatic wheels with sealed joints



A - Seat width	465
B - Total width	900
C - Backrest height	710
D - Seat Height	360
E - Seat depth	420
F - Total length	795
G - Total length with footrest	1250

CHASSIS

JOB Wheelchair	☒	JOB
----------------	-------------------------------------	-----

OTHER ACCESSORIES

	Pair of armrests	☐	JBRA
	Safety belt	☐	JCIN

	Head cushion	☐	JCUS
	Carrying bag	☐	JBOR
	Inflatable collar	☐	JCOL
	Body support	☐	JPET
	Sun canopy	☐	JSUN
	Inflatable footrest	☐	JPPG
	Third wheel kit	☐	JOB3