

 **TOPRO**

Premium rollators



EN

made in Norway



We are proud to present our revitalized and modern rollator portfolio

The updated versions the TOPRO Troja 5G, TOPRO Olympos ATR, TOPRO Troja Original and TORPO Troja Walker² which replaces the TOPRO Troja 2G, TOPRO Olympos, TOPRO Troja Classic and TOPRO Troja Walker (have a look at our walker brochure) demonstrate that TOPRO is consolidating its position as market leader in the premium rollator segment.

We are continuing our product updating programme with the new, modern designs which we began with TOPRO Pegasus in the autumn of 2019. We are enthusiastically introducing many updated practical features that our customers have requested.

TOPRO has, for more than 30 years, worked with the development, production, and sale of quality mobility aids. We are proud of the products we make, and we are happy to see the benefits these products provide for our customers.

Through a continuous focus on improvements for the benefit of our customers and the environment, we have also placed significant emphasis on sustainability in our production process. Our factory in Norway uses hydropower and environmentally friendly aluminium to minimize the carbon footprint. Our injection moulding plant uses clean materials that are recyclable to contribute

to the circular economy. We have also implemented a “zero waste” principle to ensure no plastic contaminates the natural environment for the benefit of future generations.

We hope that our products contribute to stable, safe walking and that they motivate everyone to „Keep on moving“.



Asbjørn Eskild
CEO TOPRO Industri AS

Advantages of the TOPRO rollator family

Our mobility and design specialists are constantly working to improve our products with optimized functionality and modern designs for the twenty twenties.

We believe that good design, new colour options, practical high-quality functions, give added value to our products.

Safety and functionality are our top priorities and we want you to feel safe and comfortable using a TOPRO rollator for years to come.



Ergonomic handles

The TOPRO Ergo Grip ensures correct walking posture and comfortable armrests while seated. Together with the soft tyres it provides double shock absorbing effect for stable driving comfort. (Standard on TOPRO Pegasus, Troja 5G, Olympos ATR and Odysseé)

Easy to fold

Simply pull the strap on the seat with only one hand. Practical to transport, stays firmly on its wheels and takes up little space when stowed.

Tilt function and edge guard

Enables the user to easily negotiate thresholds and kerbs, while keeping both hands on the rollator for safety.

Unique memory function

Saves the preferred preset handle height, for easy re-use. (Standard on Troja 5G, Olympos ATR and Troja Original)

Easy wheel change

The interior braking system (TOPRO IBS) in combination with the TOPRO Quick Release system for easy wheel change, allows the use of different tyre treads. (Standard on TOPRO Pegasus, Troja 5G, Olympos ATR and Odysseé)



Scan this QR code and get help by our online rollator guide to find the best TOPRO rollator that meets your needs.

TOPRO PEGASUS

CARBON FOR LIGHTWEIGHT, TOPRO FOR QUALITY



Modern, slender and lightweight carbon rollator with functionality and elegant design. This design combines to guarantee the highest level of quality, safety and stability at a weight of just 6.1 kg*.

Tried-and-tested features have been further developed:

- User friendly method for setting the height of the handles using a **numeric scale**.
- **Partly integrated brake cables** and modern reflective panels on the frame.
- Elegant detachable shopping bag, which, with its retracted position provides **increased forward view to improve user safety**.
- Easy handling and manoeuvrability.
- **Reflectors** on the frame, lock wheels and shopping bag provide **360 degree visibility for safety in the dark**.
- **Shock-absorbing tyres** and **ball bearings** in wheels and forks for a comfortable ride.

Available in colour:



Measurements, colour options and other product details are available on page 19.



* Weight including shopping bag 6.6 kg.

“TOPRO Pegasus represents a new trend in the development of mobility aids where design, functionality and safety are combined to meet users’ expectations.”



made in Norway



Features of the TOPRO Pegasus

Safety and quality

Safety, high durability, a patent-pending brake and various sophisticated functions are notable features of the TOPRO Pegasus. The modern, slender and lightweight design makes it unique.



Numeric scale enables user friendly height adjustment of the handles.

Elegant and practical shopping bag

The elegant rollator design continues with the new shopping bag. Smart and practical, with internal pocket providing a storage option for small items, such as keys and a purse. Carries up to 10 kg.



TROJA 5G

TO MAKE YOU FEEL BETTER



Upgraded level of handling comfort and stability of the ground breaking Troja Classic. The TOPRO Troja 5G offers a modern, slender design and a high level of quality, safety and driving comfort.

Innovative upgrades to provide increased comfort and ease of use:

- User friendly method for setting the height of the handles using a **numeric scale**.
- **Shock-absorbing tyres** provides a comfortable ride and smooth suspension.
- **Bell** to attract attention when necessary. Easy and safe to operate from the handle.
- Easy handling and manoeuvrability due to **ball bearings in forks and wheels**.

- **Detachable shopping bag** with reflective band, zipper and shoulder strap. Closeable and holds up to 10 kg. Two inner pockets for small items.
- Easy to fold, stays firmly on its wheels and takes up little space when stored.

Measurements, colour options and other product details are available on page 20.

Model M & S available in on-trend matt coated colours:



“I have had a TOPRO rollator for a year now and it has transformed my life. I used to have a stick before, but once I got the rollator, I found I could walk so much better.”

| Audrey from the UK



The new TOPRO Troja 5G

The slender design gives a modern and light appearance – matching it perfectly to the revitalized TOPRO rollator family.

We know how important it is to adapt the rollator to individual style and needs. With the new product line, we are proud to offer a wide range of colour options and accessories to help you find your favourite TOPRO rollator and keep on moving.



The detachable shopping bag is redesigned with an on-trend grey fabric print on the lid. Reflector band for visible safety in the dark and practical carrying strap.



Numeric scale for easy height adjustment of the handles, with improved lock clamp.



To ensure that our users are comfortable walking with a rollator, TOPRO's rollator specialists have given it an attractive design. The new frames have a more slender and lighter appearance and the new wheels a slimmer design, while retaining TOPRO's well-known stability and shock-absorbing function.

New on-trend matt colour coating and a wide range of colour* options available:



TROJA ORIGINAL

THE EASY WAY TO GET AROUND



TOPRO's basic rollator with state-of-the-art design and high level of quality and features. The slender and modern design is infused with the functionality needed for a safe walk and an active everyday life.

The functionality is based on TOPRO's high and continuous quality over many years:

- **Lightweight** at only 6.8 kg, and with a **weight bearing capacity** of 165 kg.
- Innovative **sideways folding technology** that is easy to operate with one hand by pulling the strap on the seat.
- **Ball bearings in forks and wheels** to ensure an easy ride and manoeuvrability.
- **Practical net basket** for easy access to items. **Reflective band** for safety and visibility in the dark.
- Available in three different sizes.
- Valued by users for being easy to handle and operate, having great manoeuvrability which makes it a trusty companion in everyday situations.
- Numeric scale for easy height adjustment of the handles, with improved lock clamp.

Model M, S and X available in:



Silver

Measurements, colour options and other product details are available on page 21.



Learn more
TOPRO Troja Original



made in Norway



OLYMPOS ATR

THE ALL-TERRAIN ROLLATOR



With large front wheels the Olympos ATR provides reliable and comfortable support on uneven roads and cobblestones as well as on sand or gravel foundations. The perfect choice for outdoor activities.

Innovative features to provide maximum support for outdoor use:

- **Wider design** to ensure stability in an outdoor environment.
 - **Large and shock-absorbing tyres** to enable a comfortable ride and maintain stable and smooth manoeuvrability on rough terrain.
 - **Detachable spacious water-resistant shopping bag** in on-trend grey fabric print. With reflective band and zipper. Closeable and holds up to 15 kg.
- User friendly method for setting the height of the handles with **numeric scale**.
 - Easy to ride and manoeuvre due to **ball bearings in forks and wheels**.

Measurements, colour options and other product details are available on page 22.

Model M & S available in on-trend matt coated colours:



The ergonomic handle (**TOPRO Ergo Grip**) makes it easier and safer to get up from the sitting position. Together with the soft tyres it provides double shock absorbing effect for stable driving comfort.



made in Norway



ODYSSÉ

THE TRAVEL ROLLATOR



The Odysseé takes up no more space than a sports bag when folded for transport, can easily fit into a car boot or stowed on a coach or train and easily reassembled in seconds!

Innovative features which make the TOPRO Odysseé the ideal travel companion:

- **Compact size** for transport with its patented 2-way folding system that folds the rollator both vertically and horizontally.
- Practical **transport bag** included as standard.
- Can be quickly unfolded when needed, and is **small and easy to handle** for storage.
- Fits in the ordinary luggage compartment on trains so does not block the aisles.
- Due to its compact size some airlines allow it to be stored as hand luggage in the cabin. Remember to ask the airline in advance!
- **Spacious, expandable shopping bag** in water-resistant fabric. Large cover with **reflective band**. Holds up to 10 kg.
- Easy to push and manoeuvre due to **ball bearings in forks and wheels**.

Measurements, colour options and other product details are available on page 23.

Available in colour:

-  Ultramarine
-  Wine red



“The TOPRO Odysseé is a great travel companion.”



made in Norway



Accessories



Crutch holder
can be fitted on either side.

TOPRO Pegasus Art. No. 815236
Troja Original/Troja 5G,
Olympos ATR Art. No. 815358
Olympos Art. No. 814040
Odysseé Art. No. 814457

TOPRO Pegasus, Troja 5G,
Olympos ATR, Troja Original,
Odysseé, Troja Neuro



Tray
Snaps in securely, easy to remove
and clean, carries up to 5 kg

TOPRO Pegasus, Troja Original/5G,
Odysseé Art. No. 814728
Olympos ATR M Art. No. 815212


TOPRO Pegasus, Troja 5G,
Olympos ATR, Troja Original,
Odysseé



Back support with padding
ergonomically height adjustable
according to preferred seating
position. Available in two
lengths.

short 67 cm Art. No. 815270
long 75 cm Art. No. 815272

TOPRO Pegasus, Troja 5G, Olympos
ATR, Odysseé, Troja Neuro



**LED lamp with clip on and USB
plug for safety in the dark.**
Chargeable.

Art. No. 815260

TOPRO Pegasus, Troja 5G,
Olympos ATR, Troja Original,
Odysseé, Troja Neuro



Back support with padding
ergonomically height adjustable
according to preferred seating
position. Available in two
lengths.

short 68 cm Art. No. 815271
long 76 cm Art. No. 815273

Troja Original



Umbrella with attachment
keep dry under the umbrella
in rainy weather or use it as a
parasol.

Art. No. 814068

TOPRO Pegasus, Troja 5G,
Olympos ATR, Troja Original,
Odysseé, Troja Neuro



Additional accessories can be found on our website.

Learn more
Accessories

Technical data for TOPRO Pegasus



Recommended user height	150–200 cm
Max. user weight	130 kg
Max. capacity shopping bag	10 kg / 15 l
Handle height	79–97.6 cm
Width between handles	47 cm
Handles	Ergonomic, several holding positions (TOPRO Ergo Grip)
Maximum length	65 cm
Maximum width	64 cm
Seating width	43.5 cm
Seat height	63 cm
Turning circle	84 cm
Diameter wheels	20 cm
Width wheels	3.6 cm
Grip distance driving brake	7 cm
Height when folded	80.3 cm
Width when folded	27.3 cm
Length when folded	65 cm
Transport dimensions (l × w × h)	65 × 27.3 × 80.3 cm
Product weight (excl. bag)	6.1 kg
Product weight (incl. bag)	6.6 kg
Area of use	Indoors & outdoors
Article No.	815020 + item No.

Available in colour:



White



Steel blue

Technical data for Troja 5G

Size	M (medium)	S (small)
Recommended user height	150–200 cm	135–170 cm
Max. user weight	150 kg	125 kg
Max. capacity shopping bag	10 kg / 15 l	10 kg / 15 l
Handle height	78–100 cm	71–86 cm
Width between center points of handles	47 cm	47 cm
Handles	Ergonomic, several holding positions (TOPRO Ergo Grip)	
Maximum length	65 cm	65 cm
Maximum width	60 cm	60 cm
Seating width	45.5 cm	45.5 cm
Seat height	63,5 cm	54 cm
Turning circle	84 cm	84 cm
Diameter wheels	20 cm	20 cm
Width wheels	3.6 cm	3.6 cm
Grip distance driving brake	7 cm	7 cm
Height when folded	81.5 cm	74 cm
Width when folded	24 cm	24 cm
Length when folded	65 cm	65 cm
Transport dimensions (l × w × h) (excl. wheels and shopping bag)	60 × 24 × 76 cm	60 × 24 × 72 cm
Transport weight (excl. wheels and shopping bag)	5.6 kg	5.4 kg
Product weight (incl. wheels and shopping bag)	7.8 kg	7.6 kg
Area of use	Indoors & outdoors	
Article No.	814600 + item No.	814601 + item No.



Model M & S available in on-trend matt coated colours:



Technical data for Troja Original

Size	M (medium)	S (small)	X (extra small)
Recommended user height	150–200 cm	135–170 cm	125–160 cm
Max. user weight	165 kg	130 kg	80 kg
Max. capacity basket	10 kg / 20 l	10 kg / 20 l	5 kg / 11 l
Handle height	78–100 cm	71–86 cm	67–79 cm
Width between center points of handles	47 cm	47 cm	37 cm
Handles	Anatomical		
Maximum length	65 cm	65 cm	65 cm
Maximum width	60 cm	60 cm	52 cm
Seating width	45.5 cm	45.5 cm	36 cm
Seat height	63 cm	55 cm	52 cm
Turning circle	84 cm	84 cm	79 cm
Diameter wheels	20 cm	20 cm	20 cm
Width wheels	3.6 cm	3.6 cm	3.6 cm
Grip distance driving brake	6.5 cm	6.5 cm	6.5 cm
Height when folded	80 cm	74 cm	65 cm
Width when folded	24 cm	24 cm	24 cm
Length when folded	65 cm	65 cm	65 cm
Transport dimensions (l × w × h)	65 × 24 × 80 cm	65 × 24 × 74 cm	65 × 24 × 65 cm
Product weight (incl. wheels and basket)	6.8 kg	6.6 kg	6.1 kg
Area of use	Indoors & outdoors		
Article No.	814750 -110000	814770-111000	814740-112000



Model M, S, X available in:



Technical data for Olympos ATR

Size	M (medium)	S (small)
Recommended user height	150–195 cm	140–170 cm
Max. user weight	150 kg	150 kg
Max. capacity shopping bag	15 kg / 12 l	15 kg / 10 l
Handle height	78.5–94.5 cm	70.5–86.5 cm
Width between center points of handles	49 cm	50 cm
Handles	Ergonomic, several holding positions (TOPRO Ergo Grip)	
Maximum length	77 cm	77 cm
Maximum width	68 cm	68 cm
Seating width	47 cm	47 cm
Seat height	61 cm	54.5 cm
Turning circle	97 cm	97 cm
Diameter wheels	28 cm / 23 cm	28 cm / 23 cm
Width wheels	4.2 cm	4.2 cm
Grip distance driving brake	7 cm	7 cm
Height when folded	87.5 cm	85 cm
Width when folded	27 cm	27 cm
Length when folded	77 cm	77 cm
Transport dimensions (l × w × h)	77 × 27 × 87.5 cm	77 × 27 × 85 cm
Product weight (incl. wheels and shopping bag)	8.9 kg	8.7 kg
Area of use	Outdoors	
Article No.	814300 + item No.	814305 + item No.



Model M & S available in on-trend matt coated colours:



Technical data for Odyssé

Recommended user height	150–200 cm
Max. user weight	150 kg
Volume capacity of shopping bag	13–22.5 l
Max. load in shopping bag	10 kg
Handle height	81.5–100 cm
Width between handles	47 cm
Handles	Ergonomic, several holding positions (TOPRO Ergo Grip)
Maximum length	65 cm
Maximum width	60 cm
Seating width	45.5 cm
Seat height	62 cm
Turning circle	84 cm
Diameter wheels	20 cm
Width wheels	3.6 cm
Grip distance driving brake	7 cm
Height when laterally folded	81.5 cm
Transport dimensions completely folded (l × w × h)	75 × 24.5 × 44.5 cm
Transport weight (excl. wheels and shopping bag)	6.0 kg
Product weight (incl. wheels and shopping bag)	8.2 kg
Area of use	Indoors & outdoors
Article No. Ultramarine / Wine red	815000 / 815002



Available in colour:



Premium rollators from Norway

Manufacturer:

TOPRO

TOPRO Industri AS

Rambekkevegen 5

2816 Gjøvik

Norway

Tel. (+47) 61 13 46 00

info@topromobility.com

www.topromobility.com



Learn more



topromobility.com

Norwegian design and made in Norway.

TOPRO vouches for top quality and offers a 7-year warranty (excluding assembly).

Setup and variations in fixed accessories for the rollators may differ between countries.

Quality Assurance Standards

EN ISO 9001 – quality

EN ISO 14001 – environment

EN ISO 11199-2 – rollator

EN ISO 10993-5 – biocompatibility

made in Norway



Dealer