

2026

Eco-Buggy.

ottobock.

☐ Quotation

☐ Order

Billing Address

PO no.

Customer Name

Street

County PC

Shipping address (if different from billing address)

Customer Name

Street

County PC

Name Date Signature



Base Price - £738.00

Please email all orders to
Bockuk@ottobock.com

Valid from 01.03.2026 to 31.12.26
Please ensure you are using the latest order form
v2

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals. Wheelchairs with CLT package cannot be returned or exchanged

Standard equipment Selection (unless alternative selected)

- Lightweight folding buggy
- Seat and back support upholstery
- Leg support, height-adjustable; divisible, swing-away foot plate
- Twin-roller casters with swivel lock (pin lock)
- Wheel locks at the rear wheels
- Lap belt, standard
- Folding storage bag
- Frame colour anthracite metallic
- Upholstery blue

max. load capacity 50 kg

The number codes contained in this order form are not intended to be used as article numbers for spare parts orders.
Please scan the QR code to access the spare parts catalogue for any spare parts orders.



Eco-Buggy.

Eco-Buggy

☐ HR32100001-011	Eco-Buggy, rehab folding buggy	£738.00
---	--------------------------------	---------

Accessories

☐ HR32070020-343	Padding with lateral guides, coro blue	£396.00
☐ HR32135300-343	Simple padding, coro blue	£150.00
☐ HR32102300	Five-point belt (not in combination with padding with lateral guides, coro blue)	£73.00
☐ HR32102200	Sun canopy	£114.00
☐ HR16001102	Buggy board „Techno“	£254.00
☐ HR32070060-031	Head support size 1, padded (not available in combination with padding)	£243.00


Technical data

Age	3 – 8 years	Back support height	550 mm	Lower leg length/space between seat edge and foot plate	200 – 340 mm*
Seat width	350 mm	Back support angle	90° - 110°	Max. load capacity	50 kg
Seat depth	300 mm	Wheel ø front/rear	170 / 210 mm	Weight (standard equipment)	9 kg
Folding size (L x B x H)	1150 x 420 x 290 mm	* max. length requires caster wheel adjustment			

Information



Folding storage bag



Padding with lateral guides, coro blue



Simple padding, coro blue



Five-point belt



Sun canopy



Head support size 1, padded