

2026

ottobock.

Power wheelchair *Juvo B4.*

☐ Quotation

☐ Order

Billing Address

PO no.

Customer Name

Street

County PC

Shipping address (if different from billing address)

Customer Name

Street

County PC

Name Date Signature



Rear-wheel drive

figure shows add charge options



Rear-wheel drive

figure shows add charge options



Rear-wheel drive

figure shows add charge options

Article number: 490E75=1_C

max. load capacity: xxx kg (standard equipment)

Base Price - £4305.00

Please email all orders to
Bockuk@ottobock.com

Valid from 01.03.2026 to 31.12.26

Please ensure you are using the lastet order form
v2

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals. Wheelchairs with CLT package cannot be returned or exchanged

The number codes contained in this order form are not intended to be used as article numbers for spare parts orders.
Please scan the QR code to access the spare parts catalogue for any spare parts orders.




Power wheelchair Juvo B4.

Basic data (voluntary information)

Gender ☐ male ☐ female ☐ divers

Indication



Mandatory information*

User height [cm]	
User weight [kg]	
1 Seat width [cm]	
2 Seat depth [cm]	
3 Trunk depth [cm]	
4 Arm support height [cm]	
5 Armpit height [cm]	
6 Shoulder height [cm]	
7 Head height [cm]	
8 Forearm length to middle of palm	
9 Trunk width [cm]	
10 Lower leg length [cm]	

* This information will help to support user-oriented settings and product development. Ottobock guarantees the confidential and purely internal use of these data.

Product classification

☐ EA001	Category A (Indoor, wheels: 12"/8")	£4305.00
☐ EA002	Category B (In-/Outdoor) (electric lighting mandatory in Europe)	£4305.00

Max. user weight (mandatory information; user + personal luggage, information required for chassis suspension)

☐	75 kg	Inclusive
☐	100 kg (for category A up to 12.5"/8")	Inclusive
☐	113 kg (not possible in combination with category A (indoor) see exclusions 2.1)	Inclusive
☐	125 kg (not possible in combination with category A (indoor) see exclusions 2.1)	Inclusive
☐	140 kg (not possible in combination with category A (indoor) see exclusions 2.1)	Inclusive
☐	150 kg (not possible in combination with category A (indoor) see exclusions 2.1)	Inclusive
☐	160 kg (not possible in combination with category A (indoor) see exclusions 2.1)	Inclusive

Seating systems (mandatory information)

☐ EC001	Standard seat Junior (max. user weight 75 kg)	£157.00
☐ EC002	Standard seat Small	Inclusive
☐ EC003	Standard seat Large	Inclusive
☐ EC004	Standard seat XL	£202.00
☐ EC005	VAS (Vario Adjust Seat) (possible only in combination with max. load capacity <= 125 kg) see exclusions 3.3)	£463.00

Power wheelchair Juvo B4.

Settings for Standard seat Junior see exclusions 3.1) 3.2) 4.0) 4.1) 4.4) 5.0)

Seat width SW in cm ☐ 34 ☐ 36 ☐ 38 ☐ 40

Seat depth SD in cm ☐ 34 ☐ 36 ☐ 38 ☐ 40

Back support height BH 47 cm (adjustable to 42 cm or 37 cm)

Settings for Standard seat Small see exclusions 3.1) 4.0) 4.1) 4.4) 5.0)

Seat width SW in cm ☐ 38 ☐ 40 ☐ 42

Seat depth SD in cm ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46

Back support height BH 55 cm (adjustable to 50 cm or 45 cm)

Settings for Standard seat Large see exclusions 4.0) 5.0)

Seat width SW in cm ☐ 44 ☐ 46 ☐ 48

Seat depth SD in cm ☐ 42 ☐ 44 ☐ 46 ☐ 48 ☐ 50

Back support height BH 55 cm (adjustable to 50 cm or 45 cm)

Settings for Standard seat XL see exclusions 4.0) 5.0)

Seat width SW in cm ☐ 50 ☐ 52 ☐ 54 ☐ 56

Seat depth SD in cm ☐ 42 ☐ 44 ☐ 46 ☐ 48 ☐ 50

Back support height BH 55 cm (adjustable to 50 cm or 45 cm)

Settings for VAS see exclusions 4.0) 5.0)

Seat width SW in cm ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48
☐ 50 ☐ 52 ☐ 54

Seat depth SD in cm ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48
☐ 50 ☐ 52 ☐ 54 ☐ 56 ☐ 58

Back support height BH 55 cm (adjustable to 53 cm, 51 cm or 45 cm)

☐ Lower leg length LLL cm (15 ... 47 cm) (without Seat cushion)

Overview about seating systems, seat width and seat depth

Seat width

Seat	Seat width SW [cm]											
	34	36	38	40	42	44	46	48	50	52	54	56
Standard seat Junior	☐	☐	☐	☐								
Standard seat Small			☐	☐	☐							
Standard seat Large						☐	☐	☐				
Standard seat XL									☐	☐	☐	☐
VAS			☐	☐	☐	☐	☐	☐	☐	☐	☐	

Power wheelchair *Juvo B4*.

Seat depth

Seat	Seat depth SD [cm]												
	34	36	38	40	42	44	46	48	50	52	54	56	58
Standard seat Junior	○	○	○	○									
Standard seat Small			○	○	○	○	○						
Standard seat Large					○	○	○	○	○				
Standard seat XL					○	○	○	○	○				
VAS			○	○	○	○	○	○	○	○	○	○	○

Seat cushion (only possible in combination with standard seat or VAS, mandatory information)

☐ EC011	Without seat cushion	Inclusive
☐ EC012	Seat cushion, sblack cover (functional height 5 cm)	Inclusive
☐ EC034	Livo (max. user weight 125 kg; functional height 5 cm) see exclusions 3.3) 7.6)	£88.00
☐ EC035	Evoflair (max. user weight 200 kg; functional height 8 cm) see exclusions 7.7)	£313.00
☐ EC014	Terra (up to 150 kg; functional height 6 cm) see exclusions 7.1) 7.2)	£214.00
☐ EC015	Terra Aquos (up to 150 kg; functional height 6 cm) see exclusions 7.1) 7.2)	£241.00
☐ EC016	Terra Flair (up to 150 kg; functional height 6 cm) see exclusions 7.1) 7.3)	£481.00
☐ EC018	Cloud (functional height 11 cm) see exclusions 7.5)	£589.00
☐ EC020	Contour, low contoured, Dartex® fabric cover (functional height 7 cm) see exclusions 7.8) bis 7.11)	£260.00
☐ EC021	Contour, deep contoured, Dartex® fabric cover (functional height 9 cm) see exclusions 7.8) bis 7.11)	£260.00

Please refer to our seat cushion brochure for information on additional seat cushions.

Accessories for Contour seat cushion with Dartex® fabric cover 72.0)

☐ EC039	Incontinence barrier (fitted by specialist dealer)	£27.00
--------------------------------	--	--------

Back support version (mandatory information)

☐ ED032	Without back support upholstery	Inclusive
☐ ED033	Back support upholstery, adaptable	Inclusive
☐ ED034	Back support cushion low contoured, Dartex® fabric cover (SD reduction by approx. 2 cm) exclusions see 7.8) 7.9)	£371.00
☐ ED035	Back support cushion deep contoured, Dartex® fabric cover (SD reduction by approx. 2 cm) exclusions see 7.8) 7.9)	£371.00
☐ ED050	Back support cushion low contoured AD, Dartex® fabric cover (adaptable, SD reduction by approx. 2 cm) exclusions see 9.1)	£433.00
☐ ED051	Back support cushion deep contoured AD, Dartex® fabric cover (adaptable, SD reduction by approx. 2 cm) exclusions see 9.1)	£433.00
☐ ED036	Baxx aluminium flat top back support (max. user weight 113 kg, installation accessories Quick Release is released for transportation in a vehicle, SD reduction by approx. 2 cm, back support width = chosen seat width) exclusions see 3.3) 9.1) 9.3)	£716.00

Power wheelchair *Juvo B4.*

Accessories for back version Contour 11.0)

☐	ED039	Side panels (only for back support pad deep contoured) exclusions see 11.1)	£98.00
☐	ED040	Lumbar support exclusions see 11.1)	£75.00

Accessories for arm support exclusions see 12)

☐	EL021	Padded covers for arm supports exclusions see 12.1)	£97.00
☐	EL022	Padded covers for clothing guard	£75.00

Back support height for Baxx exclusions see 14.0) 14.1)

☐	ED041	High (18"/46 cm)	Inclusive
☐	ED042	Extra high (20"/51 cm)	Inclusive

Back support contour for Baxx exclusions see 14.0) 14.1)

☐	ED043	Active (2,5"/6 cm)	Inclusive
☐	ED044	Deep (4"/10 cm)	£30.00
☐	ED045	Extra deep (6"/15 cm)	£58.00

Colour for Baxx (possible only in combination with Baxx Aluminium Flat Top back support; mandatory information) see exclusions 14.0) 14.1)

☐	ED014	Candy red	Inclusive	☐	EDM51	Sunny orange (texture)	Inclusive
☐	ED023	White RAL 9016	Inclusive	☐	EDM52	Rich lime (texture)	Inclusive
☐	EDM48	Obsidian black (texture)	Inclusive	☐	EDM55	Blueberry (texture)	Inclusive
☐	EDM50	Aubergine (texture)	Inclusive				

Seat inclination (pre-tilt)

☐	-3°	(not possible in combination with EC005) see exclusions 17.1)	Inclusive
☐	0°		Inclusive
☐	3°	(not possible in combination with EC005 möglich) see exclusions 17.1)	Inclusive
☐	6°		Inclusive
☐	9°		Inclusive

Power seat adjustment (only possible in combination with seat inclination 0°)

☐	EC025	Power seat tilt 45° with centre of gravity shifting (max. load capacity: VAS 100 kg, standard-/contour seat 140 kg) see exclusions 3.1) 18.1) 24.9)	£887.00
☐	EC026	Power seat height adjustment 35 cm (max. load capacity: VAS 100 kg, standard-/contour seat 140 kg; not possible in combination with EU060 = VR2 (90A) or EU051 = VR2 (70 A) see exclusions 3.1) 18.1) 24.9)	£2448.00
☐	EC027	Power seat height adjustment 35 cm with 45° seat tilt, with centre of gravity shifting (max. load capacity: VAS 100 kg, standard-/contour seat 140 kg; not possible in combination with EU060 or EU051) see exclusions 3.1) 18.1) 24.9)	£2860.00
☐	EC028	Power seat tilt 20° (max. load capacity: standard-/contour seat 140 kg, not possible in combination with EC005) see exclusions 3.1) 17.1) 18.1)	£619.00

Power wheelchair *Juvo B4*.

Seat adjustment accessories

☐	ES307	Stabilisation for seat height adjustment see exclusions 68.1)	£81.00
--------------------------	-------	--	--------

Caster fork (mandatory information)

☐	EF001	Suspension (not possible in combination with EF007 and EF008) exclusions see 2.1)	£265.00
☐	EF002	No suspension	Inclusive

Caster wheels (mandatory information)

☐	EF005	9" pneumatic, puncture-resistant see exclusions 2.1)	Inclusive
☐	EF006	9" puncture-proof see exclusions 2.1)	Inclusive
☐	EF007	10" pneumatic, puncture-proof (not possible in combination with EF001) see exclusions 2.1) 20.1)	Inclusive
☐	EF008	10" puncture-proof (not possible in combination with EF001) see exclusions 2.1) 20.1)	Inclusive
☐	EF004	8" pneumatic, puncture-proof (not possible in combination with EG003) see exclusions 20.2)	£77.00

Drive wheels (mandatory information)

☐	EG001	14" pneumatic, puncture-resistant see exclusions 2.1)	Inclusive
☐	EG002	14" puncture-proof see exclusions 2.1)	Inclusive
☐	EG003	12" puncture-proof (not possible in combination with EF004) see exclusions 20.2)	£171.00

Tyre colour (mandatory information)

☐	EG011	Grey	Inclusive
☐	EG012	Black (not possible in combination with EG003 and EF004) see exclusions 22.1)	Inclusive

Tread (mandatory information)

☐	EG013	Ribs (only with tyre colour grey) see exclusions 23.1)	Inclusive
☐	EG014	Lugs (not possible in combination with wheels 12" and 8") see exclusions 22.1)	Inclusive
☐	EG015	Caster wheels with ribs, drive wheels with lugs (not possible in combination with EG003 and EF004) see exclusions 22.1) 23.1)	Inclusive

Power wheelchair Juvo B4.

Seat height (at 0° seat inclination, without seat cushion, suspension fully extended/without user, mandatory information)

Seat height (at 0° seat inclination, without seat cushion, suspension fully extended, without seat, mandatory information)																								
Wheels		1.1)	1.1)	1.2)	1.3)	1.1)	1.2)	1.3)	1.1)	1.1)	1.2)	1.3)	1.2)	1.1)	1.2)	1.3)	1.3)	1.1)	1.2)	1.3)				
Seat height	cm	38	40			42,5			44	45			46	47,5			48,5	50						
	inch	14,6	15,6			16,6			17,2	17,6			18	18,6			19	19,5						
Notices and ex- clusions	24.37)	24.37)	24.33)	24.16)	24.0)	24.37)	24.31)	24.29)	24.25)	24.37)	24.37)	24.31)	24.27)	24.26)	24.37)	24.37)	24.36)	24.24)	24.21)	24.37)	24.37)	24.35)	24.23)	24.22)
	24.32)	24.17)	24.13)	24.28)	24.2)	24.17)	24.27)	24.15)	24.0)	24.34)	24.0)	24.2)	24.0)	24.0)	24.20)	24.19)	24.18)	24.0)	24.0)	24.12)	24.0)	2.1)	24.3)	24.2)
	24.4)	24.14)			24.3)	24.11)			24.9)	24.2)	24.3)	24.4)	24.7)			24.2)			24.8)					
	24.5)									24.0)				24.6)				24.6)						
	24.14)									24.10)														
	24.0)									24.1)														
	24.1)									24.6)														
	24.2)																							
	24.21)																							
Standard seat	 2.2)																							
VAS		 2.2)			 2.2)		 2.2)																	
Standard seat, tilt 20°																								
Standard seat, tilt 45°								 2.2)																
VAS, tilt 45°								 2.1)																
Standard seat, seat lift 35 cm + tilt 45°								 2.1)																
VAS, seat lift 35 cm + tilt 45°																								

Wheels

- 1.1) 8"/12,5"
 1.2) 9" caster wheels
 1.3) 9" caster wheels with suspension or 10" wheels

Weight

- 2.1) max. 75 kg
 2.2) max. 100 kg

- combination possible
 □ combination not possible

Back support angle adjustment (mandatory information)

○ ED046	Mechanical in increments (standard seat: 0/10/20/30° or -9/1/11/21°; VAS: -5/0/10/20/30°)	Inclusive
○ ED047	Power (30°)	£756.00

Mounting kit for head/neck support

○ EL001	Gooseneck	£141.00
○ EL002	Multi-axis, offset	£130.00
○ EL003	Multi-axis, straight	£125.00
○ EL031	UniLink holder for head support, multi-axis (suitable only for Advanced line head supports; not possible in combination with EL006, EL007, EL013 to EL016)	£259.00

Power wheelchair *Juvo B4*.

Accessories UniLink see exclusions 69.0)

☐	EL039	UniLink extension element, 38 mm (may be used to compensate the distance at kyphotic head posture)	£90.00
☐	EL038	UniLink extension element, 50 mm (may be used to compensate the distance at kyphotic head posture)	£66.00
☐	EL040	UniLink extension element, 76 mm (may be used to compensate the distance at kyphotic head posture)	£90.00
☐	EL031	UniLink holder for head support, multi-axis (suitable only for Advanced line head supports; not possible in combination with EL006, EL007, EL013 to EL016)	£259.00

Head/neck supports see exclusions 27.0) 27.1)

Care line

☐	EL006	Head support, small see exclusions 27.2)	£98.00
☐	EL007	Head support, large see exclusions 27.2)	£98.00
☐	EL013	Head/neck support, small see exclusions 27.2)	£93.00
☐	EL014	Head/neck support, large see exclusions 27.2)	£93.00
☐	EL015	Head/neck support X-shape, small see exclusions 27.2)	£113.00
☐	EL016	Head/neck support X-shape, large see exclusions 27.2)	£113.00

Advanced line

☐	EL033	Comfort Plus, size S	£174.00
☐	EL034	Comfort Plus, size M	£200.00
☐	EL035	Comfort Plus, size L	£200.00
☐	EL036	Combo (3-pad), size S	£344.00
☐	EL037	Combo (3-pad), size L	£358.00

Belts see exclusions 28.0)

☐	EP001	Lap belt standard	Inclusive
☐	EP002	2-point lap belt, padded, metal buckle	£164.00
☐	EP003	2-point lap belt, padded, plastic buckle	£149.00
☐	EP020	2-point lap belt with metal buckle and rewriter (buckle on side of joystick) see exclusions 7.8) 28.1)	£252.00

Chest/shoulder harnesses

☐	EP005	Dynamic harness, size S	£214.00
☐	EP006	Dynamic harness, size M	£214.00
☐	EP007	Dynamic harness, size L	£214.00
☐	EP008	Sternum harness, size S	£205.00
☐	EP009	Sternum harness, size M	£205.00
☐	EP010	Sternum harness, size L	£205.00
☐	EP011	Slimline sternum harness, size S	£198.00
☐	EP012	Slimline sternum harness, size M	£198.00
☐	EP013	Slimline sternum harness, size L	£198.00
☐	EP014	Shoulder harness, size S	£220.00
☐	EP015	Shoulder harness, size M	£220.00
☐	EP016	Shoulder harness, size L	£220.00

Power wheelchair *Juvo B4*.

VAS accessories see exclusions 9.1)

☐	ED048	Push handles	£75.00
☐	ES302	Thigh lateral supports, size M (1 pair of outer thigh lateral supports with UniLink mounting materials multi-axis adjustment) see exclusions 3.1) 30.1)	£779.00
☐	ES303	Thigh lateral supports, size S (1 pair of outer thigh lateral supports with UniLink mounting materials multi-axis adjustment) see exclusions 3.1) 30.2)	£769.00

Back extension for VAS 9.1)

☐	15 cm 71.1) 71.2) 71.3)	Inclusive
☐	15 cm with traverse 71.1) 71.2) 71.3)	Inclusive

Arm supports

☐	EE023	Standard	Inclusive
☐	EE024	Flip-up (effective seat width +/- 1 cm per side; on Standard seat Small, SW 42 cm; Standard seat Large, SW 48 cm; Standard seat XL, SW 54 cm; VAS plus 2 cm; max. user weight: 125 kg; not in combination with EC001) see exclusions 3.3) 7.8) 24.5) 67.1) bis 67.4)	£285.00
☐	EE026	Reinforced see exclusions 9.1) 30.1) 30.2) 67.5) 67.6)	£47.00

Arm pad

☐	EE060	Pads short (27 cm)	Inclusive
☐	EE061	Pads large (37 cm)	£49.00

Leg support version

☐	EB011	Standard (basic design) (information LLL required) see exclusions 7.8) 31.1) 31.6) 31.7) 31.8)	Inclusive
☐	EB001	Standard (new design) see exclusions 31.1) 31.7) 31.8)	£237.00
☐	EB002	Mechanically elevating with length adjustment (max. user weight 140 kg; not available for LLL < 21 cm; not in combination with EB012, EB009) see exclusions 3.1) 31.1) bis 31.4) 31.7) 31.8)	£658.00
☐	EB003	Power elevating with length adjustment and separate control (max. load 140 kg; not available for LLL < 21 cm; not in combination with EB012, EB009) see exclusions 3.1) 31.1) bis 31.4) 31.7) 31.8)	£1340.00
☐	EB004	Power elevating with length adjustment and simultaneous control (not available for LLL < 21 cm; not in combination with EB012) see exclusions 31.1) bis 31.4) 31.7) 31.8)	£1340.00
☐	EB005	Without leg supports see exclusions 31.5)	Inclusive

Foot plate version see exclusions 32.0) 32.1)

☐	EB006	Individual, plastic foot plates (max. load 140 kg; not available on elevating leg supports; not in combination with EB002, EB003, EB004) see exclusions 3.1) 31.6) 32.2) bis 32.4)	Inclusive
☐	EB007	Individual aluminium foot plates, new design (max. user weight 160 kg; not in combination with EB011) see exclusions 32.6) bis 32.8) 32.10)	£285.00
☐	EB009	Single-panel, aluminium foot plate (max. user weight 160 kg; not in combination with EB002, EB003) see exclusions 32.6) bis 32.9)	£166.00

Power wheelchair *Juvo B4.*

Leg support accessories

☐	EB010	Lateral leg pad (for mechanical and power elevating leg support adjustment, pair) see exclusions 33.1)	£243.00
☐	EB017	Lock for single panel aluminium foot plate (only possible in combination with EB009) see exclusions 33.2) 33.3) 33.4)	£88.00
☐	EB022	Amputation leg support, right see exclusions 3.1) 33.5)	£661.00
☐	EB023	Amputation leg support, left see exclusions 3.1) 33.5)	£661.00

Batteries

☐	EU036	AGM 56 Ah (C5) / 60,4 Ah (C20)	Inclusive
☐	EU032	Gel 63 Ah (C5) / 74 Ah (C20)	£672.00
☐	EU500	Gel 62 Ah (C5) / 79,6 Ah (C20) Basic	£312.00

Battery charger

☐	EU040	Without battery charger	Inclusive
☐	EU041	Battery charger 8 A, fanless (suitable only for batteries with 56 Ah (C5)) see exclusions 35.3) 35.4)	Inclusive
☐	EU044	Battery charger 10 A, fanless (IP54) (suitable for batteries 56 Ah (C5) and 63 Ah (C5)) see exclusions 35.2)	£295.00
☐	EU046	Battery charger 10 A, with fan (suitable for batteries 56 Ah (C5) and 63 Ah (C5))	£308.00
☐	EU042	Battery charger 12 A, incl. external charging receptacle (magnetic, installation side right) see exclusions 35.4)	£522.00

Additional charging receptacle

☐	EU043	External magnetic charging receptacle (installation side right; possible only with EU042 see exclusions 36.1)	£212.00
--------------------------	-------	--	---------

Wheelchair colour (colour for caster fork and metal strip; frame colour: KTL black; mandatory information)

☐	000	Without coloured caster wheel swing arm cover and without coloured metal strut	Inclusive	☐	050	Aubergine (texture)	£99.00
☐	014	Candy red	Inclusive	☐	051	Sunny orange (texture)	£99.00
☐	023	White RAL 9016	Inclusive	☐	052	Rich lime (texture)	£99.00
☐	048	Obsidian black (texture)	Inclusive	☐	055	Blueberry (texture)	£99.00

Chassis accessories

☐	EA010	Splash guard for drive wheels see exclusions 2.1) 39.5) 39.6)	£123.00
☐	EA011	Splash guard for caster wheels see exclusions 2.1)	£73.00

Caster wheel swivel lock

☐	EF020	Mechanical	£487.00
--------------------------	-------	-------------------	---------

Chassis safety accessories

☐	EF021	Curb climbing assist see exclusions 2.1)	£393.00
☐	ES002	Kit for vehicle transport (not possible in combination with ED036; max. load capacity 140 kg) see exclusions 3.1) 41.3)	Inclusive
☐	ES013	Ottobock installation kit for Dahl docking system (not possible in combination with EA001 and EF021; max. load capacity 140 kg) see exclusions 2.1) 3.1) 24.5) 41.2)	£123.00

Power wheelchair Juvo B4.

General accessories

☐	ES003	Pocket for mobile phone	£61.00
☐	ES004	Luggage carrier see exclusions 2.1) 3.1)	£273.00
☐	ES005	Crutch holder	£141.00
☐	ES006	Rear-view mirror (installation side across from control panel side, not possible in combination with ES012 möglich) see exclusions 42.1)	£209.00
☐	ES007	Tool kit	£73.00
☐	ES010	Beverage holder	£31.00
☐	ES011	Airman pump see exclusions 42.4)	£294.00
☐	ET001	External horn	£571.00
☐	ES012	Therapy tray (swivel mechanism across from control panel side; not possible in combination with ES006 möglich) see exclusions 42.2)	£480.00
☐	EU095	Hand warmer (single hood for handmodule) see exclusions 42.3)	£1478.00
☐	ES023	Backpack hooks (at rear)	£90.00
☐	ES305	Backpack with Ottobock logo	£65.00
☐	ES306	Catheter bag, large	£48.00
☐	ES024	Protective frame for hand module	£65.00
☐	ES300	Transfer handle, left see exclusions 30.1) 30.2) 42.5) 42.6) 42.7) 42.8)	£101.00
☐	ES301	Transfer handle, right see exclusions 30.1) 30.2) 42.5) 42.6) 42.7) 42.8)	£101.00

Receptacle for external power supply, right

☐	ET002	12 V see exclusions 43.1) 43.4)	£366.00
☐	ET003	24 V see exclusions 43.2)	£148.00
☐	ET004	USB see exclusions 43.3)	£165.00

Lighting see exclusions 44.0)

☐	ES021	Electric lighting according to German Motor Vehicle Safety Standards (StVZO)	£560.00
-----------------------	-------	--	---------

Control device

☐	EU059	VR2 (70 A), without power functions and lighting (recommended for load capacity up to 140 kg) see exclusions 46.1) 46.2) 46.4)	Inclusive
☐	EU060	VR2 (70 A), max. 2 power functions and lighting possible (recommended for load capacity up to 140 kg, not possible in combination with EC026 or EC027) see exclusions 46.2) 46.3) 46.5) 46.8)	Inclusive
☐	EU050	VR2 (90 A), without power functions and lighting (recommended for load capacity up to 160 kg) see exclusions 2.1) 46.1) 46.4)	Inclusive
☐	EU051	VR2 (90 A), max. 2 power functions and lighting possible (recommended for load capacity up to 160 kg, not possible in combination with EC026 or EC027) see exclusions 2.1) 46.3) 46.5) 46.8)	£209.00
☐	EU052	TEN° (90 A), 1 power function, choose option EU056 for more power functions, lighting possible, Colour: Shimano matt (recommended for load capacity up to 160 kg)	£966.00
☐	EU058	R-net LED-L (90 A), 1 power function, choose option EU056 for more power functions, lighting possible (recommended for load capacity up to 160 kg)	£670.00

Power wheelchair *Juvo B4.*

Control panel button (only possible in combination with EU052) see exclusions 47.0)

☐	EU067	Piko button, diameter 50 mm (red/black)	£435.00
☐	EU068	Piko button, diameter 30 mm (red/black)	£435.00

Control panel holder (mandatory information)

☐	EU070	Control panel holder standard (not possible in combination with TEN° (EU052)) see exclusions 50.1)	Inclusive
☐	EU071	Control panel holder swing-away in parallel, 6 cm height-adjustable	£140.00
☐	EU072	Control panel holder swing-away in parallel, 14 cm height-adjustable	£166.00
☐	EU440	Control panel holder, swing-away in parallel, 6 cm height-adjustable, pluggable	£380.00
☐	EU441	Control panel holder, swing-away in parallel, 14 cm height-adjustable, pluggable	£380.00

Control panel holder mounting side (mandatory information)

☐		Right	
☐		Left	

Joystick top

☐	EU080	Fork for tetraplegics, horizontal (adjustable for hand width 70 - 90 mm)	£140.00
☐	EU081	Stick S80 (special joystick handle 80 mm)	£123.00
☐	EU082	Soft ball, black cover	£209.00
☐	EU083	Ball top, small (diameter 43 mm)	£141.00
☐	EU084	Ball top, large (diameter 48 mm)	£129.00
☐	EU085	Soft ball (diameter 30 mm, black, 4 pieces)	£252.00

Modules for TEN°/R-net LED-L Colour control devices see exclusions 53.0)

☐	EU054	TEN° LCD module (when using 1-switch scan drive, programming by service required; please order button separately; adapter switch box not required)	£1974.00
--------------------------	-------	--	----------

Holder for TEN° LCD module see exclusions 54.0) 54.1)

☐	EU073	Standard (installation across from control panel side; not possible in combination with ES012) see exclusions 42.1)	£102.00
☐	EU074	With gooseneck (installation across from control panel side; not possible in combination with ES012) see exclusions 42.1)	£313.00
☐	EU075	With gooseneck for tray (included)	£235.00

Actuator module for TEN°/R-net LED-L/R-net Colour control device see exclusions 53.0) 55.0) 55.1)

☐	EU056	Actuator module (required for power seat functions on R-net LED-L, R-net colour and TEN°, when more than one power seat function is needed)	£512.00
--------------------------	-------	---	---------

Control device accessories

☐	EU057	Attendant control see exclusions 56.1)	£859.00
--------------------------	-------	--	---------

Attendant control installation side see exclusions 57.0) 57.1)

☐		Right see exclusions 57.2)	
☐		Left	

Power wheelchair *Juvo B4*.

Track stabiliser see exclusions 53.0) 59.0)

☐ EF023	Power (gyro module) (not possible in combination with EU059, EU060, EU050 or EU051; depending on terrain, the specified driving distance range can be reduced when using the gyro) see exclusions 59.1) 59.2)	£527.00
--------------------------------	--	---------

Drive unit (mandatory information)

☐ ET016	Performance (usable up to a gross weight of 320 kg)	Inclusive
☐ ET011	High Performance (only possible in combination with EU052, EU058 or EU064) see exclusions 2.1) 53.0)	£681.00

Speed

☐ ET012	6 km/h	Inclusive
☐ ET013	7,2 km/h see exclusions 2.1) 61.1)	Inclusive
☐ ET014	10 km/h see exclusions 2.1) 61.1)	Inclusive

Drive unit option see exclusions 62.0) 62.1)

☐ ET008	Individual type approval for speed > 6 km/h (required only for Germany)	£615.00
--------------------------------	---	---------

Drive unit warranty extension

☐ EW002	From 24 month to 48 month (the Mobility Solutions GmbH terms and exclusions apply)	£98.00
--------------------------------	---	--------

Additional instructions for use (additional IFU for non-national language)

☐ ES028	Qualified personnel	Inclusive
☐ ES029	User	Inclusive

Power wheelchair *Juvo B4*.



EC001 - EC004 Standard seat



EC005 VAS



EC020/ED034 Contour seat with low contour




EC021/ED035 Contour seat with deep contour




ED050/051 Contour back support cushion AD



ED036 Baxx aluminium flat top back support



ED039 Side panels



ED040 Lumbar support



EC025 Seat tilt 45°



EC026 Seat height adjustment 35 cm



ED047 power back support angle adjustment

Power wheelchair *Juvo B4.*



EG013 Ribs



EG014 Lugs



EG015 Caster wheels with ribs



EL001 Gooseneck



EL002 Multi-axis, offset



EL003 Multi-axis, straight



EL031 UniLink head support holder, multi-axis



EL006/007 head support, small/large



EL013/014 Combination head/neck support small/large



EL015/016 head/neck support X-shape, small/large



EL033/034/035 Comfort Plus



EL036/037 Combo



EE024 Flip-up arm support



EP001 Standard lap belt



EP002 Lap belt with metal buckle

Power wheelchair *Juvo B4*.



EP003 Lap belt with plastic buckle



EP020 Lap belt with rewind-er, metal buckle



EP005 - EP007 Dynamic harness



EP008 - EP010 Sternum harness




EP011 - EP013 Slimline ster-num harness



EP014 - EP016 Shoulder har-ness



EB011 Standard leg supports (basic design)



EB001 Standard leg supports (new design)



EB002 Mechanically elevating leg supports with length ad-justment



EB003 Power elevating leg support, separate control



EB004 Power elevating leg supports, simultaneous control



EB012 Individual foot plate, plastic



EB007 Individual foot plate, aluminium



EB017 Lock for single panel foot plate



EB022/023 Amputation leg support right/left

Power wheelchair *Juvo B4*.



EU043 External magnetic charging receptacle



EA010 Splash guard for drive wheels



EU041 Battery charger 8 A, fanless



EU044 Battery charger 10 A, fanless (IP54)



EU046 Battery charger 10 A, with fan



EU042 Battery charger 12 A, incl. external charging receptacle



EF021 Curb climbing assist



ES013 Mounting kit for Dahl docking system



ES003 Pocket for mobile phone



ES004 Luggage carrier




ES005 Crutch holder



ES006 Rear-view mirror



ES010 Beverage holder



ES023 Hooks for backpack



ES305 Backpack with Otto-bock logo



ES306 Catheter bag, large

Power wheelchair *Juvo B4*.



**ET002 Receptacle for 12 V
ET003 Receptacle für 24 V
ET004 USB**



ES021 Electr. lighting



EU050/059 VR2



EU051/060 VR2



EU052 TEN°



EU058 R-net LED-L



**EU080 Fork for Tetraplegics,
horizontal**



EU081 Stick S80




EU082 Soft ball, black cover



EU083 Ball top, small



EU084 Ball top, large



EU085 Soft ball



EU054 TEN° LCD module



EU057 Attendant control



EU125 P&G Wi-Fi Dongle



**EU095 Hand warmer (single
hood for handmodule)**

Power wheelchair Juvo B4.

Exclusions

- 2.1) Not possible in conjunction with product classification Category A (Indoor)
- 3.1) Possible only in combination with maximum user weight <= 140 kg
- 3.2) Possible only in combination with maximum user weight <= 75 kg
- 3.3) Possible only in combination with maximum user weight <= 125 kg
- 4.0) See table "seat depth" (490E75=1-Sitztiefe)
- 4.1) Not possible in conjunction with the combination of maximum user weight > 75 kg and standard seat, Junior or VAS
- 4.4) Possible only in combination with maximum load capacity <= 113 kg
- 5.0) See table "seat width" (490E75=1-Sitzbreite)
- 7.1) Possible only in combination with maximum load capacity <= 150 kg
- 7.2) See table "seat cushion Terra, Terra Aquos" (490E75=1-Terra_T-Aquos)
- 7.3) See table "seat cushion Terra-Flair" (490E75=1-Terra-Flair)
- 7.5) See table "seat cushion Cloud" (490E75=1-Cloud)
- 7.6) See table "seat cushion Livo" (490E75=1-Livo)
- 7.7) See table "seat cushion EvoFlair" (490E75=1-EvoFlair)
- 7.8) Cannot be combined with standard seat, Junior
- 7.9) Not possible in conjunction with the combination of VAS (Vario Adjust Seat), seat width 38 cm or 40 cm and seat depth > 46 cm, see table "Seat cushion Contour with low or deep contour" (490E75=0-Contour)
- 7.10) Not possible in conjunction with the combination of VAS (Vario Adjust Seat), seat width > 40 cm and seat depth 38 cm or 40 cm, see table "Seat cushion Contour with low or deep contour" (490E75=0-Contour)
- 7.11) Not possible in conjunction with the combination of VAS (Vario Adjust Seat) and seat depth 58 cm, see table "Seat cushion Contour with low or deep contour" (490E75=0-Contour)
- 9.1) Possible only in combination with VAS
- 9.3) Cannot be combined with seat width 54 cm
- 11.0) Possible only in combination with back version Contour with low or deep contour
- 11.1) Possible only in combination with back cushion with deep contours
- 12) Mandatory information in combination with Seat systems standard seat, small and SW 38 cm
- 12.1) Not possible in conjunction with seat systems standard seat, small and SW 38 cm and unselected Upholstery covers for arm supports
- 14.0) Possible only in combination with Baxx aluminium flat top back
- 14.1) Mandatory information in combination with Baxx aluminium flat top back
- 17.1) Cannot be combined with VAS
- 18.1) Possible only with seat inclination 0°
- 18.2) Cannot be combined with seat inclination -3° / 6° / 9°
- 20.1) Not possible in conjunction with caster wheel swing arm spring mounted
- 20.2) Not possible in conjunction with product classification Category B (Indoor/Outdoor) or product classification Category A (Indoor) (version HCD)
- 22.1) Not possible in conjunction with caster wheels 8" or drive wheels 12"
- 23.1) Possible only in conjunction with tyre colour grey
- 24.0) Seat height >= lower leg length +3 cm
- 24.1) Possible only in combination with product classification category A (Indoor)
- 24.2) Cannot be combined with seat tilt 45° or seat height adjustment 35 cm or seat height adjustment 35 cm with seat tilt 45°
- 24.3) Not possible in conjunction with seat systems VAS and caster wheel swing arm unsprung and caster wheels 9"
- 24.4) Not possible in combination with VAS and maximum user weight > 75 kg
- 24.5) Cannot be combined with seat tilt 20°
- 24.6) Possible only in combination with seat tilt 45° or seat height adjustment 35cm or seat height adjustment 35 cm with seat tilt 45°
- 24.7) Possible only in combination with caster wheel swing arm unsprung and caster wheels 9"
- 24.8) Possible only in combination with caster wheel swing arm spring mounted or caster wheels 10"
- 24.9) Not possible in conjunction with the combination of maximum user weight > 100 kg and VAS
- 24.10) Not possible in conjunction with the combination of seat height adjustment up to 35 cm with seat tilt 45° and VAS
- 24.11) Not possible in conjunction with the combination of VAS and SW > 42 cm and seat inclination 6°, 9° and caster wheel swing arm spring mounted or caster wheels 10"
- 24.12) Not possible in conjunction with the combination of VAS and caster wheel swing arm spring mounted or caster wheels 10"
- 24.13) Cannot be combined with caster wheels 10" or caster wheel swing arm spring mounted
- 24.14) Cannot be combined seat inclination -3°
- 24.15) Not possible in conjunction with the combination of seat inclination -3° and caster wheels 10" or caster wheel swing arm spring mounted
- 24.16) Not possible in conjunction with seat systems standard seat large / XL and seat tilt 20° and caster wheels 8"
- 24.17) Not possible in conjunction with seat systems standard seat large / XL and seat tilt 20° and caster wheel swing arm unsprung and caster wheels 9"
- 24.18) Not possible in conjunction with the combination of VAS and seat inclination 9° and caster wheel swing arm unsprung and caster wheels 9"
- 24.19) Not possible in conjunction with the combination of VAS and seat inclination 9° and caster wheel swing arm spring mounted or caster wheels 10"
- 24.20) Not possible in conjunction with seat systems VAS and seat inclination 9° and caster wheel swing arm unsprung and caster wheels 10" or caster wheel swing arm sprung and caster wheels 9"
- 24.21) Not possible in conjunction with seat systems VAS and caster wheels 8"
- 24.22) Not possible in conjunction with seat systems standard seat Junior / small / large / XL and seat inclination 6°, 9° and caster wheel swing arm unsprung and caster wheels 10" or caster wheel swing arm sprung and caster wheels 9"
- 24.23) Not possible in conjunction with seat systems standard seat Junior / small / large / XL and seat inclination 3°, 6°, 9° and caster wheel swing arm unsprung and caster wheels 9"
- 24.24) Not possible in conjunction with seat systems standard seat Junior / small / large / XL and seat inclination 9° and caster wheel swing arm unsprung and caster wheels 9"
- 24.25) Not possible in conjunction with the combination of caster wheels 10" or caster wheel swing arm spring mounted and seat tilt 20°
- 24.26) Not possible in conjunction with seat systems standard seat Junior / small / large / XL and seat inclination 9° and caster wheels 8"
- 24.27) Not possible in conjunction with seat systems standard seat large / XL and seat tilt 20° and caster wheel swing arm unsprung and caster wheels 10" or caster wheel swing arm sprung and caster wheels 9"

Power wheelchair Juvo B4.

- 24.28) Not possible in conjunction with the combination of caster wheel swing arm unsprung and caster wheels 9" and seat tilt 20°
- 24.29) Not possible in conjunction with the combination of VAS and SW > 42 cm and caster wheels 8"
- 24.31) Not possible in conjunction with the combination of VAS and seat inclination 9° and caster wheels 8"
- 24.32) Not possible in conjunction with seat systems standard seat large / XL
- 24.33) Not possible in conjunction with seat systems standard seat large / XL and caster wheel swing arm unsprung and caster wheels 9"
- 24.34) Not possible in combination with seat height adjustment up to 35 cm with seat tilt 45° and maximum user weight > 75 kg
- 24.35) Not possible in combination with seat tilt 20° and seat inclination 3° and caster wheel swing arm unsprung and caster wheels 9"
- 24.36) Not possible in combination with seat tilt 20° and seat inclination 3° and caster wheel swing arm unsprung and caster wheels 8"
- 24.37) See table "seat hight" (490E75=1-Sitzhoehe)
- 27.0) Possible only in combination with head/neck support mounting kit
- 27.1) Mandatory information in combination with head/neck support mounting kit
- 27.2) Possible only in combination with head/neck support mounting kit Goose Neck or Multi-axis, offset or Multi-axis, straight
- 28.0) Mandatory information in combination with category B
- 28.1) Not possible in conjunction with the combination of standard seat, small, seat depth < 42 cm and seat width > 38 cm
- 30.1) Possible only in combination with selected lower leg length
- 30.2) Not possible in conjunction with the combination of standard seat or VAS and lower leg length < 21 cm
- 31.1) Possible only in combination with selected lower leg length
- 31.2) Not possible in conjunction with the combination of standard seat or VAS and lower leg length < 21 cm
- 31.3) Not possible in conjunction with the combination of maximum user weight > 75 kg and standard seat and lower leg length < 28 cm
- 31.4) Not possible in conjunction with the combination of maximum user weight > 75 kg and VAS and lower leg length < 27 cm
- 31.5) Not possible in conjunction with selected lower leg length
- 31.7) Not possible in conjunction with the combination of seat depth 38 cm, belts lap belt with rewriter, metal buckle, accessories for VAS thigh lateral support, size M or thigh lateral support, size S and arm support, standard
- 31.8) Not possible in conjunction with the combination of seat depth < 44 cm, belts lap belt, standard, lap belt, padded, metal buckle or lap belt, padded, plastic buckle, accessories for VAS thigh lateral support, size M or thigh lateral support, size S and arm support, flip-up
- 31.6) Not possible in conjunction with the combination of seat height (SH) < 45 cm and lower leg length (LL) <= 30 cm
- 32.0) Cannot be combined with leg support version, without leg support
- 32.1) Mandatory information in combination with leg support version, standard, mechanically elevating or electrically elevating
- 32.2) Not possible in conjunction with the combination of standard seat and lower leg length < 25 cm
- 32.3) Not possible in conjunction with the combination of VAS and lower leg length < 24 cm
- 32.4) Cannot be combined with leg support version, mechanically elevating or electrically elevating
- 32.6) Cannot be combined with leg support version standard (Basic design)
- 32.7) Not possible in conjunction with the combination of standard seat or VAS and lower leg length < 17 cm and leg support version, standard
- 32.8) Not possible in conjunction with the combination of maximum user weight > 75 kg and standard seat or VAS and lower leg length < 24 cm and leg support version, standard
- 32.9) Cannot be combined with leg support version, mechanically elevating or electrically elevating, separate controllable
- 32.10) Not possible in conjunction with the combination of maximum user weight > 140 kg and leg support version, electrically elevating, simultaneously controllable
- 33.1) Possible only in combination with leg support version, mechanically or electrically elevating
- 33.2) Possible only in combination with Foot plate version Single-panel, aluminium foot plate
- 33.3) Not possible in conjunction with the combination of Leg support version Standard (New Design) and Lower leg length < 24 cm
- 33.4) Not possible in conjunction with the combination of Leg support version Electrically elevating, simult. controllable and Lower leg length < 28 cm
- 33.5) Possible only in combination with Leg support version Standard (Basic design) or Standard (New Design) or Without leg support
- 35.2) Not possible in conjunction with the combination of Argentina, Australia, Brazil, Hong Kong, Japan, Canada, Columbia, Mexico, Thailand, United Kingdom or United States of America as the country of destination
- 35.3) Cannot be combined with AGM 63 Ah (C5) or Gel 63 Ah (C5) batteries
- 35.4) Not possible in conjunction with the combination of Argentina, Brazil, Canada, Columbia, Mexico or United States of America as the country of destination
- 36.1) Possible only in conjunction with battery charger, without
- 39.5) Not possible in conjunction with the combination of SH 40 cm and caster wheel swing arm unsprung and caster wheels 9" and VAS or standard seat large / XL
- 39.6) Not possible in conjunction with the combination of SH 42,5 cm and caster wheel swing arm Spring mounted or caster wheels 10" and VAS or standard seat large / XL
- 41.2) Cannot be combined with chassis safety accessories curb climbing assist
- 41.3) Not possible in conjunction with the combination of Electric seat adjustment Seat tilt 20° Electrical, Belts 2-point Lap belt with rewriter, metal buckle and SD <44 cm
- 42.1) Cannot be combined with therapy tray
- 42.2) Cannot be combined with rearview mirror
- 42.3) Cannot be combined with Germany as the country of destination
- 42.4) Not possible in conjunction with Australia as the country of destination
- 42.5) Not possible in conjunction with the combination of standard seat, small and seat depth < 42 cm
- 42.6) Not possible in conjunction with the combination of standard seat, small and electric seat adjustment seat tilt 45°, seat height adjustment up to 35 cm or seat height adjustment up to 35 cm with seat tilt 45°
- 42.7) Not possible in conjunction with the combination of VAS (Vario Adjust Seat) and seat depth < 46 cm
- 42.8) Possible only with arm support, flip-up
- 43.1) Cannot be combined with socket for external power supply, right 24 V and socket for external power supply, right USB
- 43.2) Cannot be combined with socket for external power supply, right 12 V and socket for external power supply, right USB
- 43.3) Cannot be combined with socket for external power supply, right 12 V and socket for external power supply, right 24 V
- 43.4) Cannot be combined with batteries/battery capacity AGM 63 Ah (C5) / 75 Ah (C20)
- 44.0) Mandatory information in conjunction with the combination of product classification category B and all countries of destination except Other countries, Argentina, Australia, Brazil, France, Hong Kong, Japan, Canada, Columbia, Korea, Mexico, Russia, Switzerland, Thailand, Turkey, United States of America or United Kingdom
- 44.1) Not possible in conjunction with the combination of Germany as the country of destination and product classification category B
- 45.1) Theoretical total weight (without special requests) = sum of the selected characteristic value weights

Power wheelchair Juvo B4.

- 46.1) Cannot be combined with electric functions
- 46.2) Possible only in combination with total weight <= 260 kg
- 46.3) Possible only in combination with max. 3 electric functions
- 46.4) Cannot be combined with lighting acc. to German StVZO
- 46.5) Not possible in conjunction with the combination of country of destination Germany and more than 2 electric functions
- 46.8) Not possible in conjunction with the combination of electric seat adjustment seat height adjustment up to 35 cm, back angle adjustment, electric and and leg support version, electrically elevating, simultaneously controllable
- 50.1) Cannot be combined with control unit TEN°
- 53.0) Possible only in combination with control unit, TEN° or R-net LED-L (90 A) or R-net Colour (90 A)
- 54.0) Possible only in combination with TEN° LCD module
- 54.1) Mandatory information in combination with TEN° LCD module
- 55.0) Mandatory information in conjunction with the combination of control unit, TEN° (90A) or R-net LED-L (90 A) or R-net Colour (90 A) and > 1 electric function
- 55.1) Mandatory information in conjunction with the combination of lighting acc. to German StVZO and control unit, TEN° or R-net LED-L (90 A)
- 56.1) Cannot be combined with Control panel holder swing-away in parallel, 6 cm adjustment range, pluggable or swing-away in parallel, 14 cm adjustment range, pluggable
- 57.0) Possible only in combination with attendant control
- 57.1) Mandatory information in combination with attendant control
- 57.2) Not possible in conjunction with the combination of standard seat, Junior and actuator module for control unit, TEN°
- 58.1) Possible only in combination with control unit, VR2
- 59.0) Not possible in conjunction with Australia, Canada or the United States of America as the country of destination
- 59.1) Not possible in conjunction with the combination of batteries Gel 63 Ah (C5) / 74 Ah (C20) and General accessories external horn
- 59.2) Not possible in conjunction with the combination of batteries AGM 63 Ah (C5) / 75 Ah (C20) and general accessories external horn
- 61.1) For Germany as the country of destination, this is possible only in conjunction with the combination of rear warning sign and lighting acc. to German StVZO
- 62.0) Possible only in conjunction with the combination of Germany as the country of destination and speed > 6 km/h
- 62.1) Mandatory information in conjunction with the combination of Germany as the country of destination and speed > 6 km/h
- 63.0) Not possible in conjunction with Australia as the country of destination
- 64.0) Possible only in combination with additional instructions for use, qualified personnel
- 64.1) Mandatory information in combination with additional instructions for use, qualified personnel
- 65.0) Possible only in combination with additional instructions for use, user
- 65.1) Mandatory information in combination with additional instructions for use, user
- 67.1) Not possible in conjunction with the combination of standard seat, small and seat width 38 cm or 40 cm
- 67.2) Not possible in conjunction with the combination of standard seat, large and seat width 44 cm or 46 cm
- 67.3) Not possible in conjunction with the combination of standard seat, XL and seat width 50 cm or 52 cm
- 67.4) Not possible in conjunction with 2-point Lap belt with rewinder, metal buckle
- 67.5) Cannot be combined with seat depth < 48 cm
- 67.6) Possible only in combination with maximum user weight >= 125 kg
- 68.1) Possible only in combination with electric seat adjustment seat height adjustment up to 35 cm or seat height adjustment up to 35 cm with seat tilt 45°
- 69.0) Possible only in combination with head/neck support mounting kit UniLink, multi-axis
- 71.1) Not possible in conjunction with the combination of back version back upholstery, adaptable and seat width < 46 cm
- 71.2) Cannot be combined with back version back cushion Contour with deep contours, AD or back cushion Contour with flat contours, AD or Baxx aluminium flat top back
- 71.3) Cannot be combined with accessories for VAS push handles
- 72.0) Possible only in combination with seat cushion Contour with flat or deep contours