2021/22

ottobock.

Kimba Family



Kimba

Products for combinations



Kimba seating unit



Kimba outdoor mobility base



Kimba Cross



Kimba Inline



Kimba Home



Kimba indoor mobility base

What c	lo you want to order?	
	Kimba	see page 3
	Kimba seating unit - Please ensure you order 2 seat units if ordering the Kimba Inlne	see page 12
	Kimba outdoor mobility base	see page 18
	Kimba Cross	see page 21
	Kimba Inline - Please ensure you order 2 seat units if ordering this frame	see page 23
	Kimba Home	see page 24
	Kimba indoor mobility base	see page 26

	Dimensions kid	Dimensions Kimba size 1	Dimensions Kimba size 2	
① Back length [mm]		410 - 610	560 - 700	
② Seat depth [mm]		190 - 310 (140 shortened)	240 - 400	
Seat width [mm]		200 - 300	260 - 350	
③ Lower leg length [mm]		190 - 310 (122 shortened)	200 - 370	
Load capacity			nit: 40 kg* lity base: 55 kg	
a) User weight				
b) Medical Equipment				
* With option "footrest angle adjustab	le for short lower leg length" the	maximum load decreases to 35 kg	g.	

		Order no.:	Kimba Frame and Seat
			Article number: 470G71=5_AA04_C
4Valid from 01.04.22 to 31.03.23 Please ensure youa re using the latest o	order form		
Provide estimate only	□Order		
Attachment for custom fabr	rication enclos	ed	
Delivery address			
Sales person	Name		Figure includes accessorie
	 Date Si	ignature	Standard equipment Seat and back support Seat width, seat depth, back support height
By placing the order, you agree to comp nuals for accessories and products of of with an Ottobock product, in addition t of the respective Ottobock products. A and accepted, our general terms and co placing an order, to conduct your own of cordance with national requirements if attachments that do not meet the spec service manuals.	ther manufacturer to the instructions lso you confirm to poditions. In additi conformity assessr and when you imp	s when combining them for use/service manuals have read, understood on, you agree, upon nent or approval in ac- plement modifications or	 and lower leg length are continuously adjustable Seat inclination and knee angle continuously adjustable Back support angle adjustable Optimized quick-release mechanism to remove the seat Position the seat in/against the direction of travel Leg support assembly with flip-up foot plate Lateral side supports

Adjustable pelvic lateralsFrame colour: blackPadding colour grey

• Optimized rear axle suspension

platform for respirators)Automatic folding lock

• EVA tyres with frictionbrake

wheelchair accessible vehicle

grade aluminium

locking pin

limitation

swivel lock

Folding outdoor mobility base made of high-

Push bar with ratchet joint, height-adjustableSeat tilt with tilt activation at the frame incl.

Integrated storage space (also usable as a

• Optimized seat adapter with automatic tilt

• Front wheels 7", swivelling with integrated

• Integrated anchor points for mounting to

Seating unit						
KC01		Size 1 (from approx. 1	- 6 years)	£945.00		
KC02		Size 2 (from approx. 4	ze 2 (from approx. 4 - 10 years)			
Head support		[Γ			
KD01		T	Without head support	Inclusive		
KD02		The second se	With integrated head support with slight lateral support	£196.00		
KD03B KD03			Head support (laterally adjustable) B = enclosed mounted	£156.00 £161.00		
KL01B			Head support concave, S, enclosed (incl. holder; is delivered in combination with a short back)	£233.00		
KL02B			Head support concave, M, enclosed (incl. holder; is delivered in combination with a short back)	£237.00		
KL03B			Head support lateral, S, enclosed (incl. holder; is delivered in combination with a short back)	£342.000		
KL06B			Head support lateral adjustable, M, enclosed (incl. holder; is delivered in combination with a short back)	£342.00		
KL05B			Head support holder, enclosed (is delivered in combination with a short back)	£143.00		
Padding						
KD00			Without (for individual seat and upholstery adjustment)	Inclusive		
KD41B			Padding, enclosed	Inclusive		
KD41			Padding, mounted	£13.00		
Inlays	-					
KD33B			Inlay grey, enclosed	£59.00		
KD34B			Inlay turquoise, enclosed	£59.00		
KD35B			Inlay pink, enclosed	£59.00		
KD36B			Inlay green, enclosed	£59.00		

Pelvic belts		1			
KP02B KP02		Lap belt, padded	□ B = enclosed	mounted	£59.00 £68.00
KP18B KP18		4-Point paediactric lap belt, with plastic buckles (width: 2,5 cm; length: 100 cm)	□ B = enclosed	mounted	£99.00 £109.00
KP06B KP06	633	Groin strap	☐ B = enclosed	mounted	£61.00 £71.00
KP15B KP15		Pediatric lap belt, with buckles (width: 2,5 cm; length: 100 cm)	□ B = enclosed	□ mounted	£94.00 £99.00
Chest and should	er supports				
KP16B KP16		Shoulder harness Size S	B = enclosed	□ mounted	£158.00 £163.00
KP17B KP17		Shoulder harness Size M	□ B = enclosed	□ mounted	£158.00 £163.00
KP19B KP19	EB .	Chest and shoulder pad Size S	☐ B = enclosed	 mounted	£145.00 £151.00
KP20B KP20	1241	Chest and shoulder pad Size M	☐ B = enclosed	□ mounted	£145.00 £151.00
KP21B KP21	2	Slimline sternum harness	□ B = enclosed	mounted	£132.00 £138.00
KP22B KP22		Dynamic harness Size XS	B = enclosed	□ mounted	£137.00 £141.00
KP23B KP23		Dynamic harness Size S	B = enclosed	mounted	£142.00 £147.00
KP24B KP24		Dynamic harness Size M	B = enclosed	 mounted	£142.00 £147.00
KP03B KP03		4-point belt, padded	□ B = enclosed	mounted	£101.00 £106.00
KP05B KP05		5-Point belt, small buckle (recommended for seating unit size 1)	☐ B = enclosed	□ mounted	£101.00 £112.00
KP04B KP04		5-Point belt, large buckle	□ B = enclosed	 mounted	£101.00 £112.00
KP10B KP10		Positioning vest	☐ B = enclosed	 mounted	£81.00 £91.00

Preassembly				
KP14			Mounting kit strap/vest, mounted	£27.00
Trunk support				
KD04B KD04			Thoracic support $B = enclosed$ mounted	£81.00 £91.00
KD31B			Swing away thoracic supports (laterals), enclosed	£352.00
KD10B			Lumbar support padding, enclosed	£49.00
KD11B		4	Neck support enclosed	£58.00
КС10			Seat depth reducer (seat depth possible up to 14 cm)	£101.00
KC13B			Inlay small, enclosed (spacer fabric grey, only for size 1)	£101.00
Arm supports				
KE01		F	Arm support, Standard, mounted	£96.00
KE02			Arm support hight/low, adjustable, mounted	£135.00
Tray	T	Γ		
KS01B		T	Tray, grey, enclosed (only possible in direction of travel and in combination with KE 02 "arm support high/ low" in lowest setting only)	£194.00
KS02B		17	Tray, transparent, enclosed (only possible in direction of travel and in combination with KE 02 "arm support high/ low" in lowest setting only))	£314.00

Leg and foot su	pport		F			
КВОЗ			Leg support, Standard, mounted (from LLL 19 cm)			Inclusive
KB05		A REAL PROPERTY OF	Leg support, angle-adjustable, mounted (from LLL 19 cm)			£61.00
KB10			Leg support, angle-adjustable for short lower leg length, mounted (only for size 1, from LLL 14 cm)			£80.00
Accessories leg	and fo	oot support				
KB08B KB08			Slip resistant footrest cover (not mounted in combination with KB10)	□ B = enclosed	mounted	£21.00 £25.00
КВО9В		Ő	Padding, enclosed			£58.00
KB12B		P	Foot support – sandals, enclosed			£91.00
КВ07В			Velcro footstraps, instep padded			£25.00
КС06В КС06			Abductor size 1, only for seat size 1, (L 9 cm × W 6 cm × H 9 cm)	B = enclosed	□ mounted	£91.00 £97.00
КС07В КС07			Abductor size 2 (L 10 cm × W 7,5 cm × H 10 cm)	□ B = enclosed	 mounted	£96.00 £100.00
КС08В КС08			Abductor size 3, only for seat size 2 (L 12 cm × W 9 cm × H 10 cm)	B = enclosed	□ mounted	£99.00 £104.00
Canopy						
KW13			Canopy in grey-melange, inside grey, incl. rain cover and mounted adapter			£243.00

Slip sack				
KW01B			Summer/winter slipsack Size S, enclosed	£247.00
KW02B		15	Summer/winter slipsack Size M, enclosed	£247.00
КѠ03В			Summer/winter slipsack Size L, enclosed	£247.00
Accessories sea	t			
KQ04B			Grip bar with padding, enclosed	£51.00
KD20			Back angle limitation, mounted	£21.00
KW04B		C.S.	Rain cape, grey, enclosed	£45.00
KD12B			Incontinence cover, enclosed	£121.00

Outdoor mobility base size					
KA01		Size 1 (from approx. 1	Size 1 (from approx. 1 to 6 years) £1326.0		
KA05		Size 2 (from approx. 4	Size 2 (from approx. 4 to 10 years) £1427.0		
Frame colour					
013		Jet black RAL 9005 ma	Jet black RAL 9005 matte		
023		White RAL 9016	White RAL 9016		
Reflector sticker					
KA13			Classic	£25.00	
KA14		Sporty £42.00			



Classic

Colour reflector sticker Grey \Box Green Turquoise

Tyres	-	-		
KF03			EVA	Inclusive
KF04			Pneumatic	£45.00
Rear wheel type				
KG01		Wheel lock (actuated	by the foot via wheel lock bar)	Inclusive
KG02		With frictionbrake (ma	nually actuated via brake levers at the push bar)	£266.00
Accessories				
КА07В			Frame padding, pair, enclosed	£30.00
KS03B		A good	Buggy board, enclosed	£118.00
KA08			Tip-assist, left	£27.00
KA12			Tip-assist, right	£27.00
KS14B		\sim	Tether band, enclosed	£17.00
Accessories general				
KS07B			Cupholder, enclosed	£21.00
KS13B			Diaper bag, grey-melange, enclosed (Contains extensive accessories such as: changing mat, dummy bag, isolation pocket for bottles and a plastic uten- sil bag, as well as numerous compartments and zipped pockets)	£121.00

Products for combinations

Seating unit



Kimba seating unit

Outdoor mobility bases





Kimba outdoor mobility base

Kimba Cross

Indoor mobility bases



Kimba Home



Kimba indoor mobility base



Kimba Inline

		Order no.:			
Valid from 01.04.22 to 31.03.23 Please ensure youa re using the latest orc	der form				
Provide estimate only	□ Order				
Attachment for custom fabric	ation enclose	d			
Billing address/customer number					
Delivery address					
Sales person	Name				
	 Date Sig	gnature			
By placing the order, you agree to comply with the instructions for use/service ma- nuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in ac- cordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/ service manuals.					

Kimba Seating unit

Maximum load capacity: 40 kg Article number: 470G71=5_AA02_C



Figure includes accessories

Standard Equipment

- Seat and back support
- Seat width, seat depth, back support height and lower leg length are continuously adjustable
- Seat inclination and knee angle continuously adjustable
- Back support angle adjustable
- Optimized quick-release mechanism to remove the seat
- Position the seat in/against the direction of travel
- Leg support assembly with flip-up foot plate
- Lateral side supports
- Adjustable pelvic laterals
- Frame colour: black
- Padding colour grey

Notice:

This product can be combined with outdoor mobility bases Kimba 470G71=5_AA03_C, Kimba Cross and Kimba Inline and also with indoor mobility bases Kimba Home and Kimba. Please fill out the corresponding order sheets separately.

Please complete two of these Seat Unit formsif you are ordering a Kimba Inline Frame

Size		
KC01	Size 1 (from approx. 1 to 6 years)	£945.00
KC02	Size 2 (from approx. 4 to 10 years)	£1046.00

Head support			
KD01		Without head support Incl	usive
KD02	III III	With integrated head support with slight lateral support £19	96.00
KD03B KD03			56.00 61.00
KL01B		Head support concave, S,enclosed£2:(incl. holder; is delivered in combination with a short back)	33.00
KL02B		Head support concave, M, enclosed£2.(incl. holder; is delivered in combination with a short back)	37.00
KL03B		Head support lateral, S, enclosed £34 (incl. holder; is delivered in combination with a short back)	42.00
KL06B		Head support lateral adjustable, M, enclosed £34 (incl. holder; is delivered in combination with a short back)	42.00
KL05B		Head support holder, enclosed £14 (is delivered in combination with a short back)	43.00
Padding			
KD00		Without (for individual seat and upholstery adjustment) Incl	usive
KD41B		Padding, Incl enclosed	usive
KD41		Padding, £	13.00
Inlays			
KD33B		Inlay grey, enclosed £	59.00
KD34B		Inlay turquoise, enclosed	59.00
KD35B	G	Inlay pink, £	59.00
KD36B		Inlay green, £	59.00

Pelvic belts					
KP02B KP02		Lap belt, padded	☐ B = enclosed	mounted	£59.00 £68.00
KP18B KP18		4-Point paediactric lap belt, with plastic buckles (width: 2,5 cm; length: 100 cm)	□ B = enclosed	mounted	£99.00 £109.00
KP06B KP06		Groin strap	□ B = enclosed	 mounted	£61.00 £71.00
KP15B KP15		Pediatric lap belt, with buckles (width: 2,5 cm; length: 100 cm)	□ B = enclosed	 mounted	£94.00 £99.00
Chest and shoulder	supports				
KP16B KP16	NR()	Shoulder harness Size S	B = enclosed	 mounted	£158.00 £163.00
KP17B KP17		Shoulder harness Size M	B = enclosed	mounted	£158.00 £163.00
KP19B KP19	AB	Chest and shoulder pad Size S	B = enclosed	 mounted	£145.00 £151.00
KP20B KP20		Chest and shoulder pad Size M	B = enclosed	 mounted	£145.00 £151.00
KP21B KP21	2	Slimline sternum harness	☐ B = enclosed	 mounted	£132.00 £138.00
KP22B KP22	4	Dynamic harness Size XS	□ B = enclosed	 mounted	£137.00 £141.00
KP23B KP23	M.	Dynamic harness Size S	B = enclosed	 mounted	£142.00 £147.00
KP24B KP24		Dynamic harness Size M	B = enclosed	mounted	£142.00 £147.00
KP03B KP03		4-point belt, padded	☐ B = enclosed	 mounted	£101.00 £106.00
KP05B KP05	M	5-Point belt, small buckle (recommended for seating unit size 1)	☐ B = enclosed	 mounted	£101.00 £112.00
KP04B KP04		5-Point belt, large buckle	B = enclosed	 mounted	£101.00 £112.00
KP10B KP10		Positioning vest	☐ B = enclosed	mounted	£81.00 £91.00

Preassembly							
KP14			Mounting kit strap/vest, mounted	£27.00			
Trunk support			-				
KD04B KD04			Thoracic support $B = enclosed$ mounted	£81.00 £91.00			
KD31B			Swing away thoracic supports (laterals), enclosed	£352.00			
KD10B			Lumbar support padding, enclosed	£49.00			
KD11B		4	Neck support, enclosed	£58.00			
КС10			Seat depth reducer (seat depth possible up to 14 cm)	£101.00			
KC13B			Inlay small, enclosed (spacer fabric grey, only for size 1)	£101.00			
Arm supports							
KE01		F	Arm support, Standard, mounted	£96.00			
KE02			Arm support hight/low, adjustable, mounted	£135.00			
Tray							
KS01B		T	Tray, grey, enclosed (only possible in direction of travel and in combination with KE 02 "arm support high/ low" in lowest setting only)	£194.00			
KS02B		TT	Tray, transparent, enclosed (only possible in direction of travel and in combination with KE 02 "arm support high/ low" in lowest setting only))	£314.00			

Leg and foot	support	Γ	1			
КВОЗ			Leg support, Standard, mounted (from LLL 19 cm)			Inclusive
KB05		- And	Leg support, angle-adjustable, mounted (from LLL 19 cm)			£61.00
КВ10			Leg support, angle-adjustable for short lower leg length, mounted (only for size 1, from LLL 14 cm)	£80.00		
Accessories le	g and fo	pot support	-			
KB08B KB08			Slip resistant footrest cover (not mounted in combination with KB10)	B = enclosed	mounted	£21.00 £25.00
КВО9В		Ĩ	Padding, enclosed			£58.00
KB12B		P	Foot support – sandals, enclosed			£91.00
KB07B			Velcro footstraps, instep padded			£25.00
KC06B KC06	·		Abductor size 1, only for seat size 1, (L 9 cm × W 6 cm × H 9 cm)	□ B = enclosed	mounted	£91.00 £97.00
КС07В КС07	КС07 КС08В		Abductor size 2 (L 10 cm × W 7,5 cm × H 10 cm)	☐ B = enclosed	mounted	£96.00 £104.00
KC08B KC08			Abductor size 3, only for seat size 2 (L 12 cm × W 9 cm × H 10 cm)	B = enclosed	mounted	£99.00 £104.00
Canopy						
KW13			Canopy in grey-melange, inside grey, incl. rain cover and mounted adapter			£243.00

Slip sack				
KW01B			Summer/winter slipsack Size S, enclosed	£247.00
KW02B			Summer/winter slipsack Size M, enclosed	£247.00
KW03B			Summer/winter slipsack Size L, enclosed	£247.00
Accessories sea	at			
KQ04B			Grip bar with padding, enclosed	£51.00
KD20			Back angle limitation, mounted	£21.00
KW04B		C.S.	Rain cape, grey, enclosed	£45.00
KD12B			Incontinence cover, enclosed	£121.00
KC16B			Retrofit kit for Kimba Neo street frame (without tilt; not possible in conbimation with costum fabrication platform for respirators up to 90 kg/ 198,4 lb)	£61.00

				ottobock
				Kimba Outdoor
			Auftrags-Nr.:	mobility base
				Maximum load capacity: 55 kg (including seating unit)
Valid from 01.04.22 to 31.03 Please ensure youa re using		test order form		Article number: 470G71=5_AA03_C
				a a
Provide estimate o	nly	□ Order		
Attachment for cus	stom	fabrication enclose	ed	
Billing address/custome	r num	ıber		Contraction of the second seco
Delivery address				6
				 Standard equipment Folding outdoor mobility base made of high-Folding outdoor mobility base made of
Sales person		Name		 high-grade aluminium Optimized rear axle suspension Push bar with ratchet joint, height-adjustable
		Dete C		 Seat tilt with tilt activation at the frame incl. locking pin
By placing the order, you ag nuals for accessories and pr with an Ottobock product, i of the respective Ottobock and accepted, our general t placing an order, to conduc cordance with national requ attachments that do not me service manuals.	oducts n addi produc erms a t your uireme	comply with the instruct of other manufacturers tion to the instructions f ts. Also you confirm to h nd conditions. In additic own conformity assessm nts if and when you imp	when combining them for use/service manuals have read, understood on, you agree, upon hent or approval in ac- lement modifications or	 Integrated storage space (also usable as a platform for respirators) Automatic folding lock Optimized seat adapter with automatic tilt limitation EVA tyres with frictionbrake Front wheels 7", swivelling with integrated swivel lock Integrated anchor points for mounting to wheelchair accessible vehicle
Size		-		
KA01		Size 1 (from approx.	1 to 6 years)	£1326.00
KA05		Size 2 (from approx.	4 to 10 years)	£1427.00
Adapter Sitzschale	1	1		
КСОЗ			Seating shell interface	£125.00
KB01			Leg support assistance f	or seating shell interface £382.00
KC14			Universal adapter (seat s (only for country of dest	+8100
KC15B	1 🗆		Retrofit-kit: Kimba Neo s	eat £41.00

Frame colour					
013	Jet black RAL 9005 ma	black RAL 9005 matte			Inclusive
023	White RAL 9016	Vhite RAL 9016			Inclusive
Reflector sticker					
KA13		Classic			£25.00
KA14		Sporty			£42.00





Classic

Colour reflector sticker					
			Grey		
			Green		
			Turquoise		
Tyres					
KF03		I	EVA	Inclusive	
KF04			Pneumatic	£45.00	
Rear wheel type	÷				
KG01		Wheel lock (actuated by the foot via wheel lock bar) Inclus		Inclusive	
KG02		With frictionbrake (ma	With frictionbrake (manually actuated via brake levers at the push bar) £266.00		

Accessories			
КА07В		Frame padding, pair, enclosed	£25.00
KS03B	M	Buggy board, enclosed	£118.00
KA08	1	Tip-assist, left	£27.00
KA12		Tip-assist, right	£27.00
KS14B	\sim	Tether band, enclosed	£17.00
Accessories general			
KS07B		Cupholder, enclosed	£21.00
KS13B		Diaper bag, grey-melange, enclosed (Contains extensive accessories such as: changing mat, dummy bag, isolation pocket for bottles and a plastic uten- sil bag, as well as numerous compartments and zipped pockets)	£121.00

		Order no.:	Kimba Cro Maximum load capacity of
Valid from 01.04.22 to 31.03.23 Please ensure youa re using the latest o	order form		base: 50 kg/110 lbs (includir Article number: HR3401=00
Provide estimate only	□Order		
Attachment for custom fabr Billing address/customer number		ed	
Delivery address			 Standard equipment Three-wheeled outdoor telescoping frame, for me neo seating unit, size 1 of custom seating shell Can be folded/disassemb
Sales person	Name Date Si	gnature	 Elastomer spring system Drum brake 18" rear wheels with pne Rigidly mounted 16" fron pneumatic tyre, includes Height-adjustable push b Frame colour: black
By placing the order, you agree to comp nuals for accessories and products of of with an Ottobock product, in addition to of the respective Ottobock products. A and accepted, our general terms and ccc placing an order, to conduct your own of cordance with national requirements if attachments that do not meet the spect service manuals.	ply with the instruct ther manufacturers to the instructions f Iso you confirm to 1 onditions. In additic conformity assessm and when you imp	tions for use/service ma- s when combining them for use/service manuals have read, understood on, you agree, upon hent or approval in ac- ilement modifications or	For children from approx. 1

DSS

the outdoor mobility ng seating unit) 000_K



- mobility base with ounting a Kimba r size 2, or a
- bled
- umatic tyres
- nt wheel with splash guard
- bar

to 10 years

Outdoor mobility base				
		Kimba Cross outdoor mobility base	£1686.00	
For use in conjunction with				
		Kimba seating unit, size 1		
		Kimba seating unit, size 2		
Seat adapter for Kimba se	ating	unit		
HR34011300		Size 1	Inclusive	
HR34021300		Size 2	Inclusive	
HR32x48300		Seating shell interface Kimba	£143.00	

Accessories for outdoor n	nobilit	ty base	
HR34012100		12" swivelling wheel with quick-lock attachment	£202.00
HR34017000		Anti-tipper, complete	£140.00
HR34018000		Clothing guard (plug-on wheel covers)	£202.00
HR34018100		Clothing guard extension around the wheel (not available in combination with spoke protector)	£80.00
HR34019000		Spoke protector, transparent	£100.00
HR34015000		Storage bag (e.g. for the swivelling front wheel)	£62.00
HR34016100		Bicycle draw-bar, complete (includes coupling from Weber company, type E, safety strap)*	£251.00
HR34016200		Equipment package required according to German Road Traffic Licensing Regulations (StVZO) (lighting/pennant/roll bar)	£141.00
HR16001102		Buggy board	£126.00
HR32xx2700-		Foot support assembly for seating shell interface Kimba	£330.00
* According to a leaflet for bi	cvcle t	- railers issued by the German Federal Ministry of Transport: (www.pdeleuw.de/fahrrad/fakten) – Recommen-	

* According to a leaflet for bicycle trailers issued by the German Federal Ministry of Transport: (www.pdeleuw.de/fahrrad/fakten) – Recommended overall weight of the bicycle trailer incl. load: 40 kg. Recommended in combination with four-point belt, groin strap, spoke protectors.

HR34016200 – Equipment package required according to German Road Traffic Licensing Regulations (StVZO)

According to a leaflet for bicycle trailers issued by the German Ministry of Transport: (www.pdeleuw.de/fahrrad/fakten)

Recommended overall weight of the bicycle trailer, including load: 40 kg. Recommended in combination with four-point belt, groin strap,

spoke protectors. Seating shell interface and footrest assembly in combination with Squiggles and Mygo not possible.

HR34018000 – Clothing guard is a plug-on wheel cover

HR34018100 - Clothing guard extension around the wheel cannot be used in combination with spoke protector

Accessories only for corresponding size Footrest assembly only available in combination with seating shell interface

		Order no.:	Kimba Inline Maximum load capacity of the outdoor mobility base: 50 kg/110 lbs (including seating units)
Valid from 01.04.22 to 31.03.23 Please ensure youa re using the latest o	order form		Article number: HR3255=0000_K
Provide estimate only	□Order		
Attachment for custom fabi		ed	
Delivery address			Standard equipment Kimba Inline twin outdoor mobility base for two Kimba seating units, size 1 • Pneumatic tyres provide suspension comfort • Rear wheels with drum brake • Frame colour: black
Sales person	Name		Twin buggy for children from approx. 1–6 years
By placing the order, you agree to comp nuals for accessories and products of or with an Ottobock product, in addition to of the respective Ottobock products. A and accepted, our general terms and cc placing an order, to conduct your own cordance with national requirements if attachments that do not meet the spec	ply with the instruc- ther manufacturer to the instructions Iso you confirm to onditions. In additi conformity assess and when you imp	s when combining them for use/service manuals have read, understood on, you agree, upon nent or approval in ac- plement modifications or	— * Notice: Seating unit, size 2, is available in custom fabrication upon request.

Please complete two Kimba Seat Unit forms when ordering this frame

service manuals.

Inline outdoor mobility base				
HR32551000-		Inline twin outdoor mobility base for up to two seating systems positioned one behind the other, with pneumatic front wheel 200 and pneumatic rear wheel 300, complete with drum brake	£1434.00	
Accessories				
HR16001102		Buggy board	£133.00	
HR32548300		Front seating shell interface	£128.00	
HR32548301		Rear seating shell interface	£128.00	
HR32582700		Footrest assembly for front seating shell interface	£299.00	
HR32582701		Footrest assembly for rear seating shell interface	£299.00	

	Order no.:	
		S
Valid from 01.04.22 to 31.03.23 Please ensure youa re using the latest orc	der form	
Provide estimate only	Order	
Billing address/customer number		
Delivery address		
Sales person	Name	•
	 Date Signature	•
By placing the order, you agree to comply nuals for accessories and products of othe with an Ottobock product, in addition to of the respective Ottobock products. Also and accepted, our general terms and com- placing an order, to conduct your own con cordance with national requirements if ar attachments that do not meet the specific service manuals.	er manufacturers when combining them the instructions for use/service manuals o you confirm to have read, understood ditions. In addition, you agree, upon nformity assessment or approval in ac- nd when you implement modifications or	

ottobock. Kimba Home

Max. load capacity (including seating unit): Size 1: 50 kg Size 2: 75 kg



- Indoor chassis for Kimba seating units size 1 and size 2
- Tilt adjustment of the seat unit via gas pressure spring 500 N (tilt adjustment range -5° till +35°)
- Height adjustment via gas pressure spring 1200 N (height of seat unit adapter min. 19 cm, max. 56 cm)
- 4 castors (75 mm), of these 2 with locking device
- Frame colour: Kimba white
- Push bar incl. release lever for seat angle adjustment

Basic model				
470G71=90001		Kimba Home indoor chassis size 1	£1032.00	
470G71=90002		Kimba Home indoor chassis size 2	£1032.00	
Accessories				
471A71=SK910		Anti-tip mechanism	£402.00	
471C71=SK730		Kimba seating shell interface, size 1	£192.00	
471C71=SK830		Kimba seating shell interface, size 2	£199.00	
HR32692700-000		Foot support for seating shell interface, complete, size 2	£192.00	
471A71=SK921		Gas spring (2 x 300 N) as support for seat lift function, recommended recommended from 30 kg load	£243.00	
471A71=SK922		Gas spring (2 x 500 N) as support for seat lift function, recommended recommended from 40 kg load	£243.00	

Information



Castor with locking device



Foot pedal for height adjustment



Push bar with release lever for seat angle adjustment



Seat angle adjustment



Seat heght adjustment



Adapter for Kimba seat unit



Anit-tip mechanism

ottobock. Kimba Indoor mobility base

		Order no.:	Article number: HR2233=0000_K
Valid from 01.04.22 to 31.03.23 Please ensure youa re using the latest orc	der form		
Provide estimate only	□ Order		
Billing address/customer number			000
Delivery address			 Height-adjustable steel constru Caster with swivel lock, diamet
Sales person	Name Date Sig	Inature	 Angle-adjustable seat interface Adjustable wheel base Frame colour: anthracite metal
By placing the order, you agree to comply nuals for accessories and products of othe with an Ottobock product, in addition to of the respective Ottobock products. Also and accepted, our general terms and com- placing an order, to conduct your own con cordance with national requirements if an attachments that do not meet the specific service manuals.	er manufacturers the instructions fo you confirm to h ditions. In addition nformity assessmend when you impl	when combining them or use/service manuals ave read, understood n, you agree, upon ent or approval in ac- ement modifications or	



• Height-adjustable steel construction

- Caster with swivel lock, diameter 75 mm
- Angle-adjustable seat interface
- Adjustable wheel base
- Frame colour: anthracite metallic

Standard model					
HR22330000-011		Size 1	£572.00		
HR22430000-011		Size 2	£572.00		
For use in conjunction wit	For use in conjunction with				
		Kimba Neo seating unit, size 1			
		Kimba Neo seating unit, size 2			
Caster with swivel lock					
HR22321150		Diameter 75 mm	Inclusive		
HR22321160		Diameter 125 mm	£175.00		
Seat interface					
HR32046x00		Parallel	£128.00		
HR32x48300		Kimba	£128.00		
Accessories					
HR22327100		Push handle	£210.00		
HR32044x00-		Footrest assembly for seating shell interface	£299.00		