

# Kimba Cross

The Kimba Cross with its three wheel mobility base is versatile, practical and also looks good with its aluminium oval tubes. On rugged terrain, the seat suspension absorbs the impact allowing the child to remain comfortable at all times.

The large pneumatic tires, drum brakes and swivelling front wheel allows the buggy to be easy to manoeuvre. The front wheel can be easily removed by the quick release mechanism, which then allows the bicycle bar to be attached and then Kimba Cross can then be towed by a push bike for more adventures.

## What is special about me?

- Swivelling wheel with quick lock attachment for confined spaces
- Cycle trailer package (draw bar with quick lock, cycle coupling, pennant, light)

- 16" front wheel, pneumatic tyre, rigid quick lock attachment
- 18" rear wheels, quick release axle
- Drum brake with lock option
- Height adjustable pushbar
- Can be folded
- Attractive design

- Removable, washable covers
- Adjustable seat depth, width and back height
- Lateral support pads, height adjustable
- Knee angle and height of footrest assembly can be adjusted

## Standard Model

### Outdoor Mobility Base

- Telescoping frame fits various sized seating units
- Spring system for the seat

### Kimba Multifunctional Seating Unit

- Available in 2 sizes, size 1: 1-6years and size 2: 4-10 years approx.
- Seating unit made of aluminium
- Powder coated

## Options

- Anti-tipper, clothing protectors
- Spoke protectors
- Storage bag
- Choice of material colours

Quality for life



## Technical Data

	Outdoor Mobility Base	Size 1	Size 2
Push bar height	90cm – 115cm		
Width	66cm		
Length	140cm – 170cm		
Height	75cm – 120cm		
Seat depth		19cm – 31cm	22cm – 40cm
Seat width		20cm – 30cm	26cm – 35cm
Lower leg length		19cm – 31cm	20cm – 37cm
Backrest height/angle		46cm – 66cm / 0° to + 90°	56cm – 70cm / 0° to +90°
Wheel diameter front/back	40cm / 45cm (swivelling wheel 30cm / 45cm)		
Seat tilt	-5° to +30°		
Load capacity	50cm	35cm	35cm
Weight	15.8kg (incl. splash guard)	7cm	8.8cm
Folding size (LxWxH) cm	94 x 60 x 32	72 x 40 x 22	82 x 44 x 27