

2026

ottobock.

Rehab folding buggy *Lisa*.

☐ Quotation

☐ Order

Billing Address

PO no.

Customer Name

Street

County PC

Shipping address (if different from billing address)

Customer Name

Street

County PC

Name _____ Date _____ Signature _____



Standard equipment Selection (unless alternative selected)

- Compact folding buggy
- Lightweight aluminium design
- Seat upholstery and reinforced back support upholstery
- Back support adjustment for one-hand operation
- Adjustable and divisible leg support
- Front and rear wheels with ball bearings
- Front wheels with integrated swivel lock
- Wheel lock
- Screw connections guided by sleeves
- Standard lap belt
- Pre-mounted tip-assist, right
- Upholstery colour blue
- Frame colour: silver

Standard tyres:

- EVA wheels front 180 mm
- EVA wheels rear 250 mm

Base Price -

Please email all orders to
Bockuk@ottobock.com

Valid from 01.03.2026 to 31.12.26
Please ensure you are using the latest order form3
v2

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals. Wheelchairs with CLT package cannot be returned or exchanged

The number codes contained in this order form are not intended to be used as article numbers for spare parts orders.
Please scan the QR code to access the spare parts catalogue for any spare parts orders.



Rehab folding buggy *Lisa*.

Standard model

☐	HR32130001-027	Lisa rehab folding buggy, size 1 with seat upholstery, seat depth presetting: 30 cm/12"	£1655.00
☐	HR32130002-027	Lisa rehab folding buggy, size 1 with seat upholstery, seat depth presetting: 35 cm/14"	£1655.00
☐	HR32140001-027	Lisa rehab folding buggy, size 2 with seat upholstery, seat depth presetting: 35 cm/14"	£1655.00
☐	HR32140002-027	Lisa rehab folding buggy, size 2 with seat upholstery, seat depth presetting: 40 cm/16"	£1655.00

Accessories

☐	HR32137100-027	Arm support, piece (please order two pieces)	£82.00
☐	HR16001102	Buggy board	£254.00
☐	HR32135300-343	Simplified padding, size 1 (coro-blue)	£150.00
☐	HR32145300-343	Simplified padding, size 2 (coro-blue)	£127.00
☐	HR32135200-343	Seat pad with multifunctional back support, size 1 (coro-blue)	£378.00
☐	HR32145200-343	Seat pad with multifunctional back support, size 2 (coro-blue)	£378.00
☐	HR32137700-031	Head support size 1, matching multifunctional back support size 1	£243.00
☐	HR32147700-031	Head support size 2, matching multifunctional back support size 2	£243.00
☐	HR32070060-031	Head support size 1, padded (not available in combination with seat pad with multifunctional back support)	£243.00
☐	HR32086500-031	Head support size 2, padded (not available in combination with seat pad with multifunctional back support)	£243.00
☐	HR32136100	Five-point belt with neck pad	£228.00
☐	HR32135500-031	Lateral guides with variable seat width reduction (4 – 6 – 8 cm), size 1 (only available in combination with pad for multifunctional back support)	£200.00
☐	HR32145500-031	Lateral guides with variable seat width reduction (4 – 6 – 8 cm), size 2 (only available in combination with pad for multifunctional back support)	£157.00
☐	HR32137200-031	Flip-down pommel size 1	£323.00
☐	HR32147200-031	Flip-down pommel size 2	£331.00
☐	HR32075800-031	Groin strap for multifunctional back support size 1	£203.00
☐	HR32085800-031	Groin strap for multifunctional back support size 2	£191.00
☐	HR32137800	Grab rail size 1	£130.00
☐	HR32147800	Grab rail size 2	£144.00
☐	HR32135700-031	Padding for grab rail size 1	£66.00
☐	HR32145700-031	Padding for grab rail size 2	£66.00
☐	HR32010080	Hook-and-loop foot straps with instep pad size 1	£40.00
☐	HR32020080	Hook-and-loop foot straps with instep pad size 2	£39.00
☐	HR32200120-343	Winter slip sack coro-blue with detachable front part, with slits for straps and accessories size 2 (seat depth up to approx. 34 cm)	£692.00
☐	HR32400120-343	Winter slip sack coro-blue with detachable front part, with slits for straps and accessories size 3 (seat depth up to approx. 39 cm)	£907.00
☐	HR32100100-344	Rain cape size 1, grey	£90.00
☐	HR32100100-344	Rain cape size 2, grey	£90.00
☐	481S64=ST005	Backpack with Ottobock logo	£75.00
☐	HR32135800	Folding storage bag, size 1	£110.00
☐	HR32145800	Folding storage bag, size 2	£110.00
☐	HR32136200-094	Canopy in upholstery colour with transparent rain cover	£260.00

Rehab folding buggy *Lisa*.

Technical data

	Size 1	Size 2
Target group	Children from approx. 3 – 8 years	Children from approx. 6 – 16 years
Seat width (with/without padding)	32/36 cm	37/41 cm
Seat depth	30 oder 35 cm	35 oder 40 cm
Back support height	65 cm	70 cm
Back support adjustment	30° (in steps 4 x 7,5°)	30° (in steps 4 x 7,5°)
Lower leg length	27 – 40 cm	27 – 40 cm
Max. load capacity	50 kg	60 kg
Weight	13,6 kg	14,6 kg