

2026

ottobock.

Adaptive wheelchair *Motus HEMI.*

☐ **Quotation**

☐ **Order**

Billing Address

PO no.	
Customer Name	
Street	
County PC	

Shipping address (if different from billing address)

Customer Name	
Street	
County PC	

Name _____ Date _____ Signature _____



Base Price - £1457.00

Please email all orders to
Bockuk@ottobock.com

Valid from 01.03.2026 to 31.12.26
Please ensure you are using the latest order form
v2

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.
Wheelchairs with CLT package cannot be returned or exchanged

Article number: 480F61=2_AA04_C

Standard equipment Selection (unless alternative selected)

- Aluminium tube frame with cross brace
- Swing-away or integrated leg support holders
- Individual foot plate, angle-adjustable or single-panel foot plate, angle-adjustable
- Back support tubes with different height setting options
- Back support upholstery, adaptable
- Aluminium side panels with cold protection
- Drive wheels with hollow rims and quick-release axles, 1" tyres
- Knee lever wheel lock

Foldable


Max. user weight: 125 kg/ 275 lbs or 140 kg/ 309 lbs
Crash-tested according to ISO 7176-19

The number codes contained in this order form are not intended to be used as article numbers for spare parts orders.
Please scan the QR code to access the spare parts catalogue for any spare parts orders.



Adaptive wheelchair *Motus HEMI*.

Filling aid – Wheelchair measurements




SW= cm

Seat Width
Outside edge of frame –
Outside edge of frame



SD= cm

Seat Depth
Front of seat upholstery –
Back support tube



ASH= cm

Anterior Seat Height
Ground –
Top of upholstery



PSH= cm

Posterior Seat height
Ground –
Top of seat upholstery



BH= cm

Back Support Height
Top of seat tube –
Top of back support upholstery



LL= cm

Lower Leg Length
Top of seat upholstery –
Top of foot plate (middle)

Wheelchair Type

- ☐ Standard Inclusive
- ☐ MB50 **Foot-propelled wheelchair** (reduced clearance*) Inclusive

* Risk due to getting caught on obstacles. Please observe safety information in the instructions for use.

Foot Space Enlargement

- ☐ MA63 ☐ Yes £88.00
- ☐ No

Maximum Payload

- ☐ up to 125 kg / 276 lb
- ☐ up to 140 kg / 309 lb

Seat Width [SW] in cm

- ☐ 35.5 ☐ 38 ☐ 40.5 ☐ 43 ☐ 45.5 ☐ 48 ☐ 50.5 ☐ 53 ☐ 55.5

Seat Depth (SD) in cm

- ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48 ☐ 50 ☐ 52 ☐ 54

Anterior Seat Height [ASH] in cm (see table in annex)

22" drive wheel

- ☐ 38 ☐ 39 ☐ 40*

24" drive wheel

- ☐ 41 ☐ 42 ☐ 43 ☐ 44 ☐ 45 ☐ 46 ☐ 47 ☐ 48 ☐ 49 ☐ 50 ☐ 51
- ☐ 52 ☐ 53 ☐ 54 ☐ 55

* if PSH=40, then 24"

Lower Leg Length [LL]

LL= cm Top of seat upholstery – Top of foot plate (center)

Adaptive wheelchair *Motus HEMI.*

Posterior Seat Height [PSH] in cm

22" drive wheel

☐ 36 ☐ 37 ☐ 38 ☐ 39

24" drive wheel

☐ 40 ☐ 41 ☐ 42 ☐ 43 ☐ 44 ☐ 45 ☐ 46 ☐ 47 ☐ 48 ☐ 49 ☐ 50
☐ 51 ☐ 52

Back Height [BH] in cm

☐ 30 ☐ 32.5 ☐ 35 ☐ 37.5 ☐ 40 ☐ 42.5 ☐ 43 ☐ 45 ☐ 47.5 ☐ 48 ☐ 50

Frame Accessories

☐ MA51	Tip assist, right	£49.00
☐ MA52	Tip assist, left	£49.00
☐ MA55	Anti-tipper, swing-away, right	£99.00
☐ MA56	Anti-tipper, swing-away, left	£99.00
☐ MA57	Crutch holder with tip assist, right	£96.00
☐ MA58	Crutch holder with tip assist, left	£96.00

Caster Wheels

☐ MF30	4" Solid rubber, black	Inclusive
☐ MF31	5" Solid rubber, black	Inclusive
☐ MF04	5.5" Soft rubber	Inclusive
☐ MF32	6" Solid rubber, black	Inclusive
☐ MF08	7" Solid rubber, black	Inclusive

Horizontal Drive Wheel Position (see table in annex)

☐ A ☐ B ☐ C ☐ D

Drive Wheel Camber

☐ 0° ☐ 2° ☐ 3° ☐ 4°

Drive Wheel Rim

☐ MG06	Hollow rims	Inclusive
☐ MG15	Infinity Ultralight	£143.00
☐ MG80	With drum brake for attendant *	£399.00
☐ MG83	One-arm drive (double push ring left) *	£637.00
☐ MG83	One-arm drive (double push ring right) *	£637.00

* no combination with drive wheel camber

Tires

☐ MG31	Schwalbe RightRun 1"	Inclusive
☐ MG30	Pneumatic tire 1 3/8"	Inclusive
☐ MG34	PU 1 3/8" *	Inclusive
☐ MG341	PU 1", black *	£66.00
☐ MG37	Schwalbe Marathon Plus Evolution 1"	£76.00
☐ MG200	Schwalbe One 1"	£190.00

* no combination with MG15

Adaptive wheelchair *Motus HEMI*.

Handrims

☐	MG52	Aluminium, silver anodized, high	Inclusive
☐	MG53	Aluminium, silver anodized, low	£74.00
☐	MG54	Rubber-coated, high	£110.00
☐	MG56	Stainless steel, polished, high	£57.00
☐	MG65	Gekko (ergonomic shape)	£594.00

Wheelbase Extension

☐	MG82	With second fitting (incl. wheel lock adapter)	£362.00
-----------------------	------	--	---------

Rear Wheel Accessories

☐	MG87	High-pressure air pump, 14 bar	£135.00
☐	MG88	Airman-pump (battery-operated)	£258.00

Wheel Locks

☐	MH01	Push-to-lock	Inclusive
☐	MH11	Pull-to-lock *	£28.00
☐	MH05	Scissor wheel lock *	£206.00
☐	MH06	For user and attendant	£413.00
☐	MH08	One-hand activated, right	£147.00
☐	MH09	One-hand activated, left	£147.00

* to be preferred with short SD / DV frame

Wheel Lock Accessories

☐	MH10	Wheel lock lever extension	£40.00
-----------------------	------	----------------------------	--------

Back Support System

☐	MD01	Back support, straight	Inclusive
☐	MD09	Back support with lumbar curve 8°	Inclusive
☐	MD10	Back support, angle adjustment 30° *	£375.00
☐	MD11	Foldable	£198.00

*combination with BH= 43 or 48 only / max. cap. 125kg/ 275lb

Back Design

☐	MD03	Adaptable, padded	Inclusive
☐	MD05	Adaptable, breathable	£109.00
☐	MD93	Baxx Back (Aluminium)	£716.00
☐	MD04	Standard *	Inclusive
☐	MD06	Without	Inclusive

* chose with MD10

Widening Back Pad

☐	MD80	+ 2 cm	£17.00
☐	MD81	+ 4 cm	£17.00
☐	MD82	Without	Inclusive

Height for Baxx Back

☐	MD91	Low (10"/25 cm)	Inclusive
☐	MD94	Medium (13"/33 cm)	Inclusive
☐	MD95	Standard (16"/41 cm)	Inclusive

Adaptive wheelchair *Motus HEMI.*

Back Contour for Baxx Back

☐	MD84	Active (2.5"/6 cm)	Inclusive
☐	MD85	Deep (4"/10 cm)	Inclusive

Push Handles

☐	MD50	Short	Inclusive
☐	MD51	Long	Inclusive
☐	MD53	Telescoping, short	£99.00
☐	MD52	Telescoping, long	£99.00
☐	MD55	Height adjustable/ detachable	£150.00
☐	MD45	Height adjustable/ detachable for Baxx	£157.00
☐	MD56	Folding *	£93.00

* max. cap. 125kg/ 275lb

Back Height Range

☐	30 - 37.5 cm (if MD01 / +MD50/51/56)	Inclusive
☐	37.5 - 40 cm (if MD01 / +MD50/51/56)	Inclusive
☐	37.5 - 45 cm (if MD01 / +MD50/51/56)	Inclusive
☐	40 - 45 cm (if MD01 / +MD50/51/56)	Inclusive
☐	30 - 32.5 cm (if MD01 / +MD53/52/55)	Inclusive
☐	32.5 - 40 cm (if MD01 / +MD53/52/55)	Inclusive
☐	42.5 - 50 cm (if MD01 / +MD53/52/55)	Inclusive
☐	35 - 40 cm (if MD09 / +MD53/52/55)	Inclusive
☐	45 - 50 cm (if MD09)	Inclusive
☐	47.5 - 50 cm (if MD09 / +MD50/51/56 or MD11/ + MD51)	Inclusive

Back Support Accessories

☐	MD58	Stabiliser Bar	£119.00
-----------------------	------	----------------	---------

Thoracic Supports

☐	MD23	Lateral Trunk Support STS, medium, right	£132.00
☐	MD24	Lateral Trunk Support STS, medium, left	£132.00
☐	MD62	Lateral Trunk Support STS, large, left	£166.00
☐	MD63	Lateral Trunk Support STS, large, left	£166.00

Side Panels

☐	ME01	Aluminium with cold protection	Inclusive
☐	ME02	Aluminium with splash- and cold protection	Inclusive
☐	ME11	Plug-on side panel, plastic	Inclusive
☐	ME04	Carbon with splash protection	£472.00
☐	ME05	Desk side panel with arm support, short	£293.00
☐	ME06	Desk side panel with arm support, long	£300.00
☐	ME17	Desk side panel with arm support, depth adjustable	£340.00

Side Panel Accessories

☐	ME50	Arm support padded, swing away *	£290.00
		☐ short ☐ long ☐ long, extra high	
☐	ME60	Tray	£227.00
☐	ME67	Half-sided folding tray, left	£372.00
☐	ME68	Half-sided folding tray, right	£372.00
☐	ME61	Swivel unit for channel forearm support, right	£251.00
☐	ME62	Swivel unit for channel forearm support, left	£251.00

* no combination with ME05/06/17

Adaptive wheelchair *Motus HEMI.*

Channel Forearm Support

☐	ME70	One-piece channel forearm support, left	£71.00
☐	ME71	One-piece channel forearm support, right	£71.00

Seat Upholstery

☐	MC13	Standard	Inclusive
☐	MC14	Adaptable	£106.00
☐	MC06	Without	Inclusive

Seat Cushions

☐	MC01	Standard (functional height 5 cm)	£43.00
☐	MC02	Cubic Foam (functional height 5 cm)	£151.00
☐	MC05	Evolight (functional height 6 cm)	£215.00
☐	MC20	Terra (functional height 6 cm)	£200.00
☐	MC21	Terra Aquos (functional height 6 cm)	£241.00
☐	MC22	Terra Flair (functional height 6 cm)	£603.00
☐	MC23	Cloud (functional height 11 cm)	£589.00
☐	MC29	Livo (functional height 5 cm)	£90.00
☐	MC30	EvoFlair (functional height 8 cm)	£323.00

Leg Support

☐	MB10	Plastic, Individual, angle adjustable	Inclusive
☐	MB03	Aluminium, single panel, angle-adjustable	Inclusive
☐	MB09	Aluminium, individual, angle adjustable	£119.00
☐	MB20	Single panel, short lower leg length	£184.00
☐	MB21	Amputation leg support, right *	£223.00
☐	MB22	Amputation leg support, left *	£223.00
☐	MB23	Amputation leg support, pair	£443.00
☐	MB85	Stealth Amputation leg support, Left	£492.00
☐	MB86	Stealth Amputation leg support, Right	£492.00
☐	MB87	Stealth Amputation leg support, Pair	£983.00
☐	MB26	Elevating leg support, right *	£151.00
☐	MB27	Elevating leg support, left *	£151.00
☐	MB25	Elevating leg support, pair	£304.00
☐	MB99	Without	Inclusive

* MB10 on opposite side

Version Foot Plate Bar

☐	MB200	Standard	Inclusive
☐	MB201	Swing away, left side	£82.00
☐	MB202	Swing away, right side	£82.00

Length Foot Plate Bar without seat cushion (see table in annex)

☐	MB40	Short	Inclusive
☐	MB41	Medium	Inclusive
☐	MB42	Long	Inclusive

Adaptive wheelchair *Motus HEMI.*

Other Accessories

☐	MZ51	Lap belt with metal buckle	£77.00
☐	MZ52	Tool kit	£57.00

Spoke Protector

☐	MG90	Transparent	£136.00
-----------------------	------	-------------	---------

Color: Frame

☐	002	Anthracite		Inclusive
☐	016	Silver metallic		£109.00
☐	023	White RAL 9016		£109.00
☐	048	Obsidian Black (texture)		£109.00
☐	051	Sunny Orange (texture)		Inclusive
☐	058	Mint Green		Inclusive
☐	059	Raspberry Red		Inclusive
☐	060	Azure Blue		Inclusive
☐	056	Pearl White		£59.00

Color: Side Panel

☐		Obsidian Black (texture)	
☐	ME20	Frame color *	£55.00

* no combination with ME04/11

Color: Baxx Back

☐	MD89	White RAL 9016	Inclusive
☐	MD90	Obsidian Black (texture)	Inclusive
☐	MD88	Frame color	Inclusive

Special Accessories

☐	MZ58	Bag for rear wheels (2 wheels / max.: 28")	£116.00
☐	MZ59	Pocket for Mobile phone	£29.00
☐	MZ89	Protector for quick-release-axle	£11.00
☐	MZ92	WheelBlades (ski blades for caster wheels)	£477.00
☐	MZ93	Backpack	£66.00

Mounting Set Head Support

☐	MZ54	Mounting set	£128.00
☐	MZ61	Clinical mounting set	£128.00

Head/Neck Supports for UniLink

☐	MZ105	Comfort Plus, size M	£480.00
☐	MZ106	Combo (3-pad), size L	£643.00

UniLink Accessories


☐	MZ107	UniLink extension 38 mm	£68.00
☐	MZ108	UniLink extension 50 mm	£66.00
☐	MZ109	UniLink extension 76 mm	£66.00

Adaptive wheelchair *Motus HEMI*.

Anterior Seat Height [ASH] → Wheel Diameter / Fork Length in cm

Caster fork	Wheel diameter				
	4"	5"	5.5"	6"	7"
short	38 - 46	40 - 46	42 - 47	42 - 47	44 - 49
long	46 - 53	47 - 54	44 - 55	44 - 55	45 - 55

Drive Wheel Position → Distance from tipping point in mm

Dimension (horizontal) from centre axis to apex point Back support tube (front edge)/upholstery bar (top edge)		
Drive Wheel Position	Distance from tipping point	
D (active)	83 mm	
C	63 mm	
B	43 mm	
A	23 mm	

Length Foot Plate Bar in mm

	Motus CV			Motus CS
	MB03/09	MB10	MB25/26/27	MB03/09
short	34 - 41	34 - 38	34 - 39	35 - 47
medium	37 - 47	38 - 50	34 - 50	45 - 55
long	42 - 53	45 - 55	42 - 55	52 - 55

Adaptive wheelchair *Motus HEMI*.

Information




MA51/52
Tip-assist



MA55/56
Anti-tipper



MA57/58
Crutch holder




MF04
Soft rubber caster wheel



MG65
Gekko push ring



MH01
Push-to-lock wheel lock



MH05
Scissor wheel lock



MH08/09
One-hand wheel lock



MD10
Back support, angle
adjustment 30°



MD11
Foldable back



MD05
Breathable back support



MD93
Baxx Back



MD55
Push handles, height
adjustable/ detachable



MD56
Push handles, folding



MD58
Stabiliser Bar



MD23/24/62/63
Lateral Trunk Support STS

Adaptive wheelchair *Motus HEMI*.

Information



ME02

Aluminium with splash- and cold protection



ME04

Carbon with splash protection



ME05/06/17

Desk side panel with arm support




ME11

Plug-on side panel plastic support




ME50

Arm support padded, swing away



ME60

Tray




ME67/68

Half-sided folding tray




ME70/71

One-piece channel forearm support




MC20

Terra



MC21

Terra Aquos




MC22

Terra Flair



MC23

Cloud



MB20

Leg support for short LL




MB21/22/23

Amputation leg support



MB85/86/87

Stealth Amputation leg support



MB26/27/25

Leg support, elevating

Adaptive wheelchair *Motus HEMI.*

Information



MZ58

Bag for rear wheels



MZ59

Pocket for Mobile phone



MZ89

Protector for quick-release-axle




MZ93

Backpack



MZ54

Mounting set head support



MZ61

Clinical mounting set head support



MZ105

Comfort Plus head support



MZ106

Combo head support