

2026

StartM6 Junior

ottobock.

☐ Quotation

☐ Order

Billing Address

PO no.

Customer Name

Street

County PC

Shipping address (if different from billing address)

Customer Name

Street

County PC

Name Date Signature



Base Price - £887.00

Please email all orders to
Bockuk@ottobock.com

Valid from 01.03.2026 to 31.12.26
Please ensure you are using the latest order form

v2

Article number: 480F53=60000_K

Standard equipment Selection (unless alternative selected)

- Aluminium tube frame with crossbrace
- Tip-assist integrated into the frame
- Individual leg supports, angle-adjustable, plastic foot plates
- Swing-away leg support holders
- Side panel with arm support, height-adjustable
- Nylon seat and back support upholstery, black
- Modular design
- Drive wheels with quick-release axles
- Knee lever wheel lock

Folding, max. user weight 90 kg/198 lbs

Crash-tested up to 75 kg according to ISO 7176-19

The number codes contained in this order form are not intended to be used as article numbers for spare parts orders.
Please scan the QR code to access the spare parts catalogue for any spare parts orders.



StartM6 Junior

Wheelchair dimensions

Seat width (SW)	cm	Anterior seat height (ASH)	cm
Seat depth (SD)	cm	Posterior seat height (PSH)	cm
Back support height (BH)	cm	Lower leg length (LL) (frame measurement)	cm

Seat measurements in cm

Seat width (SW)	☐ 28	☐ 30.5	☐ 33	☐ 35.5	☐ 38
Seat depth (SD)	☐ 31	☐ 34	☐ 37		

Seat height measurements in cm

Anterior seat height (ASH)	☐ 40	☐ 41	☐ 42	☐ 43	☐ 44				
	☐ 45	☐ 46	☐ 47	☐ 48	☐ 49				
Posterior seat height	☐ 37	☐ 38	☐ 39	☐ 40	☐ 41	☐ 42	☐ 43	☐ 44	
	☐ 45	☐ 46	☐ 47	☐ 48	☐ 49				

Seat height table

Anterior seat height in cm

Caster fork, short

Caster wheel diameter 5"
Caster wheel diameter 6"


Adjustment range

40 – 44
43 – 45

Caster fork, long

Caster wheel diameter 5"
Caster wheel diameter 6"


44 – 47
45 – 49



Seat height table

Posterior seat height in cm

Drive wheel size	Adjustment range
20"	37 – 43
22"	38 – 45
24"	41 – 49



The values indicated are rounded values which have been theoretically determined. Max. deviation of 1 cm.

The anterior seat height depends on the selected wheel size, caster fork and mounting position. Please refer to the seat height table!


Specified without seat cushion at 0° seat inclination.

The maximum possible seat inclination is 3.5 cm.

StartM6 Junior.

Frame accessories

☐	MA 55	Plug-on anti-tipper, right	£27.00
☐	MA 56	Plug-on anti-tipper, left	£27.00
☐	MA 57	Crutch holder, right	£32.00
☐	MA 58	Crutch holder, left	£32.00
☐	MA 61	Swinging anti-tipper, right	£53.00
☐	MA 62	Swinging anti-tipper, left	£53.00



MA 55 / MA 56

Caster wheels

☐	MF 03	5" solid rubber	Inclusive
☐	MF 05	6" pneumatic (150 x 30)	Inclusive
☐	MF 06	6" solid rubber (150 x 27)	Inclusive

Caster forks

☐	MF 80	Caster fork, long	Inclusive
☐	MF 81	Caster fork, short	Inclusive

Drive wheel size

☐	MG 04	20"	Inclusive
☐	MG 01	22"	Inclusive
☐	MG 02	24"	Inclusive

Wheelbase

☐	short/active	☐	long/passive	Inclusive
-----------------------	--------------	-----------------------	--------------	-----------

Drive wheel rim

☐	MG 05	Standard	Inclusive
☐	MG 80	Drive wheel with drum brake/attendant	£397.00

Tyres

☐	MG 30	Pneumatic tyres 1 3/8"	Inclusive
☐	MG 35	PU tyres 1 3/8"	Inclusive

Handrims

☐	MG 52	Handrim, aluminium, silver anodised, high	Inclusive
☐	MG 56	Handrim, polished stainless steel, high	£62.00

Drive wheel accessories

☐	MG 87	High-pressure air pump, 14 bar	£121.00
☐	MG 88	Airman pump, battery-operated	£185.00

StartM6 Junior.

Wheel lock

☐	MH 01	Knee lever wheel lock, standard	Inclusive
☐	MH 06	Knee lever wheel lock for user and attendant	£258.00

Wheel lock accessories

☐	MH 10	Wheel lock lever extension	£32.00
-----------------------	-------	----------------------------	--------

Back support system

☐	MD 09	Back support with lumbar curve, standard	Inclusive
☐	MD 01	Straight back support, standard	Inclusive
☐	MD 10	Back support angle adjustment 30° (only available with BH 37.5 cm)	£250.00
☐	MD 11	Back support, folding (only available with BH 37.5 cm)	£98.00



MD 10



MD 11

Back support height in cm

☐ 30
 ☐ 32.5
 ☐ 35
 ☐ 37.5

Back support upholstery

☐	MD 04	Standard	Inclusive
☐	MD 05	Adaptable	£116.00
☐	MD 06	Without	Inclusive

Push handles

☐	MD 51	Standard	Inclusive
☐	MD 52	Telescoping (only available in combination with 30° back support angle adjustment)	£215.00
☐	MD 54	Height-adjustable, angled (only available in combination with MD01 straight back support, standard)	£124.00

Back support accessories

☐	MD 23	Side thoracic support, right	£102.00
☐	MD 24	Side thoracic support, left	£102.00
☐	MD 58	Stabiliser bar for push handles	£51.00

Side panels

☐	ME 01	Side panel, standard, 300 mm arm support, short mounted	Inclusive
☐	ME 01	Side panel, standard, 300 mm arm support, long mounted	Inclusive
☐	ME 05	Side panel, height-adjustable, 300 mm armrest	£142.00


StartM6 Junior.



MD 54



MD 23 / MD 24



ME 05

Side panel accessories

☐	ME 60	Tray	£200.00
-----------------------	-------	------	---------

Seat upholstery

☐	MC 13	Upholstery, standard	Inclusive
☐	MC 06	Without	Inclusive

Seat cushion

☐	MC 01	Seat cushion Standard (functional height 5 cm)	£43.00
☐	MC 02	Cubic Foam seat cushion (functional height 5 cm)	£151.00
☐	MC 05	Evolight seat cushion (functional height 6 cm)	£215.00

Lower leg length in cm (without seat cushion)

Lower leg length (LL) in cm (available adjustment range 16 cm – 45 cm)

Leg support bar options

☐	MB 40	Short (27 – 40 cm)	Inclusive
☐	MB 41	Medium (32 – 45 cm)	Inclusive

Please note that the lower leg lengths vary in case of elevating leg supports.

Leg supports

☐	MB 10	Individual leg supports, angle-adjustable, plastic foot plates	Inclusive
☐	MB 25	Leg supports, elevating, without automatic length adjustment, pair	£258.00
☐	MB 26	Leg support, elevating, without automatic length adjustment, right (plus standard leg support, left)	£130.00
☐	MB 27	Leg support, elevating, without automatic length adjustment, left (plus standard leg support, right)	£130.00
☐	MB 20	Single-panel leg support for short lower leg length	£220.00
☐	MB 19	Individual leg supports for short lower leg length	£81.00

Leg support accessories

☐	MB 72	Heel strap for individual leg supports	£33.00
-----------------------	-------	--	--------

StartM6 Junior.



MB 20



MB 72


Other accessories

☐	MZ 50	Lap belt with plastic buckle	£49.00
☐	MZ 52	Tool kit	£65.00
☐	MZ 53	Head support	£124.00
☐	MZ 54	Clinical installation kit for head support	£126.00
☐	MZ 55	Head support installation kit, multi-axis, offset	£111.00
☐	MZ 70	Anchor point kit for transportation in a motor vehicle (max. loading capacity: 75 kg)	£388.00


Note: Please also use the order forms for wheelchair accessories 647F661.



MZ 53




MZ 54



MZ 55

Spoke protector

☐	MG 90	Spoke protector, transparent	£155.00
☐	MG 95	Spoke protector with decor	£248.00



MG 95

StartM6 Junior.

Frame colour

☐	009	Vinho sparkle*	£191.00
☐	017	Traffic yellow RAL 1023*	Inclusive
☐	018	Signal red RAL 3001*	Inclusive
☐	021	Sparkle granny-smith*	£191.00
☐	003	Jet black**	Inclusive
☐	051	Sunny orange (texture)**	£70.00
☐	055	Blueberry (texture)**	Inclusive
☐	058	Mint green**	£70.00
☐	059	Raspberry red**	Inclusive

* available until 06/26

** available by 07/26