

2026

ottobock.

Rigid frame adaptive wheelchair *Ventus*.

☐ Quotation

☐ Order

Billing Address

PO no.

Customer Name

Street

County PC

Shipping address (if different from billing address)

Customer Name

Street

County PC

Name Date Signature



Base Price - £1710.00

Please email all orders to
Bockuk@ottobock.com

Valid from 01.03.2026 to 31.12.26
Please ensure you are using the latest order form
v2

Standard equipment Selection (unless alternative selected)

Article number: 480A25=00000_K

- Aluminium tube frame
- Rigid axle
- Front frame layout fully configurable
- Foot plate with different angle and height setting options
- Back support folding and angle adjustable
- Aluminium side panels
- Wheel camber selectable in 3° increments
- Three different drive wheel diameters
- Knee lever wheel lock
- Push handles by choice

rigid, max. user weight 125 kg

Crash tested according to ISO 7176-19


By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals. Wheelchairs with CLT package cannot be returned or exchanged

The number codes contained in this order form are not intended to be used as article numbers for spare parts orders.
Please scan the QR code to access the spare parts catalogue for any spare parts orders.



Rigid frame adaptive wheelchair Ventus.

Filling aid – Wheelchair measurements




SW= cm

Seat Width
Outside edge of frame –
Outside edge of frame



SD= cm

Seat Depth
Front of seat upholstery –
Back support tube




ASH= cm

Anterior Seat Height
Ground –
Top of upholstery



PSH= cm

Posterior Seat height
Ground –
Top of seat upholstery



BH= cm

Back Support Height
Top of seat tube –
Top of back support upholstery



LL= cm

Lower Leg Length
Top of seat upholstery –
Top of foot plate (middle)

Seat width (SB) in cm

☐ 28 ☐ 30 ☐ 32 ☐ 34 ☐ 36 ☐ 38
☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48 ☐ 50

Seat depth (SD) in cm

☐ 30 ☐ 32 ☐ 34 ☐ 36 ☐ 38
☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48 ☐ 50

Anterior seat height (ASH) in cm

☐ 42 ☐ 43 ☐ 44 ☐ 45 ☐ 46 ☐ 47 ☐ 48 ☐ 49
☐ 50 ☐ 51 ☐ 52 ☐ 53 ☐ 54

Posterior seat height (PSH) in cm

☐ 33 ☐ 34 ☐ 35 ☐ 36 ☐ 37 ☐ 38 ☐ 39 ☐ 40 ☐ 41 ☐ 42 ☐ 43
☐ 44 ☐ 45 ☐ 46 ☐ 47 ☐ 48 ☐ 49 ☐ 50 ☐ 51 ☐ 52 ☐ 53

Front frame angle

☐ MA20 75° Inclusive
☐ MA21 85° Inclusive

Abduction/offset

☐ MA06 4 cm abduction (up to seat width 36 cm) Inclusive
☐ MA14 0 cm offset Inclusive
☐ MA15 2 cm offset (from seat width 32 cm) Inclusive
☐ MA16 4 cm offset (from seat width 36 cm) Inclusive

Rigid frame adaptive wheelchair *Ventus*.

Frame accessories

☐ MA51	Tip-assist right	£106.00
☐ MA52	Tip-assist left	£106.00
☐ MA55	Anti-tipper right	£130.00
☐ MA56	Anti-tipper left	£130.00
☐ MA60	Frame pad	£87.00

Caster wheels

☐ MF30	4" solid rubber, black	Inclusive
☐ MF35	4" soft caster, aluminium	£166.00
☐ MF31	5" solid rubber, black	Inclusive
☐ MF36	5" soft caster, aluminium	£166.00
☐ MF04	5,5" soft caster, grey	Inclusive
☐ MF05	6" pneumatic	Inclusive
☐ MF32	6" solid rubber, black	Inclusive
☐ MF37	6" soft caster, aluminium	£166.00

Caster forks

☐ MF90	Standard	Inclusive
☐ MF94	Design	£183.00
☐ MX30	Froglegs, with suspension	£468.00

Fork length (see p.11, A Anterior Seat Height (ASH) – Fork length)

☐ MF80	Long	Inclusive
☐ MF81	Short	Inclusive
☐ MF82	Extra short	Inclusive

Wheelbase broadening

☐ MF99	10 mm	£43.00
☐ MF99	15 mm	£43.00

Drive wheel size (see p.11, B Posterior seat Height (PSH) – Drive wheel / slider)

☐ MG01	22"	Inclusive
☐ MG02	24"	Inclusive
☐ MG00	25"	Inclusive

Axle with camber (see p 11, C Wheel camber – overall width)

☐ MG40	0°	Inclusive
☐ MG41	3°	Inclusive
☐ MG42	6°	Inclusive
☐ MG43	9°	Inclusive

Rigid frame adaptive wheelchair *Ventus*.

Drive wheel rim		22"	24"	25"	
☐	MG06 Hollow rim	☒	☒	☒	Inclusive
☐	MG10 Spinergy SPOX		☒		£1037.00
☐	MG11 Spinergy LX		☒		£1335.00
☐	MG15 Infinity Ultralight	☒	☒	☒	£121.00
☐	MG197 SmootherWheel dampening drive wheels, max. load 75 kg		☒		£1659.00
☐	MG198 SmootherWheel dampening drive wheels, max. load 95 kg		☒		£1659.00
☐	MG199 SmootherWheel dampening drive wheels, max. load 115 kg		☒		£1659.00
Tyres		22"	24"	25"	
☐	MG31 Schwalbe RightRun, 1"	☒	☒	☒	Inclusive
☐	MG34 PU 1 "		☒		Inclusive
☐	MG37 Schwalbe Marathon Plus Evo. 1 "	☒	☒	☒	£76.00
☐	MG341 PU 1 ", black		☒		inclusive
Drive wheel attachment device / Slider (see p.11, B Posterior seat height (PSH) – Drive wheel / slider)					
☐	MG75 Short				Inclusive
☐	MG76 Medium				Inclusive
☐	MG77 Long				Inclusive
Handrims		22"	24"	25"	
☐	MG52 Aluminium, silver anodised, high	☒	☒	☒	Inclusive
☐	MG53 Aluminium, silver anodised, low	☒	☒		£16.00
☐	MG54 Rubber-coated, high	☒	☒	☒	£166.00
☐	MG56 Stainless steel, polished, high	☒	☒	☒	£64.00
☐	MG57 Stainless steel, polished, low	☒	☒		£64.00
☐	MG63 Curve L		☒	☒	£548.00
☐	MG64 Curve L Grip		☒	☒	£604.00
☐	MG65 Gekko		☒	☒	£594.00
Handrim mounting					
☐	MG70 Mounting position narrow: 15 mm				Inclusive
☐	MG71 Mounting position wide: 25 mm				Inclusive
Drive wheel accessories					
☐	MG87 High-pressure air pump, 14 bar				£125.00
Wheelbase extension					
☐	MG85 With second axle (+80 mm)				£473.00
☐	MG86 Single (+40 mm, from Drive wheel position 6 possible)				£310.00
☐	MG46 Without				Inclusive
Drive wheel position (see p.11, D Drive wheel position – Tilt point distance)					
☐	1 (active)	☐	☐	☐	☐
☐	2	☐	☐	☐	☐
☐	3	☐	☐	☐	☐
☐	4	☐	☐	☐	☐
☐	5	☐	☐	☐	☐
☐	6	☐	☐	☐	☐
☐	7	☐	☐	☐	☐
☐	8	☐	☐	☐	☐
☐	9	☐	☐	☐	☐

Rigid frame adaptive wheelchair *Ventus*.

Wheel locks

☐ MH01	Knee lever wheel lock Standard	Inclusive
☐ MH13	Knee lever wheel lock, enhanced	£396.00
☐ MH05	Scissor wheel lock	£207.00

Wheel lock accessories

☐ MH10	Wheel lock lever extension	£43.00
----------------------------	----------------------------	--------

Back support height (BH) in cm

☐ 22,5	☐ 25	☐ 27,5	☐ 30	☐ 32,5	☐ 35	☐ 37,5	☐ 40	☐ 42,5	☐ 45	☐ 47,5	☐ 50
----------------------------	--------------------------	----------------------------	--------------------------	----------------------------	--------------------------	----------------------------	--------------------------	----------------------------	--------------------------	----------------------------	--------------------------

Back support system

☐ MD10	Back support tube straight, folding onto the seat surface	Inclusive
☐ MD12	Back support tube straight, fixed	£299.00

Back support upholstery

☐ MD03	Adaptable, padded	Inclusive
☐ MD04	Standard	Inclusive
☐ MD25	Adaptable/height-adjustable, including push handles, long (adjustment range = 10 cm; ≥ BH 27.5 cm)	£254.00
☐ MD26	Adaptable/height-adjustable, including push handles, short (adjustment range = 10 cm; ≥ BH 27.5 cm)	£264.00
☐ MD40	Fixed back support plate (≥ BH 30 cm)	£107.00

Push handles

☐ MD50	Standard short (≥ BH 25 cm)	Inclusive
☐ MD51	Standard long (≥ BH 25 cm)	Inclusive
☐ MD53	Telescoping short (≥ BH 25 cm)	£178.00
☐ MD52	Telescoping long (≥ BH 25 cm)	£178.00
☐ MD55	Height-adjustable / removable (≥ BH 30 cm)	£215.00
☐ MD56	Folding (≥ BH 30 cm)	£152.00
☐ MD59	Without	Inclusive

Side panels

☐ ME01	Aluminium	Inclusive
☐ ME02	Clothing guard, aluminium	£176.00
☐ ME03	Plastic clothing guard with different setting options	£166.00
☐ ME04	Carbon clothing guard with different setting options	£658.00
☐ ME05	with forearm support, height-adjustable (see p.12, E ATR / (hSH), exclusions ME05)	£288.00

Side panel accessories

☐ ME50	Arm support, padded, swing-away, different height setting options	£240.00
☐ ME51	Neoprene insert for clothing guard	£28.00

Seat upholstery

☐ MC13	Standard	Inclusive
☐ MC14	Adaptable	£119.00
☐ MC19	Fixed seat plate	£199.00

Rigid frame adaptive wheelchair *Ventus*.

Seat cushion

☐ MC01	Standard (functional height 5 cm)	£43.00
☐ MC02	Cubic Foam (functional height 5 cm)	£151.00
☐ MC05	Evolight (functional height 6 cm)	£215.00
☐ MC20	Terra (functional height 6 cm)	£200.00
☐ MC21	Terra Aquos (functional height 6 cm)	£241.00
☐ MC22	Terra Flair (functional height 6 cm)	£603.00
☐ MC23	Cloud (functional height 11 cm)	£389.00
☐ MC29	Livo (functional height 5 cm)	£88.00
☐ MC30	EvoFlair (functional height 8 cm)	£323.00

Lower leg length

LL  in cm	Upper seat cover – Upper foot plate (middle)
--	--

Leg supports

☐ MB03	Footrest shop adjustable (with footboard)	Inclusive
☐ MB06	Foot support, light, fixed	£124.00
☐ MB20	Leg support for short LL (20 – 39 cm)	£189.00

Leg support bar options (see p. 12, **F** Front frame offset / LL – foot stirrup length)

☐ MB40	Short (39 cm – 45 cm)	Inclusive
☐ MB42	Long (46 cm – 50 cm)	Inclusive

Leg support accessories

☐ MB71	Leg positioning calf band	£28.00
☐ MB73	Lateral heel blocks	£83.00

Other accessories

☐ MZ51	Lap belt with metal buckle	£109.00
☐ MZ52	Tool kit	£55.00

Head support

☐ MD22	Attachment for head support mounting kit	£145.00
☐ MZ61	Mounting kit for head support, horizontal	£196.00
☐ MZ62	Mounting kit for head support, horizontal and vertical	£333.00
☐ MZ63	Head support, complete incl. mounting kit, horizontal and vertical	£385.00

Note: Please also use the order form for wheelchair accessories 647F661.

Spoke protector

☐ MG90	Transparent	£162.00
☐ MG100	Motiv Graphics	£235.00
☐ MG101	Motiv Miracle	£235.00
☐ MG102	Motiv Network	£235.00
☐ MG103	Motiv Robots	£235.00
☐ MG104	Motiv Spectrum	£235.00

Rigid frame adaptive wheelchair *Ventus*.

Frame colours

☐ 003	Jet Black (glossy) *	Inclusive
☐ 023	White (glossy)	Inclusive
☐ 001	Brilliant silver (glossy) *	Inclusive
☐ 002	Anthracite metallic (glossy)	Inclusive
☐ 004	Ice blue metallic (glossy) *	Inclusive
☐ 005	Sparkle Light Orange (glossy) *	Inclusive
☐ 009	Vinho Sparkle (glossy) *	Inclusive
☐ 010	Marine blue (glossy) *	Inclusive
☐ 012	Sky Blue (glossy) *	Inclusive
☐ 014	Candy red (glossy) *	Inclusive
☐ 016	Silver metallic **	Inclusive
☐ 021	Granny Smith (glossy) *	Inclusive
☐ 027	Shimano matt metallic *	Inclusive
☐ 048	Obsidian Black (textured) **	£108.00
☐ 013	Jet Black (matte) *	Inclusive
☐ 051	Sunny orange (textured) **	£108.00
☐ 056	Pearl White metallic **	£108.00
☐ 058	Mint Green (matte) **	£108.00
☐ 059	Raspberry Red (glossy) **	£108.00
☐ 060	Azure Blue metallic **	£108.00

* available until end 06/26

** available by beginning 07/26

Special accessories

☐ MZ58	Bag for drive wheels (for 2 drive wheels up to 28")	£122.00
☐ MZ59	Pocket for mobile phone	£37.00
☐ MZ70	Kit for transport in car	£295.00
☐ MZ89	Protection for quick-release-axle	£12.00
☐ MZ93	Backpack	£72.00
☐ MZ92	Wheel Blades (mini-skis for caster wheels)	£475.00

Rigid frame adaptive wheelchair *Ventus*.



**MA51/ MA52
Tip-assist**




**MA55/ MA56
Anti-tipper**




**MF30/31/32
Solid rubber**



**MF04
Soft roll wheel, grey**



**MF05
Pneumatic**




**MF35/36/37
Soft caster,
aluminium**



**MF90
Caster fork Standard**



**MF94
Caster fork Design**



**MG06
Hollow rim**



**MG15
Infinity Ultralight**



**MG10
Spinergy SPOX**



**MG11
Spinergy LX**

Rigid frame adaptive wheelchair *Ventus*.



MG31
Right Run 1"




MG37
Schwalbe Marathon Plus
Evolution 1"




MG63/ MG64
Curve L Grip



MG65
Gekko




MG85
Wheelbase extension with
second axle



MG86
Wheelbase
easy extension



MH13
Knee lever wheel lock,
enhanced




MH05
Scissor wheel lock



ME01
Aluminium side panel



ME02
Clothing guard,
aluminium



ME03 / ME04
Clothing guard
Plastic/carbon



ME05
with forearm support,
height-adjustable



ME50
Arm support, padded, height-ad-
justable, swing-away




MB03
Single-panel leg support,
angle and depth-adjustable




MB06
Foot stirrup, fixed


Rigid frame adaptive wheelchair *Ventus*.




MC02
Cubic Foam




MC05
Evolight



MC20
Terra



MC21
Terra Aquos




MC22
Terra Flair



MC23
Cloud



MC29
Livo



MC30
Evoflair



MB51
Lap belt with metal buckle



MB52
Tool kit



MZ62
Mounting kit for head sup-
port , horizontal + vertical



MB63
Head support, compl. incl.
mounting kit



MZ58
Bag for drive wheels



MZ59
Pocket for mobile phone



MZ92
Wheel Blades



MZ93
Backpack

Rigid frame adaptive wheelchair Ventus.


A) Anterior seat height (ASH) → Fork length in cm			
Caster wheel size	Caster fork Standard		
	extra short	short	long
4"	42 - 43	44 - 48	49 - 51
5"	–	45 - 49	49 - 53
5,5"	–	47 - 49	50 - 53
6"	–	47 - 49	50 - 53
Caster wheel size	Caster fork Design		
	extra short	short	long
4"	43	46 - 49	50 - 53
5"	–	46 - 50	50 - 53
5,5"	–	48 - 51	51 - 54
6"	–	50 - 51	51 - 54
Caster wheel size	Caster fork Frogleg		
	extra short	short	long
4"	–	48 - 50	50 - 52
5"	–	50 - 51	52 - 53

B) Posterior seat height (PSH) Drive → wheel/slidderr in cm				
Drive wheel size	Axle with camber	Slider		
		short	medium	long
22"	0°/3°	35 (38*) – 40	35 (38*) – 44	40 (44*) – 49
	6°	34 (37*) – 39	34 (37*) – 43	39 (43*) – 48
	9°	33 – 38**	33 – 42**	38 – 47**
24"	0°/3°	38 (41*) – 42	38 (41*) – 46	43 (47*) – 52
	6°	37 (40*) – 42	37 (40*) – 46	42 (46*) – 51
	9°	36 – 41**	36 – 43**	41 – 50**
25"	0°/3°	39 (42*) – 43	39 (42*) – 47	44 (48*) – 53
	6°	38 (41*) – 43	38 (41*) – 47	43 (47*) – 52
	9°	37 – 42**	37 – 46**	42 – 51**

* For wheelbase extension with 2nd axle, ** Cannot be selected for wheelbase extension with 2nd axle

C) Wheel camber → Overall width in cm			
Drive wheel camber	22" drive wheel	24" drive wheel	25" drive wheel
3°	< 6	6	> 6
6°	< 12	12	> 12
9°	< 18	18	> 18

D) Drive wheel position → Distance from tipping point (in mm)	
Dimension (horizontal) from centre axis to apex point Back support tube (front edge)/upholstery bar (top edge)	
Drive wheel position	Distance from tipping point
1 (active)	130 mm
2	120 mm
3	110 mm
4	100 mm
5	90 mm
6	80 mm
7	70 mm
8	60 mm
9	50 mm



Rigid frame adaptive wheelchair *Ventus*.

E) ATR / Posterior seat height (PSH), exclusions ME05						
Posterior seat height [cm]	Drive wheel without wheelbase extension			Drive wheel with wheelbase extension		
	22"	24"	25"	22"	24"	25"
33	●	●	●	●	●	●
34	●	●	●	●	●	●
35	●	●	●	●	●	●
36	●	●	●	●	●	●
37	●	●	●	●	●	●
38	●	●	●	●	●	●
39	●	●	●	●	●	●
40	●	●	●	●	●	●
41	●	●	●	●	●	●
42	●	●	●	●	●	●
43	●	●	●	●	●	●
44	●	●	●	●	●	●
45	●	●	●	●	●	●
46	●	●	●	●	●	●
47	●	●	●	●	●	●
48	●	●	●	●	●	●
49	●	●	●	●	●	●
50	●	●	●	●	●	●
51	●	●	●	●	●	●
52	●	●	●	●	●	●
53	●	●	●	●	●	●
●	Combination is possible					
●	Combination is not possible					

F) Anterior frame offset / LL → Foot stirrup length						
Front frame angle	0 cm offset		2 cm, 4 cm offset		4 cm abduction	
	75°	85°	75°	85°	75°	85°
Leg support bar short	–	–	LLL 39 - 45 cm	LLL 39 - 48 cm	LLL 39 - 47 cm	LLL 39 - 48 cm
Leg support bar long	LLL 40 - 50 cm	LLL 40 - 50 cm	LLL >= 44 cm	LLL >= 46 cm	LLL >= 47 cm	LLL >= 47 cm