

2026

ottobock.

Power wheelchair *Wingus*.☐ Quotation☐ Order

Billing Address

PO no. Cust Name Street County PC

Shipping address (if different from billing address)

Cust Name Street County PC Name Date Signature

Base Price - £3026.00

Please email all orders to
Bockuk@ottobock.comValid from 01.03.2026 to 31.12.2026
Please ensure you are using the latest order form
v2

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.

Wheelchairs with CLT package cannot be returned or exchanged


Article number: 490E163=00011
Standard equipment Selection
(unless alternative selected)

- nVR2 control unit (40 A)
- Max. speed: [6 km/h or 7,2 km/h]
- Standard control holder, mounting side right
- Gel batteries: 27Ah (C5)/35Ah (C20) or 28,75Ah (C5)/35Ah (C20)
- Max. distance range 25 km
- Battery charger: [5 A or 10 A]
- Caster wheels 8" (grey)
- Drive wheels 12,5" (grey)
- Puncture-proof wheels
- Seat incl. seat cushion (W x D x H: 45 x 51 x 5 cm)
- Back angle adjustment mechanical (0°/10°/20°)
- Flip-up arm support
- Height-adjustable leg support with swing-up foot plate,

The number codes contained in this order form are not intended to be used as article numbers for spare parts orders. Please scan the QR code to access the spare parts catalogue for any spare parts orders.



- angle adjustable (-5°/0°/+5°)
- Seat width 40 – 46 cm, adjustable via arm support (seat plate width is always 44 cm)
- Seat depth 38 – 48 cm
- Seat height 50 cm
- LLL adjustable: 38 – 48 cm
- Lap belt with clip closure
- Frame colour: Black (CDP)
- Empty weight: 57 kg
- **max. load capacity: 120 kg**

Wingus

	Quantity	
☐ 490E163=00011 Wingus		£3026.00

Options and Wingus Accessories.

Options and accessories			Quantity	
☐	491S163=SK001	Electrical light for footpaths (mounting side right or left)		£65.18
☐	491S163=SK006	Electrical light		£545.00
☐	493U163=SK040	Control panel holder, swing-away in parallel		£110.00
☐	493U163=SK046	Holder for attendant control, cpl.		£95.19
☐	491S163=SK015	Splash guard for rear wheels		£80.00
☐	491B163=ST015	Foot positioning belt (please order 2 pieces)		£33.00
☐	471S71=SK035	Cupholder		£41.66
☐	481S00=ST102	Case for mobile phone		£44.43
☐	491S163=SK040	Crutch holder		£55.00
☐	491S163=SK030	Storage bag		£316.33
☐	491D16=SK045	Adapter for head support mounting kit		£80.00
☐	476C00=SM1384	Seat cushion SW x SD: 45 cm x 51 cm		£83.00
☐	476C00=SK1456	Seat cushion Livo SW x SD: 45 x 42 cm		£104.00
☐	476C00=SK1457	Seat cushion Livo SW x SD: 45 x 44 cm		£104.00
☐	476C00=SK1458	Seat cushion Livo SW x SD: 45 x 46 cm		£104.00
☐	476C00=SK1489	Seat cushion Evoflair SW x SD: 45 x 42 cm		£298.00
☐	476C00=SK1490	Seat cushion Evoflair SW x SD: 45 x 44 cm		£298.00
☐	476C00=SK1491	Seat cushion Evoflair SW x SD: 45 x 46 cm		£298.00
☐	476P52=SK015	Chest and shoulder belt Size S *		£217.00
☐	476P52=SK016	Chest and shoulder belt Size M *		£216.00
☐	476P52=SK017	Chest and shoulder belt Size L *		£215.00

* only possible in combination with 491S163=SK045

Options and Wingus Accessories.



Electrical light for footpaths



Electrical light front



Electrical light rear



**Control panel holder,
swing-away in parallel**



**Holder for attendant
control**



**Splash guard for rear
wheels**



Foot positioning belts



Cupholder



Case for mobile phone



Crutch holder top




Crutch holder button




Storage bag



**Adapter for head support
mounting kit**



Seat cushion Livo



Seat cushion EvoFlair



Chest and shoulder belt