

2022/23

ottobock.

Base Price - £2695.00

Order no.:

Valid from 01.04.22 to 31.03.23
Please ensure you are using the latest order form☐ Provide estimate only☐ Order☐ Attachment for costumer fabrication enclosed

Billing address/customer number

Delivery address

Sales person

Name

Date Signature

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.

Zenit

Active wheelchair

Article number: 480F160=1_AA02_C



Standard equipment

- Open aluminium tube frame, foldable
- Choice of front frame offset
- Foldable leg support
- Back support upholstery adjustable
- Aluminium side panels
- Modular design
- Drive wheels with quick release axles
- Knee lever wheel lock 4.0
- Calf strap

Folding, max. user weight 120 kg

Crash-tested according to DIN ISO 7176-19

Wheelchair measurements*

1 Anterior seat height (ASH)

= ground to top of seat upholstery

2 Lower leg length (LL)

= top of seat upholstery to top of foot plate
(without seat cushion)

3 Posterior seat height (PSH)

= ground to top of seat upholstery

4 Back support height (BH)

= top of seat tube (crossbrace) to top of back support upholstery


5 Seat depth (SD)

= front of seat upholstery to back support tube

6 Seat width (SW)

= outside edge of frame tube to outside edge of frame tube

A manufacturing tolerance of 10 mm applies to all measurements and of 10% to specified weights.



The number codes contained in this order form are not intended to be used as article numbers for spare parts orders.
Please use the spare parts catalogue for spare parts orders.

Equipment			
MA 110	☒	Standard	£2695.00
Frame tube material			
MA 08	☐	Aluminium	Inclusive
MA 12	☐	Carbon (only with 85° and 2 cm offset/abduction) (max. load capacity 100 kg)	£441.00
Spring support folding mechanism			
MA 120	☒	Standard	Inclusive
Front frame angle			
MA 20	☐	75°	Inclusive
MA 21	☐	85°	Inclusive
Front frame offset			
MA 14	☐	0 cm offset	Inclusive
MA 15	☐	2 cm offset	Inclusive
MA 30	☐	2 cm abduction (possible up to SW 42 cm)	Inclusive
Front frame extension			
MA 23	☐	3 cm	£87.00
Seat width SW in cm			
☐ 32 ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46			
Seat depth SD in cm			
☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48 ☐ 50			
Anterior seat height in cm			
☐ 38 ☐ 39 ☐ 40 ☐ 41 ☐ 42 ☐ 43 ☐ 44 ☐ 45 ☐ 46 ☐ 47 ☐ 48 ☐ 49 ☐ 50 ☐ 51 ☐ 52 ☐ 53 ☐ 54			
Frame accessories			
MA 51	☐	Tip-assist, right	+ 270 g £101.00
MA 52	☐	Tip-assist, left	+ 270 g £101.00
MA 55	☐	Anti-tipper, right	+ 517 g £178.00
MA 56	☐	Anti-tipper, left	+ 517 g £178.00
MA 60	☐	Frame pad	+ 42 g £65.00
Caster wheels			
MF 30	☐	4" solid rubber, black	Inclusive
MF 35	☐	4" soft caster, aluminium	£132.00
MF 31	☐	5" solid rubber, black	Inclusive
MF 36	☐	5" soft caster, aluminium	£132.00
MF 04	☐	5.5" soft caster, grey	Inclusive
MF 05	☐	6" pneumatic tyre	Inclusive
MF 32	☐	6" solid rubber, black	Inclusive
MF 37	☐	6" soft caster, aluminium	£132.00
Caster wheel axle			
MF 90	☒	Caster wheel forks with threaded axles	Inclusive

Posterior seat height in cm							
☐ 37	☐ 38	☐ 39	☐ 40	☐ 41	☐ 42	☐ 43	☐ 44
☐ 44	☐ 45	☐ 46	☐ 47	☐ 48	☐ 49	☐ 50	

Drive wheel size			
MG 02	☐	24"	Inclusive
MG 00	☐	25"	Inclusive
Drive wheel camber ☐ 0° ☐ 2° ☐ 3° ☐ 4°			

Drive wheel rim		24"	25"		
MG 06	☐	Hollow rims	•	•	Inclusive
MG 10	☐	Spinergy SPOX - 360 g	•		£851.00
MG 11	☐	Spinergy LX - 400 g	•		£893.00
MG 12	☐	Spinergy LXL Blade, black spokes - 660 g	•		£780.00
MG 12	☐	Spinergy LXL Blade, Blade, white spokes - 660 g	•		£780.00
MG 15	☐	Infinity Ultralight (not possible in combination with MG 34) - 720 g	•	•	£95.00
MG 150	☐	Infinity Ultralight, coloured (not possible in comb. with MG 34) - 720 g	•		£162.00
		Colour set Infinity Ultralight/caster wheel fork	24"	25"	
		☐ Sparkle light orange	•		
		☐ Marine blue	•		
		☐ Candy red	•		
		☐ Sparkle granny-smith	•		

Tyres		24"	25"		
MG 31	☐	Schwalbe RightRun 1"	•	•	Inclusive
MG 34	☐	PU 1" (not possible in combination with MG 15/MG 150) + 344 g	•		Inclusive
MG 341	☐	PU 1", black (not possible in combination with MG 15/MG 150) + 344 g	•		Inclusive
MG 37	☐	Schwalbe Marathon Plus Evolution 1" + 374 g	•	•	£52.00
MG 371	☐	Schwalbe Marathon Plus Evolution 1.1" + 379 g	•		£61.00
MG 200	☐	Schwalbe ONE 1" - 226 g	•		£84.00

Handrims		24"	25"		
MG 50	☐	Aluminium, black coated, high	•	•	Inclusive
MG 52	☐	Aluminium, silver anodised, high - 40 g	•	•	£12.00
MG 53	☐	Aluminium, silver anodised, low - 120 g	•		£12.00
MG 54	☐	Rubber-coated + 500 g	•	•	£132.00
MG 56	☐	Stainless steel, polished, high + 620 g	•	•	£51.00
MG 57	☐	Stainless steel, polished, low + 580 g	•		£51.00
MG 63	☐	CURVE L (anodised) + 200 g	•	•	£335.00
MG 64	☐	CURVE L Grip for users with limited gripping function + 840 g	•	•	£362.00
MG 65	☐	GEKKO, ergonomic handrim + 580 g	•	•	£297.00

Handrim mounting			
MG 70	☐	Mounting position narrow (15 mm)	Inclusive
MG 71	☐	Mounting position wide (25 mm)	Inclusive
Tipping point in cm (measured parallel to the ground)			
☐ 5	☐ 6	☐ 7	☐ 8
☐ 9	☐ 10	☐ 11	☐ 12
☐ 13	☐ 14		
Alternative measuring method (specification of the horizontal drive wheel position of the Avantgarde)			
☐ A	☐ B	☐ C	☐ D
☐ E	☐ F	☐ G	☐ H
(Tick the drive wheel position of your Avantgarde, we calculate the tipping point.)			
Quick-release axle system			
MG 27	☐	Quick-release axle	Inclusive
MG 28	☐	Quick-release axle with Tetra release	+ 158 g £164.00
Drive wheel accessories			
MG 87	☐	High-pressure air pump, 14 bar	£99.00
MG 88	☐	Airman pump	£200.00
Wheel locks			
MH 01	☐	Knee lever wheel lock 4.0	Inclusive
MH 30	☐	Knee lever wheel lock 4.0, long lever	Inclusive
MH 05	☐	Scissor wheel lock	- 142 g £123.00
Wheel lock accessories			
MH 10	☐	Wheel lock lever extension, fixed	+ 70 g £33.00
MH 21	☐	Wheel lock lever extension, pluggable	+ 34 g £60.00
Back support height in cm			
☐ 25	☐ 27.5	☐ 30	☐ 32.5
☐ 35	☐ 37.5	☐ 40	☐ 42.5
☐ 45	☐ 47.5	☐ 50	
Back support system			
MD 01	☐	Back support, straight (measured 90° to seat tube)	Inclusive
MD 11	☐	Angle-adjustable, fold-down onto seat tubes (-10°/-5°/0°/+5°/+10°)	+ 240 g £557.00
MD 13	☐	Back support angle according to customer request (measured 90° to the ground)	£89.00
		☐ 0° ☐ -5° ☐ -10°	
Back support upholstery			
MD 03	☐	Adaptable, padded	£47.00
MD 05	☐	Adaptable, breathable	Inclusive
MD 18	☐	Ultra light	£45.00
MD 06	☐	Without	Inclusive
Back support pad widening			
MD 80	☐	+ 2 cm	£17.00
MD 81	☐	+ 4 cm	£17.00
MD 82	☐	Without	Inclusive
Use of BAXX back support system (include order form, model-specific restrictions are possible)			
MD 44	☐	BAXX Aluminium back support (not mounted)	£588.00


Push handles				
MD 50	☐	Short		Inclusive
MD 51	☐	Long	+ 128 g	Inclusive
MD 55	☐	Height-adjustable/removable	+ 580 g	£141.00
MD 56	☐	Fold-down	+ 108 g	£170.00
MD 59	☐	Without	- 258 g	Inclusive

Side panels				
ME 02	☐	Aluminium, with splash guard and protection against cold		Inclusive
ME 04	☐	Carbon, with splash guard	- 256 g	£537.00
ME 05	☐	With arm support short	+ 2490 g	£282.00
ME 06	☐	With arm support long	+ 2180 g	£282.00

Seat upholstery				
MC 15	☐	Seat upholstery with storage compartment		£35.00
MC 16	☐	Upholstery readjustable		Inclusive
MC 18	☐	Ultra light		£36.00
MC 06	☐	Without		Inclusive

Seat cushions				
MC 01	☐	Standard (thickness of cushion 5 cm)		£85.00
MC 02	☐	Cubic Foam (thickness of cushion 5 cm)		£123.00
MC 05	☐	EvoLight (thickness of cushion 6 cm)		£175.00
MC 20	☐	Terra (thickness of cushion 10 cm)		£164.00
MC 21	☐	Terra Aquos (thickness of cushion 10 cm)		£197.00
MC 22	☐	Terra Flair (thickness of cushion 10 cm)		£494.00
MC 23	☐	Cloud (thickness of cushion 10 cm)		£483.00
MC 29	☐	Livo (thickness of cushion 5 cm)		£73.00
MC 30	☐	EvoFlair (thickness of cushion 8 cm)		£264.00

Lower leg length in cm (without seat cushion)				
Lower leg length (LL) in cm (available setting range 28 cm – 51 cm)				
MB 130	☐	Adjustable		Inclusive



Leg supports (folding)				
MB 160	☐	Single-panel, aluminium		Inclusive
MB 180	☐	Single-panel, carbon	- 177 g	£131.00

Leg support accessories				
MB 70	☐	Lock for single-panel foot plate		£75.00

Other accessories				
MZ 51	☐	Lap belt with metal buckle		£68.00
MZ 59	☐	Pocket for mobile phone		£27.00

Spoke protectors			
MG 90	☐	Spoke protector, transparent	£127.00
MG 100	☐	Spoke protector, Graphics decor	£183.00
MG 101	☐	Spoke protector, Miracle decor	£183.00
MG 102	☐	Spoke protector, Network decor	£183.00
MG 103	☐	Spoke protector, Robots decor	£183.00
MG 104	☐	Spoke protector, Spectrum decor	£183.00
Frame colours			
001	☐	Brilliant silver	£384.00
002	☐	Anthracite metallic	£384.00
003	☐	Jet black RAL 9005	£384.00
004	☐	Ice blue	£384.00
005	☐	Sparkle light orange	£384.00
009	☐	Vinho sparkle	£384.00
010	☐	Marine blue	£384.00
012	☐	Sky blue	£384.00
048	☐	Jet black matte (structure)	£374.00
014	☐	Candy red	£384.00
021	☐	Sparkle granny smith	£384.00
023	☐	White RAL 9016	£384.00
027	☐	Shimano matte	£384.00
039	☐	Black, anodised, matte	Inclusive
040	☐	Anthracite, anodised, matte	Inclusive
041	☐	Gold, anodised, matte	Inclusive
042	☐	Blue, anodised, matte	Inclusive
043	☐	Orange, anodised, matte	Inclusive
044	☐	Green, anodised, matte	Inclusive
045	☐	Purple, anodised, matte	Inclusive
046	☐	Red, anodised, matte	Inclusive
Special accessories			
MZ 52	☐	Tool kit	£55.00
MZ 58	☐	Bag for drive wheels (for 2 drive wheels up to 28")	£96.00
MZ 70	☐	Kit for transportation in a motor vehicle	£286.00
MZ 89	☐	Quick-release axle protection	£8.00
MZ 92	☐	Wheelblades (Mini-skis for caster wheels)	£375.00
MZ 93	☐	Backpack	£54.00

Information



MG 06 Hollow rim



MG 10 Spinergy SPOX



MG 11 Spinergy LX




MG 15 Infinity Ultralight




MF 35/36/37 4/5/6" soft caster



MG 100 Graphic



MG 103 Robots



MG 104 Spectrum