

2024

ottobock.

Base Price - £3084.00 Valid from 01.01.24 to 31.12.24 Please ensure you are using the latest order form		Order no.: -----
☐ Provide estimate only	☐ Order	
☐ Attachment for custom fabrication enclosed		
Billing address/customer number		
Delivery address		
Sales person	Name Date Signature	
By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products, and you confirm to have read, understood and accepted our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.		

Zenit R

Active wheelchair

Article number: 480F160=1_AA04_C



Standard equipment


- Open aluminium tube frame
- Choice of front frame offset
- Back support upholstery, adaptable
- Aluminium side panels
- Modular design
- Drive wheels with quick-release axle
- Knee lever wheel lock 4.0
- Calf strap

Max. user weight 120 kg
 Crash-tested according to DIN ISO 7176-19

Wheelchair measurements*

- 1 Anterior seat height (ASH)
= ground to top of seat upholstery
- 2 Lower leg length (LL)
= top of seat upholstery right/left to top of foot plate centre (without seat cushion)
- 3 Posterior seat height (PSH)
= ground to top of seat upholstery
- 4 Back support height (BH)
= top of seat tube (crossbrace) to top of back support upholstery
- 5 Seat depth (SD)
= front of seat upholstery to back support tube
- 6 Seat width (SW)
= outside edge of frame tube to outside edge of frame tube

A manufacturing tolerance of 10 mm applies to all measurements and of 10% to specified weights.



The number codes contained in this order form are not intended to be used as article numbers for spare parts orders.
 Please use the spare parts catalogue for spare parts orders.

Equipment			
MA 110	☒	Standard	£3084.00
Frame tube material			
MA 08	☐	Aluminium	Inclusive
MA 12	☐	Carbon (available only with 85° anterior frame angle and 2 cm offset/abduction) (max. load capacity 100 kg)	£505.00
Anterior frame angle			
MA 20	☐	75°	Inclusive
MA 21	☐	85°	Inclusive
Front frame offset			
MA 14	☐	0 cm offset	Inclusive
MA 15	☐	2 cm offset	Inclusive
MA 30	☐	2 cm abduction (possible only up to SW 42 cm)	Inclusive
Front frame extension			
MA 23	☐	3 cm	£100.00
Seat width SW in cm			
☐ 32 ☐ 34 ☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46			
Seat depth SD in cm			
☐ 36 ☐ 38 ☐ 40 ☐ 42 ☐ 44 ☐ 46 ☐ 48 ☐ 50			
Anterior seat height (ASH) in cm			
☐ 38 ☐ 39 ☐ 40 ☐ 41 ☐ 42 ☐ 43 ☐ 44 ☐ 45 ☐ 46 ☐ 47 ☐ 48 ☐ 49 ☐ 50 ☐ 51 ☐ 52 ☐ 53 ☐ 54			
Frame accessories			
MA 51	☐	Tip-assist, right	+ 270 g £116.00
MA 52	☐	Tip-assist, left	+ 270 g £116.00
MA 55	☐	Anti-tipper, right	+ 517 g £204.00
MA 56	☐	Anti-tipper, left	+ 42 g £204.00
MA 60	☐	Frame pad	+ 42 g £75.00
Caster wheels			
MF 30	☐	4" solid rubber, black	Inclusive
MF 35	☐	4" soft caster, aluminium	£152.00
MF 31	☐	5" solid rubber, black	Inclusive
MF 36	☐	5" soft caster, aluminium	£152.00
MF 04	☐	5.5" soft caster, grey	Inclusive
MF 05	☐	6" pneumatic tyre	Inclusive
MF 32	☐	6" solid rubber, black	Inclusive
MF 37	☐	6" soft caster, aluminium	£152.00
Caster wheel axle			
MF 90	☒	Forks with threaded axle	Inclusive
Posterior seat height (PSH) in cm			
☐ 37 ☐ 38 ☐ 39 ☐ 40 ☐ 41 ☐ 42 ☐ 43 ☐ 44 ☐ 45 ☐ 46 ☐ 47 ☐ 48 ☐ 49 ☐ 50			

Drive wheel size				
MG 02	☐	24"		Inclusive
MG 00	☐	25"		Inclusive
Drive wheel camber ☐ 0° ☐ 2° ☐ 3° ☐ 4°				

Drive wheel rim			24"	25"	
MG 06	☐	Hollow rim	•	•	Inclusive
MG 10	☐	Spinergy SPOX - 360 g	•		£718.00
MG 11	☐	Spinergy LX - 400 g	•		£792.00
MG 12	☐	Spinergy LXL Blade, black spokes - 660 g	•		£893.00
MG 12	☐	Spinergy LXL Blade, white spokes - 660 g	•		£893.00
MG 15	☐	Infinity Ultralight (not possible in combination with MG 34) - 720 g	•	•	£109.00
MG 150	☐	Infinity Ultralight, coloured (not possible in combination with MG 34) - 720 g	•		£186.00
Infinity Ultralight/caster wheel fork colour package			24"	25"	
☐ Sparkle light orange			•		
☐ Marine blue			•		
☐ Candy red			•		
☐ Sparkle granny-smith			•		

Tyres			24"	25"	
MG 31	☐	Schwalbe RightRun 1"	•	•	Inclusive
MG 34	☐	PU 1" (not possible in combination with MG 15/MG 150) + 344 g	•		Inclusive
MG 341	☐	PU 1", black (not possible in combination with MG 15/MG 150) + 344 g	•		Inclusive
MG 37	☐	Schwalbe Marathon Plus Evolution 1" + 374 g	•	•	£52.00
MG 371	☐	Schwalbe Marathon Plus Evolution 1.1" + 379 g	•		£52.00
MG 200	☐	Schwalbe ONE 1" - 226 g	•		£85.00

Handrims			24"	25"	
MG 50	☐	Aluminium, black coated, high	•	•	Inclusive
MG 52	☐	Aluminium, silver anodised, high - 40 g	•	•	£14.00
MG 53	☐	Aluminium, silver anodised, low - 120 g	•		£14.00
MG 54	☐	Rubber-coated, high + 500 g	•	•	£152.00
MG 56	☐	Stainless steel, polished, high + 620 g	•	•	£59.00
MG 57	☐	Stainless steel, polished, low + 580 g	•		£59.00
MG 63	☐	CURVE L (anodised) + 200 g	•	•	£384.00
MG 64	☐	CURVE L Grip for users with limited gripping function + 840 g	•	•	£415.00
MG 65	☐	GEKKO, ergonomic handrim + 580 g	•	•	£340.00

Handrim mounting					
MG 70	☐	Mounting position narrow (15 mm)			Inclusive
MG 71	☐	Mounting position wide (25 mm)			Inclusive


Tipping point in cm (measured parallel to the ground)									
☐ 5	☐ 6	☐ 7	☐ 8	☐ 9	☐ 10	☐ 11	☐ 12	☐ 13	☐ 14
Alternative measuring method (specification of the horizontal drive wheel position of the Avantgarde)									
☐ A	☐ B	☐ C	☐ D	☐ E	☐ F	☐ G	☐ H		
(Tick the drive wheel position of your Avantgarde, we calculate the tipping point)									
Quick-release axle system									
MG 27	☐	Quick-release axle	Inclusive						
MG 28	☐	Quick-release axle with Tetra release	+ 158 g £161.00						
Drive wheel accessories									
MG 87	☐	High-pressure air pump, 14 bar	£114.00						
MG 88	☐	Airman pump	£173.00						
Wheel locks									
MH 01	☐	Knee lever wheel lock 4.0	Inclusive						
MH 30	☐	Knee lever wheel lock 4.0, long wheel lock lever	Inclusive						
MH 05	☐	Scissor wheel lock	- 142 g £141.00						
Wheel lock accessories									
MH 10	☐	Wheel lock lever extension, fixed	+ 70 g £38.00						
MH 21	☐	Wheel lock lever extension, pluggable	+ 34 g £69.00						
Back support height in cm									
☐ 25	☐ 27.5	☐ 30	☐ 32.5	☐ 35	☐ 37.5				
☐ 40	☐ 42.5	☐ 45	☐ 47.5	☐ 50					
Back support system									
MD 01	☐	Back support, straight (measured 90° to the seat tube)	Inclusive						
MD 11	☐	Angle-adjustable, fold-down onto seat tubes (-10°/-5°/0°/+5°/+10°)	+ 240 g £638.00						
MD 13	☐	Back support angle according to customer request (measured 90° to the ground) ☐ 0° ☐ -5° ☐ -10°	£102.00						
Back support upholstery									
MD 03	☐	Adaptable, padded	£54.00						
MD 05	☐	Adaptable, breathable	Inclusive						
MD 18	☐	Ultra light	£52.00						
MD 06	☐	Without	Inclusive						
Back support pad widening									
MD 80	☐	+ 2 cm	£20.00						
MD 81	☐	+ 4 cm	£20.00						
MD 82	☐	Without	Inclusive						
Use of BAXX back support system (include order form, model-specific restrictions are possible)									
MD 44	☐	BAXX aluminium back support (not mounted)	£673.00						

Push handles			
MD 50	☐	Short	Inclusive
MD 51	☐	Long	+ 128 g Inclusive
MD 55	☐	Height-adjustable/removable	+ 580 g £162.00
MD 56	☐	Fold-down	+ 108 g £195.00
MD 59	☐	Without	- 258 g Inclusive
Side panels			
ME 02	☐	Aluminium, with splash guard and protection against cold	Inclusive
ME 04	☐	Carbon, with splash guard	£615.00
ME 05	☐	With arm support, short	£323.00
ME 06	☐	With arm support, long	£323.00

Seat upholstery			
MC 15	☐	Upholstery with storage compartment	£41.00
MC 16	☐	Upholstery, readjustable	Inclusive
MC 18	☐	Ultra light	£42.00
MC 06	☐	Without	Inclusive

Seat cushions			
MC 01	☐	Standard (thickness of cushion 5 cm)	£98.00
MC 02	☐	Cubic Foam (thickness of cushion 5 cm)	£141.00
MC 05	☐	Evolight (thickness of cushion 6 cm)	£201.00
MC 20	☐	Terra seat (thickness of cushion 10 cm)	£188.00
MC 21	☐	Terra Aquos (thickness of cushion 10 cm)	£226.00
MC 22	☐	Terra Flair (thickness of cushion 10 cm)	£566.00
MC 23	☐	Cloud (thickness of cushion 10 cm)	£553.00
MC 29	☐	Livo (thickness of cushion 5 cm)	£84.00
MC 30	☐	Evoflair (thickness of cushion 8 cm)	£303.00

Lower leg length in cm (without seat cushion)			
Lower leg length (LL) in cm (available setting range 28 cm – 51 cm)			
MB 130	☒	Adjustable	Inclusive



Leg supports			
MB 160	☐	Single-panel, aluminium* (LL 28 – 51 cm)	Inclusive
MB 180	☐	Single-panel, carbon* (LL 28 – 51 cm)	£250.00
MB 06	☐	Aluminium light, fixed (LL 31 – 51 cm)	£153.00
* LL can be selected in dependence of wheel size, seat width and anterior seat height.			

Other accessories			
MZ 51	☐	Lap belt with metal buckle	£78.00
MZ 59	☐	Pocket for mobile phone	£31.00

Spoke protectors			
MG 90	☐	Transparent	£146.00
MG 100	☐	Graphics decor	£210.00
MG 101	☐	Miracle decor	£210.00
MG 102	☐	Network decor	£210.00
MG 103	☐	Robots decor	£210.00
MG 104	☐	Spectrum decor	£210.00
Frame colours			
001	☐	Brilliant silver	£440.00
002	☐	Anthracite metallic	£440.00
003	☐	Jet black RAL 9005	£440.00
004	☐	Ice blue	£440.00
005	☐	Sparkle light orange	£440.00
009	☐	Vinho sparkle	£440.00
010	☐	Marine blue	£440.00
012	☐	Sky blue	£440.00
048	☐	Jet black matte (texture)	£440.00
014	☐	Candy red	£440.00
021	☐	Sparkle granny-smith	£440.00
023	☐	White RAL 9016	£440.00
027	☐	Shimano matte	£440.00
039	☐	Black, matte anodised	Inclusive
040	☐	Anthracite, matte anodised	Inclusive
041	☐	Gold, matte anodised	Inclusive
042	☐	Blue, matte anodised	Inclusive
043	☐	Orange, matte anodised	Inclusive
044	☐	Green, matte anodised	Inclusive
045	☐	Violet, matte anodised	Inclusive
046	☐	Red, matte anodised	Inclusive
Special accessories			
MZ 52	☐	Tool kit	£63.00
MZ 58	☐	Bag for drive wheels (for 2 drive wheels up to 28")	£110.00
MZ 70	☐	Kit for transportation in a motor vehicle	£328.00
MZ 89	☐	Quick-release axle protection	£10.00
MZ 92	☐	Wheelblades (Mini-skis for caster wheels)	£429.00
MZ 93	☐	Backpack	£62.00

Information



MG 06 Hollow rim



MG 10 Spinerger SPOX



MG 11 Spinerger LX



MG 15 Infinity Ultralight




MF 35/36/37 4/5/6"
soft caster



MG 100 Graphic



MG 103 Robots



MG 104 Spectrum



MF 90 Caster fork with
threaded axle



MH 01 Knee lever wheel
lock 4.0



MD 55 Push handles,
height-adjustable/re-
movable



ME 04 Carbon side pan-
el, with splash guard



MD 03 Back support up-
holstery, adaptable,
padded



MB 06 Aluminium light leg
support, fixed



Drive wheel adapter, per-
manently welded