

2026

ottobock.

Active wheelchair Zenit R.

☐ Quotation

☐ Order

Billing Address

PO no.	
Customer Name	
Street	
County, PC	

Shipping address (if different from billing address)

Customer Name	
Street	
County, PC	

Name _____ Date _____ Signature _____



Base Price - £3279.00

Please email all orders to
Bockuk@ottobock.com

Valid from 01.03.2026 to 31.12.26
Please ensure you are using the latest order form
v2

Standard equipment Selection (unless alternative selected)

Article number: 480F160=1_AA04_C

- Choice of front frame offset
- Adaptable back support upholstery
- Tension adjustable seat upholstery
- Aluminium side panels
- Modular design
- Drive wheels with quick-release axle
- Knee lever wheel lock 4.0
- Calf strap

max. user weight 120 kg
crash-tested according ISO 7176-19

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals. Wheelchairs with CLT package cannot be returned or exchanged

The number codes contained in this order form are not intended to be used as article numbers for spare parts orders.
Please scan the QR code to access the spare parts catalogue for any spare parts orders.




Active wheelchair *Zenit R*.

Wheelchair measurements

- 1 Anterior seat height (aSH) = ground to top of seat upholstery
- 2 Lower leg length (LLL) = top of seat upholstery to top of foot plate (without seat cushion). Please consider the dependencies LLL and aSH.
- 3 Posterior seat height (pSH) = ground to top of seat upholstery
- 4 Back support height (BH) = top of seat tube to top of back support upholstery
- 5 Seat depth (SD) = front of seat upholstery to back support tube
- 6 Seat width (SW) = outside edge of frame tube to outside edge of frame tube

A manufacturing tolerance of 10 mm applies to all measurements and of 10 % to specified weights.



Active wheelchair Zenit R.

Equipment

☒ MA110	Standard	£3279.00
---	----------	----------

Material

☐ MA08	Aluminium	Inclusive
☐ MA12	Carbon (with 85° anterior frame angle and 2 cm offset/abduction)	£538.00

Anterior frame angle (please consider LLL on page 6)

☐ MA20	75° long frame	Inclusive
☐ MA21	85° short frame	Inclusive

Front fram offset

☐ MA14	0 cm offset	Inclusive
☐ MA15	2 cm offset	Inclusive
☐ MA30	2 cm abduction (possible only up to SW 42 cm)	Inclusive

Front frame extension

☐ MA23	3 cm	£108.00
----------------------------	------	---------

Seat width (SW) in cm

☐ 32	☐ 34	☐ 36	☐ 38	☐ 40	☐ 42	☐ 44	☐ 46
--------------------------	--------------------------	--------------------------	--------------------------	--------------------------	--------------------------	--------------------------	--------------------------

Seat depth (SD) in cm

☐ 36	☐ 38	☐ 40	☐ 42	☐ 44	☐ 46	☐ 48	☐ 50
--------------------------	--------------------------	--------------------------	--------------------------	--------------------------	--------------------------	--------------------------	--------------------------

Anterior seat height (aSH) in cm (please consider LLL on page 6)

☐ 38	☐ 39	☐ 40	☐ 41	☐ 42	☐ 43	☐ 44	☐ 45	☐ 46
☐ 47	☐ 48	☐ 49	☐ 50	☐ 51	☐ 52	☐ 53	☐ 54	

Frame accessories

☐ MA51	Tip-assist, right	£124.00
☐ MA52	Tip-assist, left	£124.00
☐ MA55	Anti-tipper, right	£218.00
☐ MA56	Anti-tipper, left	£218.00
☐ MA60	Frame pads	£81.00

Caster wheels

☐ MF30	4" solid rubber, black	Inclusive
☐ MF35	4" soft caster, aluminium	£162.00
☐ MF31	5" solid rubber, black	Inclusive
☐ MF36	5" soft caster, aluminium	£162.00
☐ MF04	5.5" soft caster, grey	Inclusive
☐ MF05	6" pneumatic, grey	Inclusive
☐ MF32	6" solid rubber, black	Inclusive
☐ MF37	6" soft caster, aluminium	£162.00

Caster axle

☒ MF90	Caster forks with threaded axles	Inclusive
--	----------------------------------	-----------

Active wheelchair Zenit R.

Posterior seat height (pSH) in cm

☐ 37	☐ 38	☐ 39	☐ 40	☐ 41	☐ 42	☐ 43
☐ 44	☐ 45	☐ 46	☐ 47	☐ 48	☐ 49	☐ 50

Drive wheel size

☐ MG02	24"	Inclusive
☐ MG00	25"	Inclusive

Drive wheel camber

☐ 0°	☐ 2°	☐ 3°	☐ 4°
--------------------------	--------------------------	--------------------------	--------------------------

Drive wheel rim

		24"	25"	
☐ MG06	Hollow rims	●	●	Inclusive
☐ MG197	SmootherWheel dampening drive wheels, max. load 75 kg	●		£1659.00
☐ MG198	SmootherWheel dampening drive wheels, max. load 95 kg	●		£1659.00
☐ MG199	SmootherWheel dampening drive wheels, max. load 115 kg	●		£1659.00
☐ MG10	Spinergy SPOX	●		£764.00
☐ MG11	Spinergy LX	●		£843.00
☐ MG12	Spinergy LXL Blade: black spokes	●		£950.00
☐ MG12	Spinergy LXL Blade: white spokes	●		£950.00
☐ MG15	Infinity Ultralight (not possible in combination with MG 34)	●	●	£117.00
☐ MG150	Infinity Ultralight, coloured (not possible in comb. with MG 34)	●		£198.00
	Colour set Infinity Ultralight/caster fork	24"	25"	
	☐ Sparkle light orange	●		Inclusive
	☐ Marine blue	●		Inclusive
	☐ Candy red	●		Inclusive
	☐ Sparkle granny-smith	●		Inclusive

Tyres

		24"	25"	
☐ MG31	Schwalbe RightRun 1"	●	●	Inclusive
☐ MG34	PU 1" (not possible in combination with Infinity rim)	●		Inclusive
☐ MG341	PU 1" schwarz (not possible in combination with Infinity rim)	●		Inclusive
☐ MG37	Schwalbe Marathon Plus Evolution 1"	●	●	£76.00
☐ MG 200	Schwalbe ONE 1"	●		£190.00

Handrims

		24"	25"	
☐ MG50	Aluminium, black coated, high	●	●	Inclusive
☐ MG52	Aluminium, silver anodised, high	●	●	£16.00
☐ MG53	Aluminium, silver anodised, deep	●		£16.00
☐ MG54	Rubber-coated, high	●	●	£162.00
☐ MG56	Stainless steel, polished, high	●	●	£63.00
☐ MG57	Stainless steel, polished, deep	●		£63.00
☐ MG63	Curve L, anodised	●	●	£409.00
☐ MG64	Curve L, Grip for users with limited gripping function	●	●	£442.00
☐ MG65	GEKKO, ergonomic handrim	●	●	£362.00

Handrim mounting

☐ MG70	Mounting position narrow (15 mm)	Inclusive
☐ MG71	Mounting position wide (25 mm)	Inclusive

Active wheelchair Zenit R.

Tipping point in cm (measured parallel to the ground)

☐ 5 ☐ 6 ☐ 7 ☐ 8 ☐ 9 ☐ 10 ☐ 11 ☐ 12 ☐ 13 ☐ 14

Alternative measuring method (specification of the horizontal drive wheel position of the Avantgarde)

☐ B ☐ C ☐ D ☐ E ☐ F ☐ G

Tick the drive wheel position of your Avantgarde, we calculate the tipping point.

Quick-release axle system

☐ MG27	Quick-release axle	Inclusive
☐ MG28	Quick-release axle with Tetra release	£172.00

Drive wheel accessories

☐ MG87	High-pressure air pump, 14 bar	£122.00
☐ MG88	Airman pump	£185.00

Wheel locks

☐ MH01	Knee lever wheel lock 4.0	Inclusive
☐ MH30	Knee lever wheel lock 4.0, long lever	Inclusive
☐ MH05	Scissor wheel lock	£151.00

Wheel lock accessories

☐ MH10	Wheel lock lever extension, fixed	£42.00
☐ MH21	Wheel lock lever extension, pluggable	£75.00

Rückenhöhe in cm

☐ 25 ☐ 27,5 ☐ 30 ☐ 32,5 ☐ 35 ☐ 37,5 ☐ 40 ☐ 42,5 ☐ 45 ☐ 47,5 ☐ 50

Back support height in cm

☐ MD01	Back support, straight (measured 90° to the seat tube)	Inclusive
☐ MD11	Angle-adjustable, fold-down onto seat tubes (-10°; -5°; 0°; +5°; +10°)	£679.00
☐ MD13	Back support angle according to customer request (measured 90° to the ground)	£110.00
	☐ 0° ☐ -5° ☐ -10°	

Back support upholstery

☐ MD03	Adaptable, padded	£58.00
☐ MD05	Adaptable, breathable	Inclusive
☐ MD18	Ultra light	£56.00
☐ MD06	Without	Inclusive

Back support pad widening

☐ MD80	+ 2 cm	£22.00
☐ MD81	+ 4 cm	£22.00
☐ MD82	Without	Inclusive

Use of BAXX back support system (include order form, model-specific restrictions are possible, take into account the reduction of the seat depth by using the Baxx back support system)

☐ MD44	Baxx aluminium back support (enclosed)	£716.00
----------------------------	--	---------

Active wheelchair Zenit R.

Push handles

☐ MD50	Short	Inclusive
☐ MD51	Long	Inclusive
☐ MD55	Height-adjustable/removable	£174.00
☐ MD45	Height-adjustable/removable for Baxx	£218.00
☐ MD56	Fold-down	£209.00
☐ MD59	Without	Inclusive

Side panels

☐ ME02	Aluminium, with splash guard and protection against the cold	Inclusive
☐ ME04	Carbon, with splash guard	£655.00
☐ ME05	With arm support, short	£344.00
☐ ME06	With arm support, long	£344.00
☐ ME17	With arm support long, depth-adjustable	£163.00

Seat upholstery

☐ MC15	Seat upholstery with storage compartment (up to SW 38 cm = 1 compartment; from SW 40 cm = 2 compartments)	£45.00
☐ MC16	Upholstery readjustable	Inclusive
☐ MC18	Ultra light	£46.00
☐ MC06	Without	Inclusive

Seat cushions

☐ MC01	Standard (functional height 5 cm)	£43.00
☐ MC02	Cubic Foam (functional height 5 cm)	£151.00
☐ MC05	EvoLight (functional height 6 cm)	£215.00
☐ MC20	Terra (functional height 6 cm)	£200.00
☐ MC21	Terra Aquos (functional height 6 cm)	£241.00
☐ MC22	Terra Flair (functional height 6 cm)	£603.00
☐ MC23	Cloud (functional height 11 cm)	£589.00
☐ MC29	Livo (functional height 5 cm)	£90.00
☐ MC30	EvoFlair (functional height 8 cm)	£323.00


Lower leg length (without seat cushion)

Lower leg length (LLL) in cm

75° front frame: aSH must be greater or equal to LLL + 2 cm

85° front frame: aSH must be greater or equal to LLL + 4 cm

☐ 28	☐ 29	☐ 30	☐ 31	☐ 32	☐ 33
☐ 34	☐ 35	☐ 36	☐ 37	☐ 38	☐ 39
☐ 40	☐ 41	☐ 42	☐ 43	☐ 44	☐ 45
☐ 46	☐ 47	☐ 48	☐ 49	☐ 50	☐ 51



☒ MB130	Adjustable	Inclusive
---	------------	-----------

Leg support

☐ MB160	Single-panel, aluminium* (LLL 28 – 51 cm)	Inclusive
☐ MB180	Single-panel, carbon* (LLL 28 – 51 cm)	£266.00
☐ MB06	Aluminium light, fixed (LLL 31 – 51 cm)	£163.00

* LLL can be selected in dependence of wheel size, seat width and anterior seat height

Active wheelchair Zenit R.

Other accessories

☐ MZ51	Lap belt with metal buckle	£84.00
☐ MZ59	Pocket for mobile phone	£68.00

Spoke protectors

☐ MG90	Transparent	£156.00
☐ MG100	Motiv Graphics	£224.00
☐ MG101	Motiv Miracle	£224.00
☐ MG102	Motiv Network	£224.00
☐ MG103	Motiv Robots	£224.00
☐ MG104	Motiv Spectrum	£224.00

Frame colours

☐ 001	Brillant silver	£469.00
☐ 002	Anthracite metallic	£469.00
☐ 003	Jet black RAL 9005	£469.00
☐ 004	Ice blue	£469.00
☐ 005	Sparkle light orange	£469.00
☐ 009	Vinho sparkle	£469.00
☐ 010	Marine blue	£469.00
☐ 012	Sky blue	£469.00
☐ 014	Candy red	£469.00
☐ 021	Sparkle granny-smith	£469.00
☐ 023	White RAL 9016	£469.00
☐ 027	Shimano matte	£469.00
☐ 048	Jet black matte structure	£469.00
☐ 039	Black, anodised, matte	Inclusive
☐ 040	Anthracite, anodised, matte	Inclusive
☐ 041	Gold, anodised, matte	Inclusive
☐ 042	Blue, anodised, matte	Inclusive
☐ 043	Orange, anodised, matte	Inclusive
☐ 044	Green, anodised, matte	Inclusive
☐ 045	Purple, anodised, matte	Inclusive
☐ 046	Red, anodised, matte	Inclusive


Special accessories

☐ MZ52	Tool kit	£68.00
☐ MZ58	Bag for drive wheels (for 2 drive wheels up to 28")	£118.00
☐ MZ70	Kit for transportation in a motor vehicle	£350.00
☐ MZ89	Protection for quick-release-axle	£12.00
☐ MZ92	Wheelblades (Mini-skis for caster wheels)	£457.00
☐ MZ93	Backpack	£66.00

Active wheelchair *Zenit R.*




MG10 Spinergy SPOX



MG06 Hollow rim



MG11 Spinergy LX



MG15 Infinity Ultralight



MG31 Right Run 1"




MG371 Marathon Plus Evolution 1.1"



Drive wheel adapter, permanently welded



MF90 Caster forks with threaded axles



MG54 Handrim rubber-coated




MG65 Gekko



MG63 Curve L (anodised)




MG64 Curve L Grip




MF35/36/37 soft caster 4", 5", 6"



MF30/31/32 Solid rubber, black 4", 5", 6"



MF04 Soft rubber 5,5"



MF05 Pneumatic tyre 6"

Active wheelchair *Zenit R.*



MD50 Push handles, short



MD55 Push handles, height-adjustable/removable



ME04 Side panel carbon, with splash guard



MD03 Back support upholstery, adaptable, padded



MB06 Leg support, aluminium light, fixed



MH01 Knee lever wheel lock 4.0



MH05 Scissor wheel lock, aluminium



MZ70 Kit for transportation in a motor vehicle



MG80 Transparent



MG100 Graphics



MG101 Miracle



MG102 Network



MG103 Robots



MG104 Spectrum