

2026

ottobock.

Active wheelchair Zenit CLT.

☐ Quotation

☐ Order

Billing Address

PO. no.

Customer Name

Street

County, PC

Shipping address (if different from billing address)

Customer Name

Street

County, PC

Name Date Signature



Base Price - £5439.00

Please email all orders to
Bockuk@ottobock.com

Valid from 01.03.2026 to 31.12.26
Please ensure you are using the latest order form
v2

Standard equipment Selection (unless alternative selected)

Article number: 480F160=1_AA03_C

- Open aluminium tube frame, foldable
- Side support (for drive wheels) firmly welded
- Back support tubes firmly welded
- Foldable leg support
- Back support upholstery adjustable
- Aluminium side panels
- Modular design
- Drive wheels with quick-release axles
- Knee lever wheel lock 4.0
- Calf strap

Folding, max. user weight 100 kg

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals. Wheelchairs with CLT package cannot be returned or exchanged

The number codes contained in this order form are not intended to be used as article numbers for spare parts orders.
Please scan the QR code to access the spare parts catalogue for any spare parts orders.




Active wheelchair *Zenit CLT*.

Wheelchair measurements

- 1 Anterior seat height (aSH) = ground to top of seat upholstery
- 2 Lower leg length (LLL) = top of seat upholstery to top of foot plate (without seat cushion). Please consider the dependencies LLL and aSH.
- 3 Posterior seat height (pSH)= ground to top of seat upholstery
- 4 Back support height (BH)= top of seat tube to top of back support upholstery
- 5 Seat depth (SD) = front of seat upholstery to back support tube
- 6 Seat width (SW) = outside edge of frame tube to outside edge of frame tube

A manufacturing tolerance of 10 mm applies to all measurements and of 10 % to specified weights.



Active wheelchair Zenit CLT.

Frame tube material

☒ MA08	Aluminium	Inclusive
--	-----------	-----------

Spring support folding mechanism

☒ MA120	Standard	Inclusive
---	----------	-----------

Front frame angle (please consider LLL on page 6)

☐ MA20	75° long frame	Inclusive
☐ MA21	85° short frame	Inclusive

Front frame offset

☐ MA14	0 cm offset	Inclusive
☐ MA15	2 cm offset	Inclusive
☐ MA30	2 cm abduction (only possible till SW 42 cm)	Inclusive

Front frame extension

☐ MA23	3 cm	£108.00
----------------------------	------	---------

Seat width SW in cm

☐ 32	☐ 34	☐ 36	☐ 38	☐ 40	☐ 42	☐ 44	☐ 46
--------------------------	--------------------------	--------------------------	--------------------------	--------------------------	--------------------------	--------------------------	--------------------------

Seat depth SD in cm

☐ 36	☐ 38	☐ 40	☐ 42	☐ 44	☐ 46	☐ 48	☐ 50
--------------------------	--------------------------	--------------------------	--------------------------	--------------------------	--------------------------	--------------------------	--------------------------

Anterior seat height (aSH) in cm (please consider LLL on page 6)

☐ 38	☐ 39	☐ 40	☐ 41	☐ 42	☐ 43	☐ 44	☐ 45	☐ 46
☐ 47	☐ 48	☐ 49	☐ 50	☐ 51	☐ 52	☐ 53	☐ 54	

Frame accessories

☐ MA51	Tip-assist, right	£124.00
☐ MA52	Tip-assist, left	£124.00
☐ MA55	Anti-tipper, right	£218.00
☐ MA56	Anti-tipper, left	£218.00
☐ MA60	Frame pads	£81.00

Caster wheels

☐ MF30	4" solid rubber, black	Inclusive
☐ MF35	4" soft caster, aluminium	£162.00
☐ MF31	5" solid rubber, black	Inclusive
☐ MF36	5" soft caster, aluminium	£162.00
☐ MF04	5.5" soft caster, grey	Inclusive
☐ MF05	6" pneumatic, grey	Inclusive
☐ MF32	6" solid rubber, black	Inclusive
☐ MF37	6" soft caster, aluminium	£162.00

Caster axle

☒ MF90	Caster forks with threaded axles	Inclusive
--	----------------------------------	-----------

Posterior seat height (pSH) in cm

☐ 37	☐ 38	☐ 39	☐ 40	☐ 41	☐ 42	☐ 43
☐ 44	☐ 45	☐ 46	☐ 47	☐ 48	☐ 49	☐ 50

Active wheelchair Zenit CLT.

Drive wheel size

☐ MG02	24" (max. pSH = 48 cm)	Inclusive
☐ MG00	25"	Inclusive

Drive wheel camber

☐ 0° ☐ 2° ☐ 3° ☐ 4°

Drive wheel rim

		24"	25"	
☐ MG06	Hollow rims	●	●	Inclusive
☐ MG197	SmootherWheel dampening drive wheels, max. load 75 kg	●		£1659.00
☐ MG198	SmootherWheel dampening drive wheels, max. load 95 kg	●		£1659.00
☐ MG199	SmootherWheel dampening drive wheels, max. load 115 kg	●		£1659.00
☐ MG10	Spinergy SPOX	●		£764.00
☐ MG11	Spinergy LX	●		£843.00
☐ MG12	Spinergy LXL Blade: black spokes	●		£950.00
☐ MG12	Spinergy LXL Blade: white spokes	●		£950.00
☐ MG15	Infinity Ultralight (not possible in combination with MG 34)	●	●	£117.00
☐ MG150	Infinity Ultralight, coloured (not possible in comb. with MG 34)	●		£198.00
Colour set Infinity Ultralight/caster fork		24"	25"	
☐	Sparkle light orange	●		Inclusive
☐	Marine blue	●		Inclusive
☐	Candy red	●		Inclusive
☐	Sparkle granny-smith	●		Inclusive

Tyres

		24"	25"	
☐ MG31	Schwalbe RightRun 1"	●	●	Inclusive
☐ MG34	PU 1" (not possible in combination with MG 15/MG 150)	●		Inclusive
☐ MG341	PU 1" black (not possible in combination with MG 15/MG 150)	●		Inclusive
☐ MG37	Schwalbe Marathon Plus Evolution 1"	●	●	£76.00
☐ MG 200	Schwalbe ONE 1"	●		£190.00

Handrims

		24"	25"	
☐ MG50	Aluminium, black coated, high	●	●	Inclusive
☐ MG52	Aluminium, silver anodised, high	●	●	£16.00
☐ MG53	Aluminium, silver anodised, deep	●		£16.00
☐ MG54	Rubber-coated	●	●	£162.00
☐ MG56	Stainless steel, polished, high	●	●	£63.00
☐ MG57	Stainless steel, polished, deep	●		£63.00
☐ MG63	Curve L, anodised	●	●	£409.00
☐ MG64	Curve L, Grip for users with limited gripping function	●	●	£442.00
☐ MG65	GEKKO, ergonomic handrim	●	●	£362.00

Handrim mounting

☐ MG70	Mounting position narrow (15 mm)	Inclusive
☐ MG71	Mounting position wide (25 mm)	Inclusive

Tipping point in cm (measured parallel to the ground)

☐ 5 ☐ 6 ☐ 7 ☐ 8 ☐ 9 ☐ 10 ☐ 11 ☐ 12 ☐ 13 ☐ 14

Active wheelchair Zenit CLT.

Alternative measuring method (specification of the horizontal drive wheel position of the Avantgarde)

☐ B ☐ C ☐ D ☐ E ☐ F ☐ G

Tick the drive wheel position of your Avantgarde, we calculate the tipping point.

Quick-release axle system

☐ MG27	Quick-release axle	Inclusive
☐ MG28	Quick-release axle with Tetra release	£172.00

Drive wheel accessories

☐ MG87	High-pressure air pump, 14 bar	£122.00
☐ MG88	Airman-pump	£185.00

Wheel locks

☐ MH01	Knee lever wheel lock 4.0	Inclusive
☐ MH30	Knee lever wheel lock 4.0, long lever	Inclusive
☐ MH05	Scissor wheel lock	£151.00

Wheel lock accessories

☐ MH10	Wheel lock lever extension, fixed	£42.00
☐ MH21	Wheel lock lever extension, pluggable	£75.00

Back support height in cm

☐ 25 ☐ 27,5 ☐ 30 ☐ 32,5 ☐ 35 ☐ 37,5 ☐ 40 ☐ 42,5 ☐ 45 ☐ 47,5 ☐ 50

Back support system

☐ MD01	Back support, straight (measured 90° to the seat tube)	Inclusive
☐ MD13	Back support angle according to customer request (measured 90° to the ground)	£110.00
	☐ 0° ☐ -5° ☐ -10°	

Back support upholstery

☐ MD03	Adaptable, padded	£58.00
☐ MD05	Adaptable, breathable	Inclusive
☐ MD18	Ultra light	£56.00
☐ MD06	Without	Inclusive

Back support pad widening

☐ MD80	+ 2 cm	£22.00
☐ MD81	+ 4 cm	£22.00
☐ MD82	Without	Inclusive

Use of BAXX back support system (include order form, model-specific restrictions are possible, take into account the reduction of the seat depth by using the Baxx back support system)

☐ MD44	Baxx aluminium back support (enclosed)	£716.00
----------------------------	--	---------

Push handles

☐ MD50	Short	Inclusive
☐ MD51	Long	Inclusive
☐ MD55	Height-adjustable/removable	£174.00
☐ MD45	Height-adjustable/removable for Baxx	£218.00
☐ MD56	Fold-down	£209.00
☐ MD59	Without	Inclusive

Active wheelchair Zenit CLT.

Side panels

☐ ME02	Aluminium, with splash guard and protection against the cold	Inclusive
☐ ME04	Carbon, with splash guard	£655.00

Seat upholstery

☐ MC15	Seat upholstery with storage compartment (up to SW 38 cm = 1 compartment; from SW 40 cm = 2 compartments)	£45.00
☐ MC16	Upholstery readjustable	Inclusive
☐ MC18	Ultra light	£46.00
☐ MC06	Without	Inclusive

Seat cushions

☐ MC01	Standard (functional height 5 cm)	£43.00
☐ MC02	Cubic Foam (functional height 5 cm)	£151.00
☐ MC05	EvoLight (functional height 6 cm)	£215.00
☐ MC20	Terra (functional height 6 cm)	£200.00
☐ MC21	Terra Aquos (functional height 6 cm)	£241.00
☐ MC22	Terra Flair (functional height 6 cm)	£603.00
☐ MC23	Cloud (functional height 11 cm)	£589.00
☐ MC29	Livo (functional height 5 cm)	£90.00
☐ MC30	EvoFlair (functional height 8 cm)	£323.00


Lower leg length (without seat cushion)

Lower leg length (LLL) in cm

75° front frame: aSH must be greater or equal to LLL + 2 cm

85° front frame: aSH must be greater or equal to LLL + 4 cm

☐ 28	☐ 29	☐ 30	☐ 31	☐ 32	☐ 33
☐ 34	☐ 35	☐ 36	☐ 37	☐ 38	☐ 39
☐ 40	☐ 41	☐ 42	☐ 43	☐ 44	☐ 45
☐ 46	☐ 47	☐ 48	☐ 49	☐ 50	☐ 51



☐ MB130	Adjustable	Inclusive
☐ MB14	Firmly welded	£108.00

Leg supports (folding)

☐ MB160	Single-panel, aluminium	Inclusive
☐ MB180	Single-panel, carbon	£160.00

Leg support accessories

☐ MB70	Lock for single panel leg support	£92.00
----------------------------	-----------------------------------	--------

Other accessories

☐ MZ51	Lap belt with metal buckle	£84.00
☐ MZ59	Pocket for mobile phone	£33.00

Active wheelchair Zenit CLT.

Spoke protectors

☐ MG90	Transparent	£156.00
☐ MG100	Motiv Graphics	£224.00
☐ MG101	Motiv Miracle	£224.00
☐ MG102	Motiv Network	£224.00
☐ MG103	Motiv Robots	£224.00
☐ MG104	Motiv Spectrum	£224.00

Frame colours

☐ 001	Brillant silver	Inclusive
☐ 002	Anthracite metallic	Inclusive
☐ 003	Jet black RAL 9005	Inclusive
☐ 004	Ice blue	Inclusive
☐ 005	Sparkle light orange	Inclusive
☐ 009	Vinho sparkle	Inclusive
☐ 010	Marine blue	Inclusive
☐ 012	Sky blue	Inclusive
☐ 048	Jet black matte (structure)	Inclusive
☐ 014	Candy red	Inclusive
☐ 021	Sparkle granny-smith	Inclusive
☐ 023	White RAL 9016	Inclusive
☐ 027	Shimano matte	Inclusive


Special accessories

☐ MZ52	Tool kit	£68.00
☐ MZ58	Bag for drive wheels (for 2 drive wheels up to 28")	£118.00
☐ MZ89	Protection for quick-release-axle	£12.00
☐ MZ92	Wheelblades (Mini-skis for caster wheels)	£457.00
☐ MZ93	Backpack	£66.00

Active wheelchair *Zenit CLT*.




MG10 Spinergy SPOX



MG06 Hollow rim



MG11 Spinergy LX




MG15 Infinity Ultralight



MG31 Right Run 1"



MH05 Scissor wheel lock



MG54 Handrim rubber-coated




MG65 Gekko



MG63 Curve L (anodised)




MG64 Curve L Grip




**MF35/36/37 4/5/6" Soft
caster**



**MF30/31/32 4/5/6" Solid rub-
ber**



MF04 5,5" Soft caster



MF05 6" Pneumatic tyre

Active wheelchair *Zenit CLT*.



MD50 Push handles, short



MD55 Push handles, height-adjustable/removable



ME04 Side panel carbon, with splash guard




MD03 Back support upholstery, adaptable, padded



MH01 Knee lever wheel lock 4.0



MH05 Scissor wheel lock, aluminium



MG80 Transparent



MG100 Graphics



MG101 Miracle



MG102 Network



MG103 Robots



MG104 Spectrum