

Panthera S3

Panthera S3. Configure each section, then transcribe selected art. nos. and prices to the compilation page.



ADJUSTED FOR 	ADJUSTED BY 	TELEPHONE · E-MAIL 	DATE 	BASE PRICE £2,896
-----------------------------	----------------------------	-----------------------------------	---------------------	----------------------

How to use this form. Tick the box next to each selected option. Fill in article numbers where indicated. Transcribe selections and prices to the compilation page.

01 Chassis

☐ S3
Seat width 33–50 cm, seat depth 40 cm

£2,896

☐ S3 Short
Seat width 30–45 cm, seat depth 27–33 cm (seat width 30–33), 35 cm (seat width 36–45)

£2,896

☐ S3 Short Low
Seat width 33–39 cm, seat depth 35 cm

£3,164

☐ S3 Large
Seat width 39–50 cm, seat depth 45 cm

£3,392

☐ S3 Long
Seat width 33–50 cm, seat depth 45 cm

£3,621

☐ S3 0°
Seat width 33–50 cm, seat depth 35/37.5/40/42.5 cm

£3,621

☐ S3 Short Abd
Seat width 27–33 cm, seat depth 27–33 cm

£3,169

02 Chassis colour

☐ Black, Anodic Black

☐ Gray, Metallic Gray

☐ Optional colour
20-day delivery · Art. no.

£299

03 Seat width (BB)

Check allowed seat width at section 01 Chassis

☐ 27 cm

☐ 30 cm

☐ 33 cm

☐ 36 cm

☐ 39 cm

☐ 42 cm

☐ 45 cm

☐ 50 cm

04 Seat depth

Option only needed for S3 0°

☐ 35 cm

☐ 37.5 cm

☐ 40 cm

☐ 42.5 cm

☐ 45 cm

Panthera S3

05 Backrest height (HH)

- | | | |
|--|--|------|
| ☐ Standard
Art. no. 523BBHH VER2 · height 20/25/30/35/40/45 cm – fill in HH | ☐ Tapered
Art. no. 525BBHH VER2 | £417 |
| ☐ Curved
Art. no. 524BBHH VER2 · frame from solid backrest | ☐ Tapered 6 cm | £417 |
| ☐ Flared 3 cm | ☐ Curved / Tapered | £417 |
- Art. no. (backrest) _____

06 Rear axle

- | | |
|--|--|
| ☐ Carbon standard
Art. no. 33023BB | ☐ Carbon +15 mm
Art. no. 33023BC · seat width BB+1 · use when choosing arm support and/or increased seat height |
| ☐ Increase the seat height
Specify (mm) +_____ – other rear axle or seat height | |

07 Castor and fork

- | | | |
|-----------------------------|------------------------------|------------------------------|
| ☐ 90 mm | ☐ 120 mm | ☐ 150 mm |
|-----------------------------|------------------------------|------------------------------|
- Art. no. _____

08 Footrest

- | | |
|---|---|
| ☐ V-shape for S (titanium)
Art. no. 36001BB | ☐ V-shape extended, reinforced (titanium)
Art. no. 36111BB · std for S3 Large |
| ☐ Foldable footplates
Art. no. _____ | ☐ Small heel bar for S series
Art. no. 3600134 (width 33) |
| ☐ Heel strap for footplates | ☐ Plastic plate for footrest
Art. no. 36010BB |
| ☐ Carbon footplate
Art. no. 65040BB | |

09 Brake

- | | | |
|---|---|--------|
| ☐ Brake S3
Art. no. 4630000 | <i>incl.</i> ☐ One hand brake, right | £153 |
| ☐ One hand brake, left | ☐ Assistant brake
Not for short / low chassis | £1,326 |
- Art. no. _____

Panthera S3



10 Rear wheel (without pushrim, all prices per pair)

☐ 22" standard

☐ 24" Spinergy X, black spokes

☐ 26" Spinergy X, black spokes

☐ 25" Spinergy X, white spokes

☐ One arm drive

£716

£388

☐ 24" standard

☐ 25" Spinergy X, black spokes

☐ 24" Spinergy X, white spokes

☐ 26" Spinergy X, white spokes

☐ Rear wheel 24x2" MTB
Includes titanium pushrims

£999

Art. no. _____

11 Tyre

All prices for tyres are for pair · reTyre does not fit PU tyres

☐ Schwalbe Right Run 22", 24", 25", 26"

☐ Schwalbe Marathon Plus 24", 25", 26"

☐ PU 24" (massive)

☐ reTyre skins 25"

incl.

£98

£98

£201

☐ Rugged 24"

☐ Schwalbe One 24", 25"

☐ reTyre skins 24"

☐ reTyre skins 26"

£98

£201

£201

Art. no. _____

12 Pushrim

All prices for pushrims are for pair · DD = wheel diameter

☐ Aluminium
Art. no. 1902430 (x2 / wheelchair)

☐ Tetragrip
Art. no. 19052DDBT x 2

☐ Titanium
Art. no. 190DD11 (x2 / wheelchair)

incl.

£196

£299

☐ Paragrip
Art. no. 19051DDBT x 2

☐ Maxgrip 22", 24"
Art. no. 1903060 x 2 (60=24", 80=22") · contains latex

£196

£299

Panthera S3



13 Side guards · Arm support

Side guard

☐ Side guard S3 £256

Art. no. _____

Arm support

☐ Arm rest S3 £379

26" std for 24" wheels for Short Low
· 26" std for Assistant brake

Art. no. _____

Side guard Carbon

☐ Without fender

☐ Small fender

☐ Large fender £518

Art. no. _____

14 Push handles

Max user weight 100 kg

Push handles

☐ Height adjustable S3 £268

Art. no. 4252060 VER2

☐ Height adjustable S3 +10 cm £268

Art. no. 4252063 VER2

Push handle foldable

☐ Foldable S3 £155

Art. no. 4245000

Push bar

☐ Push bar S3 £304

Art. no. 42400BB VER2

15 Accessories

Anti-tip

☐ Anti-tip std seat height £308

Art. no. 4353000

☐ Anti-tip S3 Short Low £308

Art. no. 4355000

☐ Anti-tip S3 Large £308

Art. no. 4353080

☐ Other seat height

Art. no. _____

Cushion · Wedge cushion

☐ Cushion 2.5 cm £116

Art. no. 572BB40

☐ Cushion 5 cm £116

Art. no. 575BB40

☐ Wedge cushion £55

Art. no. 57710BB

Panthera S3



16 Spoke protectors

Plain – £155 each

- ☐ 24" logo
Art. no. 4500000
- ☐ 22" logo
Art. no. 4500100
- ☐ 24" no logo
Art. no. 4500011
- ☐ 22" no logo
Art. no. 4500111

Themed – £204 each

- ☐ 22" Splash
Art. no. 4500322
- ☐ 24" Splash
Art. no. 4500324
- ☐ 22" Car rim
Art. no. 4500122
- ☐ 24" Car rim
Art. no. 4500124
- ☐ 22" Wild animals
Art. no. 4500422
- ☐ 24" Wild animals
Art. no. 4500424
- ☐ 22" Smiley
Art. no. 4500222
- ☐ 24" Smiley
Art. no. 4500224

17 Other accessories

- | | | | |
|--|------|--|------|
| ☐ Tube protection, front
Art. no. 9900251 | £16 | ☐ Tube protection, rear
Art. no. 9900251 | £16 |
| ☐ Crutch holder S3
Art. no. 9900301 (x1 / crutch) | £43 | ☐ Backrest wedge
Art. no. 4710012 | £94 |
| ☐ Table S3
Art. no. 6810020 · only together with arm support | £245 | ☐ FreeWheel
Art. no. 2600102 · Height 6–14 cm · Choose art. no. matching footrest height | £901 |

18 Set-up adjustments

Backrest angle

- ☐ Note measure (A) for desired backrest angle. 14–28 mm.
Standard = 20 mm
A = _____ mm

Rear axle adjustment

- ☐ Note measure (A) for desired balance. Standard = 60 mm (51 mm at raised rear axle) · 40 mm for S3 Short/Short Low · 31 mm for Short ABD (6° camber)
A = _____ mm

Footrest height – foot plates

- ☐ Note measure (A) and (B) for desired footrest
A = _____ mm · B = _____ mm

Panthera S3 — Compilation

Transcribe selected art. nos. and prices from the previous pages.



DESCRIPTION	SELECTION	ART. NO.	PRICE
			WHEELCHAIR MODEL
Wheelchair model	Panthera S3		£2,896
			CHASSIS COLOUR
Chassis colour			
			SEAT WIDTH / DEPTH
Seat width (BB)			
Seat depth (S3 0° only)			
			BACKREST
Backrest height (HH)			
			REAR AXLE
Rear axle			
			CASTOR AND FORK
Castor and fork			
			FOOTREST
Footrest			
			BRAKE
Brake			
			REAR WHEEL
Rear wheel			
			TYRE
Tyre			
			PUSHRIM
Pushrim			
			SIDE GUARDS / ARM SUPPORT
Side guard / Arm support / Carbon			
			PUSH HANDLES
Push handles			
			ACCESSORIES
Anti-tip			
Cushion / Wedge cushion			
			SPOKE PROTECTORS
Spoke protectors			
			OTHER ACCESSORIES
Other accessories			
			SET-UP ADJUSTMENTS
Set-up adjustments			

Total · GBP, ex. VAT

£

Customer signature

By signing, the customer confirms the configuration above is correct.

Signature

Date

Dealer / Prescriber

Permobil reference, name & date.

Signature

Date