

TiLite Pilot Order Form

Effective January 7, 2022

CARBON RIGID, GROWABLE FRAME

Page 1 of 4

DEALER INFORMATION		RTS INFORMATION	
Dealer:	Account #:	RTS:	Phone:
PO:	Ordered By:	Email:	
Email:		SHIP TO INFORMATION	
Mark For:		Name:	Ship Via:
TiLite Sales Rep:		Address:	
TiLite Information:	Promo Code:	City:	State: Zip:
Ship Date:	Order #:	Phone:	Fax:

Please Email to Permobil at info@permobil.co

FRAME STYLE

TPFS1	<input type="checkbox"/> TiLite Pilot INCLUDING Growth Kit (Full Adjustability. Seat Width A grows 2" and Seat Depth B grows 3".) 120 Pound Weight Limit. Price includes one-time growth kit and Transit Tie-Down option. Additional growth kits sold separately.).....	£2,423	E1235
TPFS2	<input type="checkbox"/> TiLite Pilot OMIT Growth Kit	£2,252	E1235
TPMAXSD1	<input type="checkbox"/> Permobil Smart Drive MX2+ (Not available with 8" width, rear seat heights less than 13.5", or 20" rear wheels.).....	£4,795	E0986
TPMAXSD2	PushTracker.....	STD	
	Size: <input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large		
TPMAXSD4	<input type="checkbox"/> Carrying Bag.....	£65	
TPMAXSD5	<input type="checkbox"/> Thumb Throttle Buttons.....	£364	

TRANSIT TIE-DOWN OPTION (Refer to the Transit Tie-Down Option Brochure for important information regarding the Transit Tie-Down options.)

TPWTB	<input type="checkbox"/> Transport Brackets (Consists of four Tie-Down Brackets.).....	NCO
TPTTDO	<input type="checkbox"/> Transit Tie-Down (WC-19 Tie-Down Brackets require the TiLite back and head support. 95 lb weight limit.).....	£186

FRAME FINISHES**PAINT COLOR**

☐ TPF29 Neon Pink ☐ TPF14 Ocean Blue Metallic ☐ TPF30 Acid Green Metallic ☐ TPF25 Electric Plum.....STD



ANODIZE COLOR

☐ TPAN01 Black ☐ TPAN05 Sapphire ☐ TPAN03 Red Pepper ☐ TPAN08 Ultraviolet.....STD

FRAME MEASUREMENTS**FRONT SEAT TO FLOOR HEIGHT, REAR SEAT TO FLOOR HEIGHT, CENTER OF GRAVITY AND SEAT TO FOOTREST**

Please see the table below for available configuration adjustment ranges. Center of gravity will be set at the factory. See documentation shipped with chair for detailed information on adjustments.

SEAT DEPTH	MEASUREMENT:	20" REAR WHEEL AVAILABLE RANGE (IN):	22" REAR WHEEL AVAILABLE RANGE (IN):	24" REAR WHEEL AVAILABLE RANGE (IN):
9"	FRONT SEAT HEIGHT	12.5-15	N/A	N/A
	REAR SEAT HEIGHT	12-13.5	N/A	N/A
	CENTER OF GRAVITY	0.5-1.5	N/A	N/A
10"	FRONT SEAT HEIGHT	12-15	N/A	N/A
	REAR SEAT HEIGHT	12-13.5	N/A	N/A
	CENTER OF GRAVITY	0.5-2.25	N/A	N/A
11"	FRONT SEAT HEIGHT	12-15	13.5-16.5	14-16.5
	REAR SEAT HEIGHT	12-13.5	13-15	13.5-16
	CENTER OF GRAVITY	0.5-3	0.5-2.25	0.5-2.25
12"	FRONT SEAT HEIGHT	12-15	13.5-16	14-16.5
	REAR SEAT HEIGHT	12-13.5	13-15	13.5-16
	CENTER OF GRAVITY	0.5-3.5	0.5-3.25	0.5-3.25
13"	FRONT SEAT HEIGHT	N/A	13-16	13.5-16.5
	REAR SEAT HEIGHT	N/A	13-15	13.5-16
	CENTER OF GRAVITY	N/A	0.5-3.5	0.5-3.5
14"	FRONT SEAT HEIGHT	N/A	13-16	13.5-16.5
	REAR SEAT HEIGHT	N/A	13-15	13.5-16
	CENTER OF GRAVITY	N/A	0.5-3.5	0.5-3.5
15"	FRONT SEAT HEIGHT	N/A	13-16	13.5-16.5
	REAR SEAT HEIGHT	N/A	13-15	13.5-15.5
	CENTER OF GRAVITY	N/A	1.5-3.5	1.5-3.5



TiLite Pilot Order Form

CARBON RIGID, GROWABLE FRAME

Page 2 of 4

FRAME MEASUREMENTS

SEAT WIDTH A Measure inside width from side panel to side panel. Seat width grows 2" if selected width is between 8" and 12". The 13" width will grow a maximum of 1", and the 14" width will not grow any wider. Refer to chart for seat depth availability.

TPSW1 ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14".....STD

SEAT DEPTH RANGE: The seat depth can be grown within the selected range.

☐ 9"-12" ☐ 11"-15".....STD

SEAT DEPTH B Measure from front of seat pan to front of back support cushion.

TPSD1 ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14" ☐ 15".....STD

FRONT SEAT HEIGHT C Measure from the floor to the top edge of the front of the seat pan. Adjustable in 0.5" increments.

TPFSH1 ☐ 12" ☐ 12.5" ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5".....STD

REAR SEAT HEIGHT D Measure from the floor to the top edge of the of the seat pan at the intersection of the front of the back support pad and the seat pan. Adjustable in 0.5" increments.

TPRSH1 ☐ 12" ☐ 12.5" ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5".....STD

FOOTREST TYPE

TPFTR3 ☐ Angle Adjustable Footrest.....STD

TPFTR24 ☐ Flip-Back Depth and Angle Adjustable Footrest (Not available on 8", 9" or 14" seat width).....£160 K0040

SEAT TO FOOTREST E Measure from the top of the footplate to the front of the seat pan. For an angle adjustable footrest, this measurement must be at least 2" less than front seat height. For a Flip-Back Depth and Angle Adjustable footrest, this measurement must be at least 3" less than front seat height.

☐ 3" ☐ 3.5" ☐ 4" ☐ 4.5" ☐ 5" ☐ 5.5" ☐ 6" ☐ 6.5" ☐ 7" ☐ 7.5" ☐ 8" ☐ 8.5" ☐ 9" ☐ 9.5" ☐ 10"
☐ 10.5" ☐ 11" ☐ 11.5" ☐ 12" ☐ 12.5" ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5".....STD

SEAT BACK TYPE Measure from top of back support to seat pan. Back support shell height is 9.5" with scapular cutouts.

TPSB1 ☐ Standard Non-Folding Adjustable-Height Rigid Back with Pad (11" to 13" height range).....STD

TPSB2 ☐ Contoured Non-Folding Adjustable-Height Rigid Back with Pad (11" to 13" height range. PDAC verification pending.)...£335 E2613

TPSB3 ☐ Aluminum Canes for Aftermarket back support, no Pad (Not available with height adjustable fenders.).....NCO

HEAD SUPPORT (The Comfort head support is WC-20 certified on the Acta-Back and Acta-Embrace. The head support is WC-19 certified when used with the TiLite backrest.)

TPHDR1 Comfort BodiLink head support with GlideWear cover (Uses piano taper midline hardware. Available on Acta-Back, Acta-Embrace and TiLite rigid back supports).....£298

Width: ☐ 6" ☐ 8"

E0955,
E1028

PUSH HANDLES (Height adjustable range is 33" - 40" as measured from the floor with 24" (540) wheels.)

TPSB6 ☐ Bolt-On Height-Adjustable Removable Single Push Handle.....£134

TPSB7 ☐ Bolt-On Height-Adjustable Removable Dual Push Handles.....£164

CAMBER / CAMBER TUBES (Aluminum Camber tube standard. Anodize color will match selection above.)

☐ 2° (Minimum 1.5" rear wheel spacing) ☐ 4° (Minimum 1.25" rear wheel spacing) ☐ 6° (Minimum 1" rear wheel spacing)

☐ 8° ☐ 12°.....STD

REAR WHEEL SPACING Measure from outside of side panel to inside of rear tire at top of wheel. (Not all rear wheel spacing is available with all accessories.)

TPRWS1 ☐ 0.75" ☐ 1" ☐ 1.25" ☐ 1.5" ☐ 1.75" ☐ 2".....STD

TPCBR12 ☐ Adjustable Rear Wheel Spacing Camber Plugs.....STD

TPCBR9 ☐ Additional Spacers to Increase Rear Wheel Spacing.....NCO

FRONT WHEELS (TiLite Standard Fork style will be provided.)

TPFW14 ☐ 4" X 0.75" Plastic Wheel w/Poly Tire.....STD

TPFW3 ☐ 4" X 0.75" Light-Up Micro-Caster.....£60

TPFW25 ☐ 4" X 0.75" LiteSpeed Billet Aluminum Wheel w/Poly Tire (Black Hub Only).....£112

TPFW26 ☐ 4" X 1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire.....£70

TPFW27 ☐ 4" X 1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire (Black Hub Only).....£149

TiLite Pilot Order Form

CARBON RIGID, GROWABLE FRAME

Page 3 of 4

CUSTOM CONFIGURATION OPTIONS

REAR WHEELS

- TPRW14 TiLite SHADOW™ **STD**
☐ 20"(451) *Only available with seat depth range 9"-12".*
☐ 22"(501) *Only available with seat depth range 11"-15".*
☐ 24"(540) *Only available with seat depth range 11"-15".*
- TPRW7 Spinergy SPOX ☐ Black ☐ Red ☐ Yellow ☐ Blue ☐ White **£431**
☐ 22"(501) - Std. Hub *Only available with seat depth range 11"-15".*
☐ 24"(540) - Std. Hub *Only available with seat depth range 11"-15".*
- TPRW8 Spinergy LX ☐ Black ☐ Red ☐ Yellow ☐ Blue ☐ White ☐ Go Green ☐ Bright Orange ☐ Flamingo Pink ☐ Purple **£628**
☐ 22"(501) - Std. Hub *Only available with seat depth range 11"-15".*
☐ 24"(540) - Std. Hub *Only available with seat depth range 11"-15".*
- TPRW16 Spinergy Carbon Blade **£953**
☐ 24"(540) *Not available with solid tires or airless inserts. Only available with seat depth range 11"-15".*
- TPRW17 X-Core Mag *(Not available with Short Tabs, Quad Release Axles, or Omit Handrims Option.)* ☐ Black ☐ Red ☐ Pink ☐ White **£82**
☐ 20"(451) *Only available with seat depth range 9"-12".*
☐ 22"(501) *Only available with seat depth range 11"-15".*
☐ 24"(540) *Only available with seat depth range 11"-15".*

AXLES

- TPAXL1 ☐ Stainless Quick Release **STD**
 TPAXL3 ☐ Quad Release *(Not available with Mag Rear Wheels.)* **£115**

PNEUMATIC REAR TIRES

- TPRTR1 ☐ Treaded **STD**
 TPRTR2 ☐ Airless Inserts **£70**
 TPRTR4 ☐ Primo **NCO**
 TPRTR30 ☐ Primo Sentinel *(Available on all wheel sizes, 20" (451), 22" (501), and 24" (540)).* **£97**
 TPRTR8 ☐ Schwalbe Marathon *(Only available on 24" (540)).* **£78**

E2213

HANDRIMS

- TPRIM1 ☐ Aluminum – Silver Anodized **STD**
 TPRIM21 ☐ Aluminum – Black Anodized *(Not available with 20" (451) Rear Wheels.)* **NCO**

Handrims by OUT-FRONT®

The Natural-Fit® LT *(Reduced Weight, Smaller Profile) (Not available with 20"(451) Rear Wheels.)*

- TPRIM17 ☐ Standard Grip *(Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.)* **£286**
 TPRIM18 ☐ Super Grip *(Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.)* **£328**
 TPRIM19 ☐ Without Thumb Grip **£268**
 TPRIM11 ☐ The Surge LT *(All-In-One Oval with Gription™ Strip, Smaller Profile) (Not available with 20"(451) Rear Wheels.)* **£272**
 TPRIM20 ☐ Q-Grip *(Exclusive High-Friction, Anti-Slip Coating) (Not available with 20"(451) Rear Wheels.)* **£249**

Additional Handrim Options

- TPRIM5 ☐ Plastic Coated – Black *(Not available with Mag Rear Wheels.)* **£97**
 TPRIM22 ☐ Omit Handrims **NCO**

Tab Length

- TPRIM7 ☐ Long Tabs **STD**
 TPRIM8 ☐ Short Tabs *(Not available with Mag Rear Wheels or The Natural-Fit® Handrims with Thumb Grip.)* **£15**

WHEEL LOCKS

- TPWLK30 ☐ Compact Push to Lock **NCO**
 TPWLK31 ☐ Compact Pull to Lock **NCO**

Mounting Location:


- ☐ Standard Location **F** ☐ Forward Mount Location **G** *(Mounted on forward front of side panel.)*

SOLID SEAT PAN

- TPSP1 Solid Seat Pan *(Required Option)* **£112** **E2231**

SEAT CUSHION

- TPCSH2 ☐ Comfort Curve Cushion *(Stretch-Air cover.)* **£60** **E2601**
 SCSH11 ☐ ROHO Single Compartment Low Profile Cushion **£482** **E2624**



TiLite Pilot Order Form

CARBON RIGID, GROWABLE FRAME

Page 4 of 4

CUSTOM CONFIGURATION OPTIONS

ROHO OPTIONS: *(Made to order options have a 2-3 week lead time)*

		SEAT DEPTH (B)					
		10"	11"	12"	13"	14"	15"
SEAT WIDTH (A)	8"	Made to Order	Made to Order	Made to Order	Made to Order	Made to Order	Made to Order
	9"	Made to Order	Made to Order	Made to Order	Made to Order	Made to Order	Made to Order
	10"	Made to Order	Made to Order	Made to Order	Made to Order	Made to Order	Made to Order
	11"	Made to Order	Stocked	Stocked	Stocked	Made to Order	Made to Order
	12"	Made to Order	Stocked	Stocked	Stocked	Made to Order	Made to Order
	13"	Made to Order	Made to Order	Made to Order	Stocked	Stocked	Stocked
	14"	Made to Order	Made to Order	Made to Order	Stocked	Stocked	Stocked

SIDE GUARDS: *(Side guards required)*

TPSDG30 ☐ Built in Side Guard with Height Adjustable Fender *(Not available with 20" rear wheels, .)*.....£335
 TPSDG31 ☐ Built in Side Guard WITHOUT Fender.....STD

ARMRESTS *(Armrests attach to built in Side Guards.)*

TPARM2 ☐ Height Adjustable Desk Arm w/ 10" Pad *(Works with Fendered Side Guard Attachment Option.)*.....£249 E0973
 TPARM31 ☐ Height Adjustable Armrest.....STD

UPPER BODY POSITIONING

TPUBP1 ☐ Bodypoint® PivotFit™ Dynamic shoulder harness with cam buckle mounts, rear-pull£157 E0960
☐ Extra Small ☐ Small ☐ Medium
 TPUBP2 ☐ Bodypoint® Monoflex™ chest support.....£149 E0960
☐ Extra Small ☐ Small

PELVIC POSITIONING BELTS

TPBLT1 ☐ Velcro Adjustable, black nylon webbing seat belt.....£40 E0978
 TPBLT30 ☐ Bodypoint® 1" wide two-point padded hip belt with swivel buckle.....£134 E0978
 TPBLT15 ☐ Bodypoint® 1" wide Evoflex® belt with swivel buckle.....£137 E0978
 TPBLT16 ☐ Bodypoint® 1" wide Evoflex® belt with small push-button buckle.....£137 E0978
 TPBLT32 ☐ Bodypoint® Sub-ASIS pads.....£97
☐ 100mm long ☐ 115mm long

LOWER BODY POSITIONING

TPCLF15 ☐ Bodypoint® Ankle Huggers®.....£167 E0951
☐ Extra Small ☐ Small

ANTI-TIPS

TPTIP1 ☐ Rear Aluminum - Black.....£115 E0971
 TPTIP7 ☐ Single Anti-tip.....£93 E0971
 Mounting Location: ☐ Right Side ☐ Left Side

BACKPACK

TPBCK3 ☐ Black Permobil Backpack.....£70

SPOKE GUARDS

TPGRD30 ☐ Clear polycarbonate spoke guard with TiLite Pilot sticker pack *(Not available with Mag Rear Wheels.)*.....£97 K0065

ACCESSORIES

TPCH1 ☐ Crutch Holder.....£108 E2207

GENERAL

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy.

Get the latest order forms at permobilus.com. Prices, standard features, options and accessories are subject to change without notice. All rights reserved. Trademarks are identified by symbols™ and ®. TiLite, TiFit, the TiLite logo, Ti, the Ti logo, The Titanium Experts, and TiShaft are registered trademarks of TiSport, LLC and are licensed to Permobil. Permobil® is a registered trademarks of Permobil. Adobe®, Acrobat®, Acrobat Professional® and Reader® are registered trademarks of Adobe Systems Inc. Bodypoint® is a registered trademark of Bodypoint Inc. Kenda® is a registered trademark of Kenda Rubber Industrial Co. Ltd. Ride Designs® and Corbac® are registered trademarks of Aspen Seating, LLC. ROHO®, Low Profile® and Quadro Select® are registered trademarks of ROHO Inc. Schwalbe® is a registered trademark of Ralf Bohle GmbH. Supracor®, Stimulite® and Corbee® are registered trademarks of Supracor, Inc. The Natural-Fit® and Out-Front® are registered trademarks of Three Rivers Holdings, LLC. Varilite® is a registered trademark of Cascade Designs, Inc.

NOTES

QUOTE THIS CHAIR

If you are using Adobe Acrobat, this order form may be converted into a quote request or into an order by clicking one of the buttons (this feature is available exclusively to Permobil dealers and distributors). A Quote or an Order Confirmation will be sent to you within 24 hours.

ORDER THIS CHAIR