

GO GO^e



Matte Colours Available:



Travelling Has Never Been Easier or Smoother

With the GoGo-e, travelling has never been easier! This 4-wheel scooter has everything you need for a safe and convenient trip. With two 12 Ah batteries the GoGo-e will take you further at a good speed of 4.1 mph. All with a modest price point to suit anyones budget!

Specifications	
Scooter	
Weight Capacity	18 st 13 lbs (120 kg)
Overall Length	1000 mm (39")
Overall Width	490 mm (19")
Drive Tyres	203 mm (8") Front and Rear
Ground Clearance	78 mm (3")
Colour(s)	Matte Grey, Black Matte, Matte Red
Overall Height	867 mm (34")
Mobility	
Turning Radius*	1250 mm (49")
Reversing Width (K-Turn)*	1067 mm (42")
Max Speed	4 mph (6.4 km/h)
Drive System	Rear - Wheel Drive, Single Motor, 24v DC Motor
Dual Braking System	Regenerative and Electromechanical
Battery	
Battery Requirements	2 x 12 Ah
Battery Charger	2-Amp Off-Board
Battery Weight	2.5kg
Total Weight Without Batteries*	37 kg
*This specification is for the prototype only and could be subject to change. Speak to your local representative for more details	

Features



8" Front and Rear Solid Wheels



Basket to Hold Your Everyday Essentials



Easy to use and Steer Controls



Rev 1 - 17/02/2023

