

Date.....
User.....
End user weight (kg).....
Language user manual.....

UK (English)

Company stamp

Rigid frame wheelchair

£ 3650,00

Retail price - VAT not included

Responsible.....

Keep attention to all measures.

JOKER ENERGY is a completely welded wheelchair and no modifications can be carried out after its manufacturing.

Wheelchair measures

Seat width (LS) ☐ 33 ☐ 36 ☐ 39 ☐ 42 ☐ 45 ☐ 48
The seat width is measured from the outer sides of the frame
For seat enlargement see section enlargement of mudguard
Seat depth (PS) ☐ 35 ☐ 37,5 ☐ 40 ☐ 42,5 ☐ 45 ☐ 47,5 ☐ 50
The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.

Footplate distance (DP) (cm)

(Real distance between seat and footplate) - IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)

Setting ☒ Standard ☐ Extremely prudenti ☐ Prudential ☐ Active ☐ Extreme

Backrest Height (TS) Std.Value 33,5 (min 26-max 47) adjustment steps 1.5cm

Rear wheel camber ☒ 2° ☐ 0° ☐ 4° ☐ 8° ☐ 12° Special cambers (00989 - 135£)
See reference tables N°4


Front height STD 49

Rear height STD 42


Frame

☒ Made-to-measure aluminium welded frame
Backrest with rake adjustable from -4° to +16°☐ Made-to-measure aluminium welded frame with fixed and welded backrest
Not available: bright silver!!

08932 802,00

(.....) to the ground ☐☐ (.....) to the ground


Frame Angle

☒ Frame Angle STD 95°☐ Frame Angle :☐ Frame Angle 90° to the floor
Exclude the installation of casters 5"

Frame's depth

Front frame extension

☒ 0 cm ☐ 0,5 cm ☐ 1 cm ☐ 1,5 cm ☐ 2 cm ☐ 2,5 cm
☐ 3 cm ☐ 3,5 cm ☐ 4 cm ☐ 4,5 cm ☐ 5 cm



Footplate

☒ Tubular aluminium footplate - WELDED, not adjustable in height and angle☐ with plastic plate☐ with carbon fiber plate
09193 135,00☐ Aluminium footplate EVOLUTION - (adjustable distance +/- 1,5cm)
rake and depth adjustable☒ Internal ☐ 2/3 Internal ☐ 2/3 External ☐ External

04456 270,00

☐ One piece Carbon fibre footplate - (adjustable distance +/- 1,5cm)
Rake adjustable☒ Internal ☐ External

04111 466,00

Date Order User Address

Frame colour (without extra charge)

☆ Bright silver (Not available with welded backrest)

☐ POLAR white☐ Opaque black☐ Matt orange☐ Tekno anthracite☐ Matt Midnight blue☐ Matt red☐ Matt grey**Seat**

☆ Seat cloth with pocket (under seat)

☐ Seat with tension straps

08883 240,00

☐ Carbon rigid seat

08336 509,00

☐ Pocket (under seat) [25x11x4,5]

08281 79,00

Backrest adjustment and fixing system - Not available with welded backrest

☆ Backrest with double locking (upright and folded)

with rake adjustable from -4° to +16°

☐ Backrest with upright locking system ONLY (no locking when folded) with rake adjustable from -4° to +16°☐ Fixed back(non-folding back)

with rake adjustable from -4° to +16°

Backrest

☆ Adjustable velcro straps and padded transpiring upholstery

☆ Aluminium backrest

☐ Titanium backrest

Not available with welded backrest

☐ Carbon fibre backrest

Not available with welded backrest

06811 203,00

04389 802,00



☐ **AIR** Integrated posture backrest (NOT REMOVABLE)

Carbon shell - Lateral retaining 7,5cm

Depth adjustment ±1cm. - Angle adjustment ±10°

not available for seat width 48

1.121,00

H 28		H 34		H 40	
Backrest Height (TS) min 34		Backrest Height (TS) min 38		Backrest Height (TS) min 44	
(.....)		(.....)		(.....)	
Lateral retaining wing A=high - R=regular					
A	R	A	R	A	R
					
☐	☐	☐	☐	☐	☐

Installation on:☆ Aluminium backrest tubes
WITHOUT backrest rigidizer☐ Aluminium backrest☐ Titanium backrest

Not available with welded backrest

☐ Carbon fibre backrest

Not available with welded backrest

06811 203,00

04389 802,00

Handle tube

☆ Tubes with short handle

☐ Tubes with standard handle☐ SPORT tubes Without push handles☐ Carbon fibre tubes with short handle

Lightest configuration

04390 158,00

Not available with: Posture backrest PHYSIO - Posture backrest A I R

☐ Folding handles

Minimum backrest height 30,5cm

Not available with: Push handles for companion

02958 270,00

Colour set: includes rear wheel hub and rim and front forks

Only standard Progeo 22"/24" wheels have coloured hubs and rims

☆ Black 22" 24" 25" 26"

☐ Bright silver 22" 24"

00833 135,00

☐ Bright lime 22" 24"

00833 135,00

☐ Bright blue 22" 24"

00833 135,00

☐ Bright red 22" 24"

00833 135,00

☐ Bright fuchsia 22" 24"

00833 135,00

☐ POLAR white 22" 24"

00833 135,00

☐ Tekno anthracite 22" 24"

00833 135,00

☐ Matt Midnight blue 22" 24"

00833 135,00

☐ Matt grey 22" 24"

00833 135,00

☐ Matt red 22" 24"

00833 135,00

Date	Order	User	Address
.....			

Casters

- | | | | |
|--|---|-------|-------|
| ☆ 4" Soft | | | |
| ☐ 3" Sport | <i>Lightest configuration</i> | | |
| ☐ 5" Soft | Not available with frame angle 90° to the floor | | |
| ○ 3" Active Soft (wide) - Aluminium alloy rim | | 06824 | 95,00 |
| ○ 4" Active Sport (narrow) - Aluminium alloy rim | | 01398 | 95,00 |
| ○ 4" Active Soft (wide) - Aluminium alloy rim | | 03943 | 95,00 |
| ○ 5" Active Soft (wide) - Aluminium alloy rim | Not available with frame angle 90° to the floor | 07518 | 95,00 |
| ○ 5" Pneumatic | Not available with frame angle 90° to the floor | 10844 | 68,00 |

Rear wheels

- | | | | | |
|--|-------------------------------------|------------------------|-------|-----------------|
| ☆ 24" Ultra-light Progeo - Standard! | | | | |
| □ 22" Ultra-light Progeo - Standard! | | | | |
| □ 25" Ultra-light Progeo - Standard! | (colour: black) | | | |
| □ 26" Progeo - Standard! | (colour: black) | | | |
| ○ 24" ULTRA-LIGHT PROGEO KEY SPORT | <i>Lightest configuration</i> | | 04383 | 695,00 |
| ○ 24" Ultra-light "Spinergy LX" | | | 01685 | 695,00 |
| ○ 24" Ultra-light "Spinergy LXL" | □ Black hub and white spokes | □ Hub and spokes black | 07846 | 695,00 |
| ○ 24" Ultra-light "Spinergy Flex rim" | Integrated aluminium/rubber handrim | | 02270 | 1.150,00 |
| ○ 25" Ultra-light "Spinergy LX" | | | 02101 | 695,00 |
| ○ 25" Ultra-light "Spinergy LXL" | □ Black hub and white spokes | | 07847 | 695,00 |
| ○ 26" Ultra-light "Spinergy LX" | | | 02102 | 695,00 |
| ○ 24" Off-road (aluminium handrim black) | | | 01392 | 270,00 |

Rear wheel tyre

- | | | | | |
|---|-------------------|-------------------------------|-------|--------|
| ☆ High pressure tyre | (22"-24"-25"-26") | Lightest configuration | | |
| ☐ High pressure profile tyre | (24") | | | |
| ○ Marathon Plus | (24"-25") | with anti-puncture protection | 01669 | 95,00 |
| ○ Anti-puncture - Colour Grey | (22") | | 04181 | 105,00 |
| ○ Anti-puncture - Colour Black | (24") | | 11083 | 95,00 |
| Not available with Rear wheels : 24" ULTRA-LIGHT PROGEO KEY SPORT, 24" Ultra-light "Spinergy LXL" | | | | |

Handrim

Wheel distance ★ 3 cm

□ 4 cm

- | | | | | |
|---|-------------------|---|-------|--------|
| ☆ Aluminium | (22"-24"-25"-26") | <i>Lightest configuration</i> | | |
| ○ Aluminium black SOFT-TOUCH | (24"-25") | (24" handrim, on 25" wheel is installed with spacers) | 09170 | 160,00 |
| ○ Chromed Steel | (24"-25") | (24" handrim, on 25" wheel is installed with spacers) | 01684 | 160,00 |
| ○ Titanium | (24"-25") | (24" handrim, on 25" wheel is installed with spacers) | 00016 | 290,00 |
| ○ Ultra Grepp | (24"-25") | (24" handrim, on 25" wheel is installed with spacers) | 00017 | 290,00 |
| ○ TOP FRICTION | (24"-25") | (24" handrim, on 25" wheel is installed with spacers) | 11023 | 290,00 |
| ○ Oval countour with a rubber insert | | | | |
| ○ NOVA H | (24"-25") | (24" handrim, on 25" wheel is installed with spacers) | 11025 | 290,00 |
| ○ Ergonomic black profile | | | | |
| ○ NOVA D | (24"-25") | (24" handrim, on 25" wheel is installed with spacers) | 11026 | 290,00 |
| ○ Ergonomic profile with rubber coating (ideal for quadriplegics) | | | | |
| ○ Silicon cover for standard aluminium handrim | | Available for measures (22"- 24" - 26") | 04050 | 172,00 |

Wheel locks

- | | | | | |
|---|-------------------------------|---|------------------------------|------------------------------|
| ☒ Bent lever push to lock
☐ Straight lever push to lock
☐ Straight lever pull to lock | | 

 | | |
| ☐ Sport brake | <i>Lightest configuration</i> |  | 06428 | 135,00 |
| ☐ Brake with extension lever push to lock
☐ Brake with extension lever pull to lock | | 
 | 01949
03372 | 68,00
68,00 |

Date	Order	User	Address
------------	-------------	------------	---------------

Mudguard

☆ Carbon fibre mudguard

☐ Carbon fibre Straight sideguard ☆ Aligned wheel ☐ Under wheel cm ☐ Above wheel cm

Mudguard enlargement

☆ + 0 cm ☐ + 1 cm ☐ + 2 cm 12648

☐ Removable carbon fibre mudguard 04109 270,00

☐ Removable Carbon fibre straight sideguard ☆ Aligned wheel ☐ Under wheel cm ☐ Above wheel cm 04109 270,00

☐ **Welded straight mudguard** Available only with made-to-measure aluminium welded frame with fixed and welded back 08933 398,00

☆ Aligned wheel ☐ Under wheel cm ☐ Above wheel cm

Armrests

☐ Tubular armrests (Swing away, removable and height adjustable armrest)  00694 270,00

Accessories

☐ Pad for front frame protectors 07058 56,00

☐ Backrest-seat connection cloth 10241 25,00

☐ Removable storage bag (with quick coupling) 09811 135,00

Optional 1

☐ Set-up for post-sale optional assembly

☐ Transit wheels  03290 158,00

☐ Crutches support ☐ Right 00838 79,00

☐ Left 00837 79,00



☐ Anti-tipper with rapid hookup (inside rotation to reduce the encumbrance)


☐ Right 06821 145,00

☐ Left 06822 145,00



☐ Tipping aid ☐ Right 00025 79,00

☐ Left 01046 79,00



Optional 2

☐ Push handles for companion (Removable and height adjustable)  03645 189,00

☐ Tetraclip (Quick release axle for quadriplegic) 03947 189,00

☐ Safety pelvic belt 03972 68,00

☐ Harness 04103 145,00

☐ Lumbar support (to insert in the backrest upholstery) 00013 43,00

☐ Swing away lateral support Not available Harness PHYSIO - Harness A I R

Prevents backrest folding ☐ Right 02097 270,00

☐ Left 02098 270,00



☐ Spoke guards Black (22"-24") 01043 158,00

☐ Spoke guards Transparent (22"-24") 00835 158,00

☐ Spoke guards Progeo design (22"-24") 00019 189,00

Additional Kit

☐ 24" Off-road rear wheels Off-Road tyre - aluminium handrim (black colour) 03378 534,00

☐ 24" FAT rear wheels FAT tyre 3" - aluminium handrim (black colour) 09357 895,00

Date.....

Order.....

User.....

Address.....

PHYSIO - Posture backrest
(Std. Backrest upholstery not provided)

Removable with quick release
Suitable for tubes with diameter 18÷26mm
Depth adjustment (4cm)
Angle adjustment (±18°)

☐ ~~AIR~~ Carbon shell
not available for seat width 48

BACKREST HEIGHT		28		34		40	
Shell's shape 		Lateral retaining wing A=high - R=regular					
		A	R	A	R	A	R
							

1.121,00

☐ ~~POSTURE~~ Aluminum shell

BACKREST HEIGHT	34	40	52
Shell's shape 	Lateral retaining wing R=regular		
	R	R	R

1.085,00

Configuration DEEP with side support PAD (height adjustable)

Alluminum lateral with padded, breathable upholstery.



Small (h 13cm)



☐ Right
10272

49,00

☐ Left
10272

49,00

Large (h 20cm)



☐ Right
10273


49,00

☐ Left
10273

49,00

SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve
Alluminum lateral with padded, breathable upholstery.



☐ Large (18x12)

05647 196,00

☐ Small (14x9)

05643 196,00

☐ Lateral support extension

Quantity (number pcs).....
10271 74,00

☐ Lateral support extension

Quantity (number pcs).....

☐ Small (14x9)

05642 196,00

☐ Large (18x12)

05646 196,00

☐ Headrest adjustable - Double ball mounting system

Available only with backrest Physio Posture



☐ Large (34x14)
☐ Small (28x14)

05648 385,00

Permobil LimitedUnit 7, Polaris Centre, 41 Brownfields,Welwyn Garden City, Herts, AL7 1ANUnited Kingdom - Phone: +44 1484 722888



DateOrderUserAddress

1	A Footplate width		
DP	35÷39	40÷45	46÷50
LS			
33	21	20	19
36	24	23	22
39	26	25	24
42	29	28	27
45	32	31	30
48	35	34	33

2	Setting	C :
	Extremely prudential	7,5
	Prudential	9
	Standard	10,5
	Active	12
	Extreme	13,5

3	Lt
Camber	Maximum width
0°	LS + 18
2°	LS + 21
4°	LS + 24
8°	LS + 30
12°	LS + 36

Maximum width calculated including 3 cm. distance between tyre and handrim



Wheelchair : Progeo JOKER ENERGY - Rev. 2025.1

Prices - VAT excluded
to be added to the price list

Send this page only if any optional have been chosen

☆standard

(supplied if no requests are indicated)

□on request

(without extra charge)

○optional

(cost to be added at the basic price)

customer_support_uk@permobil.com - https://www.permobil.com/en-gb