






permobil Order Form - Progeo JOKER

Permobil Limited Unit 7, Polaris Centre, 41 Brownfields, Welwyn Garden City, Herts, AL7 1AN United Kingdom - Phone: +44 1484 722888

Date	UK (English)	
User	Company stamp	
End user weight (kg)		Rigid frame wheelchair CRASH TEST APPROVED
Language user manual	Responsible	£ 2592,00 Retail price - VAT not included

Wheelchair measures - The tolerance on all the dimensions indicated and/or required is ± 10 mm			
Seat width (LS)	☐ 33 ☐ 36 ☐ 39 ☐ 42 ☐ 45 ☐ 48		
The seat width is measured from the outer sides of the frame			
Seat depth (PS)	☐ 35 ☐ 37,5 ☐ 40 ☐ 42,5 ☐ 45 ☐ 47,5 ☐ 50		
The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.			
Footplate distance (DP)	 (cm) ± 0.5 cm - adjustment steps 1.5cm		
(Real distance between seat and footplate) - IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)			
Setting	★ Standard ☐ Prudential ☐ Active ☐ Extreme		
Backrest Height (TS)	Std. Value 33,5 (min 26-max 47) adjustment steps 1.5cm		
camber	★ 0° ☐ 2° ☐ 4°		
Front height STD 49		Rear height STD 42	


	★ Aluminum rear adjustment plate		
	★ Removable mudguard ABS		
	☐ Removable ABS fibre straight sideguard		
	☐ Removable carbon fibre mudguard	03946	337,00
	☐ Removable Carbon fibre straight sideguard	08918	337,00
	☐ FULL CARBON rear adjustment plate with integrated sideguard	07646	1.006,00
	★ Version mudguard (Available only with rear wheels 24")		

JOKER		★ Aluminium footplate, rake adjustable Minimum footplate distance 35cm. ★ Internal ☐ 2/3 Internal ☐ External ☐ 2/3 External		11523	
		☐ Aluminium footplate EVOLUTION - rake and depth adjustable Minimum footplate distance 35cm. ★ Internal ☐ 2/3 Internal ☐ External ☐ 2/3 External		07649	145,00
		☐ One piece Carbon fibre footplate, rake adjustable Minimum footplate distance 35cm. ★ Internal ☐ External		03346	319,00
JOKER ZERO		☐ Plastic footrest, rake adjustable Minimum footplate distance 15cm. ★ 2/3 Internal ☐ 2/3 External		06814	203,00
		☐ One piece aluminium footplate - Flip-up!, rake and depth adjustable Minimum footplate distance 15cm. ★ Internal ☐ 2/3 Internal ☐ External ☐ 2/3 External		06856	203,00
		☐ One piece Carbon fibre footplate rake adjustable, swing away, flip-up. Minimum footplate distance 15cm. ★ Internal ☐ External		06857	522,00


Date Order User Address

Front frame


☆ Aluminium standard frame



☆



□



□

(max 4" caster)

Frame colour**Frame colour (without extra charge)**☆ Opaque black
□ Tekno anthracite□ POLAR white
□ Matt Midnight blue□ Matt orange
□ Matt grey

□ Matt red

Frame colour on request○ Bright silver
○ Bright blue○ Bright fuchsia
○ Bright lime○ Bright blue china
○ Bright red

12579

109,00**Seat**

☆ Seat cloth with pocket (under seat)

○ Seat with tension straps

08883

240,00

○ Carbon rigid seat

08336

509,00

○ Pocket (under seat) [25x11x4,5]

08281

79,00**Rear wheels axle**

☆ Aluminum rear wheels axle

○ Carbon fiber rear wheels axle

Lightest configuration

07647


158,00**Backrest adjustment and fixing system**☆ Backrest with double locking (upright and folded)
with rake adjustable from -4° to +16°□ Backrest with upright locking system ONLY (no locking when folded)
with rake adjustable from -4° to +16°□ Fixed back(non-folding back)
with rake adjustable from -4° to +16°**Backrest**

☆ Adjustable velcro straps and padded transpiring upholstery

☆ Aluminium backrest

○ Titanium backrest

○ Carbon fibre backrest



06811

203,00

04389

802,00○ **AIR** Integrated posture backrest (NOT REMOVABLE)Carbon shell - Lateral retaining 7,5cm
Depth adjustment ±1cm. - Angle adjustment ±10°
not available for seat width 48**1.121,00**

H 28	H 34	H 40			
Backrest Height (TS) min 34	Backrest Height (TS) min 38	Backrest Height (TS) min 44			
(.....)	(.....)	(.....)			
Lateral retaining wing A=high - R=regular					
A	R	A	R	A	R
☐	☐	☐	☐	☐	☐


Installation on:

☆ Backrest tubes

□ Aluminium backrest

○ Titanium backrest

○ Carbon fibre backrest



06811

203,00

04389

802,00

Date Order User Address

Handle tube

☆ Tubes with standard handle

- ☐ Tubes with short handle
☐ SPORT tubes Without push handles

☐ Carbon fibre tubes with short handle *Lightest configuration*
Not available with: Posture backrest PHYSIO - Posture backrest A I R

04390 **158,00**

☐ Folding handles **Minimum backrest height 30,5cm**
Not available with: Push handles for companion

02958 **270,00****Casters**

☆ 4" Soft

☐ 3" Sport*Lightest configuration*☐ 5" Soft☐ 6" Soft

(Available only with 100° angle)

☐ 6" Pneumatic

(Available only with 100° angle)

☐ 3" Active Soft (wide)

(Aluminium alloy rim)

06824 **95,00**☐ 4" Active Sport (narrow)

(Aluminium alloy rim)

01398 **95,00**☐ 4" Active Soft (wide)

(Aluminium alloy rim)

03943 **95,00**☐ 5" Active Soft (wide)

(Aluminium alloy rim)

07518 **95,00**☐ 5" Pneumatic10844 **68,00****Colour set: includes rear wheel hub and rim and front forks**

Only standard Progeo 22"/24" wheels have coloured hubs and rims

☆ Black

22" 24" 25" 26"

☐ Bright silver

22" 24"

00833 **135,00**☐ Bright lime

22" 24"

00833 **135,00**☐ Bright blue

22" 24"

00833 **135,00**☐ Bright red

22" 24"

00833 **135,00**☐ Bright fuchsia

22" 24"

00833 **135,00**☒ **POLAR white****22" 24"****00833 135,00**☐ Tekno anthracite

22" 24"

00833 **135,00**☐ Matt Midnight blue

22" 24"

00833 **135,00**☐ Matt grey

22" 24"

00833 **135,00**☐ Matt red

22" 24"

00833 **135,00****Rear wheels**

☆ 24" Ultra-light Progeo - Standard!

☐ 22" Ultra-light Progeo - Standard!☐ 25" Ultra-light Progeo - Standard!

(colour: black)

☐ 26" Progeo - Standard!

(colour: black)

☐ 24" ULTRA-LIGHT PROGEO KEY SPORT*Lightest configuration*04383 **695,00**☐ 24" Ultra-light "Spinergy LX"01685 **695,00**☐ 24" Ultra-light "Spinergy LXL"☐ Black hub and white spokes☐ Hub and spokes black07846 **695,00**☐ 24" Ultra-light "Spinergy Flex rim"

Integrated aluminium/rubber handrim

02270 **1.150,00**☐ 25" Ultra-light "Spinergy LX"02101 **695,00**☐ 25" Ultra-light "Spinergy LXL"☐ Black hub and white spokes07847 **695,00**☐ 26" Ultra-light "Spinergy LX"02102 **695,00**☐ 24" Off-road (aluminium handrim black)01392 **270,00****Rear wheel tyre**

☆ High pressure tyre (22"-24"-25"-26")

Lightest configuration☐ High pressure profile tyre (24")☐ Marathon Plus

(24"-25")

with anti-puncture protection

01669 **95,00**☐ Anti-puncture - Colour Grey

(22")

04181 **105,00**☐ Anti-puncture - Colour Black

(24")

11083 **95,00****Not available with Rear wheels : 24" ULTRA-LIGHT PROGEO KEY SPORT, 24" Ultra-light "Spinergy LXL"****Handrim****Wheel distance ☆ 3 cm**☐ **4 cm**

☆ Aluminium

(22"-24"-25"-26")

Lightest configuration☐ Aluminium black SOFT-TOUCH

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

09170 **160,00**☐ Chromed Steel

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

01684 **160,00**☐ Titanium

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

00016 **290,00**☐ Ultra Grepp

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

00017 **290,00**☐ TOP FRICTION

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

11023 **290,00****Oval countour with a rubber insert**☐ NOVA H

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

11025 **290,00****Ergonomic black profile**☐ NOVA D

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

11026 **290,00****Ergonomic profile with rubber coating (ideal for quadriplegics)**☐ Silicon cover for standard aluminium handrim

Available for measures (22"- 24" - 26")

04050 **172,00**

Date Order User Address

Wheel locks

☆ Bent lever push to lock

☐ Straight lever push to lock


☐ Straight lever pull to lock

☐ Sport brake

Lightest configuration

☐ Brake with extension lever push to lock

☐ Brake with extension lever pull to lock



06428 135,00

01949 68,00

03372 68,00

Armrests

☐ Tubular armrests

(Swing away, removable and height adjustable armrest)



00694 270,00

Accessories

☐ Pad for front frame protectors

07058 56,00

☐ Backrest-seat connection cloth

10241 25,00

☐ Removable storage bag (with quick coupling)

09811 135,00

SAFETY

☐ 4 points tie-down system rings

(installation must be carried out by the manufacturer)

A head restraint suitable for transportation must be available and fitted accordingly during transportation.



08654 203,00

Optional 1**Only one optional devices for side can be chosen**

☐ Transit wheels

Exclude the installation of Anti-tippers



03290 158,00

☐ Anti-tipper with rapid hookup

(inside rotation to reduce the encumbrance)

Exclude the installation of Transit wheels

☐ Right

☐ Left



06821 145,00

06822 145,00

Optional 2

☐ Push handles for companion

(Removable and height adjustable)



03645 189,00

☐ Tetraclip

(Quick release axle for quadriplegic)

03947 189,00

☐ Safety pelvic belt

03972 68,00

☐ Harness

04103 145,00

☐ Lumbar support

(to insert in the backrest upholstery)

00013 43,00

☐ Swing away lateral support

Prevents backrest folding

Not available with Posture backrest PHYSIO - Posture backrest A I R

☐ Right

☐ Left



02097 270,00

02098 270,00

☐ Spoke guards Black

(22"-24")

01043 158,00

☐ Spoke guards Transparent

(22"-24")

00835 158,00

☐ Spoke guards Progeo design

(22"-24")

00019 189,00

Additional Kit

☐ 24" Off-road rear wheels

Off-Road tyre - aluminium handrim (black colour)

03378 534,00

☐ 24" FAT rear wheels

FAT tyre 3" - aluminium handrim (black colour)

09357 895,00

Date.....

Order.....

User.....

Address.....

PHYSIO - Posture backrest
(Std. Backrest upholstery not provided)

Removable with quick release
Suitable for tubes with diameter 18÷26mm
Depth adjustment (4cm)
Angle adjustment (±18°)

- ☐ *AIR* Carbon shell
- not available for seat width 48

BACKREST HEIGHT		28		34		40	
Shell's shape		Lateral retaining wing A=high - R=regular					
7,5 cm		A	R	A	R	A	R
							
		☐	☐	☐	☐	☐	☐

1.121,00

- ☐ *POSTURE* Aluminum shell

BACKREST HEIGHT	34	40	52
Shell's shape 7,5 cm	Lateral retaining wing R=regular		
	R	R	R

1.085,00

Configuration DEEP with side support PAD (height adjustable)

Alluminum lateral with padded, breathable upholstery.



Small (h 13cm)



- ☐ Right
10272

49,00

- ☐ Left
10272

49,00

Large (h 20cm)



- ☐ Right
10273


49,00

- ☐ Left
10273

49,00

SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve
Alluminum lateral with padded, breathable upholstery.



- ☐ Large (18x12)

05647 196,00

- ☐ Small (14x9)

05643 196,00

- ☐ Lateral support extension

Quantity (number pcs) (.....)

10271 74,00

- ☐ Lateral support extension

Quantity (number pcs) (.....)

- ☐ Small (14x9)

05642 196,00

- ☐ Large (18x12)

05646 196,00

- ☐ Headrest adjustable - Double ball mounting system

Available only with backrest Physio Posture



- ☐ Large (34x14)
- ☐ Small (28x14)

05648 385,00

Permobil Limited Unit 7, Polaris Centre, 41 Brownfields, Welwyn Garden City, Herts, AL7 1AN United Kingdom - Phone: +44 1484 722888

Date Order User Address

Note

Wheelchair : Progeo JOKER - Rev. 2025.1		Prices - VAT excluded to be added to the price list		Send this page only if any optional have been chosen	
				☒	standard (supplied if no requests are indicated)
				☐	on request (without extra charge)
				☐	optional (cost to be added at the basic price)

customer_support_uk@permobil.com - https://www.permobil.com/en-gb

REFERENCE TABLES

1	Front height			
With closed frame (95°) max 5" caster With closed frame (90°) max 5" caster With closed frame (85°) max 4" caster				
FORK	Casters			
	3"	4"	5"	6"
SPORT 	max 52	max 53	max 54	
SMALL 	max 54	max 55	max 56	max 57
MEDIUM 	max 56	max 57	max 58	max 59
The size of the fork will be defined only by Progeo in accordance with wheelchair final configuration.				

2	Rear height			
Rear wheels				
22"	24"	25"	26"	
min 36	min 38	min 39	min 40	
max 44	max 46	max 47	max 48	

3	Footplate option	Footplate distance
	Aluminium footplate rake adjustable	min 35
	One piece Carbon fibre footplate	min 35
	One piece Carbon fibre footplate rake adjustable	min 35
	Plastic footrest, rake adjustable	min 35
	One piece aluminium footplate - Flip-up!, rake and depth adjustable	min 35
	One piece Carbon fibre footplate rake adjustable, swing away, flip-up.	min 35
IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)		
Calculated with the best configuration and with the footplates positioned on inside.		
NOTE! Lower footplate distances are possible only after our valuation		

4	Rear wheels	
	Model	Features
	22" Ultra-light Progeo - Standard!	(not available in anthracite color)
	24" Ultra-light Progeo - Standard!	(available for all colours)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)
	25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes)
	26" Progeo - Standard!	(colour: black)
	26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)