

**Carbon fibre frame made-to-measure.**

**Front height** (min 44-max 54)

**Rear height** (min 37-max 48)  
calculated with rear wheels 24"

**Frame Angle**  
☐ 85°   ☒ 90°   ☐ 95°

**Seat depth (PS)**  
☐ 35   ☐ 37,5  
☐ 40   ☐ 42,5  
☐ 45   ☐ 47,5

**Backrest height (TS)**  
 Std. Value 33,5  
 (min 26-max 47)  
 adjustment steps 1,5cm

**Setting C**  
☐ Extremely prudential 7,5  
☐ Prudential 9  
☒ Standard 10,5  
☐ Active 12  
☐ Extreme 13,5

**NOIR is a wheelchair without possibilities of adjustments after its manufacturing..**


See reference tables

Take care to write all the requested measures (in centimetres)

**Backrest rake adjustable from -4° to +16°**

**Seat width (LS)**  
☐ 33   ☐ 36  
☐ 39   ☐ 42  
☐ 45

**camber**  
☒ 2°   ☐ 0°   ☐ 4°

| Footplate option                                                                                                                                  |                                                                                                                                                                                                               |                                                                                       |
|---------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| ☆ Carbon fibre frame with integrated footplate, not adjustable<br><i>Lightest configuration</i>                                                   |                                                                                                                          |                                                                                       |
| <input type="checkbox"/> Aluminium footplate adjustable distance +/- 1,5cm (rake and depth adjustable)<br><b>Minimum footplate distance 35cm.</b> | <div> <input checked="" type="checkbox"/> Internal           <input type="checkbox"/> 2/3 Internal           <input type="checkbox"/> 2/3 External           <input type="checkbox"/> External         </div> |  |

| Seat                                                  |       |        |
|-------------------------------------------------------|-------|--------|
| ☆ Seat cloth with pocket (under seat)                 |       |        |
| <input type="checkbox"/> Seat with tension straps     | 08883 | 240,00 |
| <input type="radio"/> Carbon rigid seat               | 08336 | 509,00 |
| <input type="radio"/> Pocket (under seat) [25x11x4,5] | 08281 | 79,00  |

| Date  | Order | User | Address |
|-------|-------|------|---------|
| ..... |       |      |         |

*Since that the paint is sensitive to the surface material and texture, the colour selected for the wheelchair frame can also differ considerably from what is displayed on a RAL colour chart or also compared to what is displayed on a monitor.*

**CARBON**



- ☆ Polished transparent painting

- |                               |       |        |
|-------------------------------|-------|--------|
| ○ Opaque transparent painting | 08925 | 938,00 |
|-------------------------------|-------|--------|

**ESSENCE**



**Polished transparent painting**

- |                      |       |        |
|----------------------|-------|--------|
| ○ ANTHRACITE Essence | 10351 | 938,00 |
| ○ BLUE Essence       | 08928 | 938,00 |
| ○ GREEN Essence      | 11163 | 938,00 |
| ○ LIME Essence       | 08926 | 938,00 |
| ○ ORANGE Essence     | 10349 | 938,00 |
| ○ RED Essence        | 08927 | 938,00 |
| ○ SAND Essence       | 11165 | 938,00 |
| ○ VIOLET Essence     | 11167 | 938,00 |
| ○ WHITE Essence      | 10350 | 938,00 |

**Opaque transparent painting**

- |                             |       |                 |
|-----------------------------|-------|-----------------|
| ○ ANTHRACITE Opaque Essence | 08929 | <b>1,201,00</b> |
| ○ BLUE Opaque Essence       | 10354 | <b>938,00</b>   |
| ○ GREEN Opaque Essence      | 11162 | <b>938,00</b>   |
| ○ LIME Opaque Essence       | 10352 | <b>938,00</b>   |
| ○ ORANGE Opaque Essence     | 10355 | <b>938,00</b>   |
| ○ RED Opaque Essence        | 10353 | <b>938,00</b>   |
| ○ SAND Opaque Essence       | 11164 | <b>938,00</b>   |
| ○ VIOLET Opaque Essence     | 11166 | <b>938,00</b>   |
| ○ WHITE Opaque Essence      | 10356 | <b>938,00</b>   |

- my ESSENCE colour RAL

- ☐
- Glossy

---

Note

- ☐
- Opaque

- |       |        |
|-------|--------|
| 11186 | 938,00 |
|-------|--------|

## RACE



### Polished transparent painting

- |                   |       |        |
|-------------------|-------|--------|
| ○ ANTHRACITE Race | 11168 | 938,00 |
| ○ BLUE Race       | 11170 | 938,00 |
| ○ GREEN Race      | 11172 | 938,00 |
| ○ LIME Race       | 11174 | 938,00 |
| ○ ORANGE Race     | 11176 | 938,00 |
| ○ RED Race        | 11178 | 938,00 |
| ○ SAND Race       | 11180 | 938,00 |
| ○ VIOLET Race     | 11182 | 938,00 |
| ○ WHITE Race      | 11184 | 938,00 |

### Opaque transparent painting

- |                          |       |        |
|--------------------------|-------|--------|
| ○ ANTHRACITE Opaque Race | 11169 | 938,00 |
| ○ BLUE Opaque Race       | 11171 | 938,00 |
| ○ GREEN Opaque Race      | 11173 | 938,00 |
| ○ LIME Opaque Race       | 11175 | 938,00 |
| ○ ORANGE Opaque Race     | 11177 | 938,00 |
| ○ RED Opaque Race        | 11179 | 938,00 |
| ○ SAND Opaque Race       | 11181 | 938,00 |
| ○ VIOLET Opaque Race     | 11183 | 938,00 |
| ○ WHITE Opaque Race      | 11185 | 938,00 |

- ☐ my RACE colour RAL

- ☐
- Glossy

Note

- Opaque

- |       |        |
|-------|--------|
| 11187 | 938,00 |
|-------|--------|

**FLASH**



**Polished transparent painting**

- |               |       |        |
|---------------|-------|--------|
| ○ BLUE Flash  | 11190 | 938,00 |
| ○ WHITE Flash | 11192 | 938,00 |
| ○ RED Flash   | 11188 | 938,00 |

**Opaque transparent painting**

- |                      |       |        |
|----------------------|-------|--------|
| ○ BLUE Opaque Flash  | 11191 | 938,00 |
| ○ WHITE Opaque Flash | 11193 | 938,00 |
| ○ RED Opaque Flash   | 11189 | 938,00 |

**Colour set: includes rear wheel hub and rim and front forks**

Only standard Progeo 22"/24" wheels have coloured hubs and rims

- |                      |                 |       |        |
|----------------------|-----------------|-------|--------|
| ☆ Black              | 22" 24" 25" 26" |       |        |
| ○ Bright silver      | 22" 24"         | 00833 | 135,00 |
| ○ Bright lime        | 22" 24"         | 00833 | 135,00 |
| ○ Bright blue        | 22" 24"         | 00833 | 135,00 |
| ○ Bright red         | 22" 24"         | 00833 | 135,00 |
| ○ Bright fuchsia     | 22" 24"         | 00833 | 135,00 |
| ○ POLAR white        | 22" 24"         | 00833 | 135,00 |
| ○ Tekno anthracite   | 22" 24"         | 00833 | 135,00 |
| ○ Matt Midnight blue | 22" 24"         | 00833 | 135,00 |
| ○ Matt grey          | 22" 24"         | 00833 | 135,00 |
| ○ Matt red           | 22" 24"         | 00833 | 135,00 |

Permobil Limited Unit 7, Polaris Centre, 41 Brownfields, Welwyn Garden City, Herts, AL7 1AN United Kingdom - Phone: +44 1484 722888

Date.....Order.....User.....Address.....

Backrest adjustment and fixing system

☆ Backrest with double locking (upright and folded)

with rake adjustable from -4° to +16°

☐ Backrest with upright locking system ONLY (no locking when folded)

with rake adjustable from -4° to +16°

☐ Fixed back(non-folding back)

with rake adjustable from -4° to +16°

☐ Integrated backrest (Banded, not folding back)

Exclude installation of the titanium backrest

Exclude installation of the Tubular armrests

10841

1.336,00

( ..... ) to the ground ☐

☆ 0° (90° to the floor)

( ..... ) to the ground ☐

Backrest

☆ Adjustable velcro straps and padded transpiring upholstery

☆ Carbon fibre backrest

☐ Titanium backrest

Not available with: Integrated backrest (Banded, not folding back)

☐ **AIR** Integrated posture backrest (NOT REMOVABLE)

Carbon shell - Lateral retaining 7,5cm

Depth adjustment ±1cm. - Angle adjustment ±10°

not available for seat width 48

1.121,00

H 28

Backrest Height (TS)

min 34

( ..... )

H 34

Backrest Height (TS)

min 38

( ..... )

H 40

Backrest Height (TS)

min 44

( ..... )

Lateral retaining wing A=high - R=regular

A

R

A

R

A

R

☐

☐

☐

☐

☐

☐

Installation on:

☆ Aluminium backrest tubes

WITHOUT backrest rigidizer

(Carbon tubes with integrated backrest)

☐ Carbon fibre backrest

☐ Titanium backrest

Not available with: Integrated backrest (Banded, not folding back)

Handle tube

☆ Tubes with short handle

☐ SPORT tubes Without push handles

Lightest configuration

☐ Folding handles

Minimum backrest height 30,5cm

02958

270,00

Casters

☆ 3" Sport

Lightest configuration

☐ 4" Soft

☐ 3" Active Soft (wide)

☐ 4" Active Sport (narrow)

☐ 4" Active Soft (wide)

☐ 5" Active Soft (wide)

☐ 5" Pneumatic

(Aluminium alloy rim)

(Aluminium alloy rim)

(Aluminium alloy rim)

(Aluminium alloy rim)

Rear wheels - See chart N°3

☆ 24" ULTRA-LIGHT PROGEO KEY SPORT

Lightest configuration

☐ 24" Ultra-light "Spinergy LX"

☐ 24" Ultra-light "Spinergy LXL"

☐ 24" Ultra-light "Spinergy Flex rim"

☐ 25" Ultra-light "Spinergy LX"

☐ 25" Ultra-light "Spinergy LXL"

☐ 26" Ultra-light "Spinergy LX"

☐ Black hub and white spokes

Integrated aluminium/rubber handrim

☐ Black hub and white spokes

☐ Hub and spokes black

04353

1.150,00

customer\_support\_uk@permobil.com - https://www.permobil.com/en-gb

Pag. 3 / 6

UK\_202504

Rev. 2025.1\_Apr

|            |             |            |               |
|------------|-------------|------------|---------------|
| Date ..... | Order ..... | User ..... | Address ..... |
|------------|-------------|------------|---------------|

|                                                                                                          |               |                               |             |
|----------------------------------------------------------------------------------------------------------|---------------|-------------------------------|-------------|
| <b>Rear wheel tyre</b>                                                                                   |               |                               |             |
| <input checked="" type="checkbox"/> High pressure tyre                                                   | (24"-25"-26") | <i>Lightest configuration</i> |             |
| <input type="checkbox"/> High pressure profile tyre                                                      | (24")         |                               |             |
| <input type="checkbox"/> Marathon Plus                                                                   | (24"-25")     | with anti-puncture protection | 01669 95,00 |
| <input type="checkbox"/> Anti-puncture - Colour Black                                                    | (24")         |                               | 11083 95,00 |
| <b>Not available with Rear wheels : 24" ULTRA-LIGHT PROGEO KEY SPORT, 24" Ultra-light "Spinergy LXL"</b> |               |                               |             |

|                                                                        |                                                                                              |                                                       |              |
|------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|-------------------------------------------------------|--------------|
| <b>Handrim</b>                                                         | <b>Wheel distance</b> <input checked="" type="checkbox"/> 3 cm <input type="checkbox"/> 4 cm |                                                       |              |
| <input checked="" type="checkbox"/> Aluminium                          | (24"-25"-26")                                                                                | <i>Lightest configuration</i>                         |              |
| <input type="checkbox"/> Aluminium black SOFT-TOUCH                    | (24"-25")                                                                                    | (24" handrim, on 25" wheel is installed with spacers) | 09170 160,00 |
| <input type="checkbox"/> Chromed Steel                                 | (24"-25")                                                                                    | (24" handrim, on 25" wheel is installed with spacers) | 01684 160,00 |
| <input type="checkbox"/> Titanium                                      | (24"-25")                                                                                    | (24" handrim, on 25" wheel is installed with spacers) | 00016 290,00 |
| <input type="checkbox"/> Ultra Grepp                                   | (24"-25")                                                                                    | (24" handrim, on 25" wheel is installed with spacers) | 00017 290,00 |
| <input type="checkbox"/> TOP FRICTION                                  | (24"-25")                                                                                    | (24" handrim, on 25" wheel is installed with spacers) | 11023 290,00 |
| <b>Oval countour with a rubber insert</b>                              |                                                                                              |                                                       |              |
| <input type="checkbox"/> NOVA H                                        | (24"-25")                                                                                    | (24" handrim, on 25" wheel is installed with spacers) | 11025 290,00 |
| <b>Ergonomic black profile</b>                                         |                                                                                              |                                                       |              |
| <input type="checkbox"/> NOVA D                                        | (24"-25")                                                                                    | (24" handrim, on 25" wheel is installed with spacers) | 11026 290,00 |
| <b>Ergonomic profile with rubber coating (ideal for quadriplegics)</b> |                                                                                              |                                                       |              |
| <input type="checkbox"/> Silicon cover for standard aluminium handrim  | Available for measures (22" - 24" - 26")                                                     |                                                       | 04050 172,00 |

|                                                                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Wheel locks</b>                                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <input checked="" type="checkbox"/> Sport brake                  | <br><br><br><br><br> |
| <input type="checkbox"/> Bent lever push to lock                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <input type="checkbox"/> Straight lever push to lock             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <input type="checkbox"/> Straight lever pull to lock             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <input type="checkbox"/> Brake with extension lever push to lock |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <input type="checkbox"/> Brake with extension lever pull to lock |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |

|                                                                    |                                                   |                                               |                                               |
|--------------------------------------------------------------------|---------------------------------------------------|-----------------------------------------------|-----------------------------------------------|
| <b>Mudguard</b>                                                    |                                                   |                                               |                                               |
| <input checked="" type="checkbox"/> Carbon fibre mudguard          |                                                   |                                               |                                               |
| <input type="checkbox"/> Carbon fibre Straight sideguard           | <input checked="" type="checkbox"/> Aligned wheel | <input type="checkbox"/> Under wheel ..... cm | <input type="checkbox"/> Above wheel ..... cm |
| <b>Mudguard enlargement</b>                                        |                                                   |                                               |                                               |
| <input checked="" type="checkbox"/> + 0 cm                         |                                                   |                                               |                                               |
| <input type="checkbox"/> + 1 cm                                    | <input type="checkbox"/> + 2 cm                   |                                               | 12648                                         |
| <input type="checkbox"/> Removable carbon fibre mudguard           |                                                   |                                               | 04109 270,00                                  |
| <input type="checkbox"/> Removable Carbon fibre straight sideguard | <input checked="" type="checkbox"/> Aligned wheel | <input type="checkbox"/> Under wheel ..... cm | <input type="checkbox"/> Above wheel ..... cm |
|                                                                    |                                                   |                                               | 04109 270,00                                  |

|                                                                   |                                                       |                                                                                       |
|-------------------------------------------------------------------|-------------------------------------------------------|---------------------------------------------------------------------------------------|
| <b>Armrests</b>                                                   |                                                       |                                                                                       |
| Tubular armrests                                                  | (Swing away, removable and height adjustable armrest) |  |
| Not available with Integrated backrest (Banded, not folding back) |                                                       |                                                                                       |
|                                                                   |                                                       | 00694 270,00                                                                          |

|                                                                      |  |              |
|----------------------------------------------------------------------|--|--------------|
| <b>Accessories</b>                                                   |  |              |
| <input type="checkbox"/> Pad for front frame protectors              |  | 07058 56,00  |
| <input type="checkbox"/> Backrest-seat connection cloth              |  | 10241 25,00  |
| <input type="checkbox"/> Removable storage bag (with quick coupling) |  | 09811 135,00 |

|                                                        |                                             |              |
|--------------------------------------------------------|---------------------------------------------|--------------|
| <b>Optional 1</b>                                      |                                             |              |
| <input type="checkbox"/> Transit wheels                |                                             | 03290 158,00 |
| <input type="checkbox"/> Anti-tipper with rapid hookup | (inside rotation to reduce the encumbrance) |              |
|                                                        | <input type="checkbox"/> Right              | 06821 145,00 |
|                                                        | <input type="checkbox"/> Left               | 06822 145,00 |

|                                             |                                        |              |
|---------------------------------------------|----------------------------------------|--------------|
| <b>Optional 2</b>                           |                                        |              |
| <input type="checkbox"/> Safety pelvic belt |                                        | 03972 68,00  |
| <input type="checkbox"/> Harness            |                                        | 04103 145,00 |
| <input type="checkbox"/> Lumbar support     | (to insert in the backrest upholstery) | 00013 43,00  |
| <input type="checkbox"/> Tetraclip          | (Quick release axle for quadriplegic)  | 03947 189,00 |

|                                                   |                                                  |              |
|---------------------------------------------------|--------------------------------------------------|--------------|
| <b>Additional Kit</b>                             |                                                  |              |
| <input type="checkbox"/> 24" Off-road rear wheels | Off-Road tyre - aluminium handrim (black colour) | 03378 534,00 |
| <input type="checkbox"/> 24" FAT rear wheels      | FAT tyre 3" - aluminium handrim (black colour)   | 09357 895,00 |

Date.....

Order.....




User.....

Address.....

PHYSIO - Posture backrest  
(Std. Backrest upholstery not provided)

Removable with quick release  
Suitable for tubes with diameter 18÷26mm  
Depth adjustment (4cm)  
Angle adjustment (±18°)

- ☐  Carbon shell
- not available for seat width 48

| BACKREST HEIGHT                                                                                        |  | 28                                                                                |                                                                                    | 34                                                                                  |                                                                                     | 40                                                                                  |                                                                                     |
|--------------------------------------------------------------------------------------------------------|--|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| Shell's shape<br><br> |  | Lateral retaining wing A=high - R=regular                                         |                                                                                    |                                                                                     |                                                                                     |                                                                                     |                                                                                     |
|                                                                                                        |  | A                                                                                 | R                                                                                  | A                                                                                   | R                                                                                   | A                                                                                   | R                                                                                   |
|                                                                                                        |  |  |  |  |  |  |  |

1.121,00

- ☐  Aluminum shell

| BACKREST HEIGHT                                                                                     | 34                                                                                  | 40                                                                                   | 52                                                                                   |
|-----------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| Shell's shape<br> | Lateral retaining wing R=regular                                                    |                                                                                      |                                                                                      |
|                                                                                                     | R                                                                                   | R                                                                                    | R                                                                                    |
|                                                                                                     |  |  |  |

1.085,00

Configuration DEEP with side support PAD (height adjustable)

Alluminum lateral with padded, breathable upholstery.



Small (h 13cm)



- ☐ Right  
10272

49,00

- ☐ Left  
10272

49,00

Large (h 20cm)



- ☐ Right  
10273


49,00

- ☐ Left  
10273

49,00

SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve  
Alluminum lateral with padded, breathable upholstery.



- ☐ Large (18x12)

05647 196,00

- ☐ Small (14x9)

05643 196,00

- ☐ Lateral support extension

Quantity (number pcs).....  
10271 74,00

- ☐ Lateral support extension

Quantity (number pcs).....

- ☐ Small (14x9)

05642 196,00

- ☐ Large (18x12)

05646 196,00

- ☐ Headrest adjustable - Double ball mounting system  
Available only with backrest Physio Posture



- ☐ Large (34x14)
- 05648 385,00
- ☐ Small (28x14)

Permobil LimitedUnit 7, Polaris Centre, 41 Brownfields,Welwyn Garden City, Herts, AL7 1ANUnited Kingdom - Phone: +44 1484 722888

DateOrderUserAddress

REFERENCE TABLES

**Progeo (Español)** NOIR is a wheelchair without possibilities of adjustments after its manufacturing.  
(1) The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.

| 1     | A Footplate width |       |       |
|-------|-------------------|-------|-------|
| DP\LS | 35÷39             | 40÷45 | 46÷50 |
| 33    | 21,5              | 21    | 20,5  |
| 36    | 24,5              | 24    | 23,5  |
| 39    | 27,5              | 27    | 26,5  |
| 42    | 30,5              | 30    | 29,5  |
| 45    | 33,5              | 33    | 32,5  |
| 48    | 36,5              | 36    | 35,5  |

| 2 | Setting              | C :  |
|---|----------------------|------|
|   | Extremely prudential | 7,5  |
|   | Prudential           | 9    |
|   | Standard             | 10,5 |
|   | Active               | 12   |
|   | Extreme              | 13,5 |

| 3      | Maximum width |
|--------|---------------|
| Camber |               |
| 0°     | LS + 18       |
| 2°     | LS + 21       |
| 4°     | LS + 24       |

Maximum width calculated including 3 cm. distance between tyre and handrim

Wheelchair : Progeo noir 2.0 - Rev. 2025.1

Prices - VAT excluded  
to be added to the price list

Send this page only if any optional have been chosen

☆ standard

(supplied if no requests are indicated)

□ on request

(without extra charge)

○ optional

(cost to be added at the basic price)

customer\_support\_uk@permobil.com - https://www.permobil.com/en-gb