

Permobil

Progeo Noir 2.0



Carbon



Rigid



Minimized
weight



Made-to-
measure

permobil

The ultimate in carbon wheel-chair technology

The Progeo Noir's seamless frame combines stunning design with superior engineering.

For experienced active manual wheelchair users who know the preferred seating set-up and want to maximise the weight savings of a rigid frame. Combining the minimized weight with the individuality of each Noir, this gives a fine-tuned performance for a smooth driving experience.



Progeo Noir 2.0

Seat width	33 - 36 - 39 - 42 - 45 cm
Seat depth	35 - 37.5 - 40 - 42.5 - 45 - 47.5 cm
Seat height	Front from 44 to 54 cm (made-to-measure) Rear from 37 to 48 cm (made-to-measure)
Footsupport distance	From 35 to 51 cm (made-to-measure)
Backrest height	From 24.5 to 47 cm (adjustable)
Front frame angle	85° - 90° - 95°
Total length	82 cm approx.
Total width	Seat width + 17 cm (camber 0°)
Backrest angle	From 74° to 94° (adjustable)
Camber	0° - 2° - 4°
Weight without rear wheels	Starting from 4.0 kg approx.
Wheelchair weight	Starting from 6.6 kg approx.
Maximum load	125 kg

Permobil reserves the right to change specifications and design without prior notice.

Personalise your ride

CARBON
ESSENCE

LIME
ESSENCE

RED
ESSENCE

BLUE
ESSENCE

OPAQUE
ESSENCE

Actual product colors may vary.