

Quantum®

32 Wedgwood Road, Bicester,
Oxon, OX26 4UL

Phone: 01869 324600

Fax: 01869 323070

Email: sales@quantumrehab.co.uk

Web: www.quantumrehab.co.uk



Dealer April 2021



Crash
Tested



Q6 Edge 2.0 with Synergy® Seating & NE+ or Q-Logic® 3 Controls

Dealer Name:

Order Date:

Address:

Post Code:

Phone

Email:

PO No:

Order Ref:

Delivery Address:

Post Code:

Quantum Specialist:

Standard Features

- Synergy Seating System
- 4" of Seat Width Adjustment
- 4" of Seat Depth Adjustment
- 2 x 55ah batteries
- NE+ 6 Key Joystick
- 2" x 6" Front & Rear Castors, Flat Free
- 3" x 14" Pneumatic Drive Tyres
- 12 Colour Options
- 4 Mph (6 Km/h) Motors
- ATX Suspension System
- Single Post Height Adj. Removable Armrests
- Centre Mount Foot Platform
- Maximum User Weight 21 Stone 6lbs
- Synergy Simplicity Seat Cushion
- Tie-Down System for Transport

Standard = Standard Factory Feature | FOC = Select instead of Standard Factory Feature

Base Model

Q6 Edge 2.0 with Synergy® Seating and NE+ Controls - £4,595.00

Drive Wheels

- | | |
|--|----------|
| <input type="checkbox"/> 3" x 14" Pneumatic, Black Tyres | Standard |
| <input type="checkbox"/> 3" x 14" Flat Free, Black Tyres | £150.00 |
| <input type="checkbox"/> 4" x 14" Pneumatic, Black Tyres, Knobbly, Extended Traction Tyres <i>Increases overall base width by 1.5"</i> | £220.00 |
| <input type="checkbox"/> 4" x 14" Flat Free, Black Tyres, Knobbly, Extended Traction Tyres <i>Increases overall base width by 1.5"</i> | £320.00 |

Castor Wheels

- | | |
|--|----------|
| <input type="checkbox"/> 2" x 6" Flat Free, Black Tyres, Front and Rear | Standard |
| <input type="checkbox"/> 2" x 8" Flat Free, Black Tyres, Front and Rear <i>Includes extended castor arms</i> | £364.00 |

Colours



☐ Red Corvette



☐ Back in Black



☐ Pretty in Pink



☐ Lime Light



☐ Silver Bullet



☐ Orange Crush



☐ Electric Blue



☐ White Satin



☐ Deep Purple



☐ Grey Street



☐ Midnight Blue



☐ Yellow Blaze

Seat Dimensions

Select width and depth package and enter desired seat width and depth in that package range, available in 1" increments. Maximum seat size is 22"W x 22"D.

Seat Width Packages

☐ 14"W-18"W ☐ 16"W-20"W ☐ 20"W-22"W

Set Width at:_____

Seat Depth Packages

☐ 14"D-18"D ☐ 16"D-20"D ☐ 20"D-22"D

Set Depth at:_____

Seat Cushions

☐ No Seat Cushion Required

FOC

☐ Simplicity Cushion: *General Use*

Standard

☐ TRU-Comfort™2 4-Way Stretch: *Skin Protection*

£45.00

☐ Solution® Cushion: *Skin Protection & Positioning*

£100.00

☐ Solution® 1 Cushion: *Skin Protection*

£100.00

☐ Spectrum® Gel Cushion: *Positioning*

£160.00

☐ Spectrum® Foam Cushion: *Skin Protection & Positioning*

£160.00

Refer to cushion matrix for size availability and enter your selection here:

W..... X D.....

Synergy® Backrest Specifications

☐ Synergy® Backrest Assembly *Including Visco padded backrest*

Standard

☐ Tension Adjustable Backrest Upholstery *This option may increase the delivery time*

£120.00

☐ Back Canes Only *Suitable for third party backrests. A backrest will not be supplied*

FOC

Backrest/Back Cane Height:

Set Backrest/Back Cane Height at:_____ (Available in 1" increments. Min 13", Max 25") *Please refer to Synergy Backrest Matrix for seat width and backrest height compatibility*

Back Height = Seat Pan to Top of Back | Odd heights = 4" gap at bottom, Even heights = 3" gap at bottom.

TRU-Comfort Plus® Backrest Specifications

☐ TRU-Comfort Plus® Backrest Assembly *4-way stretch fabric (Positioning). Only available with seat widths 16", 18", 20" and 22"*

£220.00

Backrest Height:

☐ 27" Tall

Standard

☐ 24" Tall

FOC

Back Height = Seat Pan to Top of Back | 3" gap at bottom.

Seat Positioning Adjustments

☐ Static Seat

Standard

☐ TRU-Balance®2 Power Tilt *0-55 degree centre of gravity (CG) power tilt seating system. Uses 1 power option*

£1100.00

☐ Tru-Balance®2 10" Power Adjustable Seat Lift and Tilt *Requires Q-Logic 3 Controls. Uses 2 power options.*

£1400.00

Seat to Floor Height

Static Seat

☐ 15.875" Seat to Floor ☐ 16.875" Seat to Floor ☐ 17.875" Seat to Floor

TB2 Power Tilt

☐ 17.125" Seat to Floor ☐ 18.125" Seat to Floor ☐ 19.125" Seat to Floor

10" Seat Lift and Tilt

☐ 17.875" Seat to Floor ☐ 18.875" Seat to Floor ☐ 19.875" Seat to Floor

Backrest Positioning Adjustments

☐ Static Backrest *Tool angle adjustable from 85° to 120° in 5° increments*

Standard

Additional Static Back Options:

☐ Flip Forward Backrest *Not suitable for occupied transportation*

£200.00

☐ 1" Endo Brackets *Reduces backrest width by 2". Not available with TRU-Comfort Plus backrest on seat widths 14"-17"*

£40.00

☐ 2" Endo Brackets *Reduces backrest width by 4". Not available with TRU-Comfort Plus backrest on seat widths 14"-19"*

£40.00

☐ TRU-Balance®2 Power Recline *Only available with 16"-22" s/widths. 168° of power recline with full drive lockout at 20°. Uses 1 power option. Limited backrest heights available. See reference chart on page**

£1790.00

Armrests

☐ Single Post, Quick Height Adjustable Removable Armrests, Pair *Height adjustment from 11" - 17" in 1/2" increments*

Standard

☐ Single Post, Tool Height Adjustable Armrests, Heavy Duty, Removable, Pair *Height adjustment from 10"-14"*

£220.00

Armrests (continued)

- | | |
|---|---------|
| <input type="checkbox"/> Paediatric, Removable, Height Adjustable <i>Height adjustment from 7.75" to 9.75" in 1/2" increments. Not available with full length armpads</i> | FOC |
| <input type="checkbox"/> Cantilever Flip-Up Armrests, Pair <i>Not available with a tension adjustable backrest. Increases overall width by 7". Height adj. 11"-17" depending on backrest height. Not available with Flip Forward Backrest. Not available with TB2 Power Recline</i> | £280.00 |
| <input type="checkbox"/> 2-Post, Flip Back, Height Adjustable, Removable <i>Height range 9"-13". Not available with a seat depth less than 16"</i> | £145.00 |
| <input type="checkbox"/> TRU-Balance®2 Flip-Up Armrests, Pair <i>Only available with TB2 Power Recline</i> | FOC |

Armrests Pad Options

- | | | |
|--|---|----------|
| <input type="checkbox"/> 14"L x 2"W Straight, Pair | <input type="checkbox"/> 10"L x 2"W Straight, Pair | Standard |
| <input type="checkbox"/> 14"L x 3"W Waterfall, Pair | <input type="checkbox"/> 10"L x 3"W Waterfall, Pair | FOC |
| <input type="checkbox"/> 14"L x 3.5"W Gel Armpad, Pair | | £160.00 |
| <input type="checkbox"/> 10"L x 3.5"W Gel Armpad, Pair | | £130.00 |

Ottobock Arm Channel (Each)

- | | |
|---|---------|
| <input type="checkbox"/> Medium (4.75" x 11.5") <i>Desk length armrest required</i> | £140.00 |
| <input type="checkbox"/> Large (4.75" x 14") <i>Full length armrest required</i> | £140.00 |
| <input type="checkbox"/> One Piece (4.75" x 18") <i>Full length armrest required. Not compatible with Hand Pad.</i> | £140.00 |
| <input type="checkbox"/> Ottobock Hand Pad <i>Not available on same side as controller</i> | £50.00 |
| <input type="checkbox"/> Medium <input type="checkbox"/> Large | |

Ottobock Arm Channel Mounting Position

- ☐ Left ☐ Right ☐ Both *Please amend price for pair*

Headrests

- | | |
|---|---------|
| <input type="checkbox"/> Quantum 10" W Headrest | £140.00 |
| <input type="checkbox"/> Stealth Small Contoured 9.75"W x 3"H | £200.00 |
| <input type="checkbox"/> Stealth Large Contoured 11"W x 6"H | £200.00 |
| <input type="checkbox"/> Stealth 10" Comfort Plus | £200.00 |
| <input type="checkbox"/> Stealth 14" Comfort Plus | £200.00 |

Please visit www.stealthproducts.com for more headrest options

Centre Mount Foot Platforms

- | | |
|---|----------|
| <input type="checkbox"/> Base Mounted Foot Platform <i>Mounted to chassis, will not move with Power Tilt or Seat Lift</i> | -£100.00 |
| <input type="checkbox"/> Adjustable Centre Mount Foot Platform <i>Mounted to Seat</i> | Standard |
| <input type="checkbox"/> Centre Mount Power Articulating Foot Platform <i>Mounted to Seat. Uses 1 power options</i> | £642.00 |

Footplates

- | | | | |
|--|----|--|----------|
| <input type="checkbox"/> 10"W x 11"D Split Footplates: | or | <input type="checkbox"/> 10"W x 8"D Split Footplates | Standard |
|--|----|--|----------|

Lower Extensions *Seat pan to Heel*

- | | | | |
|--|----|--|----------|
| <input type="checkbox"/> 10-16"L <i>Seat pan to heel</i> | or | <input type="checkbox"/> 13-19"L <i>Seat pan to heel</i> | Standard |
|--|----|--|----------|

Calf Pads, Pair

- | | | | |
|--|----|--|--------|
| <input type="checkbox"/> Small 6"H x 6"W | or | <input type="checkbox"/> Large 8"H x 8"W, tapered to 6"W | £40.00 |
|--|----|--|--------|

Calf Panel

- | | |
|--|--------|
| <input type="checkbox"/> Small 8"H x 18"W | £84.00 |
| <input type="checkbox"/> Calf Panel Extension (each) | £30.00 |

Footplate Extension Kits



- | | |
|---|--------|
| <input type="checkbox"/> 4" Width Extension Kit. <i>Not available with 2" Bell Shaped, Length Extension Kit</i> | £16.00 |
|---|--------|



- | | |
|---|--------|
| <input type="checkbox"/> 2" Bell Shaped, Length Extension Kit. <i>Not available with 4" Width Extension Kit</i> | £16.00 |
|---|--------|

Centre Mount Foot Platforms (continued)

Additional Optional Accessories:

- ☐ Wheel Bumpers, Pair
- £5.00
- ☐ Footplate Connector
- £10.00

Individual Legrests

Choose instead of centre mount foot platform options

- ☐ 60-90° Swing-Away Legrests Tool adjustable in 5° Increments. Factory set at 70°. Not available with a seat width less than 16"
- £110.00

Footplate Selection Includes heel loops

- ☐ Flip-up, Extruded Footplates, 7.25"W x 6"D, Pair
- Standard
- ☐ Angle Adjustable Footplates, 6"W x 8"D, Pair Other sizes available
- £115.00

Please provide the seat pan to heel measurement required, minimum 11", maximum 19.5":

- ☐ Heavy duty, drop in, swing-away legrests with angle adjustable footplates. Not available with a seat width less than 16"
- £130.00

Legrest Angle

☐ 70°

☐ 80°

Footplate Selection Includes heel loops

☐ 6"W x 8"D

☐ 8"W x 11"D Other sizes available

Please provide the seat pan to heel measurement required, minimum 12", maximum 18":

- ☐ 90° Paediatric Swing-Away Legrests Seat Pan to Heel range of adjustment 7"-10"L. Includes 4" x 6"D angle adjustable footplates and heel loops
- £110.00

- ☐ 70° Tapered Legrests Not available with seat widths 17" or less.
- £110.00

Footplate Selection Includes heel loops

- ☐ Flip-up, Extruded Footplates, 7.25"W x 6"D, Pair
- Standard
- ☐ Angle Adjustable Footplates, 6"W x 8"D, Pair
- £115.00

Please provide the seat pan to heel measurement required, minimum 11", maximum 19.5":

Elevating Legrests

Choose instead of centre mount foot platform options

- ☐ Manual Elevating Swing-Away Legrests Not suitable for occupied transportation. Not available with a seat width less than 16"
- £130.00

Footplate Selection Includes heel loops

- ☐ Flip-up, Extruded Footplates, 7.25" W x 6"D, Pair
- Standard
- ☐ Angle Adjustable Footplates, 6" W x 8"D, Includes Wheel Bumpers, Pair Other sizes available
- £115.00

Please provide the seat pan to heel measurement required, minimum 11", maximum 18":

- ☐ 70° Power Articulating Elevating Swing-Away Legrests 4" articulation range. Not available with a seat width less than 16" . Uses 2 power options
- £1320.00

Footplate Selection Includes heel loops

- ☐ Angle Adjustable Footplates, 6"W x 8"D, Includes Wheel Bumpers, Pair Other sizes available
- Standard

Please provide the seat pan to heel measurement required, minimum 11", maximum 18":

Accessories:

☐ Small Calf Pads, Pair 5.35"W-6"H

or

☐ Large Calf Pads, Pair 5.86"Wx8"H

£88.00

☐ Gel Padded Adductor Buttons, Pair

£30.00

Operation Method: ☐ Independent ☐ Combined

Additional Legrest Options

- ☐ Stump Board support (each) Includes a 60-90° Swing-Away Legrest hanger
- ☐ Left ☐ Right
- £300.00
- ☐ Residual Limb support (each)
- ☐ Left ☐ Right
- £180.00

If both a left and right support are required, please amend price accordingly

NE+ Non-Expandable Controls

- ☐ NE+ Joystick- No Power Options
- Standard
- ☐ NE+ Joystick - Power Option Ready Required when 1 or 2 power options are selected.
- £125.00
- ☐ Mount Joystick as Attendant Control Only Please select NE+ Joystick option above
- £100.00
- ☐ Dual Control Please select NE+ Joystick option above
- £380.00

Q-Logic 3 Expandable Controls (includes Bluetooth)

- ☐ Q-Logic 3 Joystick- 2 Power Option Ready *Required when 2 power options are selected. Not compatible with iLevel* £350.00
- ☐ Q-Logic 3 Joystick- 5 Power Option Ready *Required when more than 2 power options are selected. Required for iLevel* £550.00
- ☐ Bluetooth enabled *Please tick this box if you would like the Bluetooth function turned on.* Standard
- ☐ Mount Joystick as Attendant Control Only *Please select Q-Logic 3 Joystick option above* £100.00
- ☐ Dual Control *Please select Q-Logic 3 Joystick option above* £430.00
- ☐ SCIM Module Kit *Allows use of joystick display for specialty controls. Q-Logic 3 Joystick must be selected above* £316.00
- ☐ Q-Logic 3 Enhanced Display *Includes x2 egg switches and flexible display mount. Does not include Q-Logic 3 Joystick* £1074.00
- ☐ iAccess OLED 4 button switch module *Supplied as standard with iLevel seat lift* £300.00



☐ Push Button iAccess Module



☐ Toggle lever iAccess Module



Joystick Shroud Colour Options:

- ☐ Black, Soft Touch Standard
- ☐ Red Corvette ☐ White Satin ☐ Deep Purple ☐ Grey Street £25.00
- ☐ Midnight Blue ☐ Yellow Blaze ☐ Pretty in Pink ☐ Lime Light
- ☐ Silver Bullet ☐ Orange Crush ☐ Electric Blue ☐ Back in Black

Enhanced Display Shroud Colour Options:

- ☐ Black, Soft Touch Standard
- ☐ Red Corvette ☐ White Satin ☐ Deep Purple ☐ Grey Street £25.00
- ☐ Midnight Blue ☐ Yellow Blaze ☐ Pretty in Pink ☐ Lime Light
- ☐ Silver Bullet ☐ Orange Crush ☐ Electric Blue ☐ Back in Black

Joystick Handles

- ☐ Standard Carrot Handle Standard
- ☐ Mushroom Handle ☐ Extended I-Handle ☐ Bodypoint 3 Inch Goal Post Handle £35.00 (each)
- ☐ T-Bar ☐ Soft Ball Handle ☐ Bodypoint 4 Inch Goal Post Handle
- ☐ Chin Cup ☐ Medium Ball Handle

Joystick Mounting Brackets

- ☐ Fixed, Inline Standard
- ☐ Swing-Away £75.00
- ☐ Multi-Axis £125.00
- ☐ Retract 4 Swing-Away Joystick Arm Kit *Full 180° Movement* £150.00

Mounting Position

- ☐ Left or ☐ Right
- ☐ Inline or ☐ Offset *Not available with Fixed Joystick Bracket*

Speed Program and Lights Package

- ☐ 4 Mph program *Available with NE+ or Q-Logic 3. Lights and indicators not included* Standard
- ☐ 4 Mph program with Lights and indicators £400.00
- ☐ 6 Mph program with Lights and indicators £400.00
- ☐ Lights and indicator kit surcharge *If no power options have been selected, you must tick this box in addition to the speed program and lights package you have selected above.* £260.00
- ☐ Mud guard lights (Fender Lights) *Includes single switch button for operation. Not road legal. Does not include indicators* £190.00

Quantum Accessories

Some accessories will be shipped with the unit, not installed. Multiple accessories can be mounted along the various tracks, some interference may occur.

- ☐ USB Mobile Device Charger *Please specify fitting location i.e. left seat rail or right seat rail* £110.00
- ☐ Glove Box *Supplied uninstalled* £30.00

Quantum Accessories (Continued)

<input type="checkbox"/>		Clothing Guard (pair) <i>Only available with Flip Up armrests</i>	£100.00
<input type="checkbox"/>		Quantum Backpack	£40.00
<input type="checkbox"/>		Swing-away tray table	£350.00

Price

Sub Total: £4,595.00


Options:

Total:


Special Instructions

User Measurements

Use this section to record users measurements.



- A. Shoulder Width _____
B. Chest Width _____
C. Hip Width _____
D. External Knee Width _____



- E. Maximum Sitting Height _____
F. Shoulder Height _____
G. Axilla Height _____
H. PSIS Height _____
I. Trunk Depth _____
J. Ischial Depth _____
K. Thigh Depth _____
L. Lower Leg Length _____
M. Elbow Height _____

For additional seating options such as thoracic supports and alternative drive controls please refer to secondary order forms, entitled 'Additional Seating Options' or 'Specialty Control Input Devices'.

Quantum Notes:

1. This product is suitable for occupied transportation in accordance with ISO 7176-19. We strongly recommend that all heavy accessories such as manual ELRs are removed and secured during occupied transportation. We also advise that the use of third party accessories including life Support Systems, Chin Controls and Trays be risk assessed prior to the wheelchair being transported. Full details of suitable restraints for this wheelchair can be obtained from Quantum Technical on 01869 324600 or Email us on: technical@quantumrehab.co.uk
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3 %. Overall measurements will vary based on seating and accessory selections.