

OUTBACK WITH CONTOURED SEATING

DEALERS DETAILS

Dealers Name			
Order Date			
Address			
	Postcode		
Phone		Email	
PO No		Order Ref	

DELIVERY DETAILS

Dealers Address			
	Postcode		
Quantum Specialist			
Ⓢ Clients Weight			

BASE MODEL

£11,495

7 Colour Options
 Flip-up Arm Supports
 Centre Mount Foot Platform
 Contoured Seating
 15" Pneumatic Drive Tyres

SRS Suspension System
 Q-Logic 3e Control System
 2 x 95 ah Batteries
 6 Mph (9.65 Km/h) 4-Pole Motors
 Max User Weight: 25 st (159 kg)



Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

COLOURS

	Cadet Matte	Standard
	Tan Matte	FOC
	Desert Camo	£150
	Black Camo	£150
	White Camo	£150



Please download our free PQ-365 app to view colours and other product features in more detail.

CONTOURED SEAT PACKAGES

	16"W x 16"D	FOC		20"W x 16"D*	FOC
	16"W x 18"D*	FOC		20"W x 18"D	FOC
	16"W x 20"D*	FOC		20"W x 20"D	FOC
	18"W x 16"D	FOC		20"W x 22"D*	FOC
	18"W x 18"D	Standard		22"W x 18"D*	FOC
	18"W x 20"D	FOC		22"W x 20"D	FOC
	18"W x 22"D*	FOC		22"W x 22"D	FOC

*Extended lead time. All seat packages above include a solid seat pan, contoured seat cushion, and a standard head support.

STEALTH SEAT CUSHIONS

No Seat Cushion Required.	FOC
Contoured Seat Cushion. <i>136 Kg weight limit.</i>	Standard
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	£75
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	£100
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£115
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£115
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£135
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	£135
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190

Please note that seat cushions may not be available in all width and depth combinations.
Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

Static Seat.	Standard
Powered Seat Lift. <i>10" lift and up to 3.2 Mph drive speed while elevated. Uses 1 power option.</i>	£1995
Seat to Floor Height. <i>From seat pan. Only the lowest seat to floor height can be selected with powered seat lift.</i>	19.5" 20.5" 21.5"

ARM SUPPORTS

	Std Flip-up Height Adjustable. <i>Includes 14" L contour arm pads. Do not select arm pads below.</i>				Standard
	TB Flex Height Adjustable, Removable, Drop-in T Arm. <i>Adjustment range 9" to 14" from seat pan.</i>				£100
	Single Post Height Adjustable. <i>Adjustment range 11"-17" from seat pan.</i>				£150
	Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i>				£150
	HD Height Adjustable, Removable. <i>Includes HD arm pads. Do not select arm pads below.</i>		8"-12"H.	10"-14"H.	£220
	Bariatric Removable Height Adjustable. <i>Includes Bariatric arm pads. Do not select arm pads below.</i>				£250
Arm Support Pads.	10" Standard	STD	10" Waterfall.	<i>Only with single post or paediatric arms.</i>	FOC
	14" Standard	FOC	14" Waterfall.	<i>Only with single post or paediatric arms.</i>	FOC
	10" Falcon Gel	£40	Falcon Gel Mounting Plates. <i>Required for Single Post & Paediatric arm supports.</i>		£40
	14" Falcon Gel	£40			

CENTRE MOUNTED LEG SUPPORTS

Basic Centre Mount Foot Platform						Standard			
Lower Leg Extensions. <i>From seat pan.</i>		8"-10". <i>Only available with 16" seat width.</i>	FOC		10"-14".	Standard		14"-19".	FOC
Footplates.		10"W x 9"D Split Footplates.		Standard		10"W x 12"D Split Footplates.			FOC
Calf Support.		Individual Calf Pads. <i>6"W x 6"H.</i>		Standard		Calf Panel. <i>18"W x 8"H.</i>			£43
		Individual Large Calf Pads. <i>6"W x 8"H. Only with long lower ext.</i>		FOC					
Footplate Accessories.		Footplate Connector.							£10

INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

	70° HD Drop-in Swing-away Legrests. <i>Includes heel loops.</i>	£340
70°	Lower Leg Extensions. <i>From seat pan.</i>	8.5"-12" FOC 13"-18" Standard 16"-21" FOC
	Footplate Size.	6"W x 8"D Standard 6"W x 11"D FOC 8"W x 11"D FOC
	Accessories.	Wheel Bumpers £10

The 8"W footplates may cause interference with seat widths less than 19".

STEALTH RESIDUAL LIMB SUPPORTS

LEFT	£300	Stealth Residual Limb Support.	£300	RIGHT
-------------	------	--------------------------------	------	--------------

Each limb support includes a 60-90° Swing-away leg support hanger for mounting.

PELVIC BELTS

Standard Pelvic Belts.	Auto Buckle Lap Belt, Push Button Release.	Standard
	38".	FOC
	50".	Standard
	60".	FOC
	70".	FOC
Stealth Positioning 2-Point Padded Pelvic Belt, Single Pull, Push Button Release	Medium	£78
	Large	£78
	Extra Large	£78
Retractable Lap Belt	Push Button Release.	£156

STEALTH THIGH & KNEE SUPPORTS

UniLink Thigh Supports.	LEFT	£176	UniLink Hip Guide Long Removable.	£176	RIGHT
		£176	Unilink Hip Guide Short Removable.	£176	
		£48	4x4 Pelvic Thigh Pad.	£48	
		£48	4x8 Pelvic Thigh Pad.	£48	
		£48	4x12 Pelvic Thigh Pad.	£48	
Knee Adductors.	LEFT	£156	Fixed Cane Mounted.	£156	RIGHT
		£170	Swing-Away Track Mounted.	£170	
		£36	3" Round Knee Pad.	£36	
		£36	4" Round Knee Pad.	£36	

STEALTH ELBOW SUPPORTS

UniLink Elbow Supports. <i>Only available with Single Post, Paediatric, and HD arm supports.</i>	LEFT	£133	Swing-Away, Elbow Support Bracket. <i>ARMS-ES-L ----- ARMS-ES-R</i>	£133	RIGHT
		£38	3"W x 3"H Elbow Pad. <i>ARMS-ES-3X3</i>	£38	
		£38	4"Wx4"H Elbow Pad. <i>ARMS-ES-4X4</i>	£38	

STEALTH FOOT SUPPORTS

Stealth Foot Boxes (Pair).	Small. <i>5"Wx8"Dx7"H.</i>	£361
	Medium. <i>6"Wx10"Dx9"H.</i>	£361
	Large. <i>7"Wx12"Dx11"H.</i>	£361
Stealth Shoe Holders (Pair).	Small Shoe Holder, 3/4" Webbing. <i>IFSH4211</i>	£81
	Medium Shoe Holder, 3/4" Webbing. <i>IFSH4212</i>	£81
	Large Shoe Holder, 3/4" Webbing. <i>IFSH4213</i>	£81

EXPANDABLE CONTROLS Q-Logic 3e and 3 EX hand controls are supplied with built-in Bluetooth.

Q-Logic 3e Controls. <i>No power options.</i>	Standard
Q-Logic 3e Controls. <i>3 power option ready.</i>	£500
Q-Logic 3 EX Controls. <i>No power options.</i>	£350
Q-Logic 3 EX Controls. <i>3 power option ready.</i>	£550
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i> Standard

HAND CONTROL MOUNTS

	Fixed.						Standard		
	Swing-away.						£75		
	Advanced Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i>						£100		
Mounting Position.			Left	FOC		Right	Standard		
Hand Control Joystick Handle.		Standard Carrot	Standard		T-Bar	£25		Extended I-Handle	£25
		Mushroom	£25		Medium Ball	£25		3" Goal Post	£50
		Chin Cup	£25		Soft Ball	£25		4" Goal Post	£50

SPEED PROGRAM & LIGHTS PACKAGE

6 Mph Program with Lights and indicators.	Standard
---	----------

ACCESSORIES Supplied uninstalled.

Cup Holder.	£20
Mobile Phone Holder.	£40
Clothes Guards (Pair).	£100
Crutch/Cane Holder.	£170
Glove Box.	£30
Oxygen Holder.	£50
Transfer Handles.	£125
Stealth Tablet Holder.	£342
Large Black Weather Cover.	£125
Rear View Mirror.	£15
Quantum Back Pack.	£110
Rear Basket.	£40
Saddle Bag.	£25
Additional XLR USB Charger Port.	£17
Magnetic to XLR Charger Adapter. <i>Only with Q-Logic 3e.</i>	£49
Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if requested.</i>	£379

NOTES & SPECIAL REQUESTS

OUTBACK WITH CONTOURED SEATING



BASE PRICE	£11,495
COST OPTIONS	
TOTAL	

5



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
RETAIL UK - 28/01/25 V1

