

OUTBACK WITH TB3 SEATING

DEALERS DETAILS

Dealers Name	
Order Date	
Address	
	Postcode
Phone	Email
PO No	Order Ref

DELIVERY DETAILS

Dealers Address	
	Postcode
Quantum Specialist	
Ⓢ Clients Weight	

BASE MODEL

£15,950

7 Colour Options
Height Adjustable Arm Supports
Centre Mount Foot Platform
TRU-Balance 3 Seating
15" Pneumatic Drive Tyres



SRS Suspension System
Q-Logic 3e Control System
2 x 95 ah Batteries
6 Mph (9.65 Km/h) 4-Pole Motors
Max User Weight: 21st 6 lbs (136 kg)

Standard = Standard Factory Feature
FOC = Select Instead of Standard Factory Feature

COLOURS

	Cadet Matte	Standard
	Tan Matte	FOC
	Desert Camo	£150
	Black Camo	£150
	White Camo	£150



Please download our free PQ-365 app to view colours and other product features in more detail.

SEAT PACKAGES - SEAT WIDTH

	16"-20"W.
	20"-22"W.
Set width at:	

SEAT PACKAGES - SEAT DEPTH

	16"-18"D.
	19"-22"D.
Set depth at:	

STEALTH SEAT CUSHIONS

TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	Standard
No Seat Cushion Required	FOC
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	FOC
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£115
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£115
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£135
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	£135
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190

Please note that seat cushions may not be available in all width and depth combinations.
Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

TRU-Balance 3 Powered Seat Tilt. <i>0° to 50° COG tilt. Uses 1 power option.</i>	Standard
10" Powered HD Seat Lift & Tilt. <i>10" lift & 50° tilt. Up to 3.2 Mph drive speed while elevated. Uses 2 power options.</i>	£2170
Powered Tilt Seat to Floor Height. <i>From seat pan.</i>	19" 20" 21"
Powered HD Lift & Tilt Seat to Floor Height.. <i>From seat pan.</i>	19.5" 20.5" 21.5"

BACKREST SPECIFICATIONS

TRU-Balance 3 Back Support System. <i>Select TB3 back options below. Ignore all 'Back Cane' options.</i>	Standard
Back Canes Only System. <i>Not available with Power recline or TB3 flip-back arm supports. Move directly to 'Back Cane Only Support Options' if selected.</i>	FOC

Select one option above. Both options above include a metal active contour back support shell, TRU-Comfort 2 back cushion, and a pair of plastic lateral wedge inserts

TB3 BACK SUPPORT

Active Contour Back Support with Padded Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions</i>	Standard
Privacy Flap. <i>12.75"W. Closes gap between back support and seat pan.</i>	£50
Stealth Adapter Interface Plate. <i>SSAPT4. For mounting third party back supports to TB3 back system. Not required for Q-Back 2.</i>	£341
ADI Q-Back 2. <i>Skin protection, 4 way stretch, Coolcore cover. Mounts directly to TB3 system.</i>	POA
Enter ADI Q-Back 2 part number here for quotation.	

TB3 BACK SUPPORT POSTIONING

TRU-Balance 3 Static Back. <i>85° to 105° angle adjustment via gas strut within back cover.</i>	Standard
TRU-Balance 3 Flip-forward Back. <i>Not suitable for occupied transport. Not available with TB3 flip-back arm supports or powered recline.</i>	£150
TRU-Balance 3 Powered Recline. <i>85° to 175° with sheer reduction. Uses 1 power option.</i>	£600

BACK CANE ONLY SUPPORT

	Active Contour Back Support with Padded Back Cushion and Wedge Inserts. <i>Select size in Back Support Dimensions.</i>			Standard
	Privacy Flap. <i>12.75"W. Closes gap between back support and seat pan.</i>			£50
	ADI Back Support. <i>Skin protection, 4 way stretch, Coolcore cover. Mounts to canes with cane clamps.</i>			POA
Enter ADI Back Support part number here for quotation				
Back Cane Only Support Positioning.		Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>		Standard
		Flip-forward Back Canes. <i>Not suitable for occupied transport or with Endo mounts.</i>		£200
		2" Endo Cane Mounts. <i>1" each side. Reduces back width by 2". Minimum seat width 12".</i>		£40
		4" Endo Cane Mounts. <i>2" each side. Reduces back width by 4". Minimum seat width 14".</i>		£40
Back Cane Height. <i>From seat pan.</i>		16"H.	17"H.	20"H.
		21"H.	24"H.	25"H.

BACK SUPPORT DIMENSIONS

Back width must be within +/- 4" of seat width when TB3 Back Support has been selected, and must match the seat width when the Back Canes Only option has been selected. The only exception is when the 2" or 4" Endo mounts have been selected, in which case the back can be 2" or 4" narrower than the seat width.

Backrest Support Width Package.		12"-14"W.	14"-16"W.	16"-20"W.	18"-22"W.
		20"-24"W.			
Back Support Width.	12"W.	14"W.	15"W.	16"W.	17"W.
	18"W.	19"W.	20"W.	21"W.	22"W.
Back Support Height. <i>From seat pan.</i>	14"H.	15"H.	16"H.	17"H.	18"H.
	19"H.	20"H.	21"H.	22"H.	23"H.
	24"H.	25"H.	26"H.	27"H.	28"H.
	29"H.	30"H.	31"H.	32"H.	

19" minimum backrest height when powered recline is selected.

BACK WIDTH STATIC SEAT OR POWER TILT

	12"W	14"W	15"W	16"W	17"W	18"W	19"W	20"W	21"W	22"W
14"H	T									
15"H	T									
16"H	S	S	T	S	T	S	T	S	T	S
17"H	S	S	T	S	T	S	T	S	T	S
18"H		S	T	S	T	S	T	S	T	S
19"H		S	Q	T	S	Q	T	S	Q	T
20"H		S	Q	T	S	Q	T	S	Q	T
21"H		S	Q	T	S	Q	T	S	Q	T
22"H		S	Q	T	S	Q	T	S	Q	T
23"H		S	T	S	T	S	Q	T	S	Q
24"H		S	T	S	T	S	Q	T	S	Q
25"H		S	T	S	T	S	Q	T	S	Q
26"H		S	T	S	T	S	Q	T	S	Q
27"H		S	T	S	T	S	Q	T	S	Q
28"H		S	T	S	T	S	Q	T	S	Q
29"H		S	T	S	T	S	Q	T	S	Q

BACK WIDTH POWER RECLINE

	12"W	14"W	15"W	16"W	17"W	18"W	19"W	20"W	21"W	22"W
19"H	S	S	T	S	T	S	T	S	T	S
20"H	S	S	T	S	T	S	T	S	T	S
21"H		S	T	S	T	S	T	S	T	S
22"H		S	Q	T	S	Q	T	S	Q	T
23"H		S	Q	T	S	Q	T	S	Q	T
24"H		S	Q	T	S	Q	T	S	Q	T
25"H		S	Q	T	S	Q	T	S	Q	T
26"H		S	T	S	T	S	Q	T	S	Q
27"H		S	T	S	T	S	Q	T	S	Q
28"H		S	T	S	T	S	Q	T	S	Q
29"H		S	T	S	T	S	Q	T	S	Q
30"H		S	T	S	T	S	Q	T	S	Q
31"H		S	T	S	T	S	Q	T	S	Q
32"H		S	T	S	T	S	Q	T	S	Q

Back Type.	Tru-Comfort 2	Sport Back	ADI/Q-Back 2
------------	---------------	------------	--------------

ARM SUPPORTS

TB Flex Height Adjustable, Removable, Drop-in T Arm. <i>Adjustment range 9" to 14" from seat pan.</i>	Standard
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan.</i>	£279
Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"- 9.75" from seat pan. 10" pads only.</i>	£100
HD Height Adjustable, Removable. <i>Includes HD arm pads. Do not select arm pads below.</i>	8"-12" H. 10"-14" H. £220
TB3 Flip-back. <i>Adjustment range 8"-16" from seat pan. Includes lateral in/out adjustment.</i>	£350
Down Posts for TB3 Flip-back. <i>Only with TB3 Flip-back arm supports. Minimum seat/back width 16".</i>	£100
Arm Support Pads.	10" Standard <i>Not available with HD arms.</i> STD 10" Falcon Gel £40 10" Waterfall. <i>Not available with TB Flex or HD arms.</i> FOC
	14" Standard <i>Not available with HD arms.</i> FOC 14" Falcon Gel £40 14" Waterfall. <i>Not available with TB Flex or HD arms.</i> FOC
	17" Standard. <i>TB3 only. Minimum seat depth 16".</i> £50 Falcon Gel Mounting Plates. <i>Required for 2-post, Single Post & Paediatric arm supports.</i> £40

BASIC CENTRE MOUNTED LEG SUPPORTS

Basic Centre Mount Foot Platform.							Standard
Lower Leg Extensions. <i>From seat pan.</i>	8"-10".	FOC	10"-14".	Standard	14"-19".	FOC	
Split Footplates.	10"W x 9"D.		Standard	10"W x 12"D.	<i>Tapered footplates. Not available with 8"-10" lower extensions.</i>	FOC	
Calf Support.	Individual Calf Pads. 6"W x 6"H.	Standard	Calf Panel. 18"W x 8"H.	£43			
	Individual Large Calf Pads. 6"W x 8"H. <i>Only available with 14"-19" lower extensions.</i>	FOC					
Footplate Accessories.	Footplate Connector.	£10					

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

POWERED CENTRE MOUNTED LEG SUPPORTS

Power Elevating Articulating Centre Mount Foot Platform. <i>Uses 1 power option.</i>							£700
Lower Leg Extensions. <i>From seat pan.</i>	11"-17".	FOC	14"-20".				Standard
Split Footplates.	10"W x 6"D.	FOC	10"W x 8"D.	FOC	10"W x 11"D.		Standard
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>	Standard		Small Calf Panel. <i>18"W x 8"H.</i>	£43		
	Individual Large Calf Pads. <i>6"W x 8"H.</i>	FOC		Large Calf Panel. <i>24"W x 8"H.</i> <i>Only available with seat widths 20"-22".</i>	£85		
	Reverse Dartex Calf Pad Covers.	£42		Reverse Dartex Calf Panel Cover.	£62		
Footplate Accessories.	4" Width Extension Kit. <i>Not available with 2" length kit or 10"W x 6"D footplates.</i>						£16
	2" Length Extension Kit. <i>Not available with 4" width kit or 10"W x 6"D footplates.</i>						£16
	Footplate Connector.	£10	Outer Edge Supports.	£10			
	Wheel Bumpers.	£10	Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.				

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

	Power Elevating Articulating Swing-Away Leg Supports. <i>Includes calf pads and heel loops. Uses 2 power options.</i>						£1300
POWER ELEVATING ARTICULATING	Lower Leg Extensions. <i>From seat pan.</i>	 8.5"-14.5"	FOC	 12.5"-21"	Standard		
	Footplate Type.	 Angle & Depth Adjustable. <i>Select size below.</i>	Standard				
	Footplate Size.	 4"W x 6"D	FOC	 5"W x 6"D	FOC	 5"W x 8"D	FOC
		 6"W x 8"D	Standard	 6"W x 11"D	FOC	 8"W x 11"D	FOC
	Accessories.	 Small Calf Pads	Standard	 Wheel Bumpers	£10		
		 Large Calf Pads	FOC	 Gel Padded Knee Adductors	£60		
	Leg Support Operation Method.	 Independent	 Combined				
	70° HD Drop-in Swing-away Legrests. <i>Includes heel loops.</i>						£340
70°	Lower Leg Extensions. <i>From seat pan.</i>	 8.5"-12"	FOC	 13"-18"	Standard	 16"-21"	FOC
	Footplate Size.	 6"W x 8"D	Standard	 6"W x 11"D	FOC	 8"W x 11"D	FOC
	Accessories.	 Wheel Bumpers	£10				

The 8"W footplates may cause interference with seat widths less than 19".
Certain configurations may result in reduced ground clearance below footplates.
Contact Quantum Technical for information.

STEALTH RESIDUAL LIMB SUPPORTS

LEFT	£300	Stealth Residual Limb Support.	£300	RIGHT
Each limb support includes a 60-90° Swing-away leg support hanger for mounting.				

STEALTH TRUNK SOLUTIONS

UniLink Thoracic Lateral Supports. <i>Supplied as a pair.</i>		Lateral Thoracic Swing-away Hardware.	£425		Additional 1.5" Link (each).	£52									
		3" H x 5" L Curved Lateral Pad.	£72		Additional 2" Link (each).	£52									
		4" H x 6" L Curved Lateral Pad.	£72		Additional 3" Link (each).	£52									
4-Point Harness.		Extra Small	£137	}		Classic	STD*	}		Dynamic	STD*	}		Vinyl Face	FOC
		Small	£137			Contoured	FOC			Structured	FOC			Zipper Face**	FOC
		Medium	£137												
		Large	£137												
		Extra Large	£137												
Chest Strap.		Extra Small	£87		Small	£87	}		Side Release	Standard					
		Medium	£87		Large	£87			Velcro Release	£19					
		Extra Large	£87												

*STD = Standard option **Classic Only

STEALTH PELVIC & THIGH SOLUTIONS

UniLink Thigh Supports.	LEFT	£176	UniLink Hip Guide Long Removable.	£176	RIGHT		
		£176	Unilink Hip Guide Short Removable.	£176			
		£48	4x4 Pelvic Thigh Pad.	£48			
		£48	4x8 Pelvic Thigh Pad.	£48			
		£48	4x12 Pelvic Thigh Pad.	£48			
Knee Adductors.	LEFT	£156	Fixed Cane Mounted.	£156	RIGHT		
		£170	Swing-Away Track Mounted.	£170			
		£36	3" Round Knee Pad.	£36			
		£36	4" Round Knee Pad.	£36			
Pelvic Belts.	Auto Buckle Lap Belt, Push Button Release.	Standard					
	38".	FOC	70".	FOC			
	50".	Standard	80".	FOC			
	60".	FOC	90".	FOC			
Retractable Lap Belt.	Push Button Release.	£156					
Positioning 2-Point Padded Pelvic Belt.	Extra Small	£78	Single Pull	Standard	Side Release	Standard	
	Small	£78		FOC		Push Button	FOC
	Medium	£78		FOC			
	Large	£78					
	Extra Large	£78					
Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.	Extra Small	£120	Small	£120	Push Button	Standard	
	Medium	£120	Large	£120	Side Release	£20	
Medial Supports (Pommels).	Fixed Flip-down Hardware. <i>MTHW-660.</i>					£180	
	Small Wedge. <i>MTMD-12.</i>			£103	Medium Wedge. <i>MTMD-13.</i>	£103	
	Large Wedge. <i>MTMD-14.</i>			£103	Adjustable Wedge. <i>MTMD-42.</i>	£199	

Please visit stealthproducts.com for alternative lateral, pelvic & lower extremity solutions.

STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>	£219
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£219
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£219
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	£327
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£88

Please visit stealthproducts.com for alternative head supports.

STEALTH ELBOW SOLUTIONS

UniLink Elbow Supports.	LEFT	£133	Swing-Away, Elbow Support Bracket. <i>ARMS-ES-L ----- ARMS-ES-R</i>	£133	RIGHT
		£38	3"W x 3"H Elbow Pad. <i>ARMS-ES-3X3</i>	£38	
		£38	4"Wx4"H Elbow Pad. <i>ARMS-ES-4X4</i>	£38	

STEALTH FOOT SOLUTIONS

Stealth Foot Boxes (Pair).	Small. <i>5"Wx8"Dx7"H.</i>	£361
	Medium. <i>6"Wx10"Dx9"H.</i>	£361
	Large. <i>7"Wx12"Dx11"H.</i>	£361
Stealth Shoe Holders (Pair).	Small Shoe Holder, 3/4" Webbing. <i>IFSH4211</i>	£81
	Medium Shoe Holder, 3/4" Webbing. <i>IFSH4212</i>	£81
	Large Shoe Holder, 3/4" Webbing. <i>IFSH4213</i>	£81

EXPANDABLE CONTROLS

Q-Logic 3e and 3 EX hand controls are supplied with built-in Bluetooth. Q-Logic 3 Enhanced Display is supplied with built-in Bluetooth and Infrared.

Q-Logic 3e Controls. <i>3 power option ready.</i>	Standard
Q-Logic 3e Controls. <i>5 power option ready.</i>	£500
Q-Logic 3 EX Controls. <i>3 power option ready.</i>	£350
Q-Logic 3 EX Controls. <i>5 power option ready.</i>	£550
Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches & flexible display mount.</i>	£1074
SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i>	£316

iAccess OLED 4 Button Switch Module.
*Mounted to opposite armrest unless specified.
Supplied FOC when iLevel is selected*

£300

Push Button iAccess
Toggle Lever iAccess

If iAccess is selected choose either Push Button or Toggle lever option.

Control Method. Hand Control. *Select position in Hand Control Mounts below.*

Standard

HAND CONTROL MOUNTS

	Fixed.							Standard
	Swing-away.							£75
	Advanced Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i>							£100
Hand Control Joystick Handle.	Mounting Position.		Left	FOC		Right	Standard	
			Standard Carrot	Standard		T-Bar	£25	
			Mushroom	£25		Medium Ball	£25	
			Chin Cup	£25		Soft Ball	£25	
							Extended I-Handle	£25
							3" Goal Post	£50
							4" Goal Post	£50

SPEED PROGRAM & LIGHTS PACKAGE

6 Mph Program with Lights and indicators.	Standard
---	----------

ACCESSORIES Supplied uninstalled.

USB Mobile Device Charger. <i>Please specify mounting position. i.e arm support, backrest, seat rail</i>	£110
Cup Holder.	£20
Mobile Phone Holder.	£40
Clothes Guards (Pair).	£100
Crutch/Cane Holder.	£170
Glove Box.	£30
Oxygen Holder.	£50
Transfer Handles.	£125
Stealth Tablet Holder.	£342
Large Black Weather Cover.	£125
Back-up Camera.	£300
Rear View Mirror.	£15
Quantum Back Pack.	£110
Additional XLR USB Charger Port.	£17
Magnetic to XLR Charger Adapter. <i>Only with Q-Logic 3e.</i>	£49
Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if requested.</i>	£379

NOTES & SPECIAL REQUESTS

PRICE

OUTBACK WITH TB3 SEATING



BASE PRICE	£15,950
COST OPTIONS	
TOTAL	



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
RETAIL UK - 28/01/25 V1

