

Q300[®] M Mini SEDEO Pro



ORDER FORM

Billing Address Or Customer-ID		Shipping Address		☐ Identical To Billing Address
C.-ID:		C.-ID:		
Name:		Name:		
Street:		Street:		
Town/PC:		Town/PC:		
E-Mail:		E-Mail:		
Phone:		Phone:		

Order Information

Purchase Order No.:	
Order Date:	Day: Month: Year:

Base Price **£6,750**

Seating Specifications

All Bases Prepared For DAHL Docking Station		Seat Type	SEDEO PRO
Typical WC Weight	Without Batteries 135 kg	Seat Size Range	Width 380 - 560 mm, Depth 400 - 560 mm
	With Batteries 165 kg	Seat To Floor Height Range	410 mm - 490 mm

MWD Base Specifications

User Weight Capacity	136 kg	Drive Wheels	12" or 14" Drive Wheels (Black)
Base Size	Width Min. 520 mm	Tyres / Castors	7" Front / 7" Rear (Black)
Drive Control	VR2, Option R-Net	Charger	8 A Dual Mode Charger
Driving Range	Up To 31 km	Battery Size	41 & 56 Ah Batteries
Drive Motors	2-Pole 6 /10 kph, 4-Pole 6/8 kph	Colour Accents	7 Colour Accents
Ground Clearance	70 mm	Frame Colour	Black

Q300 M Mini Drive Bases			£
☒	QMP000074	Q300 M Mini (MWD Base Chair)	6750

Seat Cushion Type			£
☐	QMP020110	Omit Seat Cushion	0

SEDEO Seat Size			£
Seat Cushion Width			
☐	SW38	15" (380 mm)	
☐	SW42	16" (420 mm)	
☐	SW46	18" (460 mm)	
☐	SW50	20" (500 mm)	
☐	SW54	21" (540 mm)	
☐	SW56	22" (560mm)	
Seat Depth			
☐	SD40	15.5" (400 mm) (Not with QMP090128 OR QMP090129)	
☐	SD42	16" (420 mm)	
☐	SD44	17" (440 mm)	
☐	SD46	18" (460 mm)	
☐	SD48	19" (480 mm)	
☐	SD50	19.5" (500 mm)	
☐	SD52	20" (520 mm)	
☐	SD54	21" (540 mm)	
☐	SD56	22" (560 mm)	

Seat Rail			£
☐	QMP020085	Seat Rail, Short, For Seat Depth 400 - 520 mm	0
☐	QMP020069	Seat Rail, Long, For Seat Depth 460 - 560 mm (Not with seat width 380 - 420mm)	0

SEDEO Seat Cushion			£
☐	QMP020073	SEDEO, Seat Cushion, Profiled, Height 70 mm	160
☐	QMP020076	SEDEO, Seat Cushion, Support Memory, Height 90 mm	260

SEDEO Seat Cushion Covers			£
☐	QMP020518	Seat Cushion Cover, Incontinence Protection	0
☐	QMP020080	Seat Cushion Cover, Comfort	0

SEDEO Seat Boards			£
☐	QMP020082	SEDEO, Seat Board, Flat, Depth-Adjustable, Cushions Guide At Rear	130

JAY Seat Size			£
Seat Cushion Width			
☐	SW40	15.5" (400 mm)	
☐	SW42	16" (420 mm)	
☐	SW44	17" (440 mm)	
☐	SW46	18" (460 mm)	
☐	SW48	19" (480 mm)	
☐	SW50	19.5" (500 mm)	
☐	SW52	20" (520 mm)	
☐	SW54	21" (540 mm)	
☐	SW56	22" (560 mm)	
Seat Depth			
☐	SD40	15.5" (400 mm)	
☐	SD42	16" (410 mm)	
☐	SD44	17" (440 mm)	
☐	SD46	18" (460 mm)	
☐	SD48	19" (480 mm)	
☐	SD50	19.5" (500 mm)	
☐	SD52	20" (520 mm) (Only with seat width 560mm)	
☐	SD54	21" (540 mm) (Only with seat width 560mm)	

JAY Seat Cushions			£
☐	JAY020006	JAY Easy Visco With Flat Base	180
☐	JAY020011	JAY Balance, Fluid, Microclimatic Cover, Until Seat Width 500 mm	655
☐	JAY020013	JAY Balance, Fluid, Microclimatic Cover, From Seat Width 560 mm	703
☐	JAY020015	JAY Balance, Fluid, Microclimatic Cover, With Positioning Kit, Until Seat Width 500 mm	755
☐	JAY020017	JAY Balance, Fluid, Microclimatic Cover, With Positioning Kit, From Seat Width 560 mm	833

Seat Heights (Measured from front of seat pan to floor)			£
☐	SH41	410 mm (Only with QMP130011 Fixed tilt module)	0
☐	SH43	430 mm (Only with seat depth 400 - 520mm)	0
☐	SH45	450 mm	0
☐	SH47	47 mm	0
☐	SH49	490 mm	0

Powered Seating Options			£
☐	QMP130011	Fixed Tilt Module (Adjustable 0°, 3°, 6°, 9°, Factory Setting 6°)	0
☐	QMP130004	Posterior Tilt, Powered, 0° - 30° (Not with QMP110036) (In case of QMP110004 please select QMP110032)	1150
☐	QMP130012	Posterior Tilt, Powered, -5° - 25° (Only available with seat height 470mm) (Not with QMP110036) (In case of QMP110004 please select QMP110032)	1250
☐	QMP130007	Posterior Tilt, Powered, 0° - 50° (Not with QMP110036) (In case of R-Net please select QMP110032) (In case of QMP110004 please select QMP110032)	1350
☐	QMP130009	Seat Elevator (300 mm) / Tilt Module (Tilt 0° - 30°), Powered (Not with QMP110036) (In case of R-Net please select QMP110032)	2050
☐	QMP130006	Seat Elevator (300 mm) / Tilt Module (Tilt 0° - 50°), Powered (Not with QMP110036) (In case of R-Net please select QMP110032)	2250

Backrest Frame Type			£
☐	QMP030083	Backrest, Tool Angle Adjustable, 90° - 120°	0
☐	QMP030084	Backrest, Manual Recline, User Adjustable, 90° - 130°, Release Lever, Right (Not on the same side as the joystick)	0
☐	QMP030085	Backrest, Manual Recline, User Adjustable, 90° - 130°, Release Lever, Left (Not on the same side as the joystick)	0
☐	QMP030086	Backrest, Powered Recline, 90° - 130° (Not with QMP110036) (In case of QMP110004 please select QMP110032)	560
☐	QMP030087	Backrest, Powered Recline, 90° - 150° With 100 mm Mechanical Anti-Shear (Not with QMP110036) (In case of QMP110004 please select QMP110032)	600
☐	QMP030229	Backrest, Powered Recline, 90° - 150° With 100 mm Mechanical Anti-Shear with Relax function (Not with QMP110036) (In case of QMP110004 please select QMP110032) (Only with QMP130006)	600

Backrest Frame Option			£
☐	QMP030098	Backrest, Flip Down, Quick Release	185

Push Handles			£
☐	QMP030204	Push Handles, Height-Adjustable	210

SEDEO Pro Back Size			
		Backrest Width	Backrest Height
☐	BW38	15" (380 mm) - Effective width 38cm	☐ BH44 18" (440 mm) - Effective Height 44 - 50cm
☐	BW42	16" (420 mm) - Effective width 38 - 42cm	☐ BH50 19.5" (500 mm) - Effective Height 50 - 56cm
☐	BW46	18" (460 mm) - Effective width 42 - 46cm	☐ BH56 22" (560 mm) Effective Height 56 - 62cm
☐	BW50	19.5" (500 mm) - Effective width 46 - 50cm	
☐	BW54	21" (540 mm) - Effective width 50 - 54cm	
☐	BW56	22" (560 mm) - Effective width 52 - 56cm	

SEDEO Pro Backrest		£
☐	QMP030109 SEDEO PRO, Back Cushion, Profiled	140
☐	QMP030113 SEDEO PRO, Back Cushion, Profiled, Adjustable Straps	260
☐	QMP030114 SEDEO PRO, Back Cushion, Deep Contour, Passive	335
☐	QMP030116 SEDEO PRO, Back Cushion, Flexible, Deep Contour, Active, Adjustable Straps	450
☐	QMP030117 SEDEO PRO, Back Cushion, Flexible, Extra Deep Contour, Active, Adjustable Straps	475
	SEDEO PRO, Back Cushion, Adaptive (Choose B4Me Requirements Separately)	

SEDEO Pro Backrest Covers		£
☐	QMP030059 Back Cushion Cover, Incontinence Protection	0
☐	QMP030108 Back Cushion Cover, Comfort	0

SEDEO Lateral Thoracic Supports (Only with QMP030109 & QMP030113)		£
☐	QMP030607 SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Right	160
☐	QMP030608 SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Left	160

Headrest Receiver		£
☐	QMP030119 Headrest Receiver, Back Cushion Mount (Only with QMP030109 & QMP030113)	0
☐	QMP030120 Headrest Reveiver, Backrest Frame Mount	0

Headrest Pad		£
☐	QMP030430 SEDEO, Headrest Pad, Small, (LxH) 240 x 150 mm (Only use with QMP030433)	95
☐	QMP030431 SEDEO, Headrest Pad, Small, With Cover, (LxH) 240 x 150 mm (Only use with QMP030433)	110
☐	QMP030432 SEDEO, Headrest Pad, Flexible, (LxH) 400 x 110 mm (Only use with QMP030433)	170
☐	WHT030402 Whitmyer, Headrest Pad, Plush, (LxH) 254 x 127 mm (Only use with WHT030408 or WHT030409)	170
☐	WHT030404 Whitmyer, Headrest Pad, Adjust-A-Plush, Small, (LxH) 355 x 76 mm (Only use with WHT030408 or WHT030409)	280

Headrest Stem (Chose a Whitmyer headrest mount and headrest receiver in case you would like to add a head control)		£
☐	QMP030433 SEDEO, Headrest Mount	140
☐	WHT030408 Whitmyer, Headrest Mount, LINX²	190
☐	WHT030409 Whitmyer, Headrest Mount, AXYS	155

JAY J3 Back Size		£
☐	JAY030111 JAY J3 Back Size, 410 mm	0
☐	JAY030112 JAY J3 Back Size, 460 mm	0
☐	JAY030113 JAY J3 Back Size, 510 mm	0

JAY J3 Back Size, Mid Deep Contour		£
☐	JAY030042 JAY J3 Back, Mid Deep Contour, Upper Thoracic, Medium, 460 mm	870

Armrest Types & Side Gaurds		£
☐	QMP040070 Armrest, Reclining, Multi Adjustable, Flip Back	0
☐	QMP040019 Armrest, Single Post, Flip Back (Not with seat height 410 & 430mm)	0
☐	QMP040021 Armrest, Single Post, Removable (Not with seat height 410 & 430mm)	240
Side Guard for Reclining Flip Back Armrest		
☐	QMP040036 Side Guard for reclining flip back arm	0
Side Guards for Single Post Flip Back / Removable		
☐	QMP040033 Side Guard	0
☐	QMP040034 Side Guard, Integral Foam (Seat width -20mm, not with QMP040123)	0
☐	QMP040035 Side Guard, Wedge Shaped, Padded (Seat width -80mm, only with QMP040121)	275

Arm Pads		£
☐	QMP040023 Arm Pad, Integral Foam, 300 mm, Right	0
☐	QMP040024 Arm Pad, Integral Foam, 300 mm, Left	0
☐	QMP040025 Arm Pad, Integral Foam, 400 mm, Right	0
☐	QMP040026 Arm Pad, Integral Foam, 400 mm, Left	0
☐	QMP040027 Arm Pad, Integral Foam, 450 mm, Right (With QMP040019, seat width +20mm)	0
☐	QMP040028 Arm Pad, Integral Foam, 450 mm, Left (With QMP040019, seat width +20mm)	0
☐	QMP040204 Arm Pad, Comfort, Padded, PU Coated, 400 mm, Right	75
☐	QMP040203 Arm Pad, Comfort, Padded, PU Coated, 400 mm, Left	75
☐	QMP040031 Arm Pad, Hemiplegia, Integral Foam, Right (Not with QMP110052, QMP090040, QMP090046 OR QMP090110)	275
☐	QMP040032 Arm Pad, Hemiplegia, Integral Foam, Left (Not with QMP110053, QMP090040, QMP090046 OR QMP090118)	275

Thigh Supports & Pommel		£
☐	QMP020342 Knee Support, Lateral, 120 x 120 mm, Swing Away, Right (Not With QMP050015-17)	185
☐	QMP020341 Knee Support, Lateral, 120 x 120 mm, Swing Away, Left (Not With QMP050015-17)	185
☐	QMP020113 Thigh Support, Lateral, Small, 120 x 120 mm, Swing Away, Right (Not with QMP040035)	185
☐	QMP020114 Thigh Support, Lateral, Small, 120 x 120 mm, Swing Away, Left (Not with QMP040035)	185
☐	QMP020121 Thigh Support, Lateral, Large, 180 x 120 mm, Swing Away, Right (Not with QMP040035)	185
☐	QMP020122 Thigh Support, Lateral, Large, 180 x 120 mm, Swing Away, Left (Not with QMP040035)	185
☐	QMP030605 Abduction Pommel (Not with centre mount legrests)	255

Centre Mount Legrest (Single Piece Footplate)		£
☐	QMP050086 Centre Mount Legrest, Tool Adjustable, With Flip Up And Angle-Adjustable Footboard (In combination with lift or tilt choose QMP050087 or 050092)	0

Pad Options		
☐	QMP050087 Centre Calf Support, Integral Foam	0
☐	QMP050092 Calf Pad, Integral Foam, Pair	0

Powered Centre Mount (In combination with lift or tilt choose QMP050092)		£
☐	QMP050076 Centre Mount Legrest, Powered, 90° - 180°, 175 mm Extension	950
☐	QMP050071 Centre Mount Legrest, Powered, 90° - 180°, 175 mm Extension, Footplate To Floor (Not with QMP110036, QMP110037 or QMP110040)(In case of R-Net QMP110032 required)	950

Footplates		
☐	QMP050072 Footplate, Dual, Height- And Angle-Adjustable, 260 x 150 mm	0
☐	QMP050073 Footplate, Dual, Height- And Angle-Adjustable, Large, 310 x 170 mm	0

Pad Options		
☐	QMP050092 Calf Pad, Integral Foam, Pair	0

SEDEO Pro Swing Away Legrest		£
☐	QMP050015 SEDEO PRO, Legrest, 70°, Swing Away	0
☐	QMP050016 SEDEO PRO, Legrest, Manual Elevating, Swing Away	375
☐	QMP050017 SEDEO PRO, Legrest, Power Elevating, Swing Away	700
☐	QMP050183 SEDEO PRO, Legrest, Power Elevating, Swing Away, LLL 300mm	700
Hanger Receivers		
☐	QMP050170 SEDEO PRO, Hanger Receivers, High, For Seat Cushion Thickness 70mm and Higher	0
☐	QMP050171 SEDEO PRO, Hanger Receivers, Low, For Seat Cushion Thickness 50mm and Lower	0
Footplates		
☐	QMP050088 Footplate, For Seat Width 385 - 400 mm	0
☐	QMP050089 Footplate, For Seat Width 420 - 460 mm	0
☐	QMP050090 Footplate, For Seat Width 480 - 520 mm	0
☐	QMP050091 Footplate, For Seat Width 540 - 560 mm	0
☐	QMP050081 Footplates, Coupled (Choose also a footplate, only for QMP050015)	0
Calf Pad / Strap		
☐	QMP050084 Heel Loop	0
☐	QMP050085 Calf Strap (Only for QMP050015)	0
☐	QMP050082 Calf Support, Integral Foam, Depth-Adjustable, Pair	100
☐	QMP050083 Calf Support, Padded, Depth-Adjustable, Pair	135
☐	QMP050200 Knee Pad, Integral Foam	75
Amputee Support (3)		
☐	QMP050013 Amputee Support, Padded, Right (Extra on a set swing away	150
☐	QMP050014 Amputee Support, Padded, Left	150
VR2 Controls		£
☐	QMP110036 Remote Joystick, VR2	0
☐	QMP110037 Remote Joystick, VR2, 2 Actuators	0
☐	QMP110040 Remote Joystick, VR2, 2 Actuators And Lights	0
VR2 Power Module & Lighting Module		
☐	QMP110002 Power Module, VR2, 90 A	0
☐	QMP110030 Lighting Module, VR2 (If LED Lights & indicators selected this is required)	0
R-Net Controls		£
☐	QMP110015 Remote Joystick, R-Net, LED, Actuators	0
☐	QMP110016 Remote Joystick, R-Net, LED, Actuators And Lights	0
☐	QMP110017 Remote Joystick, R-Net, Colour LCD	600
☐	QMP110018 Remote Joystick, R-Net, Colour LCD, Heavy Duty	1200
☐	QMP110013 Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared Control Included	750
☐	QMP110023 Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared And Bluetooth Control Included	925
R-Net Power Module & ISM		
☐	QMP110003 Power Module, R-Net, 90 A	0
☐	QMP110004 Power Module, R-Net, 120 A	135
☐	QMP110032 Actuator Control & Lights Module, R-Net, ISM-6L	0
R-Net Gyro		
☐	QMP110009 Traction Control System For R-Net (GYRO Technology)	400
Joystick Mount		£
☐	QMP110050 Joystick Mounting, Fixed	0
☐	QMP110056 Joystick Mounting, Parallel, Swing Away, Capsulated	0

Joystick Mounting Location		£
☐	QMP110052 Joystick Mounting Location, Right	0
☐	QMP110053 Joystick Mounting Location, Left	0

Attendant Controls		£
☐	QMP110021 Dual Control	425
☐	QMP110020 Joystick Mounted As Attendant Control (Please choose joystick)	270

Attendant Mounting Location		
☐	QMP110025 Attendant Mounting Location, Right	0
☐	QMP110026 Attendant Mounting Location, Left	0

Table Top Controls		£
☐	QMP110143 Table Top Control, LCD, Centre Tray Mounted	2235

Tray Table Mounting Location		
☐	QMP110080 Tray Table Mounting Location, Left	0
☐	QMP110081 Tray Table Mounting Location, Right	0

Joystick Knobs		£
☐	QMP090505 Joystick Knob, Golf Ball, Ø 42 mm	40
☐	QMP090506 Joystick Knob, T-Bar, 82 mm Wide	40
☐	QMP090507 Joystick Knob, Soft Ball, Ø 70 mm	40
☐	QMP090508 Joystick Knob, Ball, Ø 35 mm	40
☐	QMP090045 Joystick Knob, Allround, Ball, Ø 48 mm	40
☐	QMP090509 Joystick Knob, Chin Cup, Ø 37 mm	40
☐	QMP090512 Joystick Knob, Stick, 80 x 15 mm	40
☐	QMP090511 Joystick Knob, Mushroom, Ø 45 mm	40

LED Lights & Indicators		£
☐	QMP110202 LED Fender Light (QMP110030 or QMP110032 required) (Not road approved) (Only with armrest QMP040123 or QMP040070) (Not with QMP130012)	150
☐	QMP110201 LED Light And Indicator StVZO Approved (QMP110030 or QMP110032 required)	400

Drive Wheels		£
☐	QMP070084 Drive Wheel, 12", Solid, Knobbly, Black	0
☐	QMP070034 Drive Wheel, 14", Pneumatic, Knobbly, Black	0
☐	QMP070033 Drive Wheel, 14", Solid, Knobbly, Black	0

Castor Wheels		£
☐	QMP080213 Castor Wheel, 7", Solid, Black	0

Motor Packages		£
☐	QMP180012 2-Pole Motor, 6 kph/ 4 mph, Performance	0
☐	QMP180015 2-Pole Motor, 10 kph / 6.5mph, Performance Package - Includes LED Lights & Indicators, Rear View Mirror & R-Net Gyro (Must select R-Net)	500
☐	QMP180009 4-Pole Motor, 6 kph / 4 mph	310

C-Me - Drive at 4.8kph in full Lift (Requirwe QMP130009/QMP130006)(Not with QMP030084/QMP030085)		
☐	QMP110034 C-Me Function Enabled (Max speed 4.5kph when lifted) (Not available with QMP130006 with QMP030087)	0
☐	QMP110035 C-Me Function Disabled (Max speed 1.5kph when lifted)	0

Shroud Colours			£
☐	QMP100576	Red	0
☐	QMP100578	Blue	0
☐	QMP100579	White	0
☐	QMP100577	Black	0

Colour Accents			£
☐	QMP140005	Apple Green (RAL 110 70 77)	0
☐	QMP140004	Black (RAL 9005)	0
☐	QMP140003	Blue (RAL 260 40 45)	0
☐	QMP140006	Orange (RAL 060 60 80)	0
☐	QMP140002	Red (RAL 040 40 67)	0
☐	QMP140007	Silver (RAL 9006)	0
☐	QMP140001	White (RAL 9010)	0

Personalised Shrouds			£
☐	QMP140008	Grey Burst	255
☐	QMP140009	Red Burst	255
☐	QMP140010	Blue Burst	255
☐	QMP140011	Pure White Burst	255
☐	QMP140012	Flame Aesteroid Grey	255
☐	QMP140013	Flame Hot Red	255
☐	QMP140014	Flame Sharp Blue	255
☐	QMP140015	Flame Pure White	255
☐	QMP140016	Floral White	255
☐	QMP140017	Carbon Black	255

Batteries			£
☐	QML160006	Battery, 41 Ah, Maintenance-Free	0
☒	QML160013	Battery, 56 Ah, Maintenance-Free	0

Chargers			£
☐	QMP160100	Charger, Automatic, 220 V / 8 A, UK	0
☐	QMP160103	Charger, Automatic, 220 V / 10 A, UK	230

Lap Belt			£
☐	QMP090128	Positioning Belt, Carseat Buckle, Short, 50 mm x 1700 mm	0
☐	QMP090129	Positioning Belt, Carseat Buckle, Long, 50 mm x 2700 mm	0

Accessories			£
☐	QMP090110	Rear View Mirror, Right (Not with same side tray QMP090040 OR QMP090046)	60
☐	QMP090118	Rear View Mirror, Left (Not with same side tray QMP090040 OR QMP090046)	60
☐	QMP090009	Crutch Holder (Not with QMP030086 or QMP030087 or QMP030229)	70
☐	QMP090119	Brackets For Bag On Backrest Frame (Not with QMP030086 or QMP030087 or QMP030229)	35
☐	QMP090120	Bag For Armrest (Not with QMP040123)	105
☐	QMP090122	Transfer Bar (Not with QMP020342/QMP020341 & QMP050015-17)	255

Tray Table (Mounted Opposite Joystick)			
☐	QMP090040	Tray Table, Width-Adjustable, Swing Away	335
☐	QMP090046	Tray Table, Half-Side, Swing Away (Not with QMP090126)	265

Tray Table Accessories			
☐	QMP090126	Raised Edges, Acrylic For Tray Table	25

5-Switch Box		£
☐	QMP110043 #1 Ctrl+5 Switch Box With Free Assignable Functions (If combined with attendant control, Ctrl+5 will be mounted to the armrest) Voice activation functionality included	575
#1 Ctrl+5 Switch Box Mounting Location		
☐	QMP110180 #1 Ctrl+5 Switch Box Mounting Location, Behind Joystick (Not possible with QMP110056)	0
☐	QMP110181 #1 Ctrl+5 Switch Box Mounting Location, Opposite Joystick	0
#1 Ctrl+5 Switch Type		
☐	QMP110175 #1 Ctrl+5 Switch Type With Push Buttons	0
☐	QMP110176 #1 Ctrl+5 Switch Type With Toggle Switches	0
Omni2 Input Control Module		£
☐	QMP110116 OMNI2 Input Control Module, R-Net, Bluetooth, Infrared	1170
☐	QMP110117 OMNI2 Input Control Module, R-Net, Scanning Mode, Bluetooth, Infrared	1170
Omni 2 Module Mount		£
☐	QMP110119 OMNI2 Module, Fixed Mount, Right	0
☐	QMP110120 OMNI2 Module, Fixed Mount, Left	0
☐	QMP110123 OMNI2 Module, Swing Away Mount, Right (Not with QMP090110)	130
☐	QMP110124 OMNI2 Module, Swing Away Mount, Left (Not with QMP090118)	130
Hand Control		£
☐	QMP110301 Hand Control, Switch-It MicroGuide	3145
☐	QMP110302 Hand Control, Switch-It MicroPilot	3495
☐	QMP110304 Hand Control, Switch-It VersaGuide	1120
☐	QMP110305 Hand Control, Switch-It VersaGuide EZ (Not with QMP090512/11/06)	1330
☐	QMP110306 Hand Control, Allround Joystick	1120
☐	QMP110307 Hand Control, Allround Joystick Light (Only use with default joystick knob)	1270
☐	QMP110308 Hand Control, Allround Joystick Heavy Duty	2095
☐	QMP110309 Hand Control, Allround Joystick, R-Net Bus Enabled (Omni2 module not required)	1190
☐	QMP110310 Hand Control, Allround Joystick Light, R-Net Bus Enabled (Omni2 module not required)	1400
☐	QMP110311 Hand Control, Allround Joystick Heavy Duty, R-Net Bus Enabled (Omni2 module not required)	2235
Hand Control Mounting		£
☐	QMP110342 Hand Control Mounting, Link-It To The Space Disk, Hand Tray (Only for Micro Series Joystick)	700
Hand Control Mounting		
☐	QMP110350 Hand Control Mounting, Link-It To The Armrest, Fixed	210
☐	QMP110353 Hand Control Mounting, Link-It To The Armrest, Swing Away	220
Hand Control Mounting Location		
☐	QMP110362 Hand Control Mounting Location, Right	0
☐	QMP110361 Hand Control Mounting Location, Left	0
Head Control (For switch selection see section E)(Chose a Whitmyer headrest mount)		£
Proportional Input Devices		
☐	QMP110370 Head Control, Switch-It Dual Pro, 3-Way Proportional, Facial Pads	5950
☐	QMP110371 Head Control, Switch-It Dual Pro, 3-Way Proportional, Spot Pads	5950

Chin Control			£
☐	QMP110470	Chin Control, Switch-It MicroGuide	3145
☐	QMP110471	Chin Control, Switch-It MicroPilot	3495
☐	QMP110472	Chin Control, Switch-It VersaGuide	1120
☐	QMP110473	Chin Control, Switch-It VersaGuide EZ	1330
☐	QMP110474	Chin Control, Allround Joystick	1120
☐	QMP110475	Chin Control, Allround Joystick Light	1330
☐	QMP110476	Chin Control, Allround Joystick, R-Net Bus Enabled (Onmi2 not required)	1190
☐	QMP110477	Chin Control, Allround Joystick Light, R-Net Bus Enabled (Onmi2 not required)	1400

Chin Control Mounting			£
☐	QMP110480	Chin Control Mounting, Link-It To Chin Harness, Small (If switches mounted on harness, you can only use Allround joystick)	630
☐	QMP110483	Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Right	560
☐	QMP110484	Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Left	560
☐	QMP110485	Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Right	2795
☐	QMP110486	Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Left	2795
☐	QMP110487	Chin Control Mounting, Link-It To Headrest, Swing Away, Left	490
☐	QMP110488	Chin Control Mounting, Link-It To Headrest, Swing Away, Right	490

Link-It Extension Kit			
☐	QMP110590	Link-It Extension Kit, Base Module, 100 mm Extension Tube, 5 Cable Clips	170
☐	QMP110591	Link-It Extension Kit, Base Module, 150 mm Extension Tube, 5 Cable Clips	210

Sip ,N Puff			£
☐	QMP110510	Sip 'N Puff With Link-It Hardware, Swing Away Mount To The Headrest	700

#1 Switches - Section E (Select 1 switch & mounting position, for a second switch please use #2 switch selection)			£
☐	QMP110611	#1 Single Switch, Ø 64 mm, Black	120
☐	QMP110612	#1 Single Switch, Ø 64 mm, Red	120
☐	QMP110613	#1 Single Switch, Ø 64 mm, Green	120
☐	QMP110523	#1 Egg Switch – Black	135
☐	QMP110524	#1 MicroLite Switch – Black	195
☐	QMP110651	#1 Single Switch, Ø 30 mm, Black	120
☐	QMP110654	#1 Single Switch, Ø 30 mm, Blue	120
☐	QMP110653	#1 Single Switch, Ø 30 mm, Green	120
☐	QMP110652	#1 Single Switch, Ø 30 mm, Red	120
☐	QMP110655	#1 Single Switch, Ø 30 mm, Yellow	120
☐	QMP110534	#1 Twister Switch – Black	195
☐	QMP110535	#1 Twister Switch – Green	195
☐	QMP110536	#1 Twister Switch – Red	195
☐	QMP110537	#1 Satellite Switch – Black (Only for Allround joysticks)	210
☐	QMP110538	#1 Satellite Switch – Green (Only for Allround joysticks)	210
☐	QMP110539	#1 Satellite Switch – Red (Only for Allround joysticks)	210

#1 Switch Mounting			£
☐	QMP110560	#1 Switch Mounting, Link-It To The Hand Control (Not with HD Allround Joystick or Space Disk Mount)	210
☐	QMP110566	#1 Switch Mounting, Link-It To The Armrest (Opposite Site As The Hand Control) (Armrest rail mount)	210
☐	QMP110561	#1 Switch Mounting, Link-It To The Headrest	280
☐	QMP110563	#1 Switch Mounting, Link-It To The Seat Rail (Operated by knee, if 2 switches selected mount one on the right and other on the left)	210
☐	QMP110565	#1 Switch Mounting, Link-It To Chin Boom (No left/right switch selection required)	210
☐	QMP110569	#1 Omit Switch Mounting (Switches attached lose)	0

#1 Switch Mount Location			£
☐	QMP110572	#1 Switch Mounting Location, Right	0
☐	QMP110571	#1 Switch Mounting Location, Left	0

#2 Swiches - Section F			£
☐	QMP110631	#2 Single Switch, Ø 64 mm, Black	120
☐	QMP110632	#2 Single Switch, Ø 64 mm, Red	120
☐	QMP110633	#2 Single Switch, Ø 64 mm, Green	120
☐	QMP110088	#2 Egg Switch – Black	135
☐	QMP110089	#2 MicroLite Switch – Black	185
☐	QMP110671	#2 Single Switch, Ø 30 mm, Black	120
☐	QMP110674	#2 Single Switch, Ø 30 mm, Blue	120
☐	QMP110673	#2 Single Switch, Ø 30 mm, Green	120
☐	QMP110672	#2 Single Switch, Ø 30 mm, Red	120
☐	QMP110675	#2 Single Switch, Ø 30 mm, Yellow	120
☐	QMP110095	#2 Twister Switch – Black	185
☐	QMP110096	#2 Twister Switch – Green	185
☐	QMP110097	#2 Twister Switch – Red	185
☐	QMP110111	#2 Satellite Switch – Black (Only for Allround joysticks)	200
☐	QMP110112	#2 Satellite Switch – Green (Only for Allround joysticks)	200
☐	QMP110113	#2 Satellite Switch – Red (Only for Allround joysticks)	200

#2 Switch Mounting			£
☐	QMP110098	#2 Switch Mounting, Link-It To The Hand Control (Not with HD Allround Joystick or Space Disk Mount)	210
☐	QMP110099	#2 Switch Mounting, Link-It To The Armrest (Opposite Site As The Hand Control) (Armrest rail mount)	210
☐	QMP110100	#2 Switch Mounting, Link-It To The Headrest	280
☐	QMP110104	#2 Switch Mounting, Link-It To The Seat Rail (Operated by knee, if 2 switches selected mount one on the right and other on the left)	210
☐	QMP110105	#2 Switch Mounting, Link-It To Chin Boom (No left/right switch selection required)	210
☐	QMP110106	#2 Omit Switch Mounting (Switches attached lose)	0

#2 Switch Mount Location			£
☐	QMP110108	#2 Switch Mounting Location, Right	0
☐	QMP110107	#2 Switch Mounting Location, Left	0

Multi-Function Switches (Use in combination with assignable button programming)			£
☐	QMP110550	Mini 2 Button (Stereo Plug) (Attached lose when ordered, no link-it fixation possible)	385
☐	QMP110551	Mini 2 Button (Dual Mono Plug) (Attached lose when ordered, no link-it fixation possible)	385
☐	QMP110552	Mini 4 Button InLine (DB9 Connector) (Attached lose when ordered, no link-it fixation possible)	645
☐	QMP110554	4-Way Toggle Switched Joystick	645

Accessories			£
☐	QMP110580	Mono Extension Cable (3.5 mm Jack Socket) (2 m)	140
☐	QMP110581	Stereo Extension Cable (3.5 mm Jack Socket) (2 m)	180
☐	QMP110582	Extension Cable DB9 (1 m)	170
☐	QMP110583	Y-Cable, 2 Female Mono To 1 Male Stereo	210

Safety Devices			£
☐	QMP090147	Remote Stop Switch (30 m Activation Range, Fit Via OBC Port) (Not with VR2)	980

Comments For Personalisations



Q300 M Mini with SEDEO Pro			£
Base Price			
Total Of Additional Equipment			
Packaging And Transport			
Exclusive VAT			
VAT			
Total Price			

Date			Signature
Day	Month	Year	

Send E-Mail