

Q300® R Kids



ORDER FORM

Errors and imissions excepted.


Billing Address Or Customer-ID		Shipping Address ☐ Identical To Billing Address	
C.-ID:		C.-ID:	
Name:		Name:	
Street:		Street:	
Town/PC:		Town/PC:	
E-Mail:		E-Mail:	
Phone:		Phone:	
Order Information		User Weight	
☐ Order ☐ Calculation ☐ Built-4-Me®		User Weight Capacity: 75 kg	
Purchase Order No.:		User Weight: kg	
Art. No.: 2200030R-000 Reimbursement Code:			
Order Date: Day Month Year			

Base Price

£6500

Client Measurements

a) Top Of Head To Seat		mm
b) Top Of Shoulder To Seat		mm
c) Underarm To Seat		mm
d) Elbow To Seat		mm
e) Top Of Pelvis To Seat		mm
f) Shoulder Width		mm
g) Trunk Width		mm
i) Hip Width		mm
j) Actual Sitting Depth		mm
m) Foot Length		mm
n) Lower Leg Length		mm
o) Actual Seat Width		mm



Seating Specifications

All Bases Prepared For DAHL Docking Station		Seat Type	Kids Seat
Typical WC Weight	Without Batteries 80 kg	Seat Size Range	Width 280 mm - 400 mm, Depth 260 mm - 400 mm
	With Batteries from 103 kg	Seat To Floor Height Range	470-490 mm

RWD Base Specifications

User Weight Capacity	75 kg	Drive Wheels	12" or 14" Drive Wheels (Black)
Base Size	Width Min. 550mm (12" Drive Wheels)	Tyres / Castors	8" Solid / 9" Solid/Pneumatic (Black)
Drive Control	VR2, Option R-Net	Charger	8 A Dual Mode Charger
Driving Range	Up To 32 km	Battery Size	50, 56 & 64 Ah Batteries
Drive Motors	2-Pole 6 /10 kph, 4-Pole 6/10 kph	Colour Accents	5 Shroud Options, 7 Colour Accents
Ground Clearance	90 mm	Frame Colour	Black

Base

Q300 R Kids Drive Base

£

☐	QZK000165	Q300 R (RWD Base Chair)	6500
--------------------------	-----------	--------------------------------	------

Motor Packages

£

☐	QZK180014	2-Pole Motor, 4 kph / 2.5 mph, Via Programming	0
☐	QZK180011	2-Pole Motor, 6 kph / 4 mph, Standard	0
☐	QZK180012	2-Pole Motor, 6 kph / 4 mph, Performance	0
☐	QZK180013	2-Pole Motor, 10 kph / 6.5 mph, Performance	0
☐	QZK180009	4-Pole Motor, 6 kph / 4 mph (8 kph Motor With Reduced Speed)	310
☐	QZK180007	4-Pole Motor, 10 kph / 6.5 mph	390

Tracking

£

☐	QZK110009	Traction Control System For R-Net (Gyro Technology)	400
☐	QZK110014	Omit Traction Control System For R-Net	0

Drive Wheels

£

☐	QZK070084	Drive Wheel, 12", Solid, Knobbly, Black	0
☐	QZK070034	Drive Wheel, 14", Pneumatic, Knobbly, Black	0
☐	QZK070033	Drive Wheel, 14", Solid, Knobbly, Black	0

Castor Wheels

£

☐	QZK080059	Castor Wheel, 8", Solid, Black	0
☐	QZK080210	Castor Wheel, 9", Solid, Black Not Available With 80 ° Life Style Hanger (QZK050001)	0
☐	QZK080211	Castor Wheel, 9", Pneumatic, Black Not Available With 80 ° Life Style Hanger (QZK050001)	0

Shroud Colours

£

☐	QZK100577	Black (RAL 9005) Black Natural Non-Painted finish	0
☐	QZK100576	Red (RAL 040 40 67)	0
☐	QZK100578	Blue (RAL 260 40 45)	0
☐	QZK100585	Yellow (RAL 1016)	0
☐	QZK100586	Purple (RAL 4005)	0

Accent Colours

£

☐	QZK140005	Apple Green (RAL 110 70 77)	0
☐	QZK140004	Black (RAL 9005)	0
☐	QZK140006	Orange (RAL 060 60 80)	0
☐	QZK140007	Silver (RAL 9006)	0
☐	QZK140002	Red (RAL 040 40 67)	0
☐	QZK140003	Blue (RAL 260 40 45)	0
☐	QZK140001	White (RAL 9010)	0

Kerb Climber

£

☐	QZK090038	Boot Style Kerb Climber Not with Centre Mount Legrest (QZK050112)	200
☐	QZK090037	Omit Kerb Climber	0

LED Lights & Indicators

£

☐	QZK110202	Fender Lights VR2/R-NET Lighting Module Required, Linx Compatible Power Module Required, Fender Light - Not Road Approved and Without Indicators	160
☐	QZK110201	LED Light And Indicator StVZO Approved VR2/R-NET Light Module Required, Linx Compatible Power Module Required	400
☐	QZK110209	Omit LED Light And Indicator	0

Seat Heights (fixed at 3°, tilt and lift & tilt at 0°, negative tilt at minus 5°)

	QZK130002	QZK130004	QZK130012	QZK130009
	Fixed Tilt Module	Posterior Tilt 0 - 30°	Posterior Tilt -5 - 25°	Seat Elevator / Tilt Module 0 - 30°
☐ SH47	*	*	*	
☐ SH49	*	*	*	*

Batteries

£

Batteries

☐ QZK160016	Battery, 50 Ah, Maintenance-Free	0
☐ QZK160013	Battery, 55Ah/56 Ah, Maintenance-Free	130
☐ QZK160030	Battery, 63Ah/64 Ah, Maintenance-Free Not with ISM in UK Market	130

Chargers

☐ QZK160100	Charger, Automatic, 220 V / 8 A, UK	0
☐ QZK160103	Charger, Automatic, 220 V / 10 A, UK	230

Powered Seating**Powered Seating Options**

£

☐ QZK130002	Manual tilt (3° preset tilt)	0
☐ QZK130004	Powered C.O.G tilt: 0° to 30° Not With QZK110036 Remote Joystick, VR2 or QZK110810 Linx Joystick REM110	1150
☐ QZK130012	Powered negative C.O.G tilt: -5° to 25° Not With QZK110036 Remote Joystick, VR2 or QZK110810 Linx Joystick REM110	1290
☐ QZK130009	Powered lift & tilt: 30cm lift, 0° to 30° tilt Not With QZK110036 Remote Joystick, VR2 or QZK110810 Linx Joystick REM110 In Case Of R-Net: QZK110032 Actuator Control & Lights Module, R-Net, ISM-6L Required In Case of Linx: QZK110820 2 Actuator Control Module ACT-200 Required	2050

C-Me

£

☐ QZK110034	C-Me Function, Enabled Max. Speed 3 mph/4,5km/h When Lifted, Only with Powered Lift and Tilt QZK130009	0
☐ QZK110035	C-Me Function, Disabled	0

Backrest Frame Type

£

☐ QZK030018	Mechanical fixed recline back (-4 to 16° in 3 steps)	0
☐ QZK030026	Powered recline back (0 to 30°)	560

Back Cushion Type

£

Sling Back (See Section A)

☐ QZK030150	Sling Back, With Tension Adjustable Straps	0
------------------------------------	--	---

SEDEO LITE (See Section B)

☐ QZK030151	SEDEO LITE, Back Cushion	0
------------------------------------	--------------------------	---

JAYJ3 (See Section C)

☐ QZK030054	JAY J3, Back Cushion	0
☐ QZK030000	Omit Backrest Upholstery	0

JAY Shape (See Section D)

☐ QZK030152	JAY Shape, Back Cushion	0
------------------------------------	-------------------------	---

A – Sling Back (Back Width Needs To Correspond To Seat Width) £

☐	QZK030301	Tension adjustable sling backrest (Factory Setting 410 mm height)	0
--------------------------	-----------	---	---

A – Backrest Height £

☐	BH38	380 mm	0
☐	BH41	410 mm	0
☐	BH435	435 mm	0
☐	BH46	460 mm	0
☐	BH485	485 mm	0

B – SEDEO LITE Backrest £

☐	QZK030214	SEDEO Lite Deep Contour Active Required To Achieve Seat Depth SD 26 (260 mm) / SD 28 (280 mm) / Not Available With SD 38 (380 mm) / SD 40 (400 mm)	350
--------------------------	-----------	---	-----

B – SEDEO LITE Back Cushion Cover £

☐	QZK030059	Back Cushion Cover, Incontinence Protection	0
☐	QZK030108	Back Cushion Cover, Comfort	0

C – JAY J3 Back Size, Mid Deep Contour £

☐	JAY030040	JAY J3 Back, Mid Deep Contour, Mid Thoracic, Short 300 mm	870
☐	JAY030044	JAY J3 Back, Mid Deep Contour, Mid Thoracic, Medium 340 mm	870
☐	JAY030041	JAY J3 Back, Mid Deep Contour, Mid Thoracic, Tall, 380 mm	870

D – JAY Shape Back £

☐	JAY080001	JAY Shape Back Not Available With SW 28, 30 ,32 or 34	270
☐	JAY080002	JAY Shape Back, Oversized Not Available With SW 28, 30 or 32. Supplied 25 mm Wider Than Seat Width	270

D – JAY Shape Backrest Height £

☐	BH41	410 mm	0
☐	BH435	435 mm	0
☐	BH46	460 mm	0
☐	BH485	485 mm	0

D – JAY Shape Lateral Thoracic Pads £

☐	JAY080011	JAY SHAPE Lateral Thoracic supports (Pair) - Support 100 x 120 mm (Pair)	270
--------------------------	-----------	--	-----

Push Handles £

☐	QZK030212	Standard Push Handle	0
☐	QZK030213	Straight Push Handle Not available with QZK110021 / QZK110840 Dual Control Or QZK110020 / QZK110842 Joystick Mounted As Dual Control	0

Headrest Pads £

☐	QZK030430	SEDEO, Headrest Pad, Small, (LxH) 240 x 150 mm Only With QZK030433 SEDEO HEADREST MOUNT Or QZK030434 SEDEO, Quick Release Headrest Mount	95
☐	QZK030431	SEDEO, Headrest Pad, Small, With Cover, (LxH) 240 x 150 mm Only With QZK030433 SEDEO HEADREST MOUNT Or QZK030434 SEDEO, Quick Release Headrest Mount	110
☐	WHT030401	Whitmyer, Headrest Pad, Plush, (LxH) 203 x 127 mm Only With WHT030408 Whitmyer LINX2 or WHT030409 AXYS Headrest Mounts	170
☐	WHT030402	Whitmyer, Headrest Pad, Plush, (LxH) 254 x 127 mm Only With WHT030408 Whitmyer LINX2 or WHT030409 AXYS Headrest Mounts	170
☐	WHT030404	Whitmyer, Headrest Pad, Adjust-A-Plush, Narrow, (LxH) 355 x 76 mm Only With WHT030408 Whitmyer LINX2 or WHT030409 AXYS Headrest Mounts	280

Headrest Mounts £

☐	QZK030433	SEDEO, Headrest Mount	140
☐	QZK030434	SEDEO, Quick Release Headrest Mount	155
☐	WHT030408	Whitmyer, Headrest Mount, LINX ²	190
☐	WHT030409	Whitmyer, Headrest Mount, AXYS	155

Headrest Receiver £

☐	QZK030120	Headrest Receiver, Backrest Frame Mount	0
☐	QZK030119	Headrest Receiver, Back Cushion Mount Only With QZK030214 SEDEO LITE Deep Contour Active	0

Seat & Cushion £

☐	QZK020120	2" Standard Foam Seat Cushion	0
☐	QZK020110	Omit Cushion	0

Standard Seat Size £

Seat Width			
☐	SW28	280 mm Only with QZK040064 Side Guard Short In Rig	0
☐	SW30	300 mm	0
☐	SW32	320 mm	0
☐	SW34	340 mm	0
☐	SW36	360 mm	0
☐	SW38	380 mm	0
☐	SW40	400 mm Not With QZK040064, Side Guard Short In Rig	0
Seat Depth			
☐	SD26	260 mm Only With QZK030214 SEDEO LITE Deep Contour Active	0
☐	SD28	280 mm Only With QZK030214 SEDEO LITE Deep Contour Active	0
☐	SD30	300 mm	0
☐	SD32	320 mm	0
☐	SD34	340 mm	0
☐	SD36	360 mm	0
☐	SD38	380 mm Not With QZK030214 SEDEO LITE Deep Contour Active	0
☐	SD40	400 mm Not With QZK030214 SEDEO LITE Deep Contour Active	0

Armrest Types £

☐	QZK040122	Flip-up, height and angle adjustable	0
--------------------------	-----------	--------------------------------------	---

Arm Pads £

☐	QZK040023	Arm Pad, Integral Foam, 300 mm, Right	0
☐	QZK040024	Arm Pad, Integral Foam, 300 mm, Left	0
☐	QZK040062	Arm Pad, Integral Foam, 280x65 mm, Right Not with combination of QZK110056 Joystick Mounting, Parallel, Swing Way, Capsulated and QZK110052 Joystick Mounting Location, Right	0
☐	QZK040063	Arm Pad, Integral Foam, 280x65 mm, Left Not with combination of QZK110056 Joystick Mounting, Parallel, Swing Way, Capsulated and QZK110053 Joystick Mounting Location, Left	0

Side Guards £

☐	QZK040064	Side Guard Short In Rig Reduces SW By 20 mm	0
☐	QZK040033	Side Guard Not With SW 28 (280 mm)	0

Centre Mount Legrest £

☐	QZK050112	Centre Mount Legrest (includes calf pads and flip up footplate)	0
--------------------------	-----------	---	---

Swing Away Legrest & Footplates £

☐	QZK050002	Hanger Angle 70° (Quickie Life style)	0
☐	QZK050001	Hanger Angle 80° (Quickie Life style) (3) Not With QZK130012 Powered Negative C.O.G Tilt -5 to 25° or 9" Castor QZT080210/11	0
☐	QZK050115	Legrest Uni 80° Swing Away	0
☐	QZK050041	Hanger Angle 90° (Zippie Salsa style with aluminium footplates and calf strap)	0
☐	QZK050016	Legrest Manual Elevating Swing Away (Includes Flip Up Footplates)	375

Footplates for Swing Away Legrests

☐	QZK050135	Footplate, Fixed, Pair Only With QZK050115 Legrest Uni 80° Swing Away	0
☐	QZK050137	Footplate, Angle Adjustable , Pair (for Uni legrest) Only With QZK050115 Legrest Uni 80° Swing Away	80
☐	QZK050050	Footplate, Divided , Aluminium, Angle-Adjustable, Flip-Up (Sideways), Heel Loop, Black Not With QZK130012 Powered Negative C.O.G Tilt -5 to 25° or 9" Castor QZT080210/11 Only With QZK050002/1 Hanger Angle 70°/80°	0
☐	QZK050167	Footplates, High Mount - Recommended for short LLL Not With QZK130012 Powered Negative C.O.G Tilt -5 to 25° or 9" Castor QZT080210/11 Only With QZK050002/1 Hanger Angle 70°/80°	130
☐	QZK050168	Short Clamp (90 mm)	0
☐	QZK050169	Long Clamp (121 mm)	0

Accessories £

☐	QZK090100	Transportation tie down points (Included as Standard)	0
☐	QZK090171	Adjustable Lap strap	0
☐	QZK090040	Swing-away tray table	335
☐	QZK090110	Rear View Mirror, Right Not In Combination With Same Entry Side QZK090040 Tray Table	60
☐	QZK090118	Rear View Mirror, Left Not In Combination With Same Entry Side QZK090040 Tray Table	60

Remote Joystick – VR2 £

☐	QZK110036	Remote Joystick, VR2	0
☐	QZK110037	Remote Joystick, VR2 , 2 Actuators	0
☐	QZK110040	Remote Joystick, VR2, 2 Actuators and Lights	0

Joystick Mounting £

☐	QZK110082	Joystick Mounting, Parallel, Swing Away	0
☐	QZK110056	Joystick Mounting, Parallel, Swing Away, Capsulated Not With QZK040062/3 Arm Pad, Integral Foam, 280x65 mm Right/Left	80

Joystick Mounting Location £

☐	QZK110052	Joystick Mounting Location, Right	0
☐	QZK110053	Joystick Mounting Location, Left	0

Joystick Knobs (optional) £

☐	QZK090505	Joystick Knob, Golf Ball, Ø 42 mm	40
☐	QZK090506	Joystick Knob, T-Bar, 82 mm Wide	40
☐	QZK090507	Joystick Knob, Soft Ball, Ø 70 mm	40
☐	QZK090508	Joystick Knob, Ball, Ø 35 mm	40
☐	QZK090045	Joystick Knob, Allround, Ball Ø 48 mm	40
☐	QZK090509	Joystick Knob, Chin Cup, Ø 37 mm	40
☐	QZK090512	Joystick Knob, Stick, 80 x 15 mm	40
☐	QZK090511	Joystick Knob, Mushroom, Ø 45 mm	40

Power Module – VR2		£
☐ QZK110001	Power module , VR2, 70A Not Recommended For Performance/ 4 pole motors	0
☐ QZK110002	Power module , VR2, 90A Not With 4/6kph Standard Motors	135

Lighting Module - VR2		£
☐ QZK110030	Lighting Module, VR2 If LED Light And Indicator Is Chosen, Then This Is Required	0

Attendant Controls		£
☐ QZK110021	Dual Control	425
☐ QZK110020	Joystick Mounted As Attendant Control Choose Also A Remote Joystick / Right Hand Attendant Control Only	270

Remote Joystick – R-Net		£
☐ QZK110015	Remote Joystick, R-Net, LED, Actuators	485
☐ QZK110016	Remote Joystick, R-Net, LED, Actuators And Lights	485
☐ QZK110017	Remote Joystick, R-Net, Colour LCD	630
☐ QZK110018	Remote Joystick, R-Net, Colour LCD, Heavy Duty Not with any joystick knob QZK090505-QZK090511	1200
☐ QZK110013	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared Control Included	750
☐ QZK110023	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared And Bluetooth Control Included	925
☐ QZK110019	Omit Remote Joystick, R-Net	0

Joystick Mounting		£
☐ QZK110082	Joystick Mounting, Parallel, Swing Away	0
☐ QZK110056	Joystick Mounting, Parallel, Swing Away, Capsulated Not With QZK040062/3 Arm Pad, Integral Foam, 280x65 mm Right/Left	80

Joystick Mounting Location		£
☐ QZK110052	Joystick Mounting Location, Right	0
☐ QZK110053	Joystick Mounting Location, Left	0

Power Modules – R-Net		£
☐ QZK110003	Power Module, R-Net, 90 A	0

Seating / Lighting Module Control – R-Net		£
☐ QZK110032	Actuator Control & Lights Module, R-Net, ISM-6L To Be Offered Only When R-Net With Actuators AND Lights is selected	0
☐ QZK110204	Actuator Control Module, R-Net, ISM-6 To Be Offered Only When R-Net With Actuators NO Lights is selected	0
☐ QZK110157	Light Module, R-Net, ISM-L	0

Attendant Controls		£
☐ QZK110021	Dual Control	425
☐ QZK110020	Joystick Mounted As Attendant Control Choose Also A Remote Joystick / Right Hand Attendant Control Only	270

Special Controls

OMNI2 Input Control Module – R-Net

£

☐	QZK110740	OMNI2 incl. 2x Single Switch (Ø 64 mm) and Flexi Mount	1170
--------------------------	-----------	--	------

OMNI2 Module Mount

£

☐	QZK110122	OMNI2 Flexi Mount, Right	0
☐	QZK110121	OMNI2 Flexi Mount, Left	0

Safety Devices

£

☐	QZK090147	Remote Stop Switch (30 m Activation Range, Fit Via OBC Port) Not With VR2	500
--------------------------	-----------	--	-----

LiNX Controls

£

☐	QZK110810	LiNX Remote Drive Only, REM110	0
☐	QZK110811	LiNX Remote Drive 2 Act, REM210	105
☐	QZK110812	LiNX Remote Drive 6 Act, REM211	105
☐	QZK110813	LiNX Remote Drive 2 Actuators + Lights, REM215	475
☐	QZK110814	LiNX Remote Drive 6 Actuators + Lights, REM216	475
☐	QZK110815	LiNX Remote Drive 6 Act, Lights, Touchscreen, REM400	630
☐	QZK110818	Omit LiNX Remote	0

Linx Power Module & ISM

☐	QZK110802	Power Module 75A LiNX Drive Only	0
☐	QZK110803	Power Module 75A LiNX, Drive, Lights & Actuators	105
☐	QZK 110805	Power Module 120A LiNX, Drive, Lights & Actuators	135
☐	QZK110820	2 Actuator Module LiNX, ACT-200	0
☐	QZK110822	Omit Actuator Module	0

Linx Gyro

☐	QZK110830	Traction Control System, for LiNX Gyro Technology	400
☐	QZK110831	Omit Traction Control	0

LiNX Joystick Mount

£

☐	QZK110836	Joystick Mounting, Parallel, Swing Away	0
☐	QZK110837	Joystick Mounting, Parallel, Swing Away, Capsulated Not With QZK040062/3 Arm Pad, Integral Foam, 280x65 mm Right/Left	80

LiNX Joystick Mounting Location

£

☐	QZK110838	Joystick Mounting Location, Right	0
☐	QZK110839	Joystick Mounting Location, Left	0

LiNX Attendant Controls

£

☐	QZK110840	Dual Control	475
☐	QZK110842	Joystick Mounted As Attendant Control Choose Also A Remote Joystick / Right Hand Attendant Control Only	80
☐	QZK110843	Omit Attendant Control	0

LiNX Speciality Display

£

☐	QZK110823	LiNX Speciality Display, Omni, DX-REM500	1225
☐	QZK110824	Omit Speciality Display	0

Mounting Location

☐	QZK110850	Speciality Display Flexi, Right	0
☐	QZK110851	Speciality Display Flexi, Left	0

Drive Wheels & Castor Wheels & Mudguards

[◀ BACK](#)


// QZK070033



// QZK080210

LED Lights & Chargers

[◀ BACK](#)


// QZK110201



// QZK160100 & QZK160101



// QZK160102



// QZK160103 & QZK160104

Back Cushion Type

[◀ BACK](#)


// QZK030214

Headrest Pads & Headrest Mounts

[◀ BACK](#)


// QZK030430 & QZK030431



// WHT030401-WH030402



// WHT030404



// QZK030433



// WHT030408



// WHT030409

Seat Cushion Type

[◀ BACK](#)


// QZK020120

JAY J3 Back

[◀ BACK](#)


// JAY030041

Armrest Type

[◀ BACK](#)


// QZK040122

Swing Away Legrest & Support Options

[◀ BACK](#)


// QZK050001

// QZK050050

// QZK050167

Remote Joystick – R-Net & Joystick Mounting

◀ BACK



// QZK110017



// QZK110013



// QZK110018



// QZK110056

Remote Joystick – VR2

◀ BACK



// QZK110036



// QZK110037

Joystick Knobs

◀ BACK



// QZK090505



// QZK090045



// QZK090507



// QZK090508



// QZK090509



// QZK090512



// QZK090511



// QZK090506

Electronic Programmers – R-Net[◀ BACK](#)

// QZK090066

Power Modules – R-Net / Attendant Controls[◀ BACK](#)

// QZK110021

R-Net OMNI2 Input Control Modules[◀ BACK](#)

// QZK110740

Comments For Adaptation



Q300 M Mini Kids

£

Base Price	
Total Of Additional Equipment	
Packaging And Transport	
Exclusive VAT	
VAT	
Total Price	

Date

Signature

Day

Month

Year

Send E-Mail

**Sunrise Medical GmbH**

Kahlbachring 2-4

D-69254 Malsch/HD

Phone: +49 (0) 72 53 / 9 80-0

Fax: +49 (0) 72 53 / 9 80-222

E-Mail: kundenservice@sunrisemedical.de

www.SunriseMedical.eu