

Q300® R Teens



ORDER FORM

Errors and imissions excepted.


Billing Address Or Customer-ID		Shipping Address ☐ Identical To Billing Address	
C.-ID:		C.-ID:	
Name:		Name:	
Street:		Street:	
Town/PC:		Town/PC:	
E-Mail:		E-Mail:	
Phone:		Phone:	
Order Information		User Weight	
☐ Order ☐ Calculation ☐ Built-4-Me®		User Weight Capacity: 120 kg	
Purchase Order No.:		User Weight: kg	
Art. No.: 2200030R-050 Reimbursement Code:			
Order Date: Day Month Year			

Base Price

£6500

Client Measurements

a) Top Of Head To Seat		mm
b) Top Of Shoulder To Seat		mm
c) Underarm To Seat		mm
d) Elbow To Seat		mm
e) Top Of Pelvis To Seat		mm
f) Shoulder Width		mm
g) Trunk Width		mm
i) Hip Width		mm
j) Actual Sitting Depth		mm
m) Foot Length		mm
n) Lower Leg Length		mm
o) Actual Seat Width		mm



Seating Specifications

All Bases Prepared For DAHL Docking Station		Seat Type	Teens Seat
Typical WC Weight	Without Batteries 80 kg	Seat Size Range	Width 300 mm - 400 mm, Depth 360 mm - 500 mm
	With Batteries from 103 kg	Seat To Floor Height Range	470 mm - 490 mm

RWD Base Specifications

User Weight Capacity	120 kg	Drive Wheels	12" or 14" Drive Wheels (Black)
Base Size	Width Min. 550mm (12" Drive Wheels)	Tyres / Castors	8" Solid / 9" Solid/Pneumatic (Black)
Drive Control	VR2, Option R-Net	Charger	8 A Dual Mode Charger
Driving Range	Up To 32 km	Battery Size	50, 56 & 64 Ah Batteries
Drive Motors	2-Pole 6 /10 kph, 4-Pole 6/10 kph	Colour Accents	5 Shroud Options, 7 Colour Accents
Ground Clearance	90 mm	Frame Colour	Black

Base

Q300 M Mini Teens Drive Base

£

☐	QZT000165	Q300 R (RWD Base Chair)	6500
--------------------------	-----------	--------------------------------	------

Motor Packages

£

☐	QZT180014	2-Pole Motor, 4 kph / 2.5 mph, Via Programming	0
☐	QZT180011	2-Pole Motor, 6 kph / 4 mph, Standard	0
☐	QZT180012	2-Pole Motor, 6 kph / 4 mph, Performance	0
☐	QZT180013	2-Pole Motor, 10 kph / 6.5 mph, Performance	0
☐	QZT180009	4-Pole Motor, 6 kph / 4 mph (8 kph Motor With Reduced Speed)	310
☐	QZT180007	4-Pole Motor, 10 kph / 6.5 mph	390

Tracking

£

☐	QZT110009	Traction Control System For R-Net (GYRO Technology)	400
☐	QZT110014	Omit Traction Control System For R-Net	0

Drive Wheels

£

☐	QZT070084	Drive Wheel, 12", Solid, Knobbly, Black	0
☐	QZT070034	Drive Wheel, 14", Pneumatic, Knobbly, Black	0
☐	QZT070033	Drive Wheel, 14", Solid, Knobbly, Black	0

Castor Wheels

£

☐	QZT080059	Castor Wheel, 8", Solid, Black	0
☐	QZT080210	Castor Wheel, 9", Solid, Black Not Available With Centre Mount Legrest (QZT050106, QZT050107) or 80° Life Style Hanger (QZT050001)	0
☐	QZT080211	Castor Wheel, 9", Pneumatic, Black Not Available With Centre Mount Legrest (QZT050106, QZT050107) or 80° Life Style Hanger (QZT050001)	0

Shroud Colours

£

☐	QZT100577	Black (RAL 9005) Black Natural Non-Painted finish	0
☐	QZT100576	Red (RAL 040 40 67)	0
☐	QZT100578	Blue (RAL 260 40 45)	0
☐	QZT100585	Yellow (RAL 1016)	0
☐	QZT100586	Purple (RAL 4005)	0

Accent Colours

£

☐	QZT140005	Apple Green (RAL 110 70 77)	0
☐	QZT140004	Black (RAL 9005)	0
☐	QZT140006	Orange (RAL 060 60 80)	0
☐	QZT140007	Silver (RAL 9006)	0
☐	QZT140002	Red (RAL 040 40 67)	0
☐	QZT140003	Blue (RAL 260 40 45)	0
☐	QZT140001	White (RAL 9010)	0

Kerb Climber

£

☐	QZT090038	Boot Style Kerb Cimber Not with Centre Mount Legrest (QZT050106 / QZT050107)	200
☐	QZT090037	Omit Kerb Climber	0

LED Lights & Indicators

£

☐	QZT110202	Fender Lights VR2/R-NET Lighting Module Required, Linx Compatible Power Module Required, Fender Light - Not Road Approved and Without Indicators	160
☐	QZT110201	LED Light And Indicator StVZO Approved VR2/R-NET Light Module Required, Linx Compatible Power Module Required	400
☐	QZT110209	Omit LED Light And Indicator	0

Seat Heights (tilt and lift & tilt at 0°, negative tilt at minus 5°)

	QZT130004	QZT130012	QZT130009
	Posterior Tilt 0 - 30°	Posterior Tilt -5 - 25°	Seat Elevator / Tilt Module 0 - 30°
☐ SH47	*	*	
☐ SH49	*	*	*

Batteries

£

Batteries

☐ QZT160016	Battery, 50 Ah, Maintenance-Free	0
☐ QZT160013	Battery, 55Ah/56 Ah, Maintenance-Free	130
☐ QZT160030	Battery, 63Ah/64 Ah, Maintenance-Free Not with ISM in UK Market	130
☐ QZT160009	Omit Batteries	0

Chargers

☐ QZT160100	Charger, Automatic, 220 V / 8 A, UK	0
☐ QZT160103	Charger, Automatic, 220 V / 10 A, UK	230
☐ QZT160109	Omit Charger	0

Powered Seating**Powered Seating Options**

£

☐ QZT130004	Powered C.O.G tilt: 0° to 30° Not With QZT110036 Remote Joystick, VR2 or QZT110810 Linx Joystick REM110	1150
☐ QZT130012	Powered negative C.O.G tilt: -5° to 25° Not With QZT110036 Remote Joystick, VR2 or QZT110810 Linx Joystick REM110	1290
☐ QZT130009	Powered lift & tilt: 30cm lift, 0° to 30° tilt Not With QZT110036 Remote Joystick, VR2 or QZT110810 Linx Joystick REM110 In Case Of R-Net: QZT110032 Actuator Control & Lights Module, R-Net, ISM-6L Required In Case of Linx: QZT110821 4 Actuator Control Module ACT-400 Required	2050

C-Me (7)

£

☐ QZT110034	C-Me Function, Enabled Max. Speed 3 mph/4,5km/h When Lifted. Only with Powered Lift and Tilt QZT130009	0
☐ QZT110035	C-Me Function, Disabled	0

Backrest Frame Type

£

☐ QZT030018	Mechanical fixed recline back (-4 to 16° in 3 steps)	0
------------------------------------	--	---

Back Cushion Type

£

Sling Back (See Section A)

☐ QZT030150	Sling Back, With Tension Adjustable Straps	0
------------------------------------	--	---

SEDEO LITE (See Section B)

☐ QZT030151	SEDEO LITE, Back Cushion	0
------------------------------------	--------------------------	---

JAYJ3 (See Section C)

☐ QZT030054	JAY J3, Back Cushion	0
☐ QZT030000	Omit Backrest Upholstery	0

JAY Shape (See Section D)

☐ QZT030152	JAY Shape, Back Cushion	0
------------------------------------	-------------------------	---

A – Sling Back (Back Width Needs To Correspond To Seat Width) £

☐	QZT030301	Tension adjustable sling backrest (Factory Setting 410 mm height)	0
--------------------------	-----------	---	---

A – Backrest Height £

☐	BH435	435 mm	0
☐	BH46	460 mm	0
☐	BH485	485 mm	0
☐	BH51	510 mm	0
☐	BH535	535 mm	0

B – SEDEO LITE Backrest £

☐	QZT030214	SEDEO Lite Deep Contour Active Required To Achieve Seat Depth SD 26 (260 mm) / SD 28 (280 mm) / Not Available With SD 38 (380 mm) / SD 40 (400 mm)	350
--------------------------	-----------	---	-----

B – SEDEO LITE Back Cushion Cover £

☐	QZT030059	Back Cushion Cover, Incontinence Protection	0
☐	QZT030108	Back Cushion Cover, Comfort	0

C – JAY J3 Back Size, Mid Deep Contour £

☐	JAY030040	JAY J3 Back, Mid Deep Contour, Mid Thoracic, Short 300 mm	870
☐	JAY030044	JAY J3 Back, Mid Deep Contour, Mid Thoracic, Medium 340 mm	870
☐	JAY030041	JAY J3 Back, Mid Deep Contour, Mid Thoracic, Tall, 380 mm	870
☐	JAY030045	JAY J3 Back, Mid Deep Contour, Upper Thoracic, Short, 420 mm	870
☐	JAY030042	JAY J3 Back, Mid Deep Contour, Upper Thoracic, Medium, 460 mm	870
☐	JAY030046	JAY J3 Back, Mid Deep Contour, Upper Thoracic, Tall, 500 mm	870

D – JAY Shape Back £

☐	JAY080001	JAY Shape Back Not Available With SW 30 ,32 or 34	270
☐	JAY080002	JAY Shape Back, Oversized Not Available With SW 30 or 32. Supplied 25 mm Wider Than Seat Width	270

D – JAY Shape Backrest Height £

☐	BH435	435 mm	0
☐	BH46	460 mm	0
☐	BH485	485 mm	0

D – JAY Shape Lateral Thoracic Pads £

☐	JAY080011	JAY SHAPE Lateral Thoracic supports (Pair) - Support 100 mm x 120 mm (Pair)	270
--------------------------	-----------	---	-----

Push Handles £

☐	QZT030212	Standard Push Handle	0
☐	QZT030213	Straight Push Handle Not Available With QZT030120 Headrest Receiver, Backrest Frame Mount , QZT110021 Dual Control Or QZT110020 Joystick Mounted As Dual Control	0

Headrest Pads £

☐	QZT030430	SEDEO, Headrest Pad, Small, (LxH) 240 x 150 mm Only With QZK030433 SEDEO HEADREST MOUNT Or QZK030434 SEDEO, Quick Release Headrest Mount	95
☐	QZT030431	SEDEO, Headrest Pad, Small, With Cover, (LxH) 240 x 150 mm Only With QZK030433 SEDEO HEADREST MOUNT Or QZK030434 SEDEO, Quick Release Headrest Mount	110
☐	QZT030432	SEDEO, Headrest Pad, Flexible, (LxH) 400 x 110 mm Only With WHT030408 Whitmyer LINX2 or WHT030409 AXYS Headrest Mounts	170
☐	WHT030402	Whitmyer, Headrest Pad, Plush, (LxH) 254 x 127 mm Only With WHT030408 Whitmyer LINX2 or WHT030409 AXYS Headrest Mounts	170
☐	WHT030404	Whitmyer, Headrest Pad, Adjust-A-Plush, Narrow, (LxH) 355 x 76 mm Only With WHT030408 Whitmyer LINX2 or WHT030409 AXYS Headrest Mounts	280

Headrest Mounts £

☐	QZT030433	SEDEO, Headrest Mount	140
☐	QZT030434	SEDEO, Quick Release Headrest Mount	155
☐	WHT030408	Whitmyer, Headrest Mount, LINX ²	190
☐	WHT030409	Whitmyer, Headrest Mount, AXYS	155

Headrest Receiver £

☐	QZT030120	Headrest Receiver, Backrest Frame Mount	0
☐	QZT030119	Headrest Receiver, Back Cushion Mount Only With QZT030214 SEDEO LITE Deep Contour Active	0

Seat & Cushion £

☐	QZT020120	2" Standard Foam Seat Cushion	0
☐	QZT020110	Omit Cushion	0

Standard Seat Size £

Seat Width			
☐	SW30	300 mm	0
☐	SW32	320 mm	0
☐	SW34	340 mm	0
☐	SW36	360 mm	0
☐	SW38	380 mm	0
☐	SW40	400 mm	0
Seat Depth			
☐	SD36	360 mm Only With QZT030214 SEDEO LITE Deep Contour Active	0
☐	SD38	380 mm Only With QZT030214 SEDEO LITE Deep Contour Active	0
☐	SD40	400 mm	0
☐	SD42	420 mm	0
☐	SD44	440 mm	0
☐	SD46	460 mm	0
☐	SD48	480 mm Not With QZT030214 SEDEO LITE Deep Contour Active	0
☐	SD50	500 mm Not With QZT030214 SEDEO LITE Deep Contour Active	0

Armrest Types £

☐	QZT040122	Flip-up, height and angle adjustable	0
--------------------------	-----------	--------------------------------------	---

Arm Pads £

☐	QZT040040	Arm Pad, Integral Foam , 380 mm, Right	0
☐	QZT040041	Arm Pad, Integral Foam , 380 mm, Left	0
☐	QZT040204	Arm Pad, Comfort, Padded, PU coated, 400 mm, Right Not With QZT110056 Joystick Mounting, Parallel, Swing Away, Capsulated	0
☐	QZT040203	Arm Pad, Comfort, Padded, PU coated, 400 mm, Left Not With QZT110056 Joystick Mounting, Parallel, Swing Away, Capsulated	0

Side Guards £

☐	QZT040033	Side Guard	0
--------------------------	-----------	------------	---

Centre Mount Legrest £

☐	QZT050106	Centre Mount footrest, Short (includes calf pads and flip up footplate)	0
☐	QZT050107	Centre Mount footrest, Medium (includes calf pads and flip up footplate)	0

Swing Away Legrest & Footplates £

☐	QZT050002	Hanger Angle 70° (Quickie Life style)	0
☐	QZT050001	Hanger Angle 80° (Quickie Life style) Not With QZT130012 Powered Negative C.O.G Tilt -5 to 25° or 9" Castor QZT080210/11	0
☐	QZT050115	Legrest Uni 80° Swing Away	0
☐	QZT050016	Legrest Manual Elevating Swing Away (Includes Flip Up Footplates)	375
☐	QZT050017	Powered elevating legrest	695
☐	QZT050018	Powered articulating legrest	695
☐	QZT050146	Uni Amputee Support , Right	150
☐	QZT050147	Uni Amputee Support , Left	150

Footplates for Swing Away Legrests

☐	QZT050135	Footplate, Fixed, Pair Only With QZT050115 Legrest Uni 80° Swingaway	0
☐	QZT050137	Footplate, Angle Adjustable , Pair Only With QZT050115 Legrest Uni 80° Swingaway	80
☐	QZT050050	Footplate, Divided , Aluminium, Angle-Adjustable, Flip-Up (Sideways), Heel Loop, Black Not With QZT130012 Powered Negative C.O.G Tilt -5 to 25° or 9" Castor QZT080210/11 Only With QZT050002/1 Hanger Angle 70°/80°	0
☐	QZT050167	Footplates, High Mount - Recommended for short LLL Not With QZT130012 Powered Negative C.O.G Tilt -5 to 25° or 9" Castor QZT080210/11 Only With QZT050002/1 Hanger Angle 70°/80°	130
☐	QZT050168	Short Clamp (90 mm)	0
☐	QZT050169	Long Clamp (121 mm)	0

Accessories £

☐	QZT090100	Transportation tie down points (Included as Standard)	0
☐	QZT090171	Adjustable Lap strap	0
☐	QZT090040	Swing-away tray table	335
☐	QZT090110	Rear View Mirror, Right Not In Combination With Same Entry Side QZT090040 Tray Table	60
☐	QZT090118	Rear View Mirror, Left Not In Combination With Same Entry Side QZT090040 Tray Table	60

Remote Joystick – VR2 £

☐	QZT110036	Remote Joystick, VR2	0
☐	QZT110037	Remote Joystick, VR2 , 2 Actuators	0
☐	QZT110040	Remote Joystick, VR2, 2 Actuators and Lights	0

Joystick Mounting £

☐	QZT110082	Joystick Mounting, Parallel, Swing Away	0
☐	QZT110056	Joystick Mounting, Parallel, Swing Away, Capsulated Not With QZT040204/3 Arm Pad, Comfort, Padded, PU Coated, 400 mm, Right/Left	80

Joystick Mounting Location £

☐	QZT110052	Joystick Mounting Location, Right	0
☐	QZT110053	Joystick Mounting Location, Left	0

Joystick Knobs (optional) £

☐	QZT090505	Joystick Knob, Golf Ball, Ø 42 mm	40
☐	QZT090506	Joystick Knob, T-Bar, 82 mm Wide	40
☐	QZT090507	Joystick Knob, Soft Ball, Ø 70 mm	40
☐	QZT090508	Joystick Knob, Ball, Ø 35 mm	40
☐	QZT090045	Joystick Knob, Allround, Ball Ø 48 mm	40
☐	QZT090509	Joystick Knob, Chin Cup, Ø 37 mm	40
☐	QZT090512	Joystick Knob, Stick, 80 x 15 mm	40
☐	QZT090511	Joystick Knob, Mushroom, Ø 45 mm	40

Power Module – VR2 £

☐	QZT110001	Power module , VR2, 70A Not Recommended For Performance/ 4 pole motors	0
☐	QZT110002	Power module , VR2, 90A Not With 4/6kph Standard Motors	135

Lighting Module - VR2 £

☐	QZT110030	Lighting Module, VR2 (3) If LED Light And Indicator Is Chosen, Then This Is Required	0
--------------------------	-----------	---	---

Attendant Controls £

☐	QZT110021	Dual Control	425
☐	QZT110020	Joystick Mounted As Attendant Control Choose Also A Remote Joystick / Right Hand Attendant Control Only	270

Remote Joystick – R-Net £

☐	QZT110015	Remote Joystick, R-Net, LED, Actuators	485
☐	QZT110016	Remote Joystick, R-Net, LED, Actuators And Lights	485
☐	QZT110017	Remote Joystick, R-Net, Colour LCD	630
☐	QZT110018	Remote Joystick, R-Net, Colour LCD, Heavy Duty Not with any joystick knob QZT090505-QZT090511	1200
☐	QZT110013	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared Control Included	750
☐	QZT110023	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared And Bluetooth Control Included	925
☐	QZT110019	Omit Remote Joystick, R-Net	0

Joystick Mounting £

☐	QZT110082	Joystick Mounting, Parallel, Swing Away	0
☐	QZT110056	Joystick Mounting, Parallel, Swing Away, Capsulated Not With QZT040204/3 Arm Pad, Comfort, Padded, PU Coated, 400 mm, Right/Left	80

Joystick Mounting Location £

☐	QZT110052	Joystick Mounting Location, Right	0
☐	QZT110053	Joystick Mounting Location, Left	0

Power Modules – R-Net £

☐	QZT110003	Power Module, R-Net, 90 A	0
--------------------------	-----------	---------------------------	---

Seating / Lighting Module Control – R-Net £

☐	QZT110032	Actuator Control & Lights Module, R-Net, ISM-6L To Be Offered Only When R-Net With Actuators AND Lights is selected	0
☐	QZT110204	Actuator Control Module, R-Net, ISM-6 To Be Offered Only When R-Net With Actuators NO Lights is selected	0
☐	QZT110157	Light Module, R-Net, ISM-L	0

Attendant Controls £

☐	QZT110021	Dual Control	425
☐	QZT110020	Joystick Mounted As Attendant Control Choose Also A Remote Joystick / Right Hand Attendant Control Only	270

Special Controls

OMNI2 Input Control Module – R-Net

£

☐	QZT110740	OMNI2 incl. 2x Single Switch (Ø 64 mm) and Flexi Mount	1170
--------------------------	-----------	--	------

OMNI2 Module Mount

£

☐	QZT110122	OMNI2 Flexi Mount, Right	0
☐	QZT110121	OMNI2 Flexi Mount, Left	0

Safety Devices

£

☐	QZT090147	Remote Stop Switch (30 m Activation Range, Fit Via OBC Port) Not With VR2	500
--------------------------	-----------	--	-----

LiNX Controls

£

☐	QZT110810	LiNX Remote Drive Only, REM110	0
☐	QZT110811	LiNX Remote Drive 2 Act, REM210	105
☐	QZT110812	LiNX Remote Drive 6 Act, REM211	105
☐	QZT110813	LiNX Remote Drive 2 Actuators + Lights, REM215	475
☐	QZT110814	LiNX Remote Drive 6 Actuators + Lights, REM216	475
☐	QZT110815	LiNX Remote Drive 6 Act, Lights, Touchscreen, REM400	630
☐	QZT110818	Omit LiNX Remote	0

Linx Power Module & ISM

☐	QZT110802	Power Module 75A LiNX Drive Only	0
☐	QZT110803	Power Module 75A LiNX, Drive, Lights & Actuators	105
☐	QZT110804	Power Module 120A LiNX Drive Only	105
☐	QZT110805	Power Module 120A LiNX, Drive, Lights & Actuators	105
☐	QZT110821	4 Actuator Module LiNX, ACT-400	0
☐	QZT110822	Omit Actuator Module	0

Linx Gyro

☐	QZT110830	Traction Control System, for LiNX Gyro Technology	400
☐	QZT110831	Omit Traction Control	0

LiNX Joystick Mount

£

☐	QZT110836	Joystick Mounting, Parallel, Swing Away	0
☐	QZT110837	Joystick Mounting, Parallel, Swing Away, Capsulated Not With QZT040204/3 Arm Pad, Comfort, Padded, PU Coated, 400 mm, Right/Left	80

LiNX Joystick Mounting Location

£

☐	QZT110838	Joystick Mounting Location, Right	0
☐	QZT110839	Joystick Mounting Location, Left	0

LiNX Attendant Controls

£

☐	QZT110840	Dual Control	425
☐	QZT110842	Joystick Mounted As Attendant Control Choose Also a Remote Joystick / Right Hand Attendant Control Only	270
☐	QZT110843	Omit Attendant Control	0

LiNX Speciality Display

£

☐	QZT110823	LiNX Speciality Display, Omni, DX-REM500	1225
☐	QZT110824	Omit Speciality Display	0

Mounting Location

☐	QZT110850	Speciality Display Flexi, Right	0
☐	QZT110851	Speciality Display Flexi, Left	0

Drive Wheels & Castor Wheels & Mudguards

◀ BACK



// QZT070033



// QZT080210

LED Lights & Chargers

◀ BACK



// QZT110201



// QZT160100 & QZT160101



// QZT160102



// QZT160103 & QZT160104

Back Cushion Type

◀ BACK



// QZT030214

Headrest Pads & Headrest Mounts

◀ BACK



// QZT030430 & QZT030431



// WHT030401-WHT030403



// WHT030404-WHT030406



// QZT030433



// WHT030408



// WHT030409

Seat Cushion Type

◀ BACK



// QZT020120

JAY J3 Back

◀ BACK



// JAY030041

Armrest Type

◀ BACK



// QZT040122

Swing Away Legrest & Support Options

◀ BACK



// QZT050001

// QZT050050

// QZT050167

Remote Joystick – R-Net & Joystick Mounting

◀ BACK



// QZT110017



// QZT110013



// QZT110018



// QZT110056

Remote Joystick – VR2

◀ BACK



// QZT110036



// QZT110037

Joystick Knobs

◀ BACK



// QZT090505



// QZT090045



// QZT090507



// QZT090508



// QZT090509



// QZT090512



// QZT090511



// QZT090506

Electronic Programmers – R-Net[◀ BACK](#)

// QZT090066

Power Modules – R-Net / Attendant Controls[◀ BACK](#)

// QZT110021

R-Net OMNI2 Input Control Modules[◀ BACK](#)

// QZT110172

Comments For Adaptation



Q300 M Mini Teens

£

Base Price	00,000.00
Total Of Additional Equipment	
Packaging And Transport	
Exclusive VAT	
VAT	
Total Price	

Date

Signature

Day

Month

Year

Send E-Mail

**Sunrise Medical GmbH**

Kahlbachring 2-4

D-69254 Malsch/HD

Phone: +49 (0) 72 53 / 9 80-0

Fax: +49 (0) 72 53 / 9 80-222

E-Mail: kundenservice@sunrisemedical.de

www.SunriseMedical.eu