



Q400® M

with SEDEO® LITE



ORDER FORM

Billing Address Or Customer-ID		Shipping Address		☐ Identical To Billing Address
C.-ID:		C.-ID:		
Name:		Name:		
Street:		Street:		
Town/PC:		Town/PC:		
E-Mail:		E-Mail:		
Phone:		Phone:		
Order Information		User Weight		
☐ Order ☐ Calculation ☐ Built-4-Me®		User Weight Capacity: 136 kg		
Purchase Order No.:		User Weight: kg		
Art. No.: Q400 M: 2280040M-000,				
Order Date: Day Month Year				

Base Price **£6,850**

Seating Specifications

All Bases Prepared For DAHL Docking Station		Seat Type	SEDEO LITE
Typical WC Weight	Without Batteries 140 kg	Seat Size Range	Width 385 - 510 mm, Depth 410 - 560 mm
	With Batteries 165 kg	Seat To Floor Height Range	415 mm - 490 mm

MWD Base Specifications

User Weight Capacity	136 kg	Drive Wheels	13" Drive Wheels (Black)
Base Size	Width Min. 620 mm	Tires / Castors	7" Front / 7" Rear (Black)
Drive Control	VR2, Option R-Net	Charger	8 A Dual Mode Charger
Driving Range	Up To 40 km	Battery Size	50, 60 And 80 Ah Batteries
Drive Motors	2-Pole 6 kph, Options 10 kph	Colour Accents	4 Colour Accents
Ground Clearance	70 mm	Frame Colour	Black

Q400 Drive Bases			£
☒	QPL000031	Q400 M (MWD Base Chair)	6850

Seat Cushion Type			£
☐	QPL020098	Seat Cushion, Standard, Height 50 mm	0
☐	QPL020110	Omit Seat Cushion	0

Standard Seat Size					£			
SEDEO / Std Seat Cushion Width			Seat Depth		JAY Seat Width - Only for JAY Cushions			
☐	SW385	15" (385 mm) (1)	☐	SD41	16" (410 mm) (3)	☐	SW38	15" (38 mm)
☐	SW41	16" (410 mm) (1)	☐	SD435	17" (435 mm) (3)	☐	SW40	15.75" (400 mm)
☐	SW435	17" (435 mm) (1)	☐	SD46	18" (460 mm) (3)	☐	SW42	16.5" (420 mm)
☐	SW46	18" (460 mm) (2)	☐	SD485	19" (485 mm) (3)	☐	SW44	17.3" (440 mm)
☐	SW485	19" (485 mm) (2)	☐	SD51	20" (510 mm) (3)	☐	SD46	18" (460 mm) (3)
☐	SW51	20" (510 mm) (2)	☐	SD535	21" (535 mm) (3) (4)	☐	SW48	19" (480 mm)
			☐	SD56	22" (560 mm) (3) (4)	☐	SW50	19.7" (500 mm)

Backrest Frame Size			£
☐	QPL030154	Backrest, Width-Adjustable, 385 - 485 mm	0
☐	QPL030155	Backrest, Width-Adjustable, 460 - 510 mm	0

Comfort Seat Cushion (5)			£
☐	QPL020073	Profiled Cushion, Height 70 mm	160
☐	QPL020075	Support Memory Cushion, Height 90 mm	260

Seat Boards			£
☐	QPL020082	Seat Board, Flat, Depth-Adjustable, Cushions Guide At Rear	130
☐	QPL020083	Seat Pan, 50 mm, Depth-Adjustable, Cushions Guide At Rear	150

JAY Seat Cushions (6)			£
☐	JAY020006	JAY Easy Visco With Flat Base	180
☐	JAY020011	JAY Balance, Fluid, Microclimatic Cover, Until Seat Width 500 mm	655
☐	JAY020015	JAY Balance, Fluid, Microclimatic Cover, With Positioning Kit, Until Seat Width 500 mm	755

(1) Comes with cushion width 410mm

(2) Comes with cushion width 485 mm

(3) Comes with cushion depth 50cm

(4) Not with SEDEO Lite back cushion QPL03130-133

(5) Comes with seat cushion cover, incontinence protection

(6) Not all JAY Cushions are available in all sizes - for availability check our online configurator

Seat Heights (1)			£
☐	SH415	415 mm (2) (3)	0
☐	SH44	440 mm	0
☐	SH465	465 mm	0
☐	SH49	490 mm	0

Powered Seating Options			£
☐	QPL130011	Fixed Tilt Module (Adjustable 0°, 3°, 6°, 9°, Factory Setting 6°)	0
☐	QPL130004	Posterior Tilt, Powered, 0° - 30° (4)	1150
☐	QPL130012	Posterior Tilt, Powered, -5° - 25° (4) (5)	1290
☐	QPL130007	Posterior Tilt, Powered, 0° - 50° (4) (6)	2250
☐	QPL130009	Seat Elevator (300 mm) / Tilt Module (Tilt 0° - 30°), Powered (4)	2050
☐	QPL130006	Seat Elevator (300 mm) / Tilt Module (Tilt 0° - 50°), Powered (4)	2305

Backrest Frame Type			£
☐	QPL030147	Backrest, Tool Angle Adjustable, 86°, 90°, 94°, 98°, 102°, 106°	0
☐	QPL030148	Backrest, Manual Recline, Attendant Operated Only, 85° - 120°	310
☐	QPL030149	Backrest, Powered Recline, 85° - 120° (4)	560

Backrest Frame Option			£
☐	QPL030156	Backrest, Fixed	0
☐	QPL030098	Backrest, Flip Down, Quick Release	185

Backrest Cushion Type			£
☐	QPL030150	Sling Back, With Tension Adjustable Straps (Section A)	0

Sling Back Height			£
☐	BH46	460 mm	0
☐	BH485	485 mm	0
☐	BH51	510 mm	0

JAY Shape Back			£
☐	JAY080001	JAY Shape Back	270
☐	JAY080006	JAY Shape Back, Oversized (7)	270

JAY Shape Backrest Height			£
☐	BH47	470 mm	0
☐	BH495	495 mm	0
☐	BH52	520 mm	0
☐	BH56	560 mm (8)	0
☐	BH585	585 mm (8)	0
☐	BH61	610 mm (8)	0

JAY Shape Back			£
☐	JAY080010	JAY Shape Lateral Thoracic Pad, (LxH) 150 x 150 mm, Pair	270
☐	JAY080011	JAY Shape Lateral Thoracic Pad, (LxH) 100 x 152 mm, Pair	270

(1) Measured from front of the seat pan to the floor

(2) Not with 80Ah batteries

(3) If with QPL130009 Seat Elevator (300 mm) / Tilt Module (Tilt 0° - 30°), Powered, or with QPL130006 Seat Elevator (300 mm) / Tilt Module (Tilt 0° - 50°), Powered, only with QPL030147 Backrest, Tool Angle Adjustable, 86°, 90°, 94°, 98°, 102°, 106°

(4) Not with QPL110036 Remote Joystick VR2. In case of R-Net QPL110032 Actuator control & lights module, R-Net ISM-6L Required

(5) Not with seat height 415 or 440 mm

(6) Only with QPL030147 Backrest tool angle adjustable

(7) Supplied 1" larger than seat width to allow greater immersion, not with seat width 485 / 510 mm

(8) Not with QML030148 Backrest, Manual Recline, Attendant Operated Only, 85° - 120° or QML030149 Backrest, Powered Recline, 85° - 120°

SEDEO Comfort Backrest		£
☐	QPL030130 SEDEO, Back Cushion, Profiled, Height 500 mm, Comfort Cover	210
☐	QPL030131 SEDEO, Back Cushion, Profiled, Adjustable Straps, Height 570 mm, Comfort Cover	295
☐	QPL030132 SEDEO, Back Cushion, Deep Contour Passive, Adjustable Straps, Height 570 mm, Comfort Cover	350
☐	QPL030133 SEDEO, Back Cushion, Deep Contour Active, Adjustable Straps, Height 570 mm, Comfort Cover	340

SEDEO Lateral Thoracic Supports		£
☐	QPL030607 SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Right	160
☐	QPL030608 SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Left	160

Flexi Backrest Size		£
Backrest Width		
☐	BW38 380mm, For Effective Back Width 380 mm	0
☐	BW42 420 mm, For Effective Back Width 380 mm - 420 mm	0
☐	BW46 460 mm, For Effective Back Width 420 - 460 mm	0
☐	BW50 500 mm, For Effective Back Width 460 - 500 mm	0
☐	BW54 540 mm, For Effective Back Width 500 - 540 mm	0
☐	BW56 560 mm, For Effective Back Width 520 - 560 mm	0
Backrest Height		
☐	BH44 440 mm, For Effective Back Height 440 - 500 mm	0
☐	BH50 500 mm, For Effective Back Height 500 - 560 mm	0
☐	BH56 560 mm, For Effective Back Height 560 - 620 mm	0

SEDEO Flexi Back Cushion		£
☐	QPL030115 SEDEO Flexi, Back Cushion Flexible, Deep Contour, Passive, Adjustable Straps	435
☐	QPL030116 SEDEO Flexi, Back Cushion Flexible, Deep Contour, Active, Adjustable Straps	450
☐	QPL030117 SEDEO Flexi, Back Cushion Flexible, Extra Deep Contour, Active, Adjustable Straps	475

SEDEO Flexi Back Cushion Cover		£
☐	QPL030059 SEDEO Flexi Cover, Incontinence Protection	0
☐	QPL030108 SEDEO Flexi Cover, Comfort	0

Headrest Receiver		£
☐	QPL030120 Headrest Receiver, Backrest Frame Mount	0
☐	QPL030119 Headrest Receiver, Back Cushion Mount	0

Headrest Pad		£
☐	QPL030430 SEDEO, Headrest Pad, Small, (LxH) 240 x 150 mm (1)	95
☐	QPL030431 SEDEO, Headrest Pad, Small, With Cover, (LxH) 240 x 150 mm (1)	110
☐	QPL030432 SEDEO, Headrest Pad, Flexible, (LxH) 400 x 110 mm (1)	170
☐	WHT030402 Whitmyer, Headrest Pad, Plush, (LxH) 254 x 127 mm (2)	170
☐	WHT030404 Whitmyer, Headrest Pad, Adjust-A-Plush, Small, (LxH) 355 x 76 mm (2)	280

Headrest Stem (3)		£
☐	QPL030433 SEDEO, Headrest Mount	140
☐	WHT030408 Whitmyer, Headrest Mount, LINX ²	190
☐	WHT030409 Whitmyer, Headrest Mount, AXYS	155
☐	WHT030411 Whitmyer, Headrest Mount, COBRA XTRA Flip Back	220

(1) Only use with SEDEO headrest mount

(2) Only use with Whitmyer headrest mount

(3) Choose a Whitmyer headrest mount and headrest receiver in case you would like to add a head control in the special controls section

JAY J3 Back Size			£
☐	JAY030111	JAY J3 Back Size, 410 mm	0
☐	JAY030112	JAY J3 Back Size, 460 mm	0
☐	JAY030113	JAY J3 Back Size, 510 mm	0
JAY J3 Back Size, Mid Deep Contour			£
☐	JAY030042	JAY J3 Back, Mid Deep Contour, Upper Thoracic, Medium, 460 mm	870
JAY J3 Back Size, Deep Contour			£
☐	JAY030021	JAY J3 Back, Deep Contour, Upper Thoracic, Medium, 460 mm	870
Armrest Types			£
☐	QPL040039	Armrest, Flip Back	0
☐	QPL040123	Armrest, Reclining, Flip Back	240
☐	QPL040022	Armrest, Single Post, SEDEO LITE (1)(2)	130
☐	QPL040019	Armrest, Single Post, Flip Back (1)(3)	250
Side Guards			£
☐	QPL040033	Side Guard (7)	0
Arm Pads			£
☐	QPL040040	Arm Pad, Integral Foam, 380 mm, Right	0
☐	QPL040041	Arm Pad, Integral Foam, 380 mm, Left	0
☐	QPL040025	Arm Pad, Integral Foam, 400 mm, Right (4)	0
☐	QPL040026	Arm Pad, Integral Foam, 400 mm, Left (4)	0
☐	QPL040204	Arm Pad, Comfort, Padded, PU Coated, 400 mm, Right	75
☐	QPL040203	Arm Pad, Comfort, Padded, PU Coated, 400 mm, Left	75
☐	QPL040031	Arm Pad, Hemiplegia, Integral Foam, Right (5)	275
☐	QPL040032	Arm Pad, Hemiplegia, Integral Foam, Left (6)	275
Thigh Supports & Pommel (8)			£
☐	QPL020113	Thigh Support, Lateral, Small, 120 x 120 mm, Swing Away, Right	185
☐	QPL020114	Thigh Support, Lateral, Small, 120 x 120 mm, Swing Away, Left	185
☐	QPL020121	Thigh Support, Lateral, Large, 180 x 120 mm, Swing Away, Right	185
☐	QPL020122	Thigh Support, Lateral, Large, 180 x 120 mm, Swing Away, Left	185
☐	QPL030605	Abduction Pommel (9)	255
Centre Mount Legrest (Single Piece Footplate)			£
☐	QPL050086	Centre Mount Legrest, Tool Adjustable, With Flip Up And Angle-Adjustable Footboard	0
Pad Options			
☐	QPL050087	Centre Calf Support, Integral Foam	0
☐	QPL050092	Calf Pad, Integral Foam, Pair	0
Fixed Centre Mount Legrest (Split Footplate)			£
☐	QPL050140	Fixed Centre Mount Legrest	125
☐	QPL050072	Footplate, Dual, Height- And Angle-Adjustable, 260 x 150 mm	0
☐	QPL050075	Calf Pad, Padded, Pair	0

(1) QPL040033 Side Guard Included

(2) QPL040040/41 Arm Pad, Integral Foam, 380 mm right/left included

(3) Not With QPL040040-41 Armpad 380 mm, Integral Foam, Right Or Left; Not With Seat To Floor Height 415 mm

(4) Not With QPL040022 Armrest, Single Post, SEDEO LITE

(5) Not With QPL110052 Joystick Mount, Right, QPL090040 Tray Table, Swing Away, QPL090046 Tray Table, Half-Side, Swing Away, QPL090110 Rear View Mirror, Right

(6) Not With QPL110053 Joystick Mount, Left, QPL090040 Tray Table, Swing Away,

QPL090046 Tray Table, Half-Side, Swing Away, QPL090118 Rear View Mirror, Left

(7) Only With QPL040039 Armrest, Flip Back & QPL040123 Armrest, Reclining, Flip Back

(8) Not With Seat To Floor Height 415 mm

(9) Not with any centre mount legrest

Powered Centre Mount		£
☐	QPL050076 Centre Mount Legrest, Powered, 90° - 180°, 175 mm Extension	1470
Footplates		
☐	QPL050072 Footplate, Dual, Height- And Angle-Adjustable, 260 x 150 mm	0
☐	QPL050073 Footplate, Dual, Height- And Angle-Adjustable, Large, 310 x 170 mm	0
Pad Options		
☐	QPL050075 Calf Pad, Padded, Pair	0
SEDEO Lite Swing Away Legrest		£
☐	QPL050100 SEDEO LITE, Legrest, 70°, Swing Away	0
☐	QPL050105 SEDEO LITE, Legrest, Manual Elevating, Swing Away	375
Uni Swing Away Legrest		£
☐	QPL050115 Legrest, UNI, 80°, Swing Away	0
Amputee Support		
☐	QPL050142 UNI Amputee Supports, Padded, Right	150
☐	QPL050143 UNI Amputee Supports, Padded, Left	150
SEDEO Pro Swing Away Legrest		£
☐	QPL050015 SEDEO PRO, Legrest, 70°, Swing Away	0
☐	QPL050016 SEDEO PRO, Legrest, Manual Articulating, Swing Away	375
☐	QPL050017 SEDEO PRO, Legrest, Power Articulating, Swing Away & Auto Connect	695
Footplates		
☐	QPL050088 Footplate, For Seat Width 385 - 400 mm	0
☐	QPL050089 Footplate, For Seat Width 420 - 460 mm	0
☐	QPL050090 Footplate, For Seat Width 480 - 520 mm	0
☐	QPL050091 Footplate, For Seat Width 540 - 560 mm	0
☐	QPL050081 Footplates, Coupled (1) (2)	0
Hanger Receivers		
☐	QPL050170 SEDEO PRO, Hanger Receivers, High, For Seat Cushion Thicknes 70mm And Higher	0
☐	QPL050171 SEDEO PRO, Hanger Receivers, Low, For Seat Cushion Thicknes 50mm And Lower	0
☐	QPL050080 SEDEO PRO, Hanger Receivers Only, High	0
☐	QPL050172 SEDEO PRO, Hanger Receivers Only, Low	0
Calf Pad / Strap		
☐	QPL050082 Calf Support, Integral Foam, Depth-Adjustable, Pair	100
☐	QPL050083 Calf Support, Padded, Depth-Adjustable, Pair	135
☐	QPL050084 Heel Loop	0
☐	QPL050085 Calf Strap (2)	0
☐	QPL050200 Knee Pad, Integral Foam	75

(1) Choose Also A Footplate

(2) Only For QPL050015 SEDEO PRO, Legrest, 70°, Swing Away

(3) Extra On A Set Of Swing Away Legrests, Attached Lose

VR2 Controls			£
☐	QPL110036	Remote Joystick, VR2	0
☐	QPL110037	Remote Joystick, VR2, 2 Actuators	0
☐	QPL110040	Remote Joystick, VR2, 2 Actuators And Lights	0
VR2 Power Module & Lighting Module			
☐	QPL110001	Power Module, VR2, 70 A (1)	0
☐	QPL110002	Power Module, VR2, 90 A (2)	135
☐	QPL110030	Lighting Module, VR2 (3)	0
R-Net Controls			£
☐	QPL110015	Remote Joystick, R-Net, LED, Actuators	485
☐	QPL110016	Remote Joystick, R-Net, LED, Actuators And Lights	485
☐	QPL110017	Remote Joystick, R-Net, Colour LCD	600
☐	QPL110018	Remote Joystick, R-Net, Colour LCD, Heavy Duty	1200
☐	QPL110013	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared Control Included	770
☐	QPL110023	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared And Bluetooth Control Included	925
R-Net Power Module & ISM			
☐	QPL110003	Power Module, R-Net, 90 A	0
☐	QPL110004	Power Module, R-Net, 120 A	135
☐	QPL110032	Actuator Control & Lights Module, R-Net, ISM-6L	0
R-Net Gyro			
☐	QPL110009	Traction Control System For R-Net (GYRO Technology)	400
Joystick Mount			£
☐	QPL110050	Joystick Mounting, Fixed	0
☐	QPL110082	Joystick Mounting, Parallel, Swing Away	0
☐	QPL110056	Joystick Mounting, Parallel, Swing Away, Capsulated (4)	80
Joystick Mounting Location			£
☐	QPL110052	Joystick Mounting Location, Right	0
☐	QPL110053	Joystick Mounting Location, Left	0
Attendant Controls			£
☐	QPL110021	Dual Control	425
☐	QPL110020	Joystick Mounted As Attendant Control (5)	270
Attendant Mounting Location			
☐	QPL110025	Attendant Mounting Location, Right	0
☐	QPL110026	Attendant Mounting Location, Left	0
Joystick Knobs			£
☐	QPL090505	Joystick Knob, Golf Ball, Ø 42 mm	40
☐	QPL090506	Joystick Knob, T-Bar, 82 mm Wide	40
☐	QPL090507	Joystick Knob, Soft Ball, Ø 70 mm	40
☐	QPL090508	Joystick Knob, Ball, Ø 35 mm	40
☐	QPL090045	Joystick Knob, Allround, Ball, Ø 48 mm	40
☐	QPL090509	Joystick Knob, Chin Cup, Ø 37 mm	40
☐	QPL090512	Joystick Knob, Stick, 80 x 15 mm	40
☐	QPL090511	Joystick Knob, Mushroom, Ø 45 mm	40
LED Lights & Indicators			£
☐	QPL110201	LED Light And Indicator StVZO Approved (6)	400

(1) Only With 6 kph

(2) Highly Recommended With User Weight >100 kg

(3) If LED Light And Indicator Is Chosen, Then This Is Required

(4) Not With QPL040022 Armrest, Single Post

(5) Choose Also Remote Joystick

(6) QPL110030 Lighting Module VR2 Or QPL110032 Actuator Control & Lights Module, R-Net, ISM-6L Required

Remote Joystick – LinX		£
☐ QPL110810	REM110, Drive Only	0
☐ QPL110811	REM211, Drive & Actuators	105
☐ QPL110812	REM211, 6 Actuators	105
☐ QPL110813	REM215 Actuators	475
☐ QPL110814	REM216 6 Actuators	475
☐ QPL110815	REM400, Drive 6 Act,	630
☐ QPL110823	REM500 Omni	1225
☐ QPL110818	Omit Joystick	0
☐ QPL110824	Omit Omni	0

Power Modules & ISM's – LinX		£
☐ QPL110802	Power Module, LinX, 75 A, Drive Only	0
☐ QPL110803	Power Module, LinX, 75 A , Actuators and Lights	105
☐ QPL110804	Power Module, LinX, 120 Ah, Drive Only	135
☐ QPL110805	Power Module, LinX, 120 Ah , Actuators and Lights	135

Joystick Mounting - LiNX		£
☐ QPL110835	Fixed Joystick	0
☐ QPL110836	Joystick Mounting, Parallel, Swing Away	0
☐ QPL110837	Parallel Joystick	80

Joystick Mount Location - LiNX		£
☐ QPL110838	Joystick Mounting Location, Right	0
☐ QPL110839	Joystick Mounting Location, Left	0

Attendant Controls - LiNX		£
☐ QPL110840	Dual Attendant Control, LiNX	425
☐ QPL110842	Joystick Mounted As Attendant Control	270
☐ QPL110843	Omit Attendant Control	0

Attendant Mounting Location		
☐ QPL110844	Attendant Mounting Location, Right	0
☐ QPL110845	Attendant Mounting Location, Left	0
☐ QPL110846	Omni Fixed Right	0
☐ QPL110847	Omni Fixed Left	0
☐ QPL110848	Omni SA Right	130
☐ QPL110849	Omni SA Left	130

Additional Modules- LiNX		£
☐ QPL110830	Traction Control System, for LiNX Gyro Technology	400
☐ QPL110831	Omit Gyro	0
☐ QPL110820	2 Actuator Module LiNX, ACT-200	0
☐ QPL110821	4 Actuator Box	0
☐ QPL110822	Omit Actuator Box	0

Drive Wheels & Castor Wheels		£
☐	QPL070065 Drive Wheel, 13", Pneumatic, Knobbly, Black	0
☐	QPL070069 Drive Wheel, 13", Solid, Knobbly, Black	0
☐	QPL070034 Drive Wheel, 14", Pneumatic, Knobbly, Black	0
☐	QPL070033 Drive Wheel, 14", Solid, Knobbly, Black	0
☐	QPL080213 Castor Wheel, 7", Solid, Black	0

Suspension Setting		£
☐	QPL150001 Standard Suspension Setting - Selected as default	0
☐	QPL150002 Firm Suspension Setting (6) (7)	0

Motor Packages		£
☐	QPL180012 2-Pole Motor, 6 kph/ 4 mph, Performance	0
☐	QPL180013 2-Pole Motor, 10 kph / 6.5 mph, Performance	500
☐	QPL180009 4-Pole Motor, 6 kph / 4 mph	310
☐	QPL180007 4-Pole Motor, 10 kph / 6.5 mph	390

Colour Accents		£
☐	QPL140005 Apple Green (RAL 110 70 77)	0
☐	QPL140004 Black (RAL 9005)	0
☐	QPL140006 Orange (RAL 060 60 80)	0
☐	QPL140007 Silver (RAL 9006)	0

Batteries		£
☐	QPL160003 Battery, 50 Ah, Gel, Maintenance-Free	0
☐	QPL160001 Battery, 60 Ah, Gel, Maintenance-Free	0
☐	QPL160005 Battery, 80 Ah, Gel, Maintenance-Free (2)	160

Chargers		£
☐	QPL160100 Charger, Automatic, 220 V / 8 A, UK	0
☐	QPL160103 Charger, Automatic, 220 V / 10 A, UK	240

Lap Belt & JAY Belts		£
☐	QPL090128 Positioning Belt, Carseat Buckle, Short, 50 x 1700 mm	0
☐	QPL090123 Seat Rail Mounting Brackets For Positioning Belt	0
☐	JAY090100 JAY Standard 2-Point Pelvic Belt, Side Release Buckle, Pad 190 x 40 mm	60
☐	JAY090101 JAY Standard 2-Point Pelvic Belt, Side Release Buckle, Pad 240 x 40 mm	60
☐	JAY090077 JAY Medial 4-Point Pelvic Belt, Push Button Buckle, Pad 140 x 40 mm (3)	80
☐	JAY090078 JAY Medial 4-Point Pelvic Belt, Push Button Buckle, Pad 165 x 40 mm (3)	80
☐	JAY090102 JAY V-Contour Pelvic Belt, Side Release Buckle, Pad Medium	110
☐	JAY090103 JAY V-Contour Pelvic Belt, Side Release Buckle, Pad Large	110

Accessories		£
☐	QPL090110 Rear View Mirror, Right (4)	60
☐	QPL090118 Rear View Mirror, Left (4)	60
☐	QPL090009 Crutch Holder (5)	70

Tray Table (Mounted Opposite Joystick)		
☐	QPL090124 Tray Table, Swing Away	335
☐	QPL090040 Tray Table, Width-Adjustable, Swing Away	335
☐	QPL090046 Tray Table, Half-Side, Swing Away	265

(2) Requires QPL160103 Battery Charger 220 V / 10 A

(3) (Choose Also QPL090123 Seat Rail Mounting Brackets For Positioning Belt

(4) Not In Combination With Same Entry Side QPL090040 Tray Table, Swing Away, QPL090046 Tray Table, Half-Side, Swing Away

(5) Not With QPL030148 Backrest, Manual Recline, Attendant Operated Only, 85° - 120° & QPL030149 Backrest, Powered Recline, 85° - 120°

(6) QPL150002 (Firm Suspension) recommended with QPL130009 Powered Lift/Tilt

(7) QPL150002 (Firm Suspension) recommended with QPL130004/QPL130012/QPL130007 Powered Tilt for clients over 100kg user weight

C-Me - Drive at 4.8kph in full Lift (7)		£
☐	QPL110034 C-Me Function, Enabled (8)	0
☐	QPL110035 C-Me Function, Disabled (9)	0

5-Switch Box		£
☐	QPL110043 #1 Ctrl+5 Switch Box With Free Assignable Functions Voice activation functionality included	575

#1 Ctrl+5 Switch Box Mounting Location

☐	QPL110180 #1 Ctrl+5 Switch Box Mounting Location, Behind Joystick (1)	0
☐	QPL110181 #1 Ctrl+5 Switch Box Mounting Location, Opposite Joystick	0

#1 Ctrl+5 Switch Type

☐	QPL110175 #1 Ctrl+5 Switch Type With Push Buttons	0
☐	QPL110176 #1 Ctrl+5 Switch Type With Toggle Switches	0

Omni2 Input Control Module		£
☐	QPL110116 OMNI2 Input Control Module, R-Net, Bluetooth, Infrared (2)	1170
☐	QPL110117 OMNI2 Input Control Module, R-Net, Scanning Mode, Bluetooth, Infrared	1170

Omni 2 Module Mount		£
☐	QPL110119 OMNI2 Module, Fixed Mount, Right	0
☐	QPL110120 OMNI2 Module, Fixed Mount, Left	0
☐	QPL110123 OMNI2 Module, Swing Away Mount, Right	130
☐	QPL110124 OMNI2 Module, Swing Away Mount, Left	130

Hand Control		£
☐	QPL110301 Hand Control, Switch-It MicroGuide	3145
☐	QPL110302 Hand Control, Switch-It MicroPilot	3495
☐	QPL110304 Hand Control, Switch-It VersaGuide	1120
☐	QPL110305 Hand Control, Switch-It VersaGuide EZ (3)	1330
☐	QPL110306 Hand Control, Allround Joystick	1120
☐	QPL110307 Hand Control, Allround Joystick Light (4)	1270
☐	QPL110308 Hand Control, Allround Joystick Heavy Duty	2095
☐	QPL110309 Hand Control, Allround Joystick, R-Net Bus Enabled (5)	1190
☐	QPL110310 Hand Control, Allround Joystick Light, R-Net Bus Enabled (5)	1400
☐	QPL110311 Hand Control, Allround Joystick Heavy Duty, R-Net Bus Enabled (5)	2235

Hand Control Mounting		£
☐	QPL110342 Hand Control Mounting, Link-It To The Space Disk, Hand Tray (6)	700

Hand Control Mounting

☐	QPL110350 Hand Control Mounting, Link-It To The Armrest, Fixed	210
☐	QPL110353 Hand Control Mounting, Link-It To The Armrest, Swing Away	220

Hand Control Mounting Location

☐	QPL110362 Hand Control Mounting Location, Right	0
☐	QPL110361 Hand Control Mounting Location, Left	0

(1) Not Possible With QPL110056 Joystick Mounting, Parallel Swing Away, Capsulated (Here We Can't Swing Away)

(2) Display With Built In Infrared & Armrest Mount, Required For Specialty Input Devices

(3) Not With QPL090512 Joystick Knob, Stick, 80 x 15 mm, QPL090511 Joystick Knob, Mushroom, Ø 45 mm, QPL090506 Joystick Knob, T-Bar, 82 mm Wide

(4) Only Use With Default Joystick Knob Or Joystick Knob, Allround, Ball Ø 48 mm

(5) OMNI Module Not Required

(6) Only For Micro Series Joysticks

(7) Only With QML1300009

(8) 4.5 kph Maximum Speed When Lifted (Tilt Angle <9°).

(9) 1.5 kph Maximum Speed When Lifted

Table Top Controls		£
☐	QPL110143 Table Top Control, LCD, Centre Tray Mounted (6)(7)(8)	2235
Tray Table Mounting Location		
☒	QPL110080 Tray Table Mounting Location, Left	0
☐	QPL110081 Tray Table Mounting Location, Right	0
Head Control (1)		£
Proportional Input Devices		
☒	QPL110370 Head Control, Switch-It Dual Pro, 3-Way Proportional, Facial Pads	5950
☐	QPL110371 Head Control, Switch-It Dual Pro, 3-Way Proportional, Spot Pads	5950
Chin Control		£
☐	QPL110470 Chin Control, Switch-It MicroGuide	3145
☐	QPL110471 Chin Control, Switch-It MicroPilot	3495
☐	QPL110472 Chin Control, Switch-It VersaGuide	1120
☐	QPL110473 Chin Control, Switch-It VersaGuide EZ	1330
☐	QPL110474 Chin Control, Allround Joystick	1120
☐	QPL110475 Chin Control, Allround Joystick Light	1330
☐	QPL110476 Chin Control, Allround Joystick, R-Net Bus Enabled (2)	1190
☐	QPL110477 Chin Control, Allround Joystick Light, R-Net Bus Enabled (2)	1400
Chin Control Mounting		£
☒	QPL110480 Chin Control Mounting, Link-It To Chin Harness, Small (3)	630
☐	QPL110481 Chin Control Mounting, Link-It To Chin Harness, Large (3)	700
☐	QPL110483 Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Right	560
☐	QPL110484 Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Left	560
☐	QPL110485 Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Right	2795
☐	QPL110486 Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Left	2795
☐	QPL110487 Chin Control Mounting, Link-It To Headrest, Swing Away, Left	490
☐	QPL110488 Chin Control Mounting, Link-It To Headrest, Swing Away, Right	490
Link-It Extension Kit		
☐	QPL110590 Link-It Extension Kit, Base Module, 100 mm Extension Tube, 5 Cable Clips	170
Sip ,N Puff (4)		£
☐	QPL110510 Sip 'N Puff With Link-It Hardware, Swing Away Mount To The Headrest (5)	700

(1) For switch selection see section H. Choose a Whitmyer headrest mount.

(2) Omni 2 module required

(3) If you need switches mounted on the harness, you can only use All Round joystick

(4) For switch selection see section H

(5) Headrest to be selected

(6) Not with SW40, SW41, SW42, SW48, SW485

(7)Not with QPL040022 Armrest, Single Post, Sedeo Lite, QPL110309 Hand Control, Allround Joystick, QPL110310 Hand Control, Allround Joystick Light, QPL110311 Hand Control, Allround Joystick heavy Duty

(8) Only R-Net Controls, comes with joystick

#1 Switches - Section H (1)			£
☐	QPL110611	#1 Single Switch, Ø 64 mm, Black	120
☐	QPL110612	#1 Single Switch, Ø 64 mm, Red	120
☐	QPL110613	#1 Single Switch, Ø 64 mm, Green	120
☐	QPL110523	#1 Egg Switch – Black	135
☐	QPL110524	#1 MicroLite Switch – Black	195
☐	QPL110651	#1 Single Switch, Ø 30 mm, Black	120
☐	QPL110654	#1 Single Switch, Ø 30 mm, Blue	120
☐	QPL110653	#1 Single Switch, Ø 30 mm, Green	120
☐	QPL110652	#1 Single Switch, Ø 30 mm, Red	120
☐	QPL110655	#1 Single Switch, Ø 30 mm, Yellow	120
☐	QPL110534	#1 Twister Switch – Black	195
☐	QPL110535	#1 Twister Switch – Green	195
☐	QPL110536	#1 Twister Switch – Red	195
☐	QPL110537	#1 Satellite Switch – Black (7)	210
☐	QPL110538	#1 Satellite Switch – Green (7)	210
☐	QPL110539	#1 Satellite Switch – Red (7)	210

#1 Switch Mounting			£
☐	QPL110560	#1 Switch Mounting, Link-It To The Hand Control (2)	210
☐	QPL110566	#1 Switch Mounting, Link-It To The Armrest (Opposite Site As The Hand Control) (3)	210
☐	QPL110561	#1 Switch Mounting, Link-It To The Headrest	280
☐	QPL110563	#1 Switch Mounting, Link-It To The Seat Rail (4)	210
☐	QPL110565	#1 Switch Mounting, Link-It To Chin Boom (5)	210
☐	QPL110569	#1 Omit Switch Mounting (6)	0

#1 Switch Mount Location			£
☐	QPL110572	#1 Switch Mounting Location, Right	0
☐	QPL110571	#1 Switch Mounting Location, Left	0

#2 Switches - Section H (1)			£
☐	QPL110631	#2 Single Switch, Ø 64 mm, Black	120
☐	QPL110632	#2 Single Switch, Ø 64 mm, Red	120
☐	QPL110633	#2 Single Switch, Ø 64 mm, Green	120
☐	QPL110088	#2 Egg Switch – Black	125
☐	QPL110089	#2 MicroLite Switch – Black	185
☐	QPL110671	#2 Single Switch, Ø 30 mm, Black	120
☐	QPL110674	#2 Single Switch, Ø 30 mm, Blue	120
☐	QPL110673	#2 Single Switch, Ø 30 mm, Green	120
☐	QPL110672	#2 Single Switch, Ø 30 mm, Red	120
☐	QPL110675	#2 Single Switch, Ø 30 mm, Yellow	120
☐	QPL110095	#2 Twister Switch – Black	185
☐	QPL110096	#2 Twister Switch – Green	185
☐	QPL110097	#2 Twister Switch – Red	185
☐	QPL110111	#2 Satellite Switch – Black (7)	200
☐	QPL110112	#2 Satellite Switch – Green (7)	200
☐	QPL110113	#2 Satellite Switch – Red (7)	200

(1) Select 1 Switch And Mounting Position Accordingly. For A Second Switch Please Use #2 Switch Selection Part Below. Single Function For Special Input Devices

(2) Not With HD Allround Joystick Or Space Disk Mount

(3) Armrest Rail Mount

(4) Knee, If 2 Switches Selected Only Left And Right Mount

(5) No Left / Right Switch Selection Required

(6) Switches Attached Lose

(7) Only For Allround Joystick, Maximum 2 Selectable

#2 Switch Mounting		£
☐	QPL110098 #2 Switch Mounting, Link-It To The Hand Control (1)	210
☐	QPL110099 #2 Switch Mounting, Link-It To The Armrest (Opposite Site As The Hand Control) (2)	210
☐	QPL110100 #2 Switch Mounting, Link-It To The Headrest	280
☐	QPL110104 #2 Switch Mounting, Link-It To The Seat Rail (3)	210
☐	QPL110105 #2 Switch Mounting, Link-It To Chin Boom (4)	210
☐	QPL110106 #2 Omit Switch Mounting (5)	0

#2 Switch Mount Location		£
☐	QPL110108 #2 Switch Mounting Location, Right	0
☐	QPL110107 #2 Switch Mounting Location, Left	0

Multi-Function Switches (7)		£
☐	QPL110550 Mini 2 Button (Stereo Plug) (8)	365
☐	QPL110551 Mini 2 Button (Dual Mono Plug) (8)	365
☐	QPL110552 Mini 4 Button InLine (DB9 Connector) (8)	645
☐	QPL110554 4-Way Toggle Switched Joystick	645

Accessories		£
☐	QPL110580 Mono Extension Cable (3.5 mm Jack Socket) (2 m)	140
☐	QPL110581 Stereo Extension Cable (3.5 mm Jack Socket) (2 m)	180
☐	QPL110582 Extension Cable DB9 (1 m)	170
☐	QPL110583 Y-Cable, 2 Female Mono To 1 Male Stereo	210

Safety Devices		£
☐	QPL090147 Remote Stop Switch (30 m Activation Range, Fit Via OBC Port) (9)	980

B4M - LinX Control System		£
☐	B4M-UK-P-000086 HCM Mini Joystick with LinX control system – Mount as Chin Control – SA ARM (Programmed by customer)	4532
☐	B4M-UK-P-000089 MicroPilot - Mount as Link-it to Space Disk, Hand Control (Programmed by customer)	1495
☐	B4M-UK-P-000090 MicroGuide - Mount as Link-it to Space Disk, Hand Control (Programmed by customer)	1495
☐	B4M-UK-P-000091 VersaGuide - Mount as a Capsulated Swing Away Hand Control (Programmed by customer)	675
☐	B4M-UK-P-000092 VersaGuide EZ - Mount as a Capsulated Swing Away Hand Control (Programmed by customer)	785
☐	B4M-UK-P-000093 MovIS Allround Joystick - Mount as a Capsulated Swing Away Hand Control (Programmed by customer)	620
☐	B4M-UK-P-000094 Movis Allround Light Joystick - Mount as a Capsulated Swing Away Hand Control (Programmed by customer)	730
☐	B4M-UK-P-000095 Movis Allround HD - Mount as a Capsulated Swing Away Hand Control (Programmed by customer)	975
☐	B4M-UK-P-000096 Dual Pro, Facial Pads - Mounted a head control with Whitmyer Axys Stem (Programmed by customer)	2440
☐	B4M-UK-P-000097 Dual Pro, Spot Pads - Mounted a head control with Whitmyer Axys Stem (Programmed by customer)	2440
☐	B4M-UK-P-000098 Dual Pro, Facial Pads - Mounted a head control with Whitmyer Linx2 headrest Stem (Programmed by customer)	2490
☐	B4M-UK-P-000099 Dual Pro, Spot Pads - Mounted a head control with Whitmyer Linx2 headrest Stem (Programmed by customer)	2490
☐	B4M-UK-P-000100 Dual Pro, Facial Pads - Mounted a head control with Whitmyer Cobra headrest Stem (Programmed by customer)	2510

(1) Not With HD Allround Joystick Or Space Disk Mount

(2) Armrest Rail Mount

(3) Knee, If 2 Switches Selected Only Left And Right Mount

(4) No Left / Right Switch Selection Required

(5) Switches Attached Lose

(6) Only For Allround Joystick, Maximum 2 Selectable

(7) Used In Combination With Assignable Button Programming

(8) Attached Lose, When Ordered. No Link-It Fixation Possible

(9) Not With VR2

☐ B4M-UK-P-000101	Dual Pro, Spot Pads - Mounted a head control with Whitmyer Cobra headrest Stem (Programmed by customer)	2510
☐ B4M-UK-P-000102	MicroPilot - Mounted vai Link-It to Chin Harness, Small (Programmed by customer)	1450
☐ B4M-UK-P-000103	MicroGuide - Mounted vai Link-It to Chin Harness, Small (Programmed by customer)	1450
☐ B4M-UK-P-000104	MicroPilot - Mounted via S/A Link-IT to Headrest (Programmed by customer)	1525
☐ B4M-UK-P-000105	MicroGuide - Mounted via S/A Link-IT to Headrest (Programmed by customer)	1525
☐ B4M-UK-P-000106	MicroPilot - Mounted vai Link-It to Manual S/A Chin Boom (Programmed by customer)	1875
☐ B4M-UK-P-000107	MicroGuide - Mounted vai Link-It to Manual S/A Chin Boom (Programmed by customer)	1875
☐ B4M-UK-P-000108	VersaGuide - Mounted vai Link-It to Manual S/A Chin Boom (Programmed by customer)	1245
☐ B4M-UK-P-000109	VersaGuide EZ - Mounted vai Link-It to Manual S/A Chin Boom (Programmed by customer)	1355
☐ B4M-UK-P-000110	MovIS Allround Joystick - Mounted vai Link-It to Manual S/A Chin Boom (Programmed by customer)	1190
☐ B4M-UK-P-000111	Movis Allround Light Joystick - Mounted vai Link-It to Manual S/A Chin Boom (Programmed by customer)	1300
☐ B4M-UK-P-000112	MicroPilot - Mounted vai Link-It to Power S/A Chin Boom (Programmed by customer)	2300
☐ B4M-UK-P-000113	MicroGuide - Mounted vai Link-It to Power S/A Chin Boom (Programmed by customer)	2300
☐ B4M-UK-P-000114	VersaGuide - Mounted vai Link-It to Power S/A Chin Boom (Programmed by customer)	1670
☐ B4M-UK-P-000115	VersaGuide EZ - Mounted vai Link-It to Power S/A Chin Boom (Programmed by customer)	1780
☐ B4M-UK-P-000116	MovIS Allround Joystick - Mounted vai Link-It to Power S/A Chin Boom (Programmed by customer)	1615
☐ B4M-UK-P-000117	Movis Allround Light Joystick - Mounted vai Link-It to Power S/A Chin Boom (Programmed by customer)	1725

Comments For Personalisations



Q400 M with SEDEO LITE

£

Base Price

Total Of Additional Equipment

Packaging And Transport

Exclusive VAT

VAT

Total Price

Date

Signature

Day

Month

Year

Send E-Mail