

Q700® M/F/R with SEDEO® PRO



01/2024

ORDER FORM

Errors and omissions excepted.


Billing Address Or Customer-ID		Shipping Address		☐ Identical To Billing Address
C.-ID:		C.-ID:		
Name:		Name:		
Street:		Street:		
Town/PC:		Town/PC:		
E-Mail:		E-Mail:		
Phone:		Phone:		
Order Information		User Weight		
☐ Order ☐ Calculation ☐ Built-4-Me®		User Weight Capacity: 160 kg		
Purchase Order No.:		User Weight: kg		
Art. No.:		Reimbursement Code: xxx		
Order Date:	Day	Month	Year	

Base Price

£7,950

Client Measurements

a) Top Of Head To Seat		mm
b) Top Of Shoulder To Seat		mm
c) Underarm To Seat		mm
d) Elbow To Seat		mm
e) Top Of Pelvis To Seat		mm
f) Shoulder Width		mm
g) Trunk Width		mm
i) Hip Width		mm
j) Actual Sitting Depth		mm
m) Foot Length		mm
n) Lower Leg Length		mm
o) Actual Seat Width		mm



Seating Specifications

All Bases Prepared For DAHL Docking Station [Q700 M Max. 160kg User Weight / Q700 F R Max. 136kg User Weight]

Typical WC Weight	Without Batteries 140 kg	Seat Type	SEDEO PRO
	With Batteries 165 kg	Seat Size Range	Width 380 - 560 mm, Depth 400 - 560 mm
		Seat To Floor Height Range	415 mm - 490 mm

MWD Base Specifications

User Weight Capacity	160 kg	Drive Wheels	14" Drive Wheels (Black)
Base Size	Width Min. 620 mm	Tires / Castors	6" Front / 6" Rear (Black)
Drive Control	VR2, Option R-Net	Charger	8 A Dual Mode Charger
Driving Range	Up To 40 km	Battery Size	60 And 80 Ah Batteries
Drive Motors	4-Pole 6 kph, Options 8, 10 And 12.5 kph	Colour Accents	7 Colour Accents
Ground Clearance	70 mm	Frame Colour	Black

FWD Base Specifications

User Weight Capacity	160 kg	Drive Wheels	14" Drive Wheels (Black)
Base Size	Width Min. 640 mm	Tires / Castors	9" / 10" Rear (Black)
Drive Control	VR2, Option R-Net	Charger	8 A Dual Mode Charger
Driving Range	Up To 40 km	Battery Size	60 And 80 Ah Batteries
Drive Motors	4-Pole 6 kph, Options 8, 10 And 12.5 kph	Colour Accents	7 Colour Accents
Ground Clearance	75 mm	Frame Colour	Black

For Optimum Comfort & Driving Experience Of Your FWD Powerchair, Recommended Seat Depth Relationship To User Weight Is:

75 kg User Weight: SD40 - SD 44 cm; 90 kg User Weight: SD46 - SD50 cm; 160 kg User Weight: SD50 - SD56 cm.

Where Wheelchair Configuration Allows, A Recommended Backrest Angle Of 95° And A Seat Tilt Of 5° Will Further Improve The Driving Experience.

RWD Base Specifications

User Weight Capacity	160 kg	Drive Wheels	14" Drive Wheels (black)
Base Size	Width Min. 640 mm	Tires / Castors	9" / 10" Front (Black)
Drive Control	VR2, Option R-Net	Charger	8 A Dual Mode Charger
Driving Range	Up To 40 km	Battery Size	60 And 80 Ah Batteries
Drive Motors	4-Pole 6 Kph, Options 8, 10 And 12.5 kph	Colour Accents	7 Colour Accents
Ground Clearance	75 mm	Frame Colour	Black

Base

Q700 Drive Bases

£

☐	SEP000009	Q700 M (MWD Base Chair)	7950
☐	SEP000017	Q700 F (FWD Base Chair)	9450
☐	SEP000018	Q700 R (MWD Base Chair)	9450

Motor Packages

£

☒	SEP180009	4-Pole Motor, 6 kph / 4 mph (8 kph Motor With Reduced Speed) Motors optimised only for use with QEP110003 Power Module, R-Net, 90A	0
☐	SEP180010	4-Pole Motor, 8 kph / 5 mph Motors optimised only for use with QEP110003 Power Module, R-Net, 90A	0
☐	SEP180007	4-Pole Motor, 10 kph / 6.5 mph If FWD Selected: Highly Recommended With Traction Control System Motors optimised only for use with QEP110004 Power Module, R-Net, 120A	390
☐	SEP180008	4-Pole Motor, 12.5 kph / 8 mph If FWD Selected: Only In Combination With Traction Control System Motors optimised only for use with QEP110004 Power Module, R-Net, 120A	500

Tracking

£

☒	SEP110009	Traction Control System For R-Net (GYRO Technology)	400
☐	SEP110014	Omit Traction Control System For R-Net	0

Drive Wheels

£

☐	SEP070034	Drive Wheel, 14", Pneumatic, Knobbly, Black	0
☐	SEP070033	Drive Wheel, 14", Solid, Knobbly, Black	0

Castor Wheels

£

MWD Base

☐	SEP080206	Castor Wheel, 6", Solid, Black	0
--------------------------	-----------	--------------------------------	---

FWD Base

☐	SEP080210	Castor Wheel, 9", Solid, Black	0
☐	SEP080211	Castor Wheel, 9", Pneumatic, Black	0

RWD Base

☐	SEP080210	Castor Wheel, 9", Solid, Black	0
☐	SEP080211	Castor Wheel, 9", Pneumatic, Black	0
☐	SEP080212	Castor Wheel, 10", Pneumatic, Black	0

Mudguards (FWD & RWD Only)

£

☐	SEP010007	Mudguard, Drive Wheel	125
☐	SEP010008	Mudguard, Castor Wheel, 9"	100

Kerb Climber (RWD Only)

£

☐	SEP090038	Kerb Climber	200
--------------------------	-----------	--------------	-----

Shroud Colours Q700 M

£

☐	SEP100577	Black (RAL 9005) Black Natural Non-Painted finish	0
--------------------------	-----------	--	---

Base

Colour Accents

£

☐	SEP140005	Apple Green (RAL 110 70 77)	0
☐	SEP140004	Black (RAL 9005)	0
☐	SEP140003	Blue (RAL 260 40 45)	0
☐	SEP140006	Orange (RAL 060 60 80)	0
☐	SEP140002	Red (RAL 040 40 67)	0
☒	SEP140007	Silver (RAL 9006)	0
☐	SEP140001	White (RAL 9010)	0

LED Lights & Indicators

£

☐	SEP110201	LED Light And Indicator SEP110030 Lighting Module VR2 Or SEP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required	400
☒	SEP110209	Omit Light And Indicator	0
☐	SEP110202	LED Fender Light SEP110031 Lighting Module, R-Net LM Required, Fender Light - Not Road Approved and Without Indicator	150

Seat Heights Measured From Front Of Seat Pan To Floor. Obtain Quote To Verify Seat To Floor Heights

£

☐	SH415	415 mm Only With 60 Ah Batteries; Not With SEP130009 Seat Elevator/Tilt Module (Tilt 0° - 30°), Powered or SEP130006 Seat Elevator/Tilt Module (Tilt 0° - 50°), Powered	0
☐	SH425	425 mm Only With 60 Ah Batteries; Only With SEP130009 Seat Elevator/Tilt Module (Tilt 0° - 30°), Powered or SEP130006 Seat Elevator/Tilt Module (Tilt 0° - 50°), Powered	0
☒	SH44	440 mm Not With SEP130009 Seat Elevator/Tilt Module (Tilt 0° - 30°), Powered or SEP130006 Seat Elevator/Tilt Module (Tilt 0° - 50°), Powered	0
☐	SH45	450 mm Only With SEP130009 Seat Elevator/Tilt Module (Tilt 0° - 30°), Powered or SEP130006 Seat Elevator/Tilt Module (Tilt 0° - 50°), Powered	0
☐	SH465	465 mm	0
☐	SH49	490 mm	0

Batteries & Chargers

£

Batteries

☒	SEP160001	Battery, 60 Ah, Gel, Maintenance-Free	0
☐	SEP160005	Battery, 80 Ah, Gel, Maintenance-Free Only With Battery Charger 220 V / 10 A To Reduce Reloading Time	160
☐	SEP160009	Omit Batteries	0

Chargers

☐	SEP160100	Charger, Automatic, 220 V / 8 A, UK	0
☐	SEP160103	Charger, Automatic, 220 V / 10 A, UK	230
☐	SEP160109	Omit Charger	0

Powered Seating

Power Seating Options

£

☒	SEP130011	Fixed Tilt Module (Adjustable 0°, 3°, 6°, 9°, Factory Setting 6°)	0
☐	SEP130004	Posterior Tilt, Powered, 0° - 30° Not With SEP110036 Remote Joystick, VR2. In Case Of SEP110004 Power Module, R-Net: SEP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required	1150
☐	SEP130012	Posterior Tilt, Powered, -5° - 25° Not With SEP110036 Remote Joystick, VR2. In Case Of SEP110004 Power Module, R-Net: SEP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required; Not With Seat Height 415 Or 440 mm	1250
☐	SEP130007	Posterior Tilt, Powered, 0° - 50° Not With SEP110036 Remote Joystick, VR2. In Case Of SEP110004 Power Module, R-Net: SEP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required	1350
☐	SEP130005	Posterior Tilt, Powered, -5° - 45° Not With SEP110036 Remote Joystick, VR2. In Case Of SEP110004 Power Module, R-Net: SEP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required; Not With Seat Height 415 Or 440 mm	1450
☐	SEP130009	Seat Elevator (300 mm) / Tilt Module (Tilt 0° - 30°), Powered Not With SEP110036 Remote Joystick, VR2. In Case Of SEP110004 Power Module, R-Net: SEP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required	2050
☐	SEP130006	Seat Elevator (300 mm) / Tilt Module (Tilt 0° - 50°), Powered Not With SEP110036 Remote Joystick, VR2; In Case Of SEP110004 Power Module, R-Net, 120 A: SEP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required	2250

C-Me Function

£

☒	SEP110034	C-Me Function, Enabled 4.5kph Maximum Speed When Lifted (Tilt Angle <16° And Recline Angle <40°)	0
☐	SEP110035	C-Me Function, Disabled 1.5kph Maximum Speed When Lifted	0

Backrest Frame Type

£

☒	SEP030083	Backrest, Tool Angle Adjustable, 90° - 120°	0
☐	SEP030084	Backrest, Manual Recline, User Adjustable, 90° - 130°, Release Lever, Right Not On Same Joystick Side	310
☐	SEP030085	Backrest, Manual Recline, User Adjustable, 90° - 130°, Release Lever, Left Not On Same Joystick Side	310
☐	SEP030086	Backrest, Powered Recline, 90° - 130° Not With SEP110036 Remote Joystick, VR2. In Case Of SEP110004 Power Module, R-Net: SEP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required	560
☐	SEP030087	Backrest, Powered Recline, 90° - 150° With 100 mm Mechanical Anti-Shear Not With SEP110036 Remote Joystick, VR2. In Case Of SEP110004 Power Module, R-Net: SEP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required	600
☐	SEP030229	Backrest, Powered Recline, 90 - 150 With 100 mm Mechanical Anti-Shearwith Relax Function Not With SEP110036 Remote Joystick, VR2; In Case Of SEPP110004 Power Module, R-Net, 120 A: SEPP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required; Only with SEP130006 Lift Tilt 0-50°	600

Backrest Frame Option

£

☐	SEP030098	Backrest, Flip Down, Quick Release	185
--------------------------	-----------	------------------------------------	-----

Back Cushion Type

£

SEDEO PRO (See Section A)

☒	SEP030104	SEDEO PRO, Back Cushion	0
-------------------------------------	-----------	-------------------------	---

JAY J3 (See Section B)

☐	SEP030054	JAY J3, Back Cushion	0
--------------------------	-----------	----------------------	---

Interfaces

☐	SEP030106	SEDEO PRO, Interface Only	0
☐	SEP030107	SEDEO PRO, Interface Only For Back Cushion, Flexible	0
☐	SEP030055	JAY J3, Interface Only	0
☐	SEP030209	Omit Back Cushion	0

A – SEDEO PRO Back Size

£

Back Cushion Width

☐	BW38	380 mm, For Effective Back Width 380 mm	0
☐	BW42	420 mm, For Effective Back Width 380 - 420 mm	0
☐	BW46	460 mm, For Effective Back Width 420 - 460 mm	0
☐	BW50	500 mm, For Effective Back Width 460 - 500 mm	0
☐	BW54	540 mm, For Effective Back Width 500 - 540 mm	0
☐	BW56	560 mm, For Effective Back Width 520 - 560 mm	0

Back Cushion Height

☐	BH44	440 mm, For Effective Back Height 440 - 500 mm Only With SEP030116 SEDEO PRO, Back Cushion, Flexible, Deep Contour, Active, Adjustable Straps & SEP030117 SEDEO PRO, Back Cushion, Flexible, Extra Deep Contour, Active, Adjustable Straps	0
☐	BH50	500 mm, For Effective Back Height 500 - 560 mm	0
☐	BH56	560 mm, For Effective Back Height 560 - 620 mm	0

A – SEDEO PRO Back Cushion

£

☐	SEP030109	SEDEO PRO, Back Cushion, Profiled	140
☐	SEP030113	SEDEO PRO, Back Cushion, Profiled, Adjustable Straps	260
☐	SEP030114	SEDEO PRO, Back Cushion, Deep Contour, Passive	335
☐	SEP030115	SEDEO PRO, Back Cushion, Flexible, Deep Contour, Passive, Adjustable Straps	435
☐	SEP030116	SEDEO PRO, Back Cushion, Flexible, Deep Contour, Active, Adjustable Straps	450
☐	SEP030117	SEDEO PRO, Back Cushion, Flexible, Extra Deep Contour, Active, Adjustable Straps	475

Powered Seating**A – SEDEO PRO Back Cushion Cover**

£

☒	SEP030059	Back Cushion Cover, Incontinence Protection	0
☐	SEP030108	Back Cushion Cover, Comfort	0

A – SEDEO Lateral Thoracic Supports Only With SEP030109 SEDEO PRO, Back Cushion, Profiled & SEP030113 SEDEO PRO, Back Cushion, Profiled, Adjustable Straps

£

☐	SEP030607	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Right	160
☐	SEP030608	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Left	160
☐	SEP030609	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Short Pad, 90 x 120 mm, Right	160
☐	SEP030610	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Short Pad, 90 x 120 mm, Left	160

B – JAY J3 Back Size

£

☐	JAY030111	JAY J3 Back Size, 410 mm	0
☐	JAY030112	JAY J3 Back Size, 460 mm	0
☐	JAY030113	JAY J3 Back Size, 510 mm	0

B – JAY J3 Back, Mid Contour

£

☐	JAY030013	JAY J3 Back, Mid Contour, Upper Thoracic, Short, 420 mm	870
--------------------------	-----------	---	-----

B – JAY J3 Back, Mid Deep Contour

£

☐	JAY030045	JAY J3 Back, Mid Deep Contour, Upper Thoracic, Short, 420 mm	870
☐	JAY030042	JAY J3 Back, Mid Deep Contour, Upper Thoracic, Medium, 460 mm	870
☐	JAY030046	JAY J3 Back, Mid Deep Contour, Upper Thoracic, Tall, 500 mm	870
☐	JAY030043	JAY J3 Back, Mid Deep Contour, Shoulder Height, Short, 530 mm	870
☐	JAY030047	JAY J3 Back, Mid Deep Contour, Shoulder Height, Medium, 570 mm	870
☐	JAY030048	JAY J3 Back, Mid Deep Contour, Shoulder Height, Tall, 610 mm	870
☐	JAY030050	JAY J3 Back, Mid Deep Contour, Extra Deep, Upper Thoracic, Medium, 460 mm	870

Powered Seating

Headrest Pad

£

☐	SEP030430	SEDEO, Headrest Pad, Small, (LxH) 240 x 150 mm Only Use With SEDEO Headrest Mount	95
☐	SEP030431	SEDEO, Headrest Pad, Small, With Cover, (LxH) 240 x 150 mm Only Use With SEDEO Headrest Mount	110
☐	SEP030432	SEDEO, Headrest Pad, Flexible, (LxH) 400 x 110 mm Only Use With SEDEO Headrest Mount	170
☐	WHT030401	Whitmyer, Headrest Pad, Plush, (LxH) 203 x 127 mm Only Use With Whitmyer Headrest Mount	170
☐	WHT030402	Whitmyer, Headrest Pad, Plush, (LxH) 254 x 127 mm Only Use With Whitmyer Headrest Mount	170
☐	WHT030403	Whitmyer, Headrest Pad, Plush, (LxH) 356 x 127 mm Only Use With Whitmyer Headrest Mount	170
☐	WHT030404	Whitmyer, Headrest Pad, Adjust-A-Plush, Small, (LxH) 355 x 76 mm Only Use With Whitmyer Headrest Mount	280
☐	WHT030406	Whitmyer, Headrest Pad, Adjust-A-Plush, Large, (LxH) 460 mm x 127 mm Only Use With Whitmyer Headrest Mount	280

Headrest Mount Choose A Whitmyer Headrest Mount And Headrest Receiver In Case You Would Like To Add A Head Control In The Special Control Section)

£

☐	SEP030433	SEDEO, Headrest Mount	140
☐	SEP030434	SEDEO, Headrest Mount, With Quick Release	155
☐	WHT030407	Whitmyer, Headrest Mount, LINX	190
☐	WHT030408	Whitmyer, Headrest Mount, LINX²	190
☐	WHT030409	Whitmyer, Headrest Mount, AXYS	155
☐	WHT030410	Whitmyer, Headrest Mount, AXYS, Offset	155
☐	WHT030411	Whitmyer, Headrest Mount, COBRA XTRA Flip Back	210

Headrest Receiver

£

☐	SEP030119	Headrest Receiver, Back Cushion Mount Only With SEP030109 SEDEO PRO, Back Cushion, Profiled & SEP030113 SEDEO PRO, Back Cushion, Profiled, Adjustable Straps & SEP030114 SEDEO PRO, Back Cushion, Deep Contour, Passive	0
☐	SEP030120	Headrest Receiver, Backrest Frame Mount	0

Push Handles

£

☐	SEP030204	Push Handles, Height-Adjustable	210
--------------------------	-----------	---------------------------------	-----

Seat Rail

£

☒	SEP020085	Seat Rail, Short, For Seat Depth 400 - 560 mm	0
☐	SEP020069	Seat Rail, Long, For Seat Depth 460 - 560 mm Not With Seat Width 380 - 420 mm	0

Seat Cushion Type

£

SEDEO (See Section C)

☒	SEP020027	SEDEO, Seat Cushion	0
-------------------------------------	-----------	---------------------	---

JAY (See Section D)

☐	SEP020515	JAY, Seat Cushion	0
--------------------------	-----------	-------------------	---

Omit (Select Seat Width And Depth in Section C)

☐	SEP020110	Omit Seat Cushion	0
--------------------------	-----------	-------------------	---

Powered Seating**C – SEDEO Seat Size**

£

Seat Cushion Width

☐	SW38	380 mm, For Effective Seat Width 380 mm	0
☐	SW42	420 mm, For Effective Seat Width 380 - 420 mm	0
☐	SW46	460 mm, For Effective Seat Width 420 - 460 mm	0
☐	SW50	500 mm, For Effective Seat Width 460 - 500 mm	0
☐	SW54	540 mm, For Effective Seat Width 500 - 540 mm	0
☐	SW56	560 mm, For Effective Seat Width 520 - 560 mm	0

Seat Cushion Depth

☐	SD40	400 mm, Comes With Cushion Size 460 mm Not with SEP090128 or SEP090129 Positioning belt Short or Long	0
☐	SD42	420 mm, Comes With Cushion Size 460 mm	0
☐	SD44	440 mm, Comes With Cushion Size 460 mm	0
☐	SD46	460 mm, Comes With Cushion Size 500 mm	0
☐	SD48	480 mm, Comes With Cushion Size 500 mm	0
☐	SD50	500 mm, Comes With Cushion Size 500 mm	0
☐	SD52	520 mm, Comes With Cushion Size 560 mm	0
☐	SD54	540 mm, Comes With Cushion Size 560 mm	0
☐	SD56	560 mm, Comes With Cushion Size 560 mm	0

**C– SEDEO Seat Cushions **

£

☒	SEP020073	SEDEO, Seat Cushion, Profiled, Height 70 mm	160
☐	SEP020075	SEDEO, Seat Cushion, Support, Height 90 mm	225
☐	SEP020076	SEDEO, Seat Cushion, Support Memory, Height 90 mm	260
☐	SEP020078	SEDEO, Seat Cushion, Stabilizer, Height 90 mm	280

C – SEDEO Seat Cushion Cover

£

☒	SEP020518	Seat Cushion Cover, Incontinence Protection	0
☐	SEP020080	Seat Cushion Cover, Comfort	0

C – SEDEO Seat Boards

£

☐	SEP020082	SEDEO, Seat Board, Flat, Depth-Adjustable, Cushions Guide At Rear	130
☐	SEP020083	SEDEO, Seat Pan, 50 mm, Depth-Adjustable, Cushions Guide At Rear	150

Powered Seating

D – JAY Seat Size

£

Seat Cushion Width

☐	SW40	400 mm	0
☐	SW42	420 mm	0
☐	SW44	440 mm	0
☐	SW46	460 mm	0
☐	SW48	480 mm	0
☐	SW50	500 mm, Comes With Cushion Size 500 mm	0
☐	SW52	520 mm, Comes With Cushion Size 500 mm	0
☐	SW54	540 mm, Comes With Cushion Size 560 mm	0
☐	SW56	560 mm, Comes With Cushion Size 560 mm	0

Seat Cushion Depth

☐	SD40	400 mm	0
☐	SD42	420 mm	0
☐	SD44	440 mm	0
☐	SD46	460 mm	0
☐	SD48	480 mm	0
☐	SD50	500 mm	0
☐	SD52	520 mm	0
☐	SD54	540 mm Only With Seat Width 560 mm	0
☐	SD56	560 mm Only With Seat Width 560 mm	0

D – JAY Seat Cushion

£

☐	JAY020011	JAY Balance, Fluid, Microclimatic Cover, Until Seat Width 500 mm	655
☐	JAY020013	JAY Balance, Fluid, Microclimatic Cover, From Seat Width 560 mm	703
☐	JAY020015	JAY Balance, Fluid, Microclimatic Cover, With Positioning Kit, Until Seat Width 500 mm	755
☐	JAY020017	JAY Balance, Fluid, Microclimatic Cover, With Positioning Kit, From Seat Width 560 mm	803

Knee & Thigh Supports

£

☐	SEP020342	Knee Support, Lateral, 120 x 120 mm, Swing Away, Right Not With SEP050015-17 SEDEO PRO, Legrest, Swing Away, SEP050080 SEDEO PRO, Hanger Receivers Only, SEP090122 Transfer Bar	185
☐	SEP020341	Knee Support, Lateral, 120 x 120 mm, Swing Away, Left Not With SEP050015-17 SEDEO PRO, Legrest, Swing Away, SEP050080 SEDEO PRO, Hanger Receivers Only, SEP090122 Transfer Bar	185
☐	SEP020113	Thigh Support, Lateral, Small, 120 x 120 mm, Removable, Right Not With Seat To Floor Height 415 mm, Not With SEP040035 Side Guard, Wedge Shaped, Padded	185
☐	SEP020114	Thigh Support, Lateral, Small, 120 x 120 mm, Removable, Left Not With Seat To Floor Height 415 mm, Not With SEP040035 Side Guard, Wedge Shaped, Padded	185
☐	SEP020121	Thigh Support, Lateral, Large, 180 x 120 mm, Removable, Right Not With Seat To Floor Height 415 mm, Not With SEP040035 Side Guard, Wedge Shaped, Padded	185
☐	SEP020122	Thigh Support, Lateral, Large, 180 x 120 mm, Removable, Left Not With Seat To Floor Height 415 mm, Not With SEP040035 Side Guard, Wedge Shaped, Padded	185

Powered Seating

Armrest Types

£

☐	SEP040021	Armrest, Single Post, Removable Not with Seat height 415/Seat height 425 (15)	0
☒	SEP040019	Armrest, Single Post, Flip Back Not with Seat height 415/Seat height 425 (15)	250
☐	SEP040123	Armrest, Reclining, Flip Back	240
☐	SEP040070	Armrest, Reclining, Multi-Adjustable, Flip Back	320

Arm Pads

£

☐	SEP040023	Arm Pad, Integral Foam, 300 mm, Right	0
☐	SEP040024	Arm Pad, Integral Foam, 300 mm, Left	0
☒	SEP040025	Arm Pad, Integral Foam, 400 mm, Right	0
☒	SEP040026	Arm Pad, Integral Foam, 400 mm, Left	0
☐	SEP040027	Arm Pad, Integral Foam, 450 mm, Right With SEP040019 Flip Back Armrest, Seat Width +20 mm	0
☐	SEP040028	Arm Pad, Integral Foam, 450 mm, Left With SEP040019 Flip Back Armrest, Seat Width +20 mm	0
☐	SEP040029	Omit Arm Pad, Right	0
☐	SEP040030	Omit Arm Pad, Left	0

Comfort Arm Pads

☐	SEP040204	Arm Pad, Comfort, Padded, PU Coated, 400 mm, Right	75
☐	SEP040203	Arm Pad, Comfort, Padded, PU Coated, 400 mm, Left	75

Hemiplegia Arm Pads

☐	SEP040031	Arm Pad, Hemiplegia, Integral Foam, Right Not With SEP110052 Joystick Mount, Right, SEP090040 Tray Table, Swing Away, SEP090046 Tray Table, Half-Side, Swing Away, SEP090110 Rear View Mirror, Right; not with SEP110362 Hand Control Mounting Location Right; not with SEP110119/SEP110123 OMNI2 Module Mount Right	275
☐	SEP040032	Arm Pad, Hemiplegia, Integral Foam, Left Not With SEP110053 Joystick Mount, Right, SEP090040 Tray Table, Swing Away, SEP090046 Tray Table, Half-Side, Swing Away, SEP090118 Rear View Mirror, Left; not with SEP110361 Hand Control Mounting Location Left; not with SEP110120/SEP110124 OMNI2 Module Mount Left	275

Side Guard Pads For Armrests

£

Side Guard Pads For Armrest, Single Post, Flip Back / Removable

☐	SEP040033	Side Guard Not With SEP040123 Armrest, Reclining, Flip Back	0
☐	SEP040034	Side Guard, Integral Foam Seat Width -20 mm Not With SEP040123 Armrest, Reclining, Flip Back	0
☐	SEP040035	Side Guard, Wedge Shaped, Padded Seat Width -80 mm, Only With SEP040021 Armrest, Single Post, Removable	275

Side Guard Pads For Armrest, Reclining, Flip Back

☐	SEP040036	Side Guard, Mounted On Seat Rail Bracket Only With SEP040123 Armrest, Reclining, Flip Back	0
--------------------------	-----------	---	---

Centre Mount Legrest

£

☐	SEP050086	Centre Mount Legrest, Tool Adjustable, With Flip Up And Angle-Adjustable Footboard In Combination With Lift Or Tilt, Choose SEP050087 Centre Calf Support, Integral Foam Or SEP050092 Calf Pad, Integral Foam, Pair; Not With RWD	0
--------------------------	-----------	--	---

Pad Options

☐	SEP050087	Centre Calf Support, Integral Foam	0
☐	SEP050092	Calf Pad, Integral Foam, Pair	0

Fixed Centre Mount Legrest 			£
☐	SEP050140	Fixed Centre Mount Legrest Not With RWD; In Combination With Lift Or Tilt, Choose SEP050092 Calf Pad, Integral Foam, Pair	0
Footplates			
☐	SEP050072	Footplate, Dual, Height- And Angle-Adjustable, 260 x 150 mm	0
Pad Options			
☐	SEP050092	Calf Pad, Integral Foam, Pair	0
Power Centre Mount Legrest 			£
☐	SEP050076	Centre Mount Legrest, Powered, 90° - 180°, 175 mm Extension Side Guard Included; Only With SEP110082 Joystick Mounting, Parallel, Swing Away; Arm Pad, Integral Foam, 380 mm Right/Left Included	1470
☐	SEP050071	Centre Mount Legrest, Powered, 90° - 180°, 175 mm Extension, Footplate To Floor Side Guard Included; Only With SEP110082 Joystick Mounting, Parallel, Swing Away; Not With SEP090040 Tray Table, Width-Adjustable, Swing Away & SEP090046 Tray Table, Half-Side, Swing Away; Not With Special Controls Mounted On Armrest, Not With SEP090120 Bag For Armrest; Arm Pad, Integral Foam, 380 mm Right/Left Included	1470
Footplates			
☐	SEP050072	Footplate, Dual, Height- And Angle-Adjustable, 260 x 150 mm	0
☐	SEP050073	Footplate, Dual, Height- And Angle-Adjustable, Large, 310 x 170 mm	0
Pad Options			
☐	SEP050092	Calf Pad, Integral Foam, Pair	0
SEDEO PRO Swing Away Legrest 			£
☒	SEP050015	SEDEO PRO, Legrest, 70°, Swing Away	0
☐	SEP050016	SEDEO PRO, Legrest, Manual Elevating, Swing Away	375
☐	SEP050017	SEDEO PRO, Legrest, Power Elevating, Swing Away	700
☐	SEP050183	Sedeo Pro, Legrest, Power Elevating, Swing Away, Short, LLL 300mm	700
☐	SEP050000	Omit Legrest	0
Hanger Receivers			
☐	SEP050170	SEDEO PRO, Hanger Receivers, High, For Seat Cushion Thickness 70 mm And Higher	0
☐	SEP050171	SEDEO PRO, Hanger Receivers, Low, For Seat Cushion Thickness 50 mm And Lower	0
☐	SEP050080	SEDEO PRO, Hanger Receivers Only, High	0
☐	SEP050172	SEDEO PRO, Hanger Receivers Only, Low	0
Footplates			
☐	SEP050088	Footplate, For Seat Width 385 - 400 mm	0
☒	SEP050089	Footplate, For Seat Width 420 - 460 mm	0
☐	SEP050090	Footplate, For Seat Width 480 - 520 mm	0
☐	SEP050091	Footplate, For Seat Width 540 - 560 mm	0
☐	SEP050081	Footplates, Coupled Choose Also A Footplate; Only For SEP050015 SEDEO PRO, Legrest, 70°, Swing Away	0
Support Options 			
☐	SEP050082	Calf Support, Integral Foam, Depth-Adjustable, Pair Not On RWD Base With SEP050015 SEDEO PRO, Legrest, 70°, Swing Away	100
☐	SEP050083	Calf Support, Padded, Depth-Adjustable, Pair Not On RWD Base With SEP050015 SEDEO PRO, Legrest, 70°, Swing Away	135
☐	SEP050084	Heel Loop	0
☒	SEP050085	Calf Strap Only For SEP050015 SEDEO PRO, Legrest, 70°, Swing Away	0
☐	SEP050200	Knee Pad, Integral Foam	75
Amputee Support			
☐	QOP050142	Amputee Support, Integral Foam, Right	145
☐	QOP050143	Amputee Support, Integral Foam, Left	145
☐	QOP050013	Amputee Support, Padded, Right Extra On A Set Of Swing Away Legrests, Attached Lose	145
☐	QOP050014	Amputee Support, Padded, Left Extra On A Set Of Swing Away Legrests, Attached Lose	145

Positioning Supports & Accessories		£
☐	SEP090128 Positioning Belt, Carseat Buckle, Short, 50 mm x 1700 mm	0
☐	SEP090129 Positioning Belt, Carseat Buckle, Long, 50 mm x 2700 mm	0
☐	SEP090110 Rear View Mirror, Right Not In Combination With Same Entry Side SEP090040 Tray Table, Swing Away, SEP090046 Tray Table, Half-Side, Swing Away	60
☐	SEP090118 Rear View Mirror, Left Not In Combination With Same Entry Side SEP090040 Tray Table, Swing Away, SEP090046 Tray Table, Half-Side, Swing Away	60
☐	SEP090009 Crutch Holder Not With SEP030086 Powered Recline, 90° - 130°, SEP030087 Powered Recline, 90° - 150° With 100 mm Mechanical Anti-Shear	70
☐	SEP090119 Brackets For Bag On Backrest Frame Not With SEP030086 Powered Recline, 90° - 130°, SEP030087 Powered Recline, 90° - 150° With 100 mm Mechanical Anti-Shear	35
☐	SEP090120 Bag For Armrest Not With SEP040123 Armrest, Reclining, Flip Back	105
☐	SEP090122 Transfer Bar Not With SEP020342/SEP020341 Knee Support, Lateral, Swing Away, Right/Left & SEP050015-17 SEDEO PRO, Legrest, Swing Away & SEP050080 SEDEO PRO, Hanger Receivers Only	255
☐	SEP030605 Abduction Pommel Not With (Power) Centre Mount Legrests	255

Tray Table (Mounted Opposite Joystick)

☐	SEP090040 Tray Table, Width-Adjustable, Swing Away	335
☐	SEP090046 Tray Table, Half-Side, Swing Away Not With SEP090125-127 Tray Table Accessories	265

Tray Table Accessories

☐	SEP090125 Partial Padding For Tray Table	25
☐	SEP090126 Raised Edges, Acrylic For Tray Table	25
☐	SEP090127 Elbow Pads, Padded, Pair For Tray Table	25

Joystick / Electronics

Remote Joystick – VR2

☒	SEP110036 Remote Joystick, VR2	0
☐	SEP110037 Remote Joystick, VR2, 2 Actuators	0
☐	SEP110040 Remote Joystick, VR2, 2 Actuators And Lights	0

Power Module – VR2

☒	SEP110002 Power Module, VR2, 90 A	0
-------------------------------------	-----------------------------------	---

Joystick Mounting

☒	SEP110050 Joystick Mounting, Fixed	0
☐	SEP110056 Joystick Mounting, Parallel, Swing Away, Capsulated	80
☐	SEP110082 Joystick Mounting, Parallel, Swing Away Only With SEP040022 Armrest, Single Post, SEDEO LITE	0

Joystick Mounting Location

☒	SEP110052 Joystick Mounting Location, Right	0
☐	SEP110053 Joystick Mounting Location, Left	0

Joystick Knobs (Optional) 			£
☐	SEP090505	Joystick Knob, Golf Ball, Ø 42 mm	40
☐	SEP090506	Joystick Knob, T-Bar, 82 mm Wide	40
☐	SEP090507	Joystick Knob, Soft Ball, Ø 70 mm	40
☐	SEP090508	Joystick Knob, Ball, Ø 35 mm	40
☐	SEP090045	Joystick Knob, Allround, Ball, Ø 48 mm	40
☐	SEP090509	Joystick Knob, Chin Cup, Ø 37 mm	40
☐	SEP090512	Joystick Knob, Stick, 80 x 15 mm	40
☐	SEP090511	Joystick Knob, Mushroom, Ø 45 mm	40

Lighting Module – VR2			£
☐	SEP110030	Lighting Module, VR2 If LED Light And Indicator Is Chosen, Then This Is Required	0

Attendant Controls			£
☐	SEP110021	Attendant Control, Dual	425
☐	SEP110020	Joystick Mounted As Attendant Control Choose Also Remote Joystick	270
Attendant Mounting Location			
☐	SEP110025	Attendant Mounting Location, Right	0
☐	SEP110026	Attendant Mounting Location, Left	0

Joystick / Electronics

Remote Joystick – R-Net 			£
☐	SEP110015	Remote Joystick, R-Net, LED, Actuators	485
☐	SEP110016	Remote Joystick, R-Net, LED, Actuators And Lights	485
☐	SEP110017	Remote Joystick, R-Net, Colour LCD	600
☐	SEP110018	Remote Joystick, R-Net, Colour LCD, Heavy Duty	1200
☐	SEP110013	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared Control Included	770
☐	SEP110023	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared And Bluetooth Control Included	925
☐	SEP110019	Omit Remote Joystick, R-Net	0

Power Module – R-Net 			£
☐	SEP110003	Power Module, R-Net, 90 A	0
☐	SEP110004	Power Module, R-Net, 120 A	135

Joystick Mounting 			£
☐	SEP110050	Joystick Mounting, Fixed	0
☐	SEP110056	Joystick Mounting, Parallel, Swing Away, Capsulated	80
☐	SEP110082	Joystick Mounting, Parallel, Swing Away Only With SEP040022 Armrest, Single Post, SEDEO LITE	0

Joystick Mounting Location			£
☐	SEP110052	Joystick Mounting Location, Right	0
☐	SEP110053	Joystick Mounting Location, Left	0

Seating / Lighting Module Control – R-Net 			£
☐	SEP110032	Actuator Control & Lights Module, R-Net, ISM-6L	0

Joysticks / Electronics

Attendant Controls

£

☐	SEP110021	Attendant Control, Dual	425
☐	SEP110020	Joystick Mounted As Attendant Control Choose OMNI2 Mount And Input Device in Sections Below And Choose Switches in Section E	270

Attendant Mounting Location

☐	SEP110025	Attendant Mounting Location, Right	0
☐	SEP110026	Attendant Mounting Location, Left	0

Table Top Controls

£

☐	SEP110143	Table Top Control, LCD, Centre Tray Mounted	2235
--------------------------	-----------	---	------

Tray Table Mounting Location

☐	SEP110080	Tray Table Mounting Location, Left	0
☐	SEP110081	Tray Table Mounting Location, Right	0

#1 Ctrl+5 Switch Box Options

£

☐	SEP110043	#1 Ctrl+5 Switch Box With Free Assignable Functions If combined with attendant control, Ctrl/+5 will be mounted to the armrest	865
--------------------------	-----------	---	-----

#1 Ctrl+5 Switch Box Mounting Location

☐	SEP110180	#1 Ctrl+5 Switch Box Mounting Location, Behind Joystick Not Possible With SEP110056 Joystick Mounting, Parallel Swing Away, Capsulated (Here We Can't Swing Away)	0
☐	SEP110181	#1 Ctrl+5 Switch Box Mounting Location, Opposite Joystick	0

#1 Ctrl+5 Switch Type

☐	SEP110175	#1 Ctrl+5 Switch Type With Push Buttons	0
☐	SEP110176	#1 Ctrl+5 Switch Type With Toggle Switches	0

Special Controls

OMNI2 Input Control Module – R-Net Choose OMNI2 Mount And Input Device in Sections Below And Choose Switches in Section E

£

☐	SEP110116	OMNI2 Input Control Module, R-Net, Bluetooth, Infrared Remote Full VGA Display With Infrared, Separate Input Box With Bluetooth, USB Charger And Speaker	1170
☐	SEP110117	OMNI2 Input Control Module, R-Net, Scanning Mode, Bluetooth, Infrared Remote Full VGA Display With Infrared, Separate Input Box With Bluetooth, USB Charger And Speaker	1170

OMNI2 Module Mount If Used With Standard Joystick Selection, It Must Be Mounted Opposite Side

£

☐	SEP110119	OMNI2 Module, Fixed Mount, Right	0
☐	SEP110120	OMNI2 Module, Fixed Mount, Left	0
☐	SEP110123	OMNI2 Module, Swing Away Mount, Right Not With SEP090110 Rear View Mirror Right	130
☐	SEP110124	OMNI2 Module, Swing Away Mount, Left Not With SEP090118 Rear View Mirror Left	130

Hand Control For Switch Selection See Section E

£

Proportional Input Devices

☐	SEP110301	Hand Control, Switch-It MicroGuide	3145
☐	SEP110302	Hand Control, Switch-It MicroPilot	3495
☐	SEP110304	Hand Control, Switch-It VersaGuide	1120
☐	SEP110305	Hand Control, Switch-It VersaGuide EZ Not With SEP090512 Joystick Knob, Stick, 80 x 15 mm, SEP090511 Joystick Knob, Mushroom, Ø 45 mm, SEP090506 Joystick Knob, T-Bar, 82 mm Wide	1330
☐	SEP110306	Hand Control, Allround Joystick	1120
☐	SEP110307	Hand Control, Allround Joystick Light Only Use With Default Joystick Knob Or Joystick Knob, Allround, Ball Ø 48 mm	1270
☐	SEP110308	Hand Control, Allround Joystick Heavy Duty	2095
☐	SEP110309	Hand Control, Allround Joystick, R-Net Bus Enabled OMNI2 Module Not Required, Obligated To Have Master R-Net Module	1190
☐	SEP110310	Hand Control, Allround Joystick Light, R-Net Bus Enabled OMNI2 Module Not Required, Obligated To Have Master R-Net Module	1400
☐	SEP110311	Hand Control, Allround Joystick Heavy Duty, R-Net Bus Enabled OMNI2 Module Not Required, Obligated To Have Master R-Net Module	2235

Hand Control Mounting Platforms

£

☐	SEP110342	Hand Control Mounting, Link-It To The Space Disk, Hand Tray Only For Micro Series Joysticks	700
--------------------------	-----------	--	-----

Hand Control Mounting

☐	SEP110350	Hand Control Mounting, Link-It To The Armrest, Fixed	210
☐	SEP110353	Hand Control Mounting, Link-It To The Armrest, Swing Away	220

Hand Control Mounting Location

☐	SEP110362	Hand Control Mounting Location, Right	0
☐	SEP110361	Hand Control Mounting Location, Left	0

Head Control For Switch Selection See Section E. Choose A Whitmyer Headrest Mount On Previous Pages

£

Proportional Input Devices

☐	SEP110370	Head Control, Switch-It Dual Pro, 3-Way Proportional, Facial Pads	5950
☐	SEP110371	Head Control, Switch-It Dual Pro, 3-Way Proportional Head Array, Spot Pads	5950

Chin Control

£

Proportional Input Devices

☐	SEP110470	Chin Control, Switch-It MicroGuide	3145
☐	SEP110471	Chin Control, Switch-It MicroPilot	3495
☐	SEP110472	Chin Control, Switch-It VersaGuide	1120
☐	SEP110473	Chin Control, Switch-It VersaGuide EZ Not with SEP090512 Joystick Knob, Stick, 80 x 15 mm, SEP090511 Joystick Knob, Mushroom, Ø 45mm, SEP090506 Joystick Knob, T-Bar, 82 mm wide	1330
☐	SEP110474	Chin Control, Allround Joystick	1120
☐	SEP110475	Chin Control, Allround Joystick Light Not with SEP090512 Joystick Knob, Stick, 80 x 15 mm, SEP090511 Joystick Knob, Mushroom, Ø 45mm, SEP090506 Joystick Knob, T-Bar, 82 mm wide	1330
☐	SEP110476	Chin Control, Allround Joystick, R-Net Bus Enabled OMNI2 Module Not Required, Obligated To Have Master R-Net Module	1190
☐	SEP110477	Chin Control, Allround Joystick Light, R-Net Bus Enabled OMNI2 Module Not Required, Obligated to Have Master R-Net Module; Not With SEP090512 Joystick Knob, Stick, 80 x 15 mm, SEP090511 Joystick Knob, Mushroom, Ø 45 mm, SEP090506 Joystick Knob, T-Bar, 82 mm Wide	1400

Chin Control Mounting

☐	SEP110480	Chin Control Mounting, Link-It To Chin Harness, Small If You Need Switches Mounted On The Harness, You Can Only Use	630
☐	SEP110483	Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Right	560
☐	SEP110484	Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Left	560
☐	SEP110485	Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Right	2795
☐	SEP110486	Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Left	2795
☐	SEP110487	Chin Control Mounting, Link-It To Headrest, Swing Away, Left	490
☐	SEP110488	Chin Control Mounting, Link-It To Headrest, Swing Away, Right	490

Link-It Extension Kit

☐	SEP110590	Link-It Extension Kit, Base Module, 100 mm Extension Tube, 5 Cable Clips	170
☐	SEP110591	Link-It Extension Kit, Base Module, 150 mm Extension Tube, 5 Cable Clips	210

Special Controls

E – Switches Select 1 Switch And Mounting Position Accordingly. For A Second Switch Please Use #2 Switch Selection Part. Single Function
For Special Input Devices 

£

#1 Switch Selection

☐	SEP110611	#1 Single Switch, Ø 64 mm, Black	120
☐	SEP110612	#1 Single Switch, Ø 64 mm, Red	120
☐	SEP110613	#1 Single Switch, Ø 64 mm, Green	120
☐	SEP110523	#1 Egg Switch – Black	135
☐	SEP110524	#1 MicroLite Switch – Black	195
☐	SEP110651	#1 Single Switch, Ø 30 mm, Black	120
☐	SEP110654	#1 Single Switch, Ø 30 mm, Blue	120
☐	SEP110653	#1 Single Switch, Ø 30 mm, Green	120
☐	SEP110652	#1 Single Switch, Ø 30 mm, Red	120
☐	SEP110655	#1 Single Switch, Ø 30 mm, Yellow	120
☐	SEP110534	#1 Twister Switch – Black	195
☐	SEP110535	#1 Twister Switch – Green	195
☐	SEP110536	#1 Twister Switch – Red	195

#1 Switch Mounting

☐	SEP110560	#1 Switch Mounting, Link-It To The Hand Control Not With HD Allround Joystick Or Space Disk Mount	210
☐	SEP110566	#1 Switch Mounting, Link-It To The Armrest (Opposite Site As The Hand Control) Armrest Rail Mount	210
☐	SEP110561	#1 Switch Mounting, Link-It To The Headrest	280
☐	SEP110563	#1 Switch Mounting, Link-It To The Seat Rail Operated By Knee, If Two Switches Selected, One Mounted On The Right Side, One On The Left	210
☐	SEP110565	#1 Switch Mounting, Link-It To Chin Boom No Left / Right Switch Selection Required	210
☐	SEP110569	#1 Omit Switch Mounting Switches Attached Lose	0

#1 Switch Mounting Location

☐	SEP110572	#1 Switch Mounting Location, Right	0
☐	SEP110571	#1 Switch Mounting Location, Left	0

#1 Satellite Switches Only For Allround Joystick, Maximum 2 Selectable

☐	SEP110537	#1 Satellite Switch – Black	210
☐	SEP110538	#1 Satellite Switch – Green	210
☐	SEP110539	#1 Satellite Switch – Red	210

#2 Switch Selection

☐	SEP110631	#2 Single Switch, Ø 64 mm, Black	120
☐	SEP110632	#2 Single Switch, Ø 64 mm, Red	120
☐	SEP110633	#2 Single Switch, Ø 64 mm, Green	120
☐	SEP110088	#2 Egg Switch – Black	135
☐	SEP110089	#2 MicroLite Switch – Black	185
☐	SEP110671	#2 Single Switch, Ø 30 mm, Black	120
☐	SEP110674	#2 Single Switch, Ø 30 mm, Blue	120
☐	SEP110673	#2 Single Switch, Ø 30 mm, Green	120
☐	SEP110672	#2 Single Switch, Ø 30 mm, Red	120
☐	SEP110675	#2 Single Switch, Ø 30 mm, Yellow	120
☐	SEP110095	#2 Twister Switch – Black	185
☐	SEP110096	#2 Twister Switch – Green	185
☐	SEP110097	#2 Twister Switch – Red	185

#2 Switch Mounting

☐	SEP110098	#2 Switch Mounting, Link-It To The Hand Control Not With HD Allround Joystick Or Space Disk Mount	210
☐	SEP110099	#2 Switch Mounting, Link-It To The Armrest (Opposite Site As The Hand Control) Armrest Rail Mount	210
☐	SEP110100	#2 Switch Mounting, Link-It To The Headrest	280
☐	SEP110104	#2 Switch Mounting, Link-It To The Seat Rail Knee, If 2 Switches Selected Only Left And Right Mount	210
☐	SEP110105	#2 Switch Mounting, Link-It To Chin Boom No Left / Right Switch Selection Required	210
☐	SEP110106	#2 Omit Switch Mounting Switches Attached Lose	0

#2 Switch Mounting Location

☐	SEP110108	#2 Switch Mounting Location, Right	0
☐	SEP11107	#2 Switch Mounting Location, Left	0

#2 Satellite Switches Only For Allround Joystick, Maximum 2 Selectable

☐	SEP110111	#2 Satellite Switch – Black	200
☐	SEP110112	#2 Satellite Switch – Green	200
☐	SEP110113	#2 Satellite Switch – Red	200

Multi-Function Switches Used In Combination With Assignable Button Programming

☐	SEP110550	Mini 2 Button (Stereo Plug) Attached Lose, When Ordered. No Link-It Fixation Possible	385
☐	SEP110551	Mini 2 Button (Dual Mono Plug) Attached Lose, When Ordered. No Link-It Fixation Possible	385
☐	SEP110552	Mini 4 Button InLine (DB9 Connector) Attached Lose, When Ordered. No Link-It Fixation Possible	645
☐	SEP110554	4-Way Toggle Switched Joystick	645

Accessories

☐	SEP110580	Mono Extension Cable (3.5 mm Jack Socket) (2 m)	140
☐	SEP110581	Stereo Extension Cable (3.5 mm Jack Socket) (2 m)	180
☐	SEP110582	Extension Cable DB9 (1 m)	170
☐	SEP110583	Y-Cable, 2 Female Mono To 1 Male Stereo	210

Safety Devices

£

☐	SEP090147	Remote Stop Switch (9 m Activation Range, Fit Via OBC Port) Not With VR2	980
--------------------------	-----------	---	-----

Comments For Adaption



Q700 M/F/R with SEDEO PRO

£

Base Price	00,000.00
Total Of Additional Equipment	
Packaging And Transport	
Exclusive VAT	
VAT	
Total Price	

Date

Signature

Day

Month

Year

Send E-Mail



Sunrise Medical GmbH

Kahlbachring 2-4

D-69254 Malsch/HD

Phone: +49 (0) 72 53 / 9 80-0

Fax: +49 (0) 72 53 / 9 80-222

E-Mail: kundenservice@sunrisemedical.de

www.SunriseMedical.eu