

Q700-UP M/F

with
SEDEO® PRO ADVANCED™




ORDER FORM

Billing Address Or Customer-ID		Shipping Address		☐ Identical To Billing Address
C.-ID:		C.-ID:		
Name:		Name:		
Street:		Street:		
Town/PC:		Town/PC:		
E-Mail:		E-Mail:		
Phone:		Phone:		
Order Information		User Weight		
☐ Order ☐ Calculation ☐ Built-4-Me®		User Weight Capacity: 120 kg		
Purchase Order No.:		User Weight:		kg
Art. no.: Q700-UP M: 2320070M-000, Q700-UP F: 2320070F-000				
Order Date:	Day	Month	Year	

Client Measurements

a) Top of head to seat		mm
b) Top of shoulder to seat		mm
c) Underarm to seat		mm
d) Elbow to seat		mm
e) Top of pelvis to seat		mm
f) Shoulder width		mm
g) Trunk width		mm
i) Hip width		mm
j) Actual sitting depth		mm
m) Foot length		mm
n) Lower leg length		mm
o) Actual seat width		mm



Seating Specifications

		Seat Type	SEDEO PRO ADVANCED UP
Typical WC Weight	Without Batteries 114 kg	Seat Size Range	Width 380-560 mm, Depth 400-540mm
Typical WC Weight	With Batteries 160 kg	Seat to Floor Height Range	450 - 480 mm
Power Seating	22° of CG Tilt And Power Recline With 180 mm Mechanical Shear Reduction / 200 mm Seat Elevator		

MWD Base Specifications

User Weight Capacity	120 kg	Drive Wheels	14" Drive Wheels (Black)
Base Size	Width Min. 622 mm	Tires / Castors	6" Front / 6" Rear (Black)
Drive Control	R-Net PM 90A/120A (Expandable)	Charger	8 A Dual Mode Charger
Driving Range	Up To 40 km	Battery Size	60 And 80 Ah Batteries
Drive Motors	4-Pole 6 kph, Options 8, 10 And 12.5 kph	Colour Accents	7 Colour Accents
Ground Clearance	90 mm	Frame Colour	Black

FWD Base Specifications

User Weight Capacity	120 kg	Drive Wheels	14" Drive Wheels (Black)
Base Size	Width Min. 640 mm	Tires / Castors	9" Rear (Black)
Drive Control	R-Net PM 90A/120A (Expandable)	Charger	8 A Dual Mode Charger
Driving Range	Up To 40 km	Battery Size	60 And 80 Ah Batteries
Drive Motors	4-Pole 6 Kph, Options 8, 10 And 12.5 Kph	Colour Accents	7 Colour Accents
Ground Clearance	71 mm	Frame Colour	Black

Q700-UP CHAIR

☐	QUE000008	Q700-UP M (MWD Base Chair)	19,995
☐	QUE000016	Q700-UP F (FWD Base Chair)	19,995

	Motor Package	£
☐	QUE180009 4-Pole Motor, 6 kph / 4 mph (8 kph Motor With Reduced Speed)	0
☐	QUE180007 4-Pole Motor, 10 kph / 6.5 mph	0
☐	QUE180008 4-Pole Motor, 12.5 kph / 8 mph	125

	Tracking	£
☐	QUE110009 Traction Control System For R-Net (GYRO Technology)	0
☐	QUE110014 Omit Traction Control System For R-Net	0

	Drive Wheels 	£
☐	QUE070034 Drive Wheel, 14", Pneumatic, Knobbly, Black	0
☐	QUE070033 Drive Wheel, 14", Solid, Knobbly, Black	0

	Castor Wheels 	
MWD Base		
☐	QUE080206 Castor Wheel, 6", Solid, Black	0
FWD Base		
☐	QUE080210 Castor Wheel, 9", Solid, Black	0
☐	QUE080211 Castor Wheel, 9", Pneumatic, Black	0

	Mudguards (FWD Only) 	£
☐	QUE010007 Mudguard, Drive Wheel	0
☐	QUE010008 Mudguard, Castor Wheel, 9"	0

	Colour Accents	£
☐	QUE140005 Apple Green (RAL 110 70 77)	0
☐	QUE140004 Black (RAL 9005)	0
☐	QUE140003 Blue (RAL 260 40 45)	0
☐	QUE140006 Orange (RAL 060 60 80)	0
☐	QUE140002 Red (RAL 040 40 67)	0
☐	QUE140007 Silver (RAL 9006)	0
☐	QUE140001 White (RAL 9010)	0

	LED Lights & Indicators 	£
☐	QUE110201 LED Light and Indicator QUE110031 Lighting Module, R-Net LM required Approved for road use	0
☐	QUE110202 LED Fender Light QUE110031 Lighting Module, R-Net LM required, Fender Light - Not road approved and without indicators	150
☐	QUE110209 Omit Light And Indicator	0

Seat Heights Measured from front of seat pan to floor. Obtain quote to verify seat-to-floor heights.			£
☐	SH45	450 mm Only with QUE160001 battery, 60 Ah, Gel, Maintenance-Free In case of FWD base, not with SD 520mm & SD 540mm	0
☐	SH465	465 mm In case of FWD base, not with SD540 mm	0
☐	SH48	480 mm	0

Batteries & Chargers	£
---	---

Batteries

☐	QUE160001	Battery, 60 Ah, Gel, Maintenance-Free	0
☐	QUE160005	Battery, 80 Ah, Gel, Maintenance-Free Only with battery charger 220 V / 10 A to reduce reloading time	150

Chargers

☐	QUE160100	Charger, Automatic, 220 V / 8 A, UK	0
☐	QUE160103	Charger, Automatic, 220 V / 10 A, UK	100

Power Seating Options	£
--	---

☐	QUE130010	Stand-Up Seat, Powered	0
		Posterior Tilt, Powered, 0° - 22°	
		Anterior Tilt, Powered, 0° - 80°	
		Seat Elevator, Powered, 200 mm	

Backrest Frame Type	£
----------------------------	---

☐	QUE030047	Backrest Frame Type, Powered Recline, 90° - 175° With 180 mm Mechanical Anti-Shear	0
--------------------------	-----------	--	---

Back Cushion Type	£
--------------------------	---

SEDEO PRO (See Section A)

☐	QUE030104	SEDEO PRO, Back Cushion	0
--------------------------	-----------	-------------------------	---

JAY Syncra (See Section B)

☐	QUE030216	JAY Syncra, Back Cushion	0
--------------------------	-----------	--------------------------	---

JAY J3 (See Section C)

☐	QUE030054	JAY J3, Back Cushion	0
--------------------------	-----------	----------------------	---

A SEDEO PRO Back Size			£
Back Cushion Width			
☐	BW38	380 mm, For Effective Back Width 380 mm	0
☐	BW42	420 mm, For Effective Back Width 380 - 420 mm	0
☐	BW46	460 mm, For Effective Back Width 420 - 460 mm	0
☐	BW50	500 mm, For Effective Back Width 460 - 500 mm	0
☐	BW54	540 mm, For Effective Back Width 500 - 540 mm	0
☐	BW56	560 mm, For Effective Back Width 520 - 560 mm	0
Back Cushion Height			
☐	BH44	440 mm, For Effective Back Height 510 - 570 mm	0
☐	BH50	500 mm, For Effective Back Height 570 - 630 mm	0
☐	BH56	560 mm, For Effective Back Height 630 - 700 mm	0
Installed Back Height (From seat pan to top of back)			
☐	BIH51	510mm (440 mm Back Cushion Height)	0
☐	BIH54	540mm (440 mm Back Cushion Height)	0
☐	BIH57	570mm (440 mm Back Cushion Height)	0
☐	BIH57_1	570mm (500 mm Back Cushion Height)	0
☐	BIH60	600mm (500 mm Back Cushion Height)	0
☐	BIH63	630mm (500 mm Back Cushion Height)	0
☐	BIH63_1	630mm (560 mm Back Cushion Height)	0
☐	BIH66	660mm (560 mm Back Cushion Height)	0
☐	BIH70	700mm (560 mm Back Cushion Height)	0
A SEDEO PRO Back Cushion 			£
☐	QUE030116	SEDEO PRO, Back Cushion, Flexible, Deep Contour, Active, Adjustable Straps	300
☐	QUE030117	SEDEO PRO, Back Cushion, Flexible, Extra Deep Contour, Active, Adjustable Straps	300
A SEDEO PRO Back Cushion Cover			£
☐	QUE030108	Back Cushion Cover, Comfort	0

B	JAY Syncra
----------	---

B	JAY Syncra Back Width	£
----------	------------------------------	----------

Back can be ordered a different width than seat pan. Armrest is adjusted to back width ordered.

☐	BW40	400 mm (Adjustment range 400 mm to 430 mm, Contour Depth 75 mm)	0
☐	BW46	460 mm (Adjustment range 460 mm to 480 mm, Contour Depth 75 mm)	0
☐	BW51	510 mm (Adjustment range 510 mm to 530 mm, Contour Depth 75 mm)	0
☐	BW56	560 mm (Adjustment range 560 mm to 585 mm, Contour Depth 75mm)	0

B	JAY Syncra Back Shell Options And Height Ranges	£
----------	--	----------

JAY Syncra back shells are mounted with a minimim 14 cm gap between the bottom of the shell to the seat pan to accommodate the mechanical anti-shear. Back height can be adjusted in the field within the ranges shown based on the back frame type.

JAY Syncra Shell Height

☐	BH40	400 mm Tall Shell (560 mm to 610 mm Back Height)	0
☐	BH46	460 mm Tall Shell (610 mm to 660 mm Back Height)	0
☐	BH51	510 mm TallShell (660 mm to 710 mm Back Height)	0
☐	BH56	560 mm Tall Shell (710 mm to 760 mm Back Height)	0

JAY Syncra Back Installed Back Height (From Pan to Top of Back)

☐	BIH56	560 mm (400mm Shell Height)	0
☐	BIH585	585 mm (400mm Shell Height)	0
☐	BIH61	610 mm (400mm Shell Height)	0
☐	BIH61_1	610 mm (460mm Shell Height)	0
☐	BIH635	635 mm (460mm Shell Height)	0
☐	BIH66	660 mm (460mm Shell Height)	0
☐	BIH66_1	660 mm (510mm Shell Height)	0
☐	BIH685	685 mm (510mm Shell Height)	0
☐	BIH71	710 mm (510mm Shell Height)	0
☐	BIH71_1	710 mm (560mm Shell Height)	0
☐	BIH735	735 mm (560mm Shell Height)	0
☐	BIH76	760 mm (560mm Shell Height)	0

B	JAY Syncra Back Upholstery	£
----------	-----------------------------------	----------

☐	JAY030060	Syncra 3DX Microclimatic	0
--------------------------	-----------	--------------------------	---

B/C	JAY Lateral Thoracic Swing-Away Supports (Requires Pad, Upholstery and Mount)	£
------------	--	----------

JAY Lateral Thoracic Support Pad Upholstery

☐	JAY090110	Syncra 3DX Microclimatic	0
--------------------------	-----------	--------------------------	---

JAY Fixed Lateral Support Accessory Hardware

☐	JAY090111	Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm - Right JAY Fixed Lateral Support is adjustable - Syncra Back only to increase the contour depth up to 120 mm - Not available with JAY J3 Back	0
☐	JAY090112	Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm - Left JAY Fixed Lateral Support is adjustable - Syncra Back only to increase the contour depth up to 120 mm - Not available with JAY J3 Back	0

C	JAY J3 Back Size	£
☐	JAY030111 JAY J3 Back Size, 410 mm	0
☐	JAY030112 JAY J3 Back Size, 460 mm	0
☐	JAY030113 JAY J3 Back Size, 510 mm	0

C	JAY J3 Back, Mid Deep Contour 	£
☐	JAY030042 JAY J3 Back, Mid Deep Contour, Upper Thoracic, Medium, 460 mm	721

C	JAY J3 Back, Deep Contour 	£
☐	JAY030021 JAY J3 Back, Deep Contour, Upper Thoracic, Medium, 460 mm	721

C	JAY J3 Installed Back Height (From Pan to Top of Back) 	£
☐	BIH53 530 mm (420mm Shell Height)	0
☐	BIH56 560 mm (420mm Shell Height)	0
☐	BIH585 585 mm (420mm Shell Height)	0
☐	BIH585_1 585 mm (460mm Shell Height)	0
☐	BIH61 610 mm (460mm Shell Height)	0
☐	BIH61_1 610 mm (500mm Shell Height)	0
☐	BIH635 635 mm (500mm Shell Height)	0
☐	BIH66 660 mm (500mm Shell Height)	0

	Headrest Pad 	£
☐	QUE030430 SEDEO, Headrest Pad, Small, (LxH) 240 x 150 mm Only use with SEDEO Headrest Mount	0
☐	QUE030431 SEDEO, Headrest Pad, Small, With Cover, (LxH) 240 x 150 mm Only use with SEDEO Headrest Mount	0
☐	QUE030432 SEDEO, Headrest Pad, Flexible, (LxH) 400 x 110 mm Only use with SEDEO Headrest Mount	145
☐	WHT030402 Whitmyer, Headrest Pad, Plush, (LxH) 254 x 127 mm Only use with Whitmyer Headrest Mount	140
☐	WHT030404 Whitmyer, Headrest Pad, Adjust-A-Plush, Small, (LxH) 355 x 76 mm Only use with Whitmyer Headrest Mount	225

	Headrest Mount  Choose a Whitmyer Headrest Mount and Headrest Receiver in case you would like to add a Head Control in the Special Control section.	£
☐	QUE030433 SEDEO, Headrest Mount	0
☐	WHT030408 Whitmyer, Headrest Mount, LINX²	150
☐	WHT030409 Whitmyer, Headrest Mount, AXYS	130
☐	WHT030411 Whitmyer, Headrest Mount, COBRA XTRA Flip Back	190

	Headrest Receiver	£
☐	QUE030120 Headrest Receiver, Backrest Frame Mount	0

	Push Handles	£
☐	QUE030202 Push Handles, Fixed	175
☐	QUE030200 Omit Push Handles	0

	Seat Rail	£
☐	QUE020085 Seat Rail, Short, For Seat Depth 400 - 540 mm	0
☐	QUE020069 Seat Rail, Long, For Seat Depth 460 - 540 mm Not available SW380mm-SW460mm in combination with SD460mm - SD500mm	0

	Seat Cushion Type	£
SEDEO (See Section D)		
☐	QUE020027 SEDEO, Seat Cushion	0
JAY Syncra (See Section E)		
☐	QUE020132 JAY Syncra, Seat Cushion	0
Omit (Select Seat Width And Depth in Section D)		
☐	QUE020110 Omit Seat Cushion	0

D	SEDEO Seat Size	£
Seat Cushion Width		
☐	SW38 380 mm, For Effective Seat Width 380 mm Minimum seat pan width (outside of seat rail to outside of seat rail) is 460 mm - Cushion will be supplied at selected size	0
☐	SW42 420 mm, For Effective Seat Width 380 - 420 mm Minimum seat pan width (outside of seat rail to outside of seat rail) is 460 mm - Cushion will be supplied at selected size	0
☐	SW46 460 mm, For Effective Seat Width 420 - 460 mm Minimum seat pan width (outside of seat rail to outside of seat rail) is 460 mm - Cushion will be supplied at selected size	0
☐	SW50 500 mm, For Effective Seat Width 460 - 500 mm	0
☐	SW54 540 mm, For Effective Seat Width 500 - 540 mm	0
☐	SW56 560 mm, For Effective Seat Width 520 - 560 mm	0
Seat Cushion Depth		
☐	SD40 400 mm, Comes With Cushion Size 460 mm	0
☐	SD42 420 mm, Comes With Cushion Size 460 mm	0
☐	SD44 440 mm, Comes With Cushion Size 460 mm	0
☐	SD46 460 mm, Comes With Cushion Size 500 mm	0
☐	SD48 480 mm, Comes With Cushion Size 500 mm	0
☐	SD50 500 mm, Comes With Cushion Size 500 mm	0
☐	SD52 520 mm, Comes With Cushion Size 560 mm	0
☐	SD54 540 mm, Comes With Cushion Size 560 mm	0

D	SEDEO Seat Cushions 	£
☐	QUE020073 SEDEO PRO, Seat Cushion, Profiled, Height 70 mm	130
☐	QUE020076 SEDEO PRO, Seat Cushion, Support Memory, Height 90 mm	220

D	SEDEO Seat Cushion Cover	£
☐	QUE020518 Seat Cushion Cover, Incontinence Protection	0
☐	QUE020080 Seat Cushion Cover, Comfort	0

D/E	SEDEO Seat Boards	£
☐	QUE020082 SEDEO, Seat Board, Flat, Depth-Adjustable, Cushions Guide At Rear	115
☐	QUE020083 SEDEO, Seat Pan, 50 mm, Depth-Adjustable, Cushions Guide At Rear	130
☐	QUE020084 SEDEO, Seat Pan, 75 mm, Depth-Adjustable, Cushions Guide At Rear	135

E JAY Syncra Seat Size			£
Seat Cushion Width			
☐	SW40	400 mm	0
☐	SW43	430 mm	0
☐	SW46	460 mm	0
☐	SW48	480 mm	0
☐	SW51	510 mm	0
☐	SW56	560 mm	0
Seat Cushion Depth			
☐	SD40	400 mm	0
☐	SD43	430 mm	0
☐	SD46	460 mm	0
☐	SD48	480 mm	0
☐	SD51	510 mm	0
E JAY Syncra Seat Cushion 			£
☐	JAY020030	JAY Syncra Cushion with Visco Foam	300
☐	JAY020031	JAY Syncra Cushion featuring JAY Cryo/Flow Fluid	500
E JAY Syncra Cushion Cover Upholstery			£
☐	JAY020040	Syncra 3DX Microclimatic	0
☐	JAY020041	Syncra Incontinence Protection	0
Knee & Thigh Supports Not Available With SW 380 mm / 400 mm / 420 mm			£
☐	QUE020342	Knee Support, Lateral, 120 x 120 mm, Swing Away, Right	125
☐	QUE020341	Knee Support, Lateral, 120 x 120 mm, Swing Away, Left	125
☐	QUE020113	Thigh Support, Lateral, Small, 120 x 120 mm, Swing Away, Right	150
☐	QUE020114	Thigh Support, Lateral, Small, 120 x 120 mm, Swing Away, Left	150
☐	QUE020121	Thigh Support, Lateral, Large, 180 x 120 mm, Removable, Right Not available with QUE040036 Side Guard	150
☐	QUE020122	Thigh Support, Lateral, Large, 180 x 120 mm, Removable, Left Not available with QUE040036 Side Guard	150
Armrest Type 			£
☐	QUE040070	Armrest, Reclining, Multi-Adjustable, Flip Back	0
Arm Pads 			£
☐	QUE040023	Arm Pad, Integral Foam, 300 mm, Right	0
☐	QUE040024	Arm Pad, Integral Foam, 300 mm, Left	0
☐	QUE040025	Arm Pad, Integral Foam, 400 mm, Right	0
☐	QUE040026	Arm Pad, Integral Foam, 400 mm, Left	0
Comfort Arm Pads			
☐	QUE040204	Arm Pad, Comfort, Padded, PU Coated, 400 mm, Right	100
☐	QUE040203	Arm Pad, Comfort, Padded, PU Coated, 400 mm, Left	100
Hemiplegia Arm Pads			
☐	QUE040031	Arm Pad, Hemiplegia, Integral Foam, Right Not available with Tray Table QUE090040 / QUE090046 / QUE110143	130
☐	QUE040032	Arm Pad, Hemiplegia, Integral Foam, Left Not available with Tray Table QUE090040 / QUE090046 / QUE110143	130

Side Guard Pads For Armrest			£
Side Guard Pads For Armrest, Reclining, Multi-Adjustable, Flip Back			
☐	QUE040036	Side Guard, Mounted On Seat Rail Bracket Not available with SW380mm/400mm/420mm	0
Power Centre Mount Legrest 			£
☐	QUE050071	Power Centre Mount Legrest, 90° - 180° Angle, 165 mm Extension	0
Footplates			
☐	QUE050072	Footplate, Dual, Height- And Angle-Adjustable, 260 x 150 mm	0
☐	QUE050073	Footplate, Dual, Height- And Angle-Adjustable, Large, 310 x 170 mm	0
Pad Options			
☐	QUE050074	Knee Fixation Pad, Integral Foam, Adjustable	0
☐	QUE050139	Knee Fixation Pad, Soft Foam, Small, Adjustable Shortest distance calf support to knee fixation + 40 mm adjustment	50
☐	QUE050092	Calf Pad, Padded, Pair	0
Positioning Supports & Accessories 			£
☐	QUE090047	Positioning Belt, Aircraft Buckle, 50 x 1300 mm	0
☐	QUE090048	Positioning Belt, Aircraft Buckle, Padded, 50 x 1300 mm	0
☐	QUE090052	Chest Strap, 2-Point, Padded	0
☐	QUE090053	Additional Chest Strap, 2-Point, Padded	75
☐	QUE090097	Chest Bar, Detachable, Round Pad Not available with Tray Table QUE090040 / QUE090046 / QUE110143	150
☐	QUE090098	Chest Bar, Detachable, Flat Pad, 200 x 100 mm Not available with Tray Table QUE090040 / QUE090046 / QUE110143	150
☐	QUE090110	Rear View Mirror, Right	45
☐	QUE090118	Rear View Mirror, Left	45
☐	QUE090122	Transfer Bar Not available with Knee Support	200
☐	QUE090600	Black Utility Backpack (Requires Push Handles)	75
☐	QUE090172	Ram Mount Phone	80
☐	QUE090173	Ram Mount Tablet	130
☐	QUE090174	Ram Mount Cup Holder	30
☐	QUE090175	Reflectors (Backrest shroud mounted)	50
Tray Table (Mounted Opposite Joystick)			
☐	QUE090040	Tray Table, Width-Adjustable, Swing Away	300
☐	QUE090046	Tray Table, Half-Side, Swing Away	200
Tray Table Accessories			
☐	QUE090126	Raised Edges, Acrylic For Tray Table Not available with QUE090046 Tray Table, Half-Side, Swing Away	30
Remote Joystick – R-Net 			£
☐	QUE110017	Remote Joystick, R-Net, Colour LCD	0
☐	QUE110018	Remote Joystick, R-Net, Colour LCD, Heavy Duty	900
☐	QUE110013	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared Control Included	0
☐	QUE110023	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared And Bluetooth Control Included	500
Joystick Mounting 			£
☐	QUE110050	Joystick Mounting, Fixed	0
☐	QUE110056	Joystick Mounting, Parallel Swing Away, Capsulated	0

Joystick Mounting Location			£
☐	QUE110052	Joystick Mounting Location, Right	0
☐	QUE110053	Joystick Mounting Location, Left	0
Joystick Knobs (Optional) 			£
☐	QUE090505	Joystick Knob, Golf Ball, Ø 42 mm	35
☐	QUE090506	Joystick Knob, T-Bar, 82 mm Wide	35
☐	QUE090507	Joystick Knob, Soft Ball, Ø 70 mm	35
☐	QUE090508	Joystick Knob, Ball, Ø 35 mm	35
☐	QUE090045	Joystick Knob, Allround, Ball, Ø 48 mm	35
☐	QUE090509	Joystick Knob, Chin Cup, Ø 37 mm	35
☐	QUE090512	Joystick Knob, Stick, 80 x 15 mm	35
☐	QUE090511	Joystick Knob, Mushroom, Ø 45 mm	35
Power Modules – R-Net 			£
☐	QUE110003	Power Module, R-Net, 90 A	0
☐	QUE110004	Power Module, R-Net, 120 A	200
Seating / Lighting Module Control – R-Net 			£
☐	QUE110031	Lighting Module, R-Net	0
☐	QUE110033	Actuator Control Module CxSM, R-Net	0
Attendant Controls 			£
☐	QUE110021	Attendant Control, Dual	400
☐	QUE110020	Joystick Mounted As Attendant Control Choose also Remote Joystick	250
Attendant Mounting Location			
☐	QUE110025	Attendant Mounting Location, Push Handle, Right	0
☐	QUE110026	Attendant Mounting Location, Push Handle, Left	0
Table Top Controls 			£
☐	QUE110143	Table Top Control, LCD, Centre Tray Mounted	1900
Tray Table Mounting Location			
☐	QUE110080	Tray Table Mounting Location, Left	0
☐	QUE110081	Tray Table Mounting Location, Right	0
Memory Seating and Pressure Relief Alert			£
Biometric Repositioning, includes Memory Seating for 6 positions, Stand UP, Pressure Relief Relax, Lie Down Supine, Transfer (Anterior Tilt), Memory Seated (Lifted), Seated.			
☐	QUE190001	Biometric Repositioning (Includes Memory Seating)	0
Pressure Relief Alert and Recording (PRAR) - Programmable reminder to assist users to use their Pressure Relief Memory Seat Position, in combination with CTRL+5.			
☐	QUE190004	PRAR Enabled (Factory Set To Relax Axis, Duration 1 Minute)	0

Table Top Contols 			£
☐	QUE110143	Table Top Control, LCD, Centre Tray Mounted	1900
Tray Table Mounting Location			
☐	QUE110080	Tray Table Mounting Location, Left	0
☐	QUE110081	Tray Table Mounting Location, Right	0
Memory Seating and Pressure Relief Alert			£
Biometric Repositioning, includes Memory Seating for 6 positions, Stand UP, Pressure Relief Relax, Lie Down Supine, Transfer (Anterior Tilt), Memory Seated (Lifted), Seated.			
☐	QUE190001	Biometric Repositioning (Includes Memory Seating)	0
Pressure Relief Alert and Recording (PRAR) - Programmable reminder to assist users to use their Pressure Relief Memory Seat Position, in combination with CTRL+5.			
☐	QUE190004	PRAR Enabled (Factory Set To Relax Axis, Duration 1 Minute)	0
#1 Ctrl+5 Switch Box Options 			£
☐	QUE110043	#1 Ctrl+5 Switch Box With Free Assignable Functions Voice activation functionality included	575
☐	QUE110179	#1 Ctrl+5 Switch Box With Free Assignable Functions And Switch-It Remote Seating To achieve full functionality of the Remote Seating, it is required to download the app from AppleStore (iOS) or Playstore (Adroid) and install it on a smartphone or tablet. Further information and instructions about the Remote Seating App can be found at www.SunriseMedical.com . Voice activation functionality included	575
#1 Ctrl+5 Switch Box Mounting Location			
☐	QUE110180	#1 Ctrl+5 Switch Box Mounting Location, Behind Joystick If combined with Attendant Control, Ctrl+5 will be mounted to the armrest	0
☐	QUE110181	#1 Ctrl+5 Switch Box Mounting Location, Opposite Joystick If combined with Attendant Control, Ctrl+5 will be mounted to the armrest Not available with Tray Table QUE090040 / QUE090046 / QUE110143	0
#1 Ctrl+5 Switch Type			
☐	QUE110175	#1 Ctrl+5 Switch Type With Push Buttons	0
☐	QUE110176	#1 Ctrl+5 Switch Type With Toggle Switches	0
SPECIAL CONTROLS			
OMNI2 Input Control Modules – R-Net 			£
Choose OMNI2 mount and input device in sections below, and choose switches in section F.			
☐	QUE110116	OMNI2 Input Control Module, R-Net, Bluetooth, Infrared Remote Full VGA Display, with separate input box with Bluetooth, Infrared, USB charger and speaker for specialty input devices	1050
☐	QUE110117	OMNI2 Input Control Module, R-Net, Scanning Mode, Bluetooth, Infrared Remote Full VGA Display, with separate input box with Bluetooth, Infrared, USB charger and speaker for specialty input devices	1050
OMNI2 Module Mount			£
Choose OMNI2 mount and input device in sections below, and choose switches in section F. If used with standard joystick selection, it must be mounted opposite side.			
☐	QUE110119	OMNI2 Module, Fixed Mount, Right	0
☐	QUE110120	OMNI2 Module, Fixed Mount, Left	0
☐	QUE110123	OMNI2 Module, Swing Away Mount, Right	100
☐	QUE110124	OMNI2 Module, Swing Away Mount, Left	100

Hand Control  For Switch Selection see Section F.		£
---	--	----------

Proportional Input Devices

☐	QUE110301	Hand Control, Switch-It MicroGuide	2300
☐	QUE110302	Hand Control, Switch-It MicroPilot	2300
☐	QUE110304	Hand Control, Switch-It VersaGuide	1050
☐	QUE110305	Hand Control, Switch-It VersaGuide EZ Not with QUE090512 Joystick Knob, Stick, 80 x 15 mm, QUE090511 Joystick Knob, Mushroom, Ø 45mm, QUE090506 Joystick Knob, T-Bar, 82 mm wide	1270
☐	QUE110306	Hand Control, Allround Joystick	1000
☐	QUE110307	Hand Control, Allround Joystick Light Only use with default joystick knob or joystick knob, Allround, Ball Ø 48 mm	1100
☐	QUE110308	Hand Control, Allround Joystick Heavy Duty	1650
☐	QUE110309	Hand Control, Allround Joystick, R-Net Bus Enabled OMNI2 module not required, obligated to have master R-Net module	1200
☐	QUE110310	Hand Control, Allround Joystick Light, R-Net Bus Enabled OMNI2 module not required, obligated to have master R-Net module	1450
☐	QUE110311	Hand Control, Allround Joystick Heavy Duty, R-Net Bus Enabled OMNI2 module not required, obligated to have master R-Net module	1995

Hand Control Mounting Platforms 		£
--	--	----------

☐	QUE110342	Hand Control Mounting, Link-It To The Space Disk, Hand Tray Only for Micro Series joystick	450
--------------------------	-----------	---	-----

Hand Control Mounting Position

☐	QUE110350	Hand Control Mounting, Link-It To The Armrest, Fixed Not with Control 5-Switch Box behind joystick QUE110180	0
☐	QUE110353	Hand Control Mounting, Link-It To The Armrest, Swing Away Not with Control 5-Switch Box behind joystick QUE110180	100

Hand Control Mounting Location

☐	QUE110362	Hand Control Mounting Location, Right	0
☐	QUE110361	Hand Control Mounting Location, Left	0

Head Control  For switch selection see Section F. Choose a Whitmyer Headrest Mount on previous pages.		£
---	--	----------

Proportional Input Devices

☐	QUE110370	Head Control, Switch-It Dual Pro, 3-Way Proportional, Facial Pads	4500
☐	QUE110371	Head Control, Switch-It Dual Pro, 3-Way Proportional Head Array, Spot Pads	4500

Chin Control 			
Proportional Input Devices			
☐	QUE110470	Chin Control, Switch-It MicroGuide	2300
☐	QUE110471	Chin Control, Switch-It MicroPilot	2300
☐	QUE110472	Chin Control, Switch-It VersaGuide	1050
☐	QUE110473	Chin Control, Switch-It VersaGuide EZ Not with QUE090512 Joystick Knob, Stick, 80 x 15 mm, QUE090511 Joystick Knob, Mushroom, Ø 45mm, QUE090506 Joystick Knob, T-Bar, 82 mm wide	1270
☐	QUE110474	Chin Control, Allround Joystick	1000
☐	QUE110475	Chin Control, Allround Joystick Light Only use with default joystick knob or joystick knob, Allround, Ball Ø 48 mm	1100
☐	QUE110476	Chin Control, Allround Joystick, R-Net Bus Enabled OMNI2 module not required, obligated to have master R-Net module	1200
☐	QUE110477	Chin Control, Allround Joystick Light, R-Net Bus Enabled (6) (3) Only use with default joystick knob or joystick knob, Allround, Ball Ø 48 mm OMNI2 module not required, obligated to have master R-Net module	1450
Chin Control Mounting 			
☐	QUE110480	Chin Control Mounting, Link-It To Chin Harness, Small If you need switches mounted on the harness, you can only use Allround Joystick	400
☐	QUE110481	Chin Control Mounting, Link-It To Chin Harness, Large If you need switches mounted on the harness, you can only use Allround Joystick	400
☐	QUE110483	Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Right	1300
☐	QUE110484	Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Left	1300
☐	QUE110485	Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Right	2000
☐	QUE110486	Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Left	2000
☐	QUE110487	Chin Control Mounting, Link-It To Headrest, Swing Away, Left	550
☐	QUE110488	Chin Control Mounting, Link-It To Headrest, Swing Away, Right	550
Link-It Extension Kit			
☐	QUE110590	Link-It Extension Kit, Base Module, 100 mm Extension Tube, 5 Cable Clips	100
☐	QUE110591	Link-It Extension Kit, Base Module, 150 mm Extension Tube, 5 Cable Clips	100
Sip 'N Puff For switch selection see Section F.			£
☐	QUE110510	Sip 'N Puff With Link-It Hardware, Swing Away Mount To The Headrest	880,-

F	Switches  Select 1 switch and mounting position accordingly. For a second switch please use #2 Switch Selection part. Single function for Special Input Devices.	£
----------	--	----------

#1 Switch Selection

☐	QUE110611	#1 Single Switch, Ø 64 mm, Black	100
☐	QUE110612	#1 Single Switch, Ø 64 mm, Red	100
☐	QUE110613	#1 Single Switch, Ø 64 mm, Green	100
☐	QUE110523	#1 Egg Switch – Black	100
☐	QUE110524	#1 MicroLite Switch – Black	140
☐	QUE110651	#1 Single Switch, Ø 30 mm, Black	100
☐	QUE110654	#1 Single Switch, Ø 30 mm, Blue	100
☐	QUE110653	#1 Single Switch, Ø 30 mm, Green	100
☐	QUE110652	#1 Single Switch, Ø 30 mm, Red	100
☐	QUE110655	#1 Single Switch, Ø 30 mm, Yellow	100
☐	QUE110534	#1 Twister Switch – Black	140
☐	QUE110535	#1 Twister Switch – Green	140
☐	QUE110536	#1 Twister Switch – Red	140

#1 Switch Mounting

☐	QUE110560	#1 Switch Mounting, Link-It To The Hand Control Not with HD Allround Joystick or Space Disk mount	220
☐	QUE110566	#1 Switch Mounting, Link-It To The Armrest (Opposite side as the Hand Control) Armrest rail mount	130
☐	QUE110561	#1 Switch Mounting, Link-It To The Headrest	160
☐	QUE110563	#1 Switch Mounting, Link-It To The Seat Rail Operated by knee, if two switches selected, one mounted on the right side, one on the left	160
☐	QUE110565	#1 Switch Mounting, Link-It To Chin Boom No left / right switch selection required	400
☐	QUE110569	#1 Omit Switch Mounting Switches attached loose	0

#1 Mounting Location

☐	QUE110572	#1 Switch Mounting Location, Right	0
☐	QUE110571	#1 Switch Mounting Location, Left	0

#1 Satellite Switches Only for Allround Joystick, maximum 2 selectable

☐	QUE110537	#1 Satellite Switch – Black	150
☐	QUE110538	#1 Satellite Switch – Green	150
☐	QUE110539	#1 Satellite Switch – Red	150

#2 Switch Selection

☐	QUE110631	#2 Single Switch, Ø 64 mm, Black	100
☐	QUE110632	#2 Single Switch, Ø 64 mm, Red	100
☐	QUE110633	#2 Single Switch, Ø 64 mm, Green	100
☐	QUE110088	#2 Egg Switch – Black	100
☐	QUE110089	#2 MicroLite Switch – Black	140
☐	QUE110671	#2 Single Switch, Ø 30 mm, Black	100
☐	QUE110674	#2 Single Switch, Ø 30 mm, Blue	100
☐	QUE110673	#2 Single Switch, Ø 30 mm, Green	100
☐	QUE110672	#2 Single Switch, Ø 30 mm, Red	100
☐	QUE110675	#2 Single Switch, Ø 30 mm, Yellow	100
☐	QUE110095	#2 Twister Switch – Black	140
☐	QUE110096	#2 Twister Switch – Green	140
☐	QUE110097	#2 Twister Switch – Red	140

#2 Switch Mounting

☐	QUE110098	#2 Switch Mounting, Link-It To The Hand Control Not with HD Allround Joystick or Space Disk mount	220
☐	QUE110099	#2 Switch Mounting, Link-It To The Armrest (Opposite side as the Hand Control) Armrest rail mount	130
☐	QUE110100	#2 Switch Mounting, Link-It To The Headrest	160
☐	QUE110104	#2 Switch Mounting, Link-It To The Seat Rail Knee, if 2 switches selected only left and right mount	160
☐	QUE110105	#2 Switch Mounting, Link-It To Chin Boom No left / right switch selection required	400
☐	QUE110106	#2 Omit Switch Mounting Switches attached loose	0

#2 Switch Mounting Location

☐	QUE110108	#2 Switch Mounting Location, Right	0
☐	QUE110107	#2 Switch Mounting Location, Left	0

#2 Satellite Switches Only for Allround Joystick, maximum 2 selectable

☐	QUE110111	#2 Satellite Switch – Black	150
☐	QUE110112	#2 Satellite Switch – Green	150
☐	QUE110113	#2 Satellite Switch – Red	150

Multi-Function Switches Used in combination with Assignable Button programming.

☐	QUE110550	Mini 2 Button (Stereo Plug) Attached loose, when ordered. No Link-It fixation possible	450
☐	QUE110551	Mini 2 Button (Dual Mono Plug) Attached loose, when ordered. No Link-It fixation possible	450
☐	QUE110552	Mini 4 Button InLine (DB9 Connector) Attached loose, when ordered. No Link-It fixation possible	600
☐	QUE110554	4-Way Toggle Switched Joystick	600

Accessories

☐	QUE110580	Mono Extension Cable (3.5 mm Jack Socket) (2 m)	110
☐	QUE110581	Stereo Extension Cable (3.5 mm Jack Socket) (2 m)	110
☐	QUE110582	Extension Cable DB9 (1 m)	110
☐	QUE110583	Y-Cable, 2 Female Mono To 1 Male Stereo	110

Safety Devices

☐	QUE090147	Remote Stop Switch (30 m Activation Range, Fit Via OBC Port)	750
--------------------------	-----------	--	-----

Comments for Special Build

Q700-UP M/F with SEDEO PRO ADVANCED

£

Base Price	00,000.00
Total Of Additional Equipment	
Packaging And Transport	
Exclusive VAT	
VAT	
Total Price	