



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Manual Retail Prescription Form

05 JANUARY 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Rea Azalea Max

Standard Features

- Max User Weight 180KG/28 Stone
- Crash Tested ISO 7176-Part 19
- Transit Propulsion with Attendant Hub Brakes and Black Solid Tyres
- Choice of 2 Seat Widths 550mm or 610mm (+50mm width adjustment, +100mm also with Transit)
- Choice of 2 Seat Heights 400mm or 450mm
- Depth Adjustable Seat 530mm-600mm
- Max Contour Seat, (Visco top, Wipeable), Laguna Backrest Contour Cushion (Wipeable)
- Tension and Height Adjustable Backrest (Set with space 'upwards'), Sideways Adjustable Headrest
- Removable Height Adjustable and Reinforced Armrests
- Assistant Operated Gas Piston Seat Tilt -1° - +23° and Backrest Recline to 30° with Tilt Lock
- Angle Adjustable Swingaway Elevating and Removable Legrests with Reinforced Calf Support and 1 Piece Footboard Converter
- Width Adjustable Footplate
- Solid Front Castors
- Height Adjustable Push Handles
- Pearl Silver Frame
- Trunk (Thoracic Lateral) Supports Removable and Multi Adjustable
- User Manoeuvred Tyre Brakes, Pelvic Belt (adjustable both sides)
- Anti Tippers/Tipping Lever

Standard Base price: £2,976



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code Description

BASE PRODUCT*

BSA0010	☐	REA Azalea Max	£2,976
BSA8000	☐	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE	POA

PROPULSION*				
BSA0020		Transit 16" wheels with hub brakes	(NC 0100, 0514, 0516, 0522, 0540, 0550)	◆
BSA0030		Self Propel 24"/22" wheels with user brakes	(NC 0080, 0512, 0524)	○
BSA0040		Self Propel 24"/22" wheels with assistant hub brakes	(NC 0080, 0512, 0522)	£140
SEAT WIDTH*				
BSA0050		550 mm / 22"		○
BSA0060		610 mm / 24"		○
SEAT WIDTH EXTENSION*				
BSA0070		+ 50 mm seat width pre-adjustment on armrests and legrests		○
BSA0080		+ 100 mm seat width pre-adjustment on armrests and legrests		○
BSA0085		Without seat width extension, set at 0mm adjustment but can be adjusted wider if required		○
SEAT DEPTH*				
BSA0090		530-600 mm / 21" - 24"		◆
SEAT HEIGHT*				
BSA0100		400 mm / 16"	(NC 0290, 0495, 0512, 0516)	○
BSA0110		450 mm / 18"	(NC 0480, 0514)	○
SEAT TYPE*				
BSA0120		Depth Adjustable		◆
SEAT CUSHION*				
BSA0130		Max Seat Cushion Visco top	(NC 1111)	◆
BSA0135		Without Seat Cushion, Seat Plate only	(NC 0145, 0155)	○
SEAT UPHOLSTERY*				
BSA0145		Black Plush TR35		○
BSA0155		Black Multi-stretch Black PU coated fabric		◆
BSA1111		No upholstery, select with option 0135 Without seat cushion		○
BACKREST*				
BSA0160		Backrest, Tension Adjustable - Pre-adjusted space upwards		◆
BACKREST CUSHION*				
BSA0170		Cover for tension adjustable upholstery		£178
BSA0180		Cushion, Laguna (Lateral Support)	(NC 0217)	◆
BACKREST UPHOLSTERY*				
BSA0205		Black Plush TR35		○
BSA0215		Black Multi-stretch Black PU coated fabric		◆
BSA0217		Black Spacer TR36		○
PUSH HANDLES*				
BSA0220		Push Handles braced		£44
BSA0230		High Push Bar		◆
SEAT TILT AND BACKREST REGULATION*				
BSA0240		Assistant operated gas piston seat and backrest tilt	(NC 2222)	◆
BSA0250		Electrical operated seat and backrest tilt	(NC 0450, 0605)	£616
TILT LOCK*				
BSA0605		Tilt Lock (locking of pre-set angle)		◆
BSA2222		No Tilt Lock, Select with 0250 Electric operated seat and backrest		○
LEGREST HANGER*				
BSA0260		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad and calf support Pair		◆

FOOTREST*			
BSA0280	☐	Footrest Lower part, Depth and angle adjustable with short tube and footboard converter, Pair	☐
BSA0290	☐	Footrest Lower part, Depth and angle adjustable with long tube and footboard converter, Pair	☐
BSA0291	☐	Legrest Amputee LH (NC 0293)	£137
BSA0292	☐	Legrest Amputee RH (NC 0293)	£137
HEEL STRAP			
BSA0293	☐	Heel strap	£32
ARMREST LH*			
BSA0300	☐	Reinforced, Wide Pad LH	◆
ARMREST RH*			
BSA0320	☐	Reinforced, Wide Pad RH	◆
ARMREST OPTIONS			
BSA0340	☐	Side Support Cover, PU fabric cover for the sideguards (Removable with tools) (NC 9999)	£51
BSA9999	☐	Side Support Cover with foam cushion. 2cm thick cushion supplied complete with PU fabric cover. Reduces seat width by 2cm for a comfortable fit. Consists of options 0680 and 0690 (Removable with tools)	£51
SEAT OPTIONS			
BSA0360	☐	Incontinence Cover	☐
HEADREST*			
BSA0370	☐	Headrest, Sideway Adjustable (NC 0400)	◆
BSA0390	☐	Headrest with cheek support, sideway adjustable (NC 0380, 0400)	£137
BSA0400	☐	Neckrest (NC 0380)	☐
HEADREST COVER			
BSA0380	☐	Cover for headrest	£31
TRUNK SUPPORTS LH			
BSA0413	☐	Trunk Support Multi Adjustable LH	£113
TRUNK SUPPORTS RH			
BSA0414	☐	Trunk Support Multi Adjustable RH	£107
BACKREST OPTIONS			
BSA0420	☐	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)	£42
BSA0430	☐	Lumbar Support Cushion	£31
BSA0440	☐	Wedges (One Pair)	£11
CHASSIS COLOUR*			
BSA0441	☐	Pearl Grey	◆
BSA0442	☐	Azurite Blue	☐
BSA0443	☐	Sand	☐
BSA0444	☐	Electric Green	☐
BSA0445	☐	Happy Red	☐
CHASSIS SAFETY			
BSA0450	☐	Chassis extender, castors extend 50mm forward for added stability	☐
CASTOR ATTACHMENT*			
BSA0460	☐	Fixed castor attachment	◆
FRONT CASTORS*			
BSA0480	☐	140x40 mm (5 1/2" x 1 5/8") Medium Hard Black tyre	£25
BSA0495	☐	200x45 mm (8" X 1 3/4") Black Solid tyre	◆

REAR WHEEL SIZE*			
BSA0512	☐	16" / 405 mm Transit rear wheel with fixed attachment (NC 0522, 0540, 0550)	◆
BSA0514	☐	22" / 560 mm Self Propel rear wheel with QR attachment	○
BSA0516	☐	24" / 600 mm Self Propel rear wheel with QR attachment	○
REAR WHEEL TYPE*			
BSA0522	☐	Standard self propel with colourless hub, clear anodized rim and Stainless Steel spokes	○
BSA0524	☐	Hub Brake 16" Transit or 22" or 24" Self Propel Transit with Black Plastic Hub and Rim Self Propel with colourless hub, clear anodized rim and Stainless Steel spokes	○
REAR WHEEL TYRE*			
BSA0535	☐	Solid Black tyres	◆
HANDRIM			
BSA0540	☐	Aluminium handrim	◆
SPOKEGUARD			
BSA0550	☐	Frost spokeguard	£166
BRAKE*			
BSA0560	☐	User manoeuvred tyre brakes	◆
SAFETY OPTIONS			
BSA0570	☐	Pelvic belt	◆
BSA0590	☐	Reflector Kit (wheels and chassis)	£22
TRANSPORTATION BRACKETS*			
BSA0595	☐	Transportation Brackets, mounted to the front of the wheelchair	◆
OTHER OPTIONS			
BSA0600	☐	Fixed Table Tray with elbow pads includes mounting kit	£209
BSA0620	☐	Handles Clinic Kit	£55

TOTAL PRICE			
-------------	--	--	--

Non Conforming Options:	
BSA0080	+ 100 mm seat width pre-adjustment on armrests and legrests
BSA0100	400 mm / 16"
BSA0145	Black Plush TR35
BSA0155	Black Multi-stretch Black PU coated fabric
BSA1111	No upholstery, select with option 0135 Without seat cushion
BSA0217	Black Spacer TR36
BSA0605	Tilt Lock (locking of pre-set angle)
BSA2222	No Tilt Lock, Select with 0250 Electric operated seat and backrest
BSA0290	Footrest Lower part, Depth and angle adjustable with long tube and footboard converter, Pair
BSA0293	Heel strap
BSA9999	Side Support Cover with foam cushion. 2cm thick cushion supplied complete with PU fabric cover. Reduces seat width by 2cm for a comfortable fit. Consists of options 0680 and 0690 (Removable with tools)
BSA0400	Neckrest
BSA0380	Cover for headrest
BSA0450	Chassis extender, castors extend 50mm forward for added stability
BSA0480	140x40 mm (5 1/2" x 1 5/8") Medium Hard Black tyre
BSA0495	200x45 mm (8" X 1 3/4") Black Solid tyre
BSA0512	16" / 405 mm Transit rear wheel with fixed attachment
BSA0514	22" / 560 mm Self Propel rear wheel with QR attachment
BSA0516	24" / 600 mm Self Propel rear wheel with QR attachment
BSA0522	Standard self propel with colourless hub, clear anodized rim and Stainless Steel spokes
BSA0524	Hub Brake 16" Transit or 22" or 24" Self Propel Transit with Black Plastic Hub and Rim Self Propel with colourless hub, clear anodized rim and Stainless Steel spokes
BSA0540	Aluminium handrim
BSA0550	Frost spokeguard



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):