

Customer ID: _____

ORDER FORM

ORDER ☐ QUOTE ☐

ICON 60

LIGHTWEIGHT, RIGID FRAME
WHEELCHAIR

Base model £ 1895.00

- Max user weight - 120kg
- Series 7000 aluminium rigid active wheelchair
- Multiple adjustments, seat depth, front and rear seat heights and centre of gravity
- Angle adjustable front forks
- 24" Daylight rear wheel with back anodised hand rim, QR axle and pneumatic tyre
- 4" Alu core castor
- Lockdown folding, angle, height and tension adjustable backrest with 3D breathable fabric and seat cushion.
- Low profile scissor lock brakes
- Aluminium side guards
- Calf strap and under seat zipped pouch included



SUITABLE FOR USE WITH POWERED
ASSISTED WHEELCHAIR SYSTEMS



Help to save paper and use the interactive PDF.

Please order by e-mail to ukadmin@rehasense.com. Thank's a lot.


Name:	Date:
Customer ID:	Order number:
Contact person:	Delivery address
E-Mail:	Str. & house nr:
Tel./Fax:	City & postcode:

PLEASE MARK TO SELECT AN OPTION


☐ Hereby I consent to the processing of my personal data for the purpose of the order fulfilment for a rehabilitation equipment in accordance with Article 4 Paragraph 11 of the Regulation (EU) 2016/679 of the European Parliament and of the Council of 27 April 2016 on the protection of natural persons with regard to the processing of personal data and on the free movement of such data, and repealing Directive 95/46/EC (General Data Protection Regulation)

☐ I am aware of the fact that all dimensions can be provided within the dimensional tolerance of ± 1 cm. A dimensional deviation within this tolerance range shall not constitute a reason to claim the wheelchair and deem it non-compliant with this contract.


SEAT WIDTH

tick to order	Part number	Description	Standard or cost option	Price	
☐	SW36	Seat width 36 cm	Standard	£ 0.00	
☐	SW38	Seat width 38 cm	Standard	£ 0.00	
☐	SW40	Seat width 40 cm	Standard	£ 0.00	
☐	SW42	Seat width 42 cm	Standard	£ 0.00	
☐	SW44	Seat width 44 cm	Standard	£ 0.00	
☐	SW46	Seat width 46 cm	Standard	£ 0.00	
☐	SW48	Seat width 48 cm	Standard	£ 0.00	


SEAT DEPTH

tick to order	Part number	Description	Standard or cost option	Price	Seat range adjustable between 40 cm – 49 cm 
☐	SD40	Seat depth 40 cm	Standard	£ 0.00	
☐	SD41.5	Seat depth 41.5 cm	Standard	£ 0.00	
☐	SD43	Seat depth 43 cm	Standard	£ 0.00	
☐	SD44.5	Seat depth 44.5 cm	Standard	£ 0.00	
☐	SD46	Seat depth 46 cm	Standard	£ 0.00	
☐	SD47.5	Seat depth 47.5 cm	Standard	£ 0.00	
☐	SD49	Seat depth 49 cm	Standard	£ 0.00	

FRONT SEAT HEIGHT

tick to order	Part number	Description	Standard or cost option	Price	
☐	FSH1	Seat height front 49.5 cm (4" castor wheel only)	Standard	£ 0.00	
☐	FSH2	Seat height front 50.5 cm (4" & 5" castor wheel)	Standard	£ 0.00	
☐	FSH3	Seat height front 51.5 cm (4" & 5" castor wheel)	Standard	£ 0.00	
☐	FSH4	Seat height front 52.5 cm (5" castor wheel only)	Option	£ 0.00	
☐	FSH5	Seat height front 53 cm (4" castor wheel only/jazz fork)	Option	£ 0.00	
☐	FSH6	Seat height front 54.5 cm (4" castor wheel only/jazz fork)	Option	£ 0.00	
☐	FSH7	Seat height front 56 cm (5" castor wheel only/jazz fork)	Option	£ 0.00	

REAR SEAT HEIGHT

tick to order	Part number	Description	Standard or cost option	Price	
☐	RSH1	Position 1: 38.5 cm (24" rear wheels only)	Standard	£ 0.00	
☐	RSH2	Position 2: 41 cm (24" rear wheels only)	Standard	£ 0.00	
☐	RSH3	Position 3: 44.5 cm (24" rear wheels only)	Standard	£ 0.00	
☐	RSH4	Position 4: 39.5 cm (25" rear wheels only)	Option	£ 0.00	
☐	RSH5	Position 5: 42 cm (25" rear wheels only)	Option	£ 0.00	
☐	RSH6	Position 6: 45.5 cm (25" rear wheels only)	Option	£ 0.00	

LOWER LEG LENGTH (FOOTREST POSITION) DEDICATED TO FOOTRESTS FR1 AND FR2

☐	LLL1	Position 1: 41.5 cm	Option	£ 0.00	
☐	LLL2	Position 2: 43 cm	Standard	£ 0.00	
☐	LLL3	Position 3: 44.5 cm	Option	£ 0.00	
☐	LLL4	Position 4: 46 cm	Option	£ 0.00	

LOWER LEG LENGTH (FOOTREST POSITION) ELEVATED FOOTREST FR3 REQUIRED

	LLL5	Min. value → 7 cm (Position "upside down") Required value cm enter the value Max. value ← 39 cm (Position "standard")	Option	£ 0.00	
---	------	---	--------	--------	---

SEAT POSITION (BALANCE)

☐	SPB1	Position 1 = +7.5 cm Most Active	Option	£ 0.00
☐	SPB2	Position 2 = + 6 cm	Option	£ 0.00
☐	SPB3	Position 3 = + 4.5 cm	Option	£ 0.00
☐	SPB4	Position 4 = + 3 cm	Standard	£ 0.00
☐	SPB5	Position 5 = +1.5 cm	Option	£ 0.00
☐	SPB6	Position 6 = 0	Option	£ 0.00
☐	SPB7	Position 7 = -1.5 cm	Option	£ 0.00
☐	SPB8	Position 8 = -3 cm Least Active	Option	£ 0.00

The diagram illustrates the relationship between a wheelchair and a person's lower body. A blue wheelchair is shown at the top right. A red circle highlights the seat pan area, with a red arrow pointing down to a larger, detailed view of the seat pan and footrest. This detailed view shows a blue seat pan with a row of white dots representing the seat's surface. Below the seat pan is a blue footrest with three white dots. Two vertical red lines are positioned on either side of the seat pan and footrest, and a horizontal red double-headed arrow is at the bottom, indicating the width of the seat area.


CAMBER REAR WHEELS

☐	CRW0	Camber 0°	Option	£ 0.00	
☐	CRW1	Camber 1°	Option	£ 0.00	
☐	CRW2	Camber 2°	Standard	£ 0.00	
☐	CRW3	Camber 3°	Option	£ 0.00	


BUSHING TYPE

☐	BT1	Inch 1/2"	Standard	£ 0.00	
☐	BT2	Metric 12 mm	Option	£ 0.00	

BACKREST ANGLE

☐	BRA1	Forward (from +1° to +10°) enter the value	Option	£ 0.00	
☐	BRA0	Vertical backrest (0°)	Standard	£ 0.00	
☐	BRA2	Backward (from -1° to -10°) enter the value	Option	£ 0.00	

FRAME COLOUR



☐	FC1	White	Option	£ 0.00	
☐	FC2	Black Matte	Standard	£ 0.00	
☐	FC3	Silver	Option	£ 0.00	
☐	FC4	Gold	Option	£ 40.00	
☐	FC5	Blue	Option	£ 40.00	
☐	FC6	Red	Option	£ 40.00	
☐	FC7	Green	Option	£ 40.00	

SEAT CUSHION

☐	SC1	Standard seat cushion	Standard	£ 0.00	
☐	SC0	Without seat cushion	Option	£ 0.00	

BACKREST OPTIONS

☐	BRO1	Height adjustable backrest without push handles	Standard	£ 0.00	
☐	BRO2	Height adjustable backrest fitted with height adjustable push handles	Option	£ 160.00	
☐	BRO3	Height adjustable backrest fitted with ergonomic backrest frame	Option	£ 140.00	
☐	BRO4	Height adjustable backrest fitted with ergonomic backrest frame and height adjustable push handles	Option	£ 160.00	
BACKREST AD-ONS (available with BRO1 - BRO4)					
☐	EV1	Height adjustable, standard backrest without push handles width adjustable between -7 to +7 cm	Option	£ 50.00	
☐	EV2	Height adjustable, standard backrest without push handles width adjustable between -3 to +3 cm	Option	£ 50.00	



BACKREST HEIGHT (CM)

☐ 32 standard		
☐ 33	☐ 34	☐ 35
☐ 36	☐ 37	☐ 38

BACKREST OPTIONS (not available with EV1, EV2)

☐	BRH1	Backrest with foldable pushhandles, single bended crossbar, frame height 25cm +6cm with handle (31cm)	Option	£ 100.00	 	Backrest height non-extendable, for folding handles in the backrest frame.
☐	BRH2	Backrest with foldable pushhandles, single bended crossbar, frame height 30cm +6cm with handle (36cm)	Option	£ 100.00		
☐	BRH3	Backrest with foldable pushhandles, single bended crossbar, frame height 35cm +6cm with handle (41cm)	Option	£ 100.00		
☐	BRH4	Backrest with foldable pushhandles, double bended crossbar, frame height 25cm +6cm with handle (31cm)	Option	£ 100.00	 	Backrest height non-extendable, for folding handles in the backrest frame. Crossbar curved, two radiuses
☐	BRH5	Backrest with foldable pushhandles, double bended crossbar, frame height 30cm +6cm with handle (36cm)	Option	£ 100.00		
☐	BRH6	Backrest with foldable pushhandles, double bended crossbar, frame height 35cm +6cm with handle (41cm)	Option	£ 100.00		

FOOTREST

☐	FR1	Standard footrest	Standard	£ 0.00	 FR1	 FR2	 FR3
☐	FR2	Reinforced footrest for track wheel	Option	£ 75,00			
☐	FR3	Elevated footrest	Option	£ 35,00			

FOOTREST OPTIONS

☐	FC1	Footrest cover (available for FR1 & FR3)	Option	£ 100.00			
-----------------------	-----	--	--------	----------	---	--	--

SIDE PANELS

☐	SP1	Aluminium side panel, flat, swing away	Standard	£ 0.00			
☐	SP2	Carbon side panel, flat, swing away	Option	£ 130.00			
☐	SP3	Aluminium fender with radius/removable black structured gloss	Option	£ 150.00			
☐	SP4	Aluminium fender with radius/removable black structured matte	Option	£ 150.00			
☐	SP5	Carbon fenders with radius/removable	Option	£ 220.00			
☐	SP6	Carbon fenders with radius drop-in/removable (only available with Rear Wheels 24")	Option	£ 210.00			
☐	SP7	Side panel with armrest cushion	Option	£ 150.00			

BRAKES

☐	B1	Horizontal (left, right) (only with long clamp BC2)	Option	£ 45.00	
☐	B2	Horizontal Premium Lock-P (L+R)	Standard	£ 0.00	
☐	B3	Horizontal Lock-A in aluminium	Option	£ 75.00	
☐	B4	Vertical - push to lock (only with long clamp BC2)	Option	£ 45.00	

 B1 Horizontal (left, right)	 B2 Lock-P in plastic	 B3 Lock-A in aluminium	 B4 Vertical push to lock
--	---	---	---

BRAKE MOUNTING CLAMPS

☐	BC1	Standard short clamp (only B2/B3)	Standard	£ 0.00	 BC1	 BC2
☐	BC2	Long clamp (all brakes)	Option	£ 35.00		

FORKS

☐	F1	Standard Light	Standard	£ 0.00	 F1	 F2
☐	F2	Jazz-Shock	Option	£ 330.00		

FRONT WHEELS

☐	FW1	Alucore - 4" 100x34 mm rim anodised, Silver coloured	Standard	£ 0.00	 FW1 Alucore - 4" 100x34 mm rim anodised, Silver coloured
☐	FW2	Alucore - 4" 100x34 mm rim anodised, black coloured	Standard	£ 0.00	
☐	FW3	Alucore - 5" 125x34 mm rim anodised, silver coloured	Option	£ 50.00	
☐	FW4	Alucore - 5" 125x34 mm rim anodised, black coloured	Option	£ 50.00	
☐	FW5	Tristar - 4" 100x34 mm rim anodised, silver coloured	Option	£ 50.00	
☐	FW6	Tristar - 5" 125x34 mm rim anodised, silver coloured	Option	£ 50.00	
☐	FW7	Lotus fibercore - 4" 100x34 mm rim plastic, black	Option	£ 55.00	
☐	FW8	Lotus fibercore - 5" 125x34 mm rim plastic, black	Option	£ 55.00	

 FW2 Alucore - 4", 100 X 34 mm, rim anodised, black	 FW3 Alucore - 5", 125 X 34 mm rim anodised, silver	 FW4 Alucore - 5", 125 X 34 mm rim anodised, black	 FW5 Tristar - 4" 100x34 mm rim anodised, silver	 FW6 Tristar - 5" 125x34 mm rim anodised, silver	 FW7 Lotus fibercore - 4" 100x34 mm rim plastic, black	 FW8 Lotus fibercore - 5" 125x34 mm rim plastic, black
--	--	---	---	--	---	---

REAR WHEELS

	RW1	DAYLITE - 24" x 1" wheel (Tyre: T1, T2, T3: Pushrim: PR1, PR2, PR3, PR4, PR5, PR6, PR7, PR8)	Standard	£ 0.00	 RW1 DAYLITE 930S 24" wheel/2 pcs.
	RW2	SKYLITE - 24" x 1" wheel (Tyres: T1, T2, T3: Pushrim: PR1, PR2, PR3, PR4, PR5, PR6, PR7, PR8)	Option	£ 220.00	
	RW3	TENSIO RAY 12S - 24" wheel (Tyres: T1, T2, T3: Pushrim: PR1, PR2, PR3, PR4, PR5, PR6, PR7)	Option	£ 550.00	
	RW4	HURRICANE 24" x 1" wheel (Tyres: T1, T2, T3: Pushrim: PR1, PR2, PR3, PR4, PR5, PR6, PR7, PR8)	Option	£ 410.00	
	RW5	OFF ROAD 24" x 2" buddy wheel (Tyres: T4; Pushrim: PR1, PR7)	Option	£ 550.00	
	RW6	SMARTLITE 24" wheel	Option	£ 220.00	
	RW7	CYCLONE 24" wheel (Tyres: T3: Pushrim: PR1)	Option	£ 400.00	
	RW8	DAYLITE - 25" wheel (Tyres: T3, T2, T3: Pushrim: PR1, PR5, PR8)	Option	£ 220.00	
	RW9	SKYLITE - 25" wheel (Tyres: T3, T2, T3: Pushrim: PR1, PR5, PR8)	Option	£ 240.00	



PUSH RIMS

	PR1	Anodised aluminium, Black - coloured	Standard	£ 0.00	 PR1 Anodised aluminium, Black - coloured
	PR2	Anodised aluminium, Silver - coloured	Option	£ 0.00	
	PR3	Steel	Option	£ 70.00	
	PR4	Titanium	Option	£ 175.00	
	PR5	Tetra Optimum F	Option	£ 200.00	
	PR7	Tetra silicone, Black	Option	£ 85.00	
	PR8	Nova H	Option	£ 210.00	



TYRES

☐	T1	Speedlite 24" x 1"	Standard	£ 0.00	 T1 Speedlite 24"
☐	T2	Gazelle 24" x 1"	Option	£ 50.00	
☐	T3	Trailblazer 24" x 1"	Option	£ 80.00	
☐	T4	Schwalbe Big Apple (only to be used with 24" Buddy)	Option	£ 0.00	



QR AXLES

☐	QR1	Quick release axle / 2pcs	Standard	£ 0.00	 QR1	 QR2
☐	QR2	Quick release tetra handle axle / 2pc	Option	£ 80.00		

OPTIONS

☐	O1	Anti-Tip Left	Option	£ 120.00	
☐	O2	Anti-Tip Right	Option	£ 120.00	
☐	O3	Tip-on Left	Option	*available Febr 2025	
☐	O4	Tip-on Right	Option	*available Febr 2025	

TRANSPORT CLAMPS

☐	TC1	Transport clamps for vehicle	Option	£ 120.00	
-----------------------	-----	------------------------------	--------	----------	---

TRACK WHEEL

☐	STWA	Single arm track wheel Aluminium	Option	£ 399.00	 STWA, STWC	 DTWA, DTWC
☐	DTWA	Double arm track wheel Aluminium	Option	£ 499.00		
☐	STWC	Single arm track wheel Carbon	Option	£ 499.00		
☐	DTWC	Double arm track wheel Carbon	Option	£ 599.00		

Your notes