

# Robustly Rugged Perfectly Convenient



**REVO 2.0**



Robustly rugged meets brilliant convenience. With the Revo® 2.0, there is nothing you can't handle! This tough and dependable scooter has powerful suspension for a smooth ride, yet easily disassembles for full portability.

## Specifications

Weight Capacity	26 st 11 lbs (170 kg)
Maximum Speed <sup>1</sup>	Up to 4 mph
Ground Clearance <sup>3</sup>	32 mm (1.25") at motor
Turning Radius <sup>3</sup>	1321 mm (52")
Overall Length <sup>3</sup>	1194 mm (47")
Overall Width <sup>3</sup>	527 mm (20.75")
Max. Obstacle Climb	63.5 mm (2.5")
Drive Tyres	Front & rear 228.5 mm (9") - solid or pneumatic
Suspension	Enhanced front and rear suspension
Total Weight Without Batteries <sup>3, 4</sup>	60 kg (132 lbs)
Heaviest Piece When Disassembled <sup>3</sup>	25.85 kg (56.9 lbs)
Drivetrain	24 volt DC motor
Dual Breaking System	Regenerative and electromechanical
Battery Requirements <sup>3</sup>	2 x 30 Ah batteries
Battery Weight <sup>4</sup>	11.5 kg (25 lbs)
Battery Charger	3.5 Amp off-board
Colour	Grey Street, True Blue

## Features

- Comfort Trac front and rear suspension
- Full LED lights package
- Solid alloy style wheels
- Underseat storage
- Large tiller basket
- Easily disassembled into five pieces

Available in:

- Grey Street
- True Blue



\*Solid tyres – available in blue only  
Pneumatic tyres – available in blue and grey

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.

<sup>2</sup> Includes standard seat weight.

<sup>3</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

<sup>4</sup> Battery weight may vary based on battery manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.