

Robby Bathing Aid

Simple and practical. Thanks to the low centre of gravity, the Robby Bathing Support needs comparatively little water in the bath. It can be folded up to a compact size allowing it to be stored easily. The compact size also makes it easy to take with you where needed.

Combining design with function, its four sets of ratchet joints and multiple positioning accessories support a wide range of sitting and lying positions for a safe and fun bath time.

What is special about me?

- What is special about me?
- Secure support while bathing
- Low body centre of gravity for optimum use of the bath
- Various sitting and lying positions

- Cover adaptable to the body shape with hook and loop straps
- Seat surface and back length grow with the child thanks to telescoping tubes and continuously adjustable supports

When and where am I used?

- Bathing at home or in facilities
- By the therapy pool or swimming pool as a means of transport
- Out and about at a family picnic or in the garden for temporary positioning

Technical Data

Dimensions (cm) and weight (kg)	Size 1	Size 2
Max user weight	30	60
Seat depth (depending on strap tension and expansion of telescoping tubes)	20 – 24	31 – 43
Backrest height	63 – 72	89 – 94
Total width	38	38
Total length (support surface)	115 – 133	147 – 173
Total weight (approx.)	3.8	4.3

Quality for life

