

| | |
|---|--------------------------------|
| <p>Base Price -</p> <p>Valid 01.04.20 to 31.03.21 Please ensure you are using the latest order form</p> | <p>Order no.:</p> <p>-----</p> |
|---|--------------------------------|

Robby Bathing support



Standard equipment

- Folding bathing and lying support that can be adapted to grow with the user
- Net upholstery (blue)
- Lap belt and leg strap
- Frame colour: silver

| | |
|--|--------------------------------|
| <input type="checkbox"/> Provide estimate only | <input type="checkbox"/> Order |
| Billing address/customer number | |
| Delivery address | |
| Sales person | Name |
| | |
| | Date Signature |
| <p>By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.</p> | |

| Robby bathing support | | | |
|-----------------------|--------------------------|--|---------|
| HR51040000-027 | <input type="checkbox"/> | Robby, size 1 (seat depth 20 – 24 cm, back support height incl. head section 63 – 72.5 cm) Bathing and lying support, includes washable net upholstery, lap belt, leg straps; max. load capacity 30 kg | £643.00 |
| HR51050000-027 | <input type="checkbox"/> | Robby, size 2 (seat depth 31 – 42.5 cm, back support height incl. head section 89 – 94 cm) Bathing and lying support, includes washable net upholstery, lap belt, leg straps; max. load capacity 60 kg | £682.00 |

| Accessories | | | |
|----------------|--------------------------|---|--------|
| HR51045300 | <input type="checkbox"/> | Trunk belt | £35.38 |
| HR51045200-103 | <input type="checkbox"/> | Adjustable head support pads (can also be used as trunk support pads) | £97.78 |

Information



Adjustable head support pads

Technical Data

| | Size 1 | Size 2 |
|--|----------------|----------------|
| Seat depth (depending on strap tension and expansion of telescoping tubes) | 200 – 240 mm | 310 – 425 mm |
| Back support height | 630 – 725 mm | 890 – 940 mm |
| Overall width | 360 mm | 380 mm |
| Overall length (support surface) | 1150 – 1330 mm | 1470 – 1730 mm |
| Max. load capacity | 30 kg | 60 kg |
| Weight | 3.8 kg | 4.3 kg |