



Bariatric Shower Chair Clean

Strongly engineered and functionally designed exclusively for bariatric users!

A combination chair with toilet opening. When the seat cover is removed, the chair becomes a bedside commode

We have designed the armrests to bend out like a butterfly, so the body doesn't get squeezed, getting in and out the chair. The butterfly armrests of the Shower Chair also provide much better functionality for the user, compared to other chairs with straight and short armrests where it is almost impossible for a bariatric user to get a good grip.

The height of the chair is adjustable from 44 to 59 cm. It is very important to adjust the chair to the right height to have success getting in and out of the chair. If the chair is too high for the user, the feet will lift off the ground, making it hard to get seated in the chair. If the chair is too low, the user will have trouble getting up.

It is often a problem for bariatric users to be seated comfortably. Sometimes the seat depth is too short, sometimes too deep and with no back support. We have designed the seat to be individually adjusted to provide a comfortable seating. The backrest can be moved +/- 6 cm and the bended tubing assures that the right comfort can be found.

It is easy to clean the Shower Chair Clean with its removable seat pad. No tools are needed and the bedpan is redesigned for a smoothly hygienic cleaning operation.

It is also made easy to both handle and move the chair from place to place. In the construction of the chair we have only used high-tech lightweight tubing, to provide the highest strength at the lowest weight.

The commode bucket is included when purchasing the chair.



Seat with lid



Split toilet seat

Shower Chair Clean	
Maximum weight capacity is	325 kg
Seat width	61 / 71 cm
Total length	80-101 cm
Seat depth	40-52 cm
Seat height	44-59 cm
Backrest height	37 cm
Distance from floor to seat underside	33-48 cm
Total outer width of front legs	69 / 79 cm
Total height	79-94 cm
Total weight	15 / 16 kg
Total outer width on armrests	73 / 83 cm
Seat / back / armrests	Polyurethan
Frame	Powder-lacquered chromium molybdenum steel
ISO - classification	ISO 09 33 03

MAX 325 kg

- Adjustable seat height
- Adjustable seat depth
- Long stable armrests
- Two seat widths
- Combined bedside commode and shower bench

Ideal use for...

- Bariatric users.

Model	Seat width cm	Seat depth cm	Order No
XXL-Rehab Clean	61	40-52	0142-061-000
XXL-Rehab Clean	71	40-52	0142-071-000
Accessories			
Commode bucket			0140-999-001

