

## Skippi / Skippi<sup>Plus</sup>

The Skippi power wheelchair supports a child's natural urge to move. It is individually adjustable and, with the enAble40 control unit, offers the ability to operate two power options.

The Skippi<sup>Plus</sup> offers additional functions with the enAble50, such as additional control options as well as special and environmental controls. The Skippi<sup>Plus</sup> equipped with power seat tilt and power seat height adjustment makes everyday life easier.

### What is special about me?

- Compact and manoeuvrable
- Easy to transport due to its modular design and low weight of components

### Standard Model

- Standard seat with back bag
- Seat cushion
- Telescoping side panels
- Drive wheels with aluminium rim, black
- Wheel cover for drive wheels
- 8" PU caster wheel
- Segmented telescoping footrests

- 28 Ah (C5) AGM batteries
- Battery charger
- **Skippi: enAble40 control unit** with the ability to control 2 power options (choice of: power back angle adjustment, power seat tilt and power seat height adjustment)
- **Skippi<sup>Plus</sup>: enAble50 control unit** with the ability to control 3 power options (choice same as Skippi as well as special and environmental controls)
- Frame colours: depending on frame colour combination, upholstery in red/black, blue/black

### Options

- Push handles, height-adjustable
- Attendant control
- Power back adjustment
- Power seat tilt
- Electric seat height adjustment (Skippi<sup>Plus</sup>)
- Bumper bar, detachable
- Manual horn
- Kit for transport by car

Quality for life



### Technical Data

Speed	6 km/h
Battery capacity	28 Ah (C5) AGM batteries
Range under test conditions approx.	25 km
Control unit	Skippi (enAble40) Skippi <sup>Plus</sup> (enable50)
Max user weight	50kg
Obstacle height the wheelchair is able to overcome (optional)	50mm
Climbing ability	12%
Length including footplates	850mm
Width	570mm
Turning Radius	520mm
Seat height (without cushion) at 3° seat tilt	430mm
Seat width	300 – 380mm
Seat depth	280 – 380mm
Back angle	-9 to + 30°
Seat tilt	25° mechanical / 20° powered
Back height	340 / 460mm
Armrest height	170 – 230mm
Lower leg length	150 – 380mm
Weight when empty	62kg
Charging time	12hr
Skippi <sup>Plus</sup> with electric seat height adjustment	200mm
Seat inclination	0° / 3° / 6°