



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

Spark

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

☐ Quote ☐ PO Number: _____

City: _____ Postal Code: _____

RTS: _____

Phone: _____

Marked For: _____

Email/Fax: _____

Select if chair quote has: ☐ Additional Ki Products with Different Quote/Order Number **Quote/Order Number:** _____

☐ **Spark - Transit Included** £2100

Effective 18/07/2022

165 lb weight capacity. Not available with upholstery.

☐ **Spark - Transit Not Included** £2100

165 lb weight capacity. Not available with upholstery.

Frame

Frame Width / Frame Depth

Select box corresponding to desired Frame Width and Depth.
See references at right.

		Frame Width						
		10"	11"	12"	13"	14"	15"	16"
Frame Depth	12"	☐	☐	☐	☐	☐	☐	☐
	13"	☐	☐	☐	☐	☐	☐	☐
	14"	☐	☐	☐	☐	☐	☐	☐
	15"	☐	☐	☐	☐	☐	☐	☐
	16"	☐	☐	☐	☐	☐	☐	☐
	17"	☐	☐	☐	☐	☐	☐	☐
	18"	☐	☐	☐	☐	☐	☐	☐



Frame Width

Spark is set up to allow 2" increase in frame width from initial frame width in 1" increments. Growth is achieved by expanding cross braces. 16" maximum frame width. 10" Frame Width expands only 1".



Frame Depth

Spark is set up to allow 4" increase in frame depth from initial frame depth in 1" increments. 19" maximum frame depth.

Seat Pan

☐ Removable Aluminum Seat Pan £210
Not available with Transit System. Pan Depth 1/2" shorter than Frame Depth.

Backrest

Back Post / Height Range / Back Height

Select Back Height in row corresponding to desired Back Post and Height Range. Both Back Post options are folding and angle adjustable in 4° increments.

Back Post	Price	Height Range	Back Height						
Stroller Handle Back Post	N/C	Short	☐ 20						
		Tall	☐ 24						
Height Adjustable Straight with Push Handle	£63	Short ¹	☐ 13	☐ 14	☐ 15	☐ 16	☐ 17		
		Medium	☐ 16	☐ 17	☐ 18	☐ 19	☐ 20		
		Tall	☐ 20 ²	☐ 21	☐ 22	☐ 23	☐ 24		

¹ Not available with Seating Dynamics Rocker Back

² Not available with Angle Adjustable Locking Extendable Flip Up Armrest

Back Post Options

- ☐ Removable Stroller Handle Extensions £105
Not available with Height Adjustable Straight with Push Handle Back Post
- ☐ Seating Dynamics Rocker Back £1680
Available only with T-Arms. Not available with 12" Frame Depths with Rear Frame Heights lower than 14.5".

Footrest

Hanger / Front Frame Bend

- ☐ 60° _____ N/C
Available only with Extension Tube Hangers
- ☐ 70° _____ N/C
- ☐ 80° _____ N/C
Not available with 7" or 8" Casters
- ☐ 90° _____ N/C
Available only with Front Mount Hangers. Not available with 7" or 8" Casters.
- ☐ Pediatric Pro Elevating Leg Rests (Peds PRO ELR) _____ £252
- ☐ Pro Elevating Leg Rests (PRO ELR) _____ £252
Not available with Seat Widths less than 14"
- ☐ Omit Hangers _____ N/C

PRO ELR Cover

Available only on PRO ELR

- ☐ Gel Pad Covers (Pair) _____ £264

Footrest / Extension / Seat to Footrest Length

Select Seat to Footrest Length in row corresponding to desired Footrest and Extension. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements are based on 70° Hanger / Front Frame Bend.

Footrest	Extension	Seat to Footrest Length												
Extension Tube Hangers (60°, 70°, 80°)	Short	☐ 11	☐ 11.5	☐ 12	☐ 12.5	☐ 13	☐ 13.5	☐ 14	☐ 14.5	☐ 15				
	Medium	☐ 13	☐ 13.5	☐ 14	☐ 14.5	☐ 15	☐ 15.5	☐ 16	☐ 16.5	☐ 17				
	Long	☐ 15.5	☐ 16	☐ 16.5	☐ 17	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5				
Front Mount Hangers ¹		Standard Selections										Requires 3" Extension for £46		
	Short (70°, 80°, 90°)	☐ 4 ⁵	☐ 4.5 ⁵	☐ 5	☐ 5.5	☐ 6	☐ 6.5				☐ 10	☐ 10.5	☐ 11	
		☐ 7	☐ 7.5	☐ 8	☐ 8.5	☐ 9	☐ 9.5				☐ 11.5	☐ 12	☐ 12.5	
	Medium (70°, 80°, 90°)	☐ 4 ⁵	☐ 4.5 ⁵	☐ 5	☐ 5.5	☐ 6	☐ 6.5	☐ 7	☐ 7.5		☐ 12	☐ 12.5	☐ 13	
		☐ 8	☐ 8.5	☐ 9	☐ 9.5	☐ 10	☐ 10.5	☐ 11	☐ 11.5		☐ 13.5	☐ 14	☐ 14.5	
	Long (70°, 80°)	☐ 5	☐ 5.5	☐ 6	☐ 6.5	☐ 7	☐ 7.5	☐ 8	☐ 8.5	☐ 9	☐ 14	☐ 14.5	☐ 15	
☐ 9.5		☐ 10	☐ 10.5	☐ 11	☐ 11.5	☐ 12	☐ 12.5	☐ 13	☐ 13.5	☐ 15.5	☐ 16	☐ 16.5		
Peds ² PRO ELR	Low Pivot	☐ 7.5 ³	☐ 8 ³	☐ 8.5 ³	☐ 9 ⁴	☐ 9.5 ⁴	☐ 10	☐ 10.5	☐ 11	☐ 11.5	☐ 12	☐ 12.5	☐ 13	
	High Pivot	☐ 7.5 ³	☐ 8 ³	☐ 8.5 ³	☐ 9 ⁴	☐ 9.5 ⁴	☐ 10	☐ 10.5	☐ 11	☐ 11.5				
PRO ELR Low Pivot Setting	Short	☐ 12 - 17.5												
	Long	☐ 15 - 21												
PRO ELR High Pivot Setting	Short	☐ 10.5 - 15.5												
	Long	☐ 13 - 19												

¹ See Minimum Front Seat Height chart below.

² Comes standard with 3" Extension

³ Not available on chairs 13" wide or more

⁴ Not available on chairs 15" wide or more

⁵ Available with 90° only

3" Extension

Available only with Front Mount Hangers. Provided as parts.

- ☐ 3" Extension for Front Mount Hanger _____ £42

Hanger Release



If no selection is made Classic Hanger Release will be used

- ☐ Classic _____ Std
- ☐ 4-Way - In line _____ N/C
See reference A to the right. Not available with PRO ELR or Peds PRO ELR.
- ☐ 4-Way - Out _____ N/C
See reference B to the right. Not available with PRO ELR or Peds PRO ELR.

Footplate

- ☐ Composite _____ N/C
Not available with Front Mount Hanger, Frame Widths less than 14"
- ☐ Composite Angle Adjustable _____ £158
Not available with Front Mount Hanger or Frame Widths less than 15"
- ☐ Aluminum Angle Adjustable _____ N/C
- ☐ Aluminum Locking Multi-Angle Adjustable _____ £263
Available only with Front Mount Hanger
- ☐ One Piece Flip Up Angle Adjustable _____ £163
Not available with Front Mount Hanger, PRO ELR or Frame Widths less than 14"

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13"	13.5"	14"	15"	16"	17"
Minimum Front Seat Height with 3" Extension	14.5"	16"	16.5"	18"	18.5"	20"



Footrest continued

Residual Limb Support

Residual Limb Support

- ☐ Left _____ £147 ☐ Right _____ £147
☐ Omit Left Hanger _ N/C ☐ Omit Right Hanger _ N/C

Cover

- ☐ Smooth ☐ Air Mesh ☐ Smooth ☐ Air Mesh

Pad Size

- ☐ 5 x 8 ☐ 9 x 14 ☐ 5 x 8 ☐ 9 x 14
☐ 7 x 10 ☐ 11 x 10 ☐ 7 x 10 ☐ 11 x 10
☐ 9 x 10 ☐ 11 x 14 ☐ 9 x 10 ☐ 11 x 14

Heel Loop

Not available with One Piece Flip Up Angle Adjustable Footplate

- ☐ Heel Loops _____ £32
Not available with Frame Widths less than 14"
☐ Adjustable Heel Loops _____ £53
☐ Omit Heel Loops _____ N/C

Calf Strap

- ☐ Velcro Adjustable _____ £32
☐ Padded Velcro Adjustable _____ £58
☐ Bodypoint Padded AeroMesh _____ £89

Front Frame Height and Casters

Caster Size / Front Frame Type / Front Frame Height

Select Front Frame Height in row corresponding to desired caster size. Selections in first three columns require Super Low Front Frame.

Caster Size	Front Frame Height													
	Requires Super Low Front Frame ¹													
3"	☐ 13.5	☐ 14		☐ 14.5	☐ 15	☐ 15.5	☐ 16	☐ 16.5	☐ 17	☐ 17.5	☐ 18	☐ 18.5		
4"	☐ 13.5	☐ 14	☐ 14.5	☐ 15	☐ 15.5	☐ 16	☐ 16.5	☐ 17	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20
5"	☐ 14.5	☐ 15	☐ 15.5	☐ 16	☐ 16.5	☐ 17	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20	☐ 20.5	
6"	☐ 15.5	☐ 16	☐ 16.5	☐ 17	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20	☐ 20.5			
7"	☐ 17	☐ 17.5		☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20	☐ 20.5					
8"	☐ 17.5	☐ 18	☐ 18.5	☐ 19	☐ 19.5	☐ 20	☐ 20.5							

¹ Price is £40

Caster Type

- ☐ .75" Rollerblade _____ N/C
Available only with 3" and 4" Caster Sizes
☐ .75" Lighted Rollerblade _____ £74
Available only with 3", 4" and 5" Caster Sizes
☐ 1" Poly _____ N/C
Not available with 3" Caster Size
☐ 1" Poly Aluminum _____ £53
Available only with 4" and 5" Caster Sizes
☐ 1" Pneumatic _____ N/C
Available only with 6" and 8" Caster Sizes
☐ 1.5" Poly _____ N/C
Available only with 4", 5" and 6" Caster Sizes
☐ 1.5" Soft Roll Aluminum _____ £53
Not available with 7" and 8" Caster Sizes
☐ 2" Poly _____ N/C
Available only with 6" Caster Size
☐ 2" Pneumatic w/ Foam Insert _____ £189
Available only with 8" Caster Size

Caster Housing

- ☐ Standard _____ N/C
☐ Performance _____ £32
Not available with 350 lb Weight Capacity

Maximum Seat Dump Reference

(Front Seat Height minus Rear Seat Height)

Seat Depth	12"	13"	14"	15"	16"	17"	18"	19"
Maximum Dump Standard Caster Housing	1.5"	1.5"	1.5"	2"	2"	2"	2.5"	2.5"
Maximum Dump Performance Caster Housing	2.5"	2.5"	3"	3"	3"	3.5"	3.5"	4"

Fork Stem

- ☐ Best Fit _____ Std
Ki Mobility will pick Fork and Stem combination to optimize performance based on Caster Size and Front Seat Height
☐ Standard Fork Stem _____ N/C
☐ 3/4" Longer Fork Stem _____ N/C
☐ 1 1/2" Longer Fork Stem _____ N/C

Fork Options

- ☐ Frog Legs Suspension Fork _____ £299
Available only with 3", 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 2" Poly Caster Type or 2" Pneumatic w/ Foam Insert Caster Type.

End User Weight Required: _____ lbs

- ☐ Caster Pin Locks _____ £168
Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork

Wheel Locks

Wheel Lock		
☐	Push to Lock	Std
☐	Pull to Lock	N/C
☐	Omit Wheel Locks	N/C

Wheel Lock Options		
☐	Extension Handles	£42
<i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>		
☐	Grade Aids	£126
<i>Available only with Push to Lock Wheel Locks and Pneumatic and Pneumatic with Airless Insert tires</i>		

Rear Frame Height and Wheels

Wheel Size / Tire Profile / Rear Frame Height

Select Rear Frame Height in row corresponding to Wheel Size and Tire Profile. 12"-14" Frame Depths not available with 24" wheel size.

☐ Requires Universal Axle Plate ☐ Any axle plate can be used (Universal, Amputee, One Arm Drive)

Wheel Size	Tire Profile	Rear Frame Height											
18"	Low	☐ 11	☐ 11.5 ¹	☐ 12	☐ 12.5 ¹	☐ 13	☐ 13.5	☐ 14	☐ 14.5	☐ 15	☐ 15.5	☐ 16	
	Full	☐ 11.5	☐ 12 ¹	☐ 12.5	☐ 13 ¹	☐ 13.5	☐ 14	☐ 14.5	☐ 15	☐ 15.5	☐ 16	☐ 16.5	
20"	Low	☐ 11.5	☐ 12 ¹	☐ 12.5	☐ 13 ¹	☐ 13.5	☐ 14	☐ 14.5	☐ 15	☐ 15.5	☐ 16	☐ 16.5	
	Full	☐ 12	☐ 12.5 ¹	☐ 13	☐ 13.5 ¹	☐ 14	☐ 14.5	☐ 15	☐ 15.5	☐ 16	☐ 16.5	☐ 17	
22"	Low	☐ 12.5	☐ 13 ¹	☐ 13.5	☐ 14 ¹	☐ 14.5	☐ 15	☐ 15.5	☐ 16	☐ 16.5	☐ 17	☐ 17.5	
	Full	☐ 13	☐ 13.5 ¹	☐ 14	☐ 14.5 ¹	☐ 15	☐ 15.5	☐ 16	☐ 16.5	☐ 17	☐ 17.5	☐ 18	
24"	Low	☐ 13.5	☐ 14 ¹	☐ 14.5	☐ 15 ¹	☐ 15.5	☐ 16	☐ 16.5	☐ 17	☐ 17.5	☐ 18	☐ 18.5	
	Full	☐ 14	☐ 14.5 ¹	☐ 15	☐ 15.5 ¹	☐ 16	☐ 16.5	☐ 17	☐ 17.5	☐ 18	☐ 18.5	☐ 19	

¹ Amputee Axle Plate not available with lowest two Rear Frame Height positions

One Arm Drive

One Arm Drive £735

Available only with 20", 22" and 24" spoke wheels and Quick Release Axles.
Spoke wheel style is exclusive to OAD. Not available with 4 degrees of camber.
Standard handrim will be Aluminum Anodized.

Drive Handrim Side: ☐ Left ☐ Right

Drive Handrim Type: ☐ Aluminum N/C
☐ Poly Coated £168

Axle Plate

No selection needed with One Arm Drive

☐ Universal Std
☐ Amputee £147

Rear Wheel Type

No selection needed with One Arm Drive

☐ MAXX Performance Spoke N/C
☐ Mag (Black) £105
Not available in 18" wheel size

Colored Mag £158
Available only in 20" and 22" wheel sizes

Color: ☐ Florescent Red ☐ White ☐ Blue ☐ Pink

Spinergy Spox £473
Available only in 22" and 24" wheel sizes

☐ Black Spokes N/C
☐ Colored Spokes £79
Color: ☐ Red ☐ White ☐ Blue ☐ Yellow

Rear Tire Type

Full Profile

☐ Street Pneumatic 1 3/8" Std
☐ Full Poly 1 3/8" N/C
☐ Pneumatic with Airless Insert 1 3/8" £95

Rear Tire Type continued

Low Profile

☐ Low Poly 1" N/C
☐ Sentinel Solid Tire 1" £95
Not available on 18" wheel size or 20" and 22" Mag wheels
☐ High Pressure 1" N/C
Not available on Mag wheel or 18" wheel size

Handrim

No selection needed with One Arm Drive

Standard Handrims

☐ Aluminum Anodized Std
☐ Plastic Coated £105
☐ Projection (8 Vertical Projections) £247
Available only with full profile tires in combination with MAXX Performance, 22"-24" Spoke, or 24" Mag wheel
☐ Omit Handrim N/C

Natural Fit Handrims

Not available on 18" or 20" wheel sizes or 22" Mag wheel. Not recommended with Mag wheel.

☐ Natural Fit LT Regular Grip £236
☐ Natural Fit Standard Regular Grip £236
☐ Natural Fit LT Super Grip £236
☐ Natural Fit Standard Super Grip £236
☐ Natural Fit LT Omit Grip £236
☐ Natural Fit Standard Omit Grip £236

Handrim Options

☐ 3M Micro Grip £79
Not available on Plastic Coated or Projection Handrims

Handrim Tabs:

Available only with Spinergy, Mag or Maxx Performance Spoke and Aluminum Anodized or Plastic Coated Handrims

☐ Short Tab N/C
☐ Long Tab N/C

Axles


Axle		Axle Position
No selection needed with One Arm Drive		
☐ Quick Release	Std	Rear Wheel Camber: ☐ 0° ☐ 2° ☐ 4°
☐ Quad Quick Release	£110	Center of Gravity: ☐ Rear ☐ Center ☐ Forward
<i>Not available with Mag or Spinergy wheels</i>		<i>No selection needed with Amputee Axle Plate</i>

Armrests

Armrest Please select Arm Pad option		Arm Pad
No selection needed with One Arm Drive		No selection needed with One Arm Drive
☐ Height Adjustable T-Arm Tall (11.5" - 15") New!	£223	Not available with Tubular Flip Up Armrests. Full length Arm Pads not available with Pediatric Height Adjustable T-Arms.
<i>Not available with Frame Depths less than 14"</i>		
☐ Height Adjustable T-Arm (8.5" - 13.5")	£223	Standard Arm Pad _____ N/C ☐ Desk ☐ Full
<i>Not available with Frame Depths less than 14"</i>		Waterfall Arm Pad _____ £21 ☐ Desk ☐ Full
☐ Height Adjustable T-Arm Low (6.5" - 10.5")	£223	<i>Not available with Angle Adjustable Locking Extendable Flip Up Armrest</i>
<i>Not available with Frame Depths less than 14". If requested configuration is unavailable, Standard Height T-Arm will be used.</i>		Foam Grip _____ N/C ☐ Desk ☐ Full
☐ Pediatric Height Adjustable T-Arm Desk (6.5" - 10")	£164	<i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest</i>
☐ Angle Adjustable Locking Extendable Flip Up	£223	
Armrest continued		
☐ Tubular Flip Up		N/C
<i>Not available with Arm Pads</i>		
☐ Omit Armrests		N/C

Accessories

Removable Side Guard		Other Accessories
No selection needed with One Arm Drive		☐ Tool Kit _____ N/C
No selection needed with One Arm Drive		☐ Removable Underseat Pouch _____ £42
<i>Not available with Height Adjustable T-Arm Armrests</i>		<i>Not available with 10" - 13" Seat Widths</i>
☐ Adult Composite Side Guards (14" - 16")	N/C	☐ Pediatric Ki Mobility Backpack _____ £63
☐ Pediatric Composite Side Guards (10" - 14")	N/C	☐ Ki Mobility Backpack _____ £63
Anti-Tip		☐ Black Neoprene Impact Guards _____ £74
☐ Rear Anti-Tips _____ £184		☐ Cane and Crutch Holder _____ £79
Pelvic Positioning Belt		☐ Canopy with Hardware _____ £236
☐ 1" Auto Buckle _____ £59		<i>See image below</i>
☐ 1" Padded Auto Buckle _____ £131		
☐ 1 1/2" Auto Buckle _____ £59		
☐ 1 1/2" Padded Auto Buckle _____ £131		
☐ 2" Airline Buckle _____ £95		
Bodypoint Evoflex - Pelvic Stabilizer _____ £163		
Size: ☐ Small (16-24") ☐ Medium (19-29") ☐ Large (25-37")		
Bodypoint Belt Mounting Kit _____ £89		
Style: ☐ Bracket ☐ Band Clamp 1"		
Spoke Guards		
☐ Clear Spoke Guards _____ £131		
<i>Not available with Mag wheel</i>		
☐ Black Spoke Guards _____ £131		
<i>Not available with 18" or 20" wheel sizes or Mag wheel</i>		
Designer Spoke Guards _____ £158		
<i>Available only with 20" and 22" wheel sizes. Not available with Mag wheel.</i>		
Design:		
☐ Blue Monster ☐ Fairy Garden ☐ Havin' a Hoot		
☐ Robo Fun ☐ Rockin' Reptiles ☐ Super Puppies		
☐ Tutu Cute ☐ Tuxedo Time		



Frame Color

Gloss Finish

☐ Black	N/C
☐ Blue Skies	N/C
☐ Bubblegum Pink	N/C
☐ Candy Blue	N/C
☐ Candy Purple	N/C
☐ Candy Red	N/C
☐ Copper	N/C
☐ Granny Smith	N/C
☐ Metallic Orange	N/C
☐ Yellow	N/C

Matte Finish

☐ Charcoal Gray	N/C
☐ Emerald Green	N/C
☐ Iridescent Burgundy	N/C
☐ Iridescent Navy	N/C
☐ Military Green	N/C
☐ Retro Red	N/C
☐ Royal Blue	N/C
☐ Shadow Black	N/C
☐ Smooth White	N/C

Hydrographics

Visit www.kimobility.com for current Hydrographic options. Note selection below

☐ _____ £630

Spark Front Frame Height Matrix																							
3" Caster				4" Caster				5" Caster				6" Caster				7" Caster				8" Caster			
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork
14.5	13.5	Std	4"	15	13.5	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	16.5	Std	6"	19	17.5	Std	6"
15	14			15.5	14			16.5	15.5	3/4		17.5	16.5	3/4		18.5	17			19.5	18.5	3/4	
15.5	14	3/4		15.5 ¹	14 ¹	3/4		17.5	16	1 1/2	18.5	17	1 1/2	18.5	17.5	3/4	20.5	n/a		1 1/2			
16	14.5			16	14.5			Std	17	15.5	3/4	17	15.5	Std	19		18	19		17.5	Std	7"	
16	15	1 1/2		16.5 ¹	15 ¹	1 1/2			16.5	15		Std	17.5		16	Std	19.5	18		19.5			18
16.5	15.5			17	15.5			16.5	15.5	3/4	18		16.5	5"	18		16.5	20		18.5			20
15.5	14.5	Std	15.5	14	Std	17	16	3/4	17.5		16.5	6"	18.5		17	Std	20.5	19	Std	7"			
16	15		16	14.5		17.5	16		1 1/2	18	17		3/4	19	17.5		Std	19.5			18.5		3/4
16.5	15	3/4	16	15	3/4	18	16.5	1 1/2		18.5	17.5	6"		19.5	18	Std		20			19		
17	15.5		16.5	15.5		16.5	15		Std	17	15.5		Std	18.5	17		1 1/2	20			18.5	20.5	n/a
17	16	1 1/2	17	15.5	1 1/2	17.5	16	Std		19	17.5	1 1/2		19	18	3/4		21	n/a		1 1/2		
17.5	16.5		17.5	16		17.5	16		3/4	18	16.5		6"	19.5	18		3/4	19.5	18.5			20.5	19
16.5	15.5	Std	16.5	15	Std	17.5	16	3/4		18	16.5	6"		18	16.5	Std		20	19	Std			
17	16		17	15.5		18	16.5		3/4	18.5	17		Std	19	17.5		Std	20.5	n/a			1 1/2	6" x 2" caster not available with 5" fork
17.5	16	3/4	17	16	3/4	18.5	17	1 1/2		19	17.5	6"		19.5	18	1 1/2		20	18.5		7"		
18	16.5		17.5	16.5		18	16.5		1 1/2	18	16.5		Std	18.5	17		3/4	20.5	19			1 1/2	
18	17	1 1/2	18	16.5	1 1/2	18.5	17	1 1/2		19	17.5	6"		18.5	17.5	3/4		21	n/a	7"			
18.5	17.5		18.5	17		19	17.5		Std	19	18		7"	19.5	18.5		1 1/2	20	19			1 1/2	6" x 2" caster not available with 5" fork
				18	16.5	Std	18	16.5		Std	19.5	18.5		1 1/2	20.5	19		1 1/2	6" x 2" caster not available with 5" fork				
				18.5	17		18.5	17	20		19	1 1/2	6" x 2" caster not available with 5" fork										
				18.5	17.5	3/4	19	17.5	19.5	18	1 1/2			6" x 2" caster not available with 5" fork									
				19	18		18.5	17.5	20	18.5		1 1/2			6" x 2" caster not available with 5" fork								
				19.5	18.5	1 1/2	19	18	20.5	19	1 1/2					6" x 2" caster not available with 5" fork							
				20	19		19.5	18.5	21	n/a		1 1/2	6" x 2" caster not available with 5" fork										
								19.5	18	6" x 2" caster not available with 5" fork													
								20	18.5					1 1/2	6" x 2" caster not available with 5" fork								
								20.5	19							1 1/2	6" x 2" caster not available with 5" fork						
								6" x 2" caster not available with 5" fork															
																6" x 2" caster not available with 5" fork							
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			
												6" x 2" caster not available with 5" fork											
				6" x 2" caster not available with 5" fork																			

¹ Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.

Axiom G, S and P Cushions

Select size in chart at right. Weight capacity: 350 lbs.

Axiom G, S and P Cushions

- ☐ Axiom G £79
☐ Axiom S £110
☐ Axiom P £155

Axiom G, S or P Cushion Quantity: ☐ 1 ☐ 2

If no selection is made, quantity 1 will be selected

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Depth	14"	☐								
	15"	☐	☐	☐						
	16"	☐	☐	☐	☐	☐	☐	☐		
	17"	☐	☐	☐	☐	☐	☐	☐		
	18"	☐	☐	☐	☐	☐	☐	☐		☐
	19"		☐	☐	☐	☐	☐	☐		☐
	20"		☐	☐	☐	☐	☐	☐		☐
	21"			☐		☐	☐	☐		☐
	22"			☐		☐	☐	☐		☐

Axiom SP Cushions

Select size in chart at right. Weight capacity: 350 lbs.

Axiom SP Cushions

- ☐ Axiom SP - Visco £340
☐ Axiom SP - Fluid £315

Axiom SP Cushion Quantity: ☐ 1 ☐ 2

If no selection is made, quantity 1 will be selected

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Depth	14"	☐	☐	☐						
	15"	☐	☐	☐	☐	☐				
	16"	☐	☐	☐	☐	☐	☐	☐		☐
	17"	☐	☐	☐	☐	☐	☐	☐		☐
	18"	☐	☐	☐	☐	☐	☐	☐	☐	☐
	19"	☐	☐	☐	☐	☐	☐	☐	☐	☐
	20"	☐	☐	☐	☐	☐	☐	☐	☐	☐
	21"		☐	☐	☐	☐	☐	☐	☐	☐
	22"		☐	☐	☐	☐	☐	☐	☐	☐

Axiom Cushion Accessories

Outer Covers

- ☐ Axiom G Extra Outer Cover £59
☐ Axiom S Extra Outer Cover £67
☐ Axiom P Extra Outer Cover £76
☐ Axiom SP Extra Visco Outer Cover £84
☐ Axiom SP Extra Fluid Outer Cover £84

Outer Cover Quantity: ☐ 1 ☐ 2

If no selection is made, quantity 1 will be selected

Inner Covers

- ☐ Axiom SP Extra Visco Inner Cover £84
☐ Axiom SP Extra Fluid Inner Cover £84

Inner Cover Quantity: ☐ 1 ☐ 2

If no selection is made, quantity 1 will be selected

Solid Seat Insert

Select size in chart below. Weight capacity: 350 lbs.

Solid Seat Insert £46

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Depth	14"	☐	☐	☐						
	15"	☐	☐	☐	☐					
	16"	☐	☐	☐	☐	☐	☐	☐		☐
	17"	☐	☐	☐	☐	☐	☐	☐		☐
	18"	☐	☐	☐	☐	☐	☐	☐	☐	☐
	19"	☐	☐	☐	☐	☐	☐	☐	☐	☐
	20"	☐	☐	☐	☐	☐	☐	☐	☐	☐
	21"		☐	☐	☐	☐	☐	☐	☐	☐
	22"		☐	☐	☐	☐	☐	☐	☐	☐

Solid Seat Insert Quantity: ☐ 1 ☐ 2


If no selection is made, quantity 1 will be selected

Axiom Back


Select back that is the same or up to 2" narrower than the width of the frame. Hardware is included with back. Low Back selections are not available with headrests, modifications, or posture supports. Weight capacity: 300 lbs.

Axiom Posterior Back (PB)

		Width					
		12"	14"	16"	18"	20"	22"
		Retail £499	£499	£499	£499	£499	£499
Height	Low	8"	☐ AXPB1208	☐ AXPB1408	☐ AXPB1608	☐ AXPB1808	☐ AXPB2008
		10"		☐ AXPB1410	☐ AXPB1610	☐ AXPB1810	☐ AXPB2010
		12"	☐ AXPB1212	☐ AXPB1412	☐ AXPB1612	☐ AXPB1812	☐ AXPB2012
	Mid	14"		☐ AXPB1414	☐ AXPB1614	☐ AXPB1814	☐ AXPB2014
		16"	☐ AXPB1216	☐ AXPB1416	☐ AXPB1616	☐ AXPB1816	☐ AXPB2016
		18"		☐ AXPB1418	☐ AXPB1618	☐ AXPB1818	☐ AXPB2018



Axiom Posterior Back - Low




Axiom Posterior Back - Mid

Axiom Posterior Back Measurement Reference				
		Back Height	A	B
Height	Low	8"	10"	6.5"
		10"	12"	8.5"
		12"	14"	10.5"
	Mid	14"	16"	8.5"
		16"	18"	9.5"
		18"	20"	10.5"


2.25"

Axiom Posterior Lateral Back (PL)


		Width					
		12"	14"	16"	18"	20"	22"
		Retail £530	£530	£530	£530	£530	£530
Height	Low	8"	☐ AXPL1208	☐ AXPL1408	☐ AXPL1608	☐ AXPL1808	☐ AXPL2008
		10"		☐ AXPL1410	☐ AXPL1610	☐ AXPL1810	☐ AXPL2010
		12"	☐ AXPL1212	☐ AXPL1412	☐ AXPL1612	☐ AXPL1812	☐ AXPL2012
	Mid	14"		☐ AXPL1414	☐ AXPL1614	☐ AXPL1814	☐ AXPL2014
		16"	☐ AXPL1216	☐ AXPL1416	☐ AXPL1616	☐ AXPL1816	☐ AXPL2016
		18"		☐ AXPL1418	☐ AXPL1618	☐ AXPL1818	☐ AXPL2018
	Tall	20"	☐ AXPL1220	☐ AXPL1420	☐ AXPL1620	☐ AXPL1820	☐ AXPL2020
		22"		☐ AXPL1422	☐ AXPL1622	☐ AXPL1822	☐ AXPL2022



Axiom Posterior Lateral Back - Low



Axiom Posterior Lateral Back - Mid



Axiom Posterior Lateral Back - Tall

Axiom Posterior Lateral Back Measurement Reference				
		Back Height	A	B
Height	Low	8"	10"	6.5"
		10"	12"	8.5"
		12"	14"	10.5"
	Mid	14"	16"	8.5"
		16"	18"	9.5"
		18"	20"	10.5"
	Tall	20"	22"	11.5"
		22"	24"	12.5"


3.65"

Axiom Back continued


Select back that is the same or up to 2" narrower than the width of the frame. Hardware is included with back. Weight capacity: 300 lbs.

Axiom Deep Lateral Back (DL)

			Width					
			12"	14"	16"	18"	20"	22"
			Retail	£610	£610	£610	£610	£610
Height	Mid	14"	☐ AXDL1214	☐ AXDL1414	☐ AXDL1614	☐ AXDL1814	☐ AXDL2014	☐ AXDL2214
		16"		☐ AXDL1416	☐ AXDL1616	☐ AXDL1816	☐ AXDL2016	
		18"	☐ AXDL1218	☐ AXDL1418	☐ AXDL1618	☐ AXDL1818	☐ AXDL2018	☐ AXDL2218
	Tall	20"	☐ AXDL1220	☐ AXDL1420	☐ AXDL1620	☐ AXDL1820	☐ AXDL2020	
		22"		☐ AXDL1422	☐ AXDL1622	☐ AXDL1822	☐ AXDL2022	☐ AXDL2222



Axiom Deep Lateral Back - Mid



Axiom Deep Lateral Back - Tall


Axiom Deep Lateral Back Measurement Reference									
		Back Height	A	B ¹					Shell Contour Depth
Height	Mid	14"	16"	8.25"	8.25"	8.25"	8.25"	8.25"	6"
	Tall	16"	18"	9.25"	9.25"	9.25"	9.25"	9.25"	
		18"	20"	10.25"	10.25"	10.25"	10.25"	10.25"	
		20"	22"	9.5"	9.75"	10"	10.25"	10.25"	
		22"	24"	10.75"	11"	11.25"	11.5"	11.75"	

¹ Measurement based on back width

Axiom Back Accessories

Headrest and Hardware

Choose headrest and hardware combination in chart below. Not available with Low Backs.



			Hardware					
			Standard Straight Back Mount Hardware with Ball		Multi Axis Removable Hardware with Horizontal Bar		Headrest Hardware with Ball	
			HEDHDW101 / HEDHDW117		HEDHDW201 / HEDHDW217		HEDHDW401 / HEDHDW417	
	Description	Item Code	Kio Code	Retail	Kio Code	Retail	Kio Code	Retail
Headrest Pad	Comfort Plus Headrest 8"	HEDPAD408	☐ HDR001	£233	☐ HDR009	£263	☐ HDR017	£326
	Comfort Plus Headrest 10"	HEDPAD410	☐ HDR002	£233	☐ HDR010	£263	☐ HDR018	£326
	Comfort Plus Headrest 14"	HEDPAD414	☐ HDR003	£233	☐ HDR011	£263	☐ HDR019	£326
	Small All Positioning Headrest Pad and Ring Only	HEDPAD201	☐ HDR004	£233	☐ HDR012	£263	☐ HDR020	£326
	Medium All Positioning Headrest Pad and Ring Only	HEDPAD202	☐ HDR005	£233	☐ HDR013	£263	☐ HDR021	£326
	Large All Positioning Headrest Pad and Ring Only	HEDPAD203	☐ HDR006	£233	☐ HDR014	£263	☐ HDR022	£326
	Tri-Comfort Plus Headrest 14"	HEDPAD502	☐ HDR007	£242	☐ HDR015	£272	☐ HDR023	£335
	Tri-Comfort Plus Headrest 19"	HEDPAD503	☐ HDR008	£242	☐ HDR016	£272	☐ HDR024	£335

Headrest Hardware Options

	Description	Item Code	Retail
☐	Headrest Flip Down Mount <i>Available only with HEDHDW401/HEDHDW417</i>	HDR100	£147
☐	Universal Mounting Plate <i>For use with other manufacturer headrests</i>	UMP	£37

Axiom Back Accessories continued

Lateral Support						
Pair	Left	Right	Description	Item Code	Pair	Each
Laterals for Posterior Back						
☐	☐	☐	3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back	LF3x4P	£256	£128
☐	☐	☐	4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back	LF4x5P	£256	£128
☐	☐	☐	3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back	LS3x4P	£546	£273
☐	☐	☐	4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back	LS4x5P	£546	£273
Laterals for Posterior Lateral Back						
☐	☐	☐	3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back	LF3x4L	£256	£128
☐	☐	☐	4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back	LF4x5L	£256	£128
☐	☐	☐	3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back	LS3x4L	£546	£273
☐	☐	☐	4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back	LS4x5L	£546	£273

Axiom Lateral Support Modifications

Extended lead times will apply. Modified Lateral Supports are not eligible for refund.

	Description	Item Code	Retail
☐	Modification for Chest Strap attached to Lateral Supports	AXMLATMT	N/C

Privacy Flap

	Description	Retail
☐	Privacy Flap	£38

Shoulder Guides

	Description	Item Code	Retail
☐	Shoulder Strap Guides and Mount <i>Selection of a modification from the Axiom Back Modifications required in order to mount Shoulder Guides. Not available with Low Backs.</i>	AXBSG	£37

Therafin Posture Support

Selection of a modification from the Axiom Back Modifications section required in order to mount posture supports. Not available with Low Backs.

	Description	Item Code	Retail
☐	Small Static Therafit 2-Piece Chest Strap with H & L Adj and 1" Safety Strap	60153T	£121
☐	Medium Static Therafit 2-Piece Chest Strap with H & L Adj and 1" Safety Strap	60154T	£121
☐	Large Static Therafit 2-Piece Chest Strap with H & L Adj and 1" Safety Strap	60155T	£121
☐	Small Dynamic Therafit Chest Strap with Dual D-Ring Adjustment	60110T	£95
☐	Medium Dynamic Therafit Chest Strap with Dual D-Ring Adjustment	60111T	£95
☐	Large Dynamic Therafit Chest Strap with Dual D-Ring Adjustment	60112T	£95
☐	Small Dynamic Therafit Comfortfit Full Vest	30589T	£140
☐	Medium Dynamic Therafit Comfortfit Full Vest	30590T	£140
☐	Large Dynamic Therafit Comfortfit Full Vest	30591T	£140
☐	XLarge Dynamic Therafit Comfortfit Full Vest	30592T	£140
☐	XXLarge Dynamic Therafit Comfortfit Full Vest	30593T	£140
☐	Small Dynamic Therafit Comfortfit Trim Vest	30618T	£142
☐	Medium Dynamic Therafit Comfortfit Trim Vest	30619T	£142
☐	Large Dynamic Therafit Comfortfit Trim Vest	30620T	£142
☐	XLarge Dynamic Therafit Comfortfit Trim Vest	30621T	£142
☐	XXLarge Dynamic Therafit Comfortfit Trim Vest	30622T	£142

Axiom Back Modifications

Extended lead times will apply for all modifications. Makes back non-refundable.

Axiom Back Modifications

Choose modification to add holes to the back shell in order to mount posture supports. Not available with Low Backs.

	Description	Item Code	Retail
☐	Modification for 4 Point Harness	AXM4PT	N/C
☐	Modification for 2 Point Chest Strap	AXM2PT	N/C
☐	Modification for Shoulder Guides	AXMSG	N/C

Painted Back Shell Modification

If no selection is made, standard Shadow Black will be used with no additional lead time

	Gloss Finish	Item Code	Retail
☐	Black	AXMBLK	N/C
☐	Blue Skies	AXMSKY	N/C
☐	Bubblegum Pink	AXMPNK	N/C
☐	Candy Blue	AXMCBL	N/C
☐	Candy Purple	AXMPUR	N/C
☐	Candy Red	AXMRED	N/C
☐	Copper	AXMCOP	N/C
☐	Granny Smith	AXMGSM	N/C
☐	Metallic Orange	AXMMOR	N/C
☐	Silver	AXMSIL	N/C
☐	Teal	AXMTEA	N/C
☐	Yellow	AXMYLW	N/C

	Matte Finish	Item Code	Retail
☐	Charcoal Gray	AXMCGM	N/C
☐	Emerald Green	AXMEGM	N/C
☐	Iridescent Burgundy	AXMIBM	N/C
☐	Iridescent Navy	AXMINM	N/C
☐	Military Green	AXMAGM	N/C
☐	Retro Red	AXMRRM	N/C
☐	Royal Blue	AXMRBM	N/C
☐	Smooth White	AXMSWM	N/C