2024

ottobock.

Base Price - £834.00	Order no.:					
Valid from 01.01.24 to 31.12.24 Please ensure you are using the latest or	rder form					
☐ Provide estimate only	□Order					
Billing address/customer number						
Delivery address						
Sales person	Name					
	Date Si	gnature				
By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock						

product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for

Start^{M6} Junior Lightweight wheelchair

Article number: 480F53=60000_K



Standard equipment

- Aluminium tube frame with crossbrace
- Tip-assist integrated into the frame
- Individual leg supports, angle-adjustable, plastic foot plates
- Swing-away leg support holders
- Side panel with arm support, height-adjust-
- · Nylon seat and back support upholstery, black
- Modular design
- Drive wheels with quick-release axles
- Knee lever wheel lock

Folding, max. user weight 90 kg/198 lbs Crash-tested up to 75 kg according to DIN ISO 7176-19

Wheelchair dimensions								
Seat width (SW)			cm	Anterior sea	at height (<i>A</i>	ASH)		 cm
Seat depth (SD)			cm	Posterior se	at height (PSH)		 cm
Back support height (BH)			cm	Lower leg le	ength (LL)	(frame meası	urement)	 cm
Seat measurements in cm								
Seat width (SW)	<u>28</u>	□30.5	□33	35.5	□38			
Seat depth (SD)	□31	□34	37					
Seat height measurements in cm (table on the next page)								
Anterior seat height (ASH)	4 0	4 1	<u>42</u>	4 3	<u>44</u>			
	☐ 45	□ 46	□ 47	□48	49			
Posterior seat height 37	□38	□39	40	41	42	43	44	
	☐ 45	□ 46	47	□48	□ 49			

The number codes contained in this order form are not intended to be used as article numbers for spare parts orders. Please use the spare parts catalogue for spare parts orders.

Information


Seat height table Anterior seat height in cm

Caster fork, short Adjustment range

Caster wheel diameter 5" 40 – 44 Caster wheel diameter 6" 43 – 45

Caster fork, long

Caster wheel diameter 5" 44 – 47 Caster wheel diameter 6" 45 – 49



Seat height table

Posterior seat height in cm

Drive wheel size	Adjustment range
20"	37 – 43
22"	38 – 45
24"	41 – 49



The values indicated are rounded values which have been theoretically determined. Max. deviation of 1 cm. The anterior seat height depends on the selected wheel size, caster fork and mounting position. Please refer to the seat height table!

Specified without seat cushion at 0° seat inclination.

The maximum possible seat inclination is 3.5 cm.

ottobock.

Frame accessories					
MA 55	Тп	Plug-on anti-tipper, right	£25.00		
MA 56	1_	Plug-on anti-tipper, left	£25.00		
MA 57		Crutch holder, right	£30.00		
MA 58		Crutch holder, left	£30.00		
MA 61		Swinging anti-tipper, right	£49.00		
MA 62		Swinging anti-tipper, left	£49.00		
Caster wheels	Caster wheels				
MF 03		5" solid rubber	Inclusive		
MF 05		6" pneumatic (150 x 30)	Inclusive		
MF 06		6" solid rubber (150 x 27)	Inclusive		
Caster forks					
MF 80		Caster fork, long	Inclusive		
MF 81		Caster fork, short	Inclusive		
Drive wheel size					
MG 04		20"	Inclusive		
MG 01		22"	Inclusive		
MG 02		24"	Inclusive		
Wheelbase					
short/active]long/passive			
Drive wheel rim					
MG 05		Standard	Inclusive		
MG 80		Drive wheel with drum brake/attendant	£373.00		
Tyres					
MG 30		Pneumatic tyres 1 3/8"	Inclusive		
MG 35		PU tyres 1 3/8"	Inclusive		
Handrims					
MG 52		Handrim, aluminium, silver anodised, high	Inclusive		
MG 56		Handrim, polished stainless steel, high	£58.00		
Drive wheel accessories					
MG 87		High-pressure air pump, 14 bar	£113.00		
MG 88		Airman pump, battery-operated	£173.00		
Wheel lock					
MH 01		Knee lever wheel lock, standard	Inclusive		
MH 06		Knee lever wheel lock for user and attendant	£242.00		
Wheel lock accessories					
MH 10		Wheel lock lever extension	£30.00		

Information



MA 55 Plug-on anti-tipper



MA 62 Swinging anti-tipper



MD 10 Back support angle adjustment 30°



MD 11 Folding back support



MD 23 Thoracic support



MD 54 Angled push handles



ME 05 Side panel, height-adjustable



MB 20 Leg support for short lower leg length



MB 72 Heel strap



MZ 53 Head support



MZ 54 Clinical installation set



MZ 55 Multi-axis installation set



MG 95 Spoke protector with decor

ottobock.

Back support system	1					
MD 09		Back support with lumbar curve, standard	Inclusive			
MD 01		Straight back support, standard	Inclusive			
MD 10		Back support angle adjustment 30° (only available with BH 37.5 cm)	£234.00			
MD 11		Back support, folding (only available with BH 37.5 cm)	£92.00			
Back support height	in cn	n				
□30 □32.5		35 37.5				
Back support uphols	tery					
MD 04		Standard	Inclusive			
MD 05		Adaptable	£108.00			
MD 06		Without	Inclusive			
Push handles						
MD 51		Standard	Inclusive			
MD 52		Telescoping (only available in combination with 30° back support angle adjustment)	£201.00			
MD 54		Height-adjustable, angled (only possible in connection with MD01)	£116.00			
Back support access	ories					
MD 23		Side thoracic support, right	£95.00			
MD 24		Side thoracic support, left	£95.00			
MD 58		Stabiliser bar for push handles	£47.00			
Side panels						
ME 01		Side panel, standard, 300 mm arm support, short mounted	Inclusive			
ME 01		Side panel, standard, 300 mm arm support, long mounted	Inclusive			
ME 05		Side panel, height-adjustable, 300 mm armrest	£132.00			
Side panel accessori						
ME 60		Tray	£188.00			
Seat upholstery						
MC 13		Upholstery, standard	Inclusive			
MC 06		Without	Inclusive			
Seat cushion						
MC 01		Seat cushion Standard (thickness 5 cm)	£108.00			
MC 02		Cubic Foam seat cushion (thickness 5 cm)	£141.00			
MC 05		Evolight seat cushion (thickness 6 cm)	£201.00			
Lower leg length in cm (without seat cushion)						
		m (available adjustment range 16 cm – 45 cm)				
Leg support bar options						
MB 40		Short (27 – 40 cm)	Inclusive			
MB 41		Medium (32 – 45 cm)	Inclusive			
	wer le	g lengths vary in case of elevating leg supports.				
. Tests that are larger regularity in case of creating leg supports.						

Leg supports					
MB 10		Individual leg supports, angle-adjustable, plastic foot plates	Inclusive		
MB 25		Leg supports, elevating, without automatic length adjustment, pair	£242.00		
MB 26		Leg support, elevating, without automatic length adjustment, right (plus standard leg support, left)	£122.00		
MB 27		Leg support, elevating, without automatic length adjustment, left (plus standard leg support, right)	£122.00		
MB 20		Single-panel leg support for short lower leg length	£206.00		
MB 19		Individual leg supports for short lower leg length	£75.00		
Leg support accesso	ries				
MB 72		Heel strap for individual leg supports	£31.00		
Other accessories					
MZ 50		Lap belt with plastic buckle	£45.00		
MZ 52		Tool kit	£61.00		
MZ 53		Head support	£116.00		
MZ 54		Clinical installation kit for head support	£118.00		
MZ 55		Head support installation kit, multi-axis, offset	£103.00		
MZ 70		Anchor point kit for transportation in a motor vehicle (max. loading capacity: 75 kg)	£364.00		
Note: Please also use the order forms for wheelchair accessories 647F661.					
Spoke protector					
MG 90		Spoke protector, transparent	£145.00		
MG 95		Spoke protector with decor	£232.00		
Frame colour					
009		Vinho sparkle	£179.00		
017		Traffic yellow RAL 1023	Inclusive		
018		Signal red RAL 3001	Inclusive		
021		Sparkle granny-smith	£179.00		