

STELLAR LEAP



The *Stellar LEAP* has a tilt range of 45° and **combines posterior and anterior tilt** in a single wheelchair. The LEAP offers enhanced user function and positioning, as well as **standing transfer support** through anterior tilt.



Hand Propulsion



Foot Propulsion



Self-Tilt



Desk/Table Accessible



Transit Tie Downs

Stellar LEAP *A unique chair with anterior and posterior tilt*

visit pdgmobility.com/leap for more information

Available options See order form

- ✓ 20-30° of anterior plus 15-25° of posterior tilt (45° total range)
- ✓ Transit tie downs WC-19 & ISO-7176-19
- ✓ Power tilt (anterior/posterior)



DYNAMIC TILT
45° Range



WEIGHT CAPACITY
250 lbs



SEAT WIDTH
14" to 22"



SEAT DEPTH
16" to 22"



SEAT HEIGHT
13" to 18"



BACK SEAT ANGLE
90° to 122°



OVERALL WIDTH
+9.5" to 10.5"



CHAIR WEIGHT
58 lbs

