




STIMULITE® TENSION ADJUSTABLE BACK

Supracor's new Stimulite Tension Adjustable Back features Stimulite honeycomb to circulate air and evaporate moisture while equalizing pressure for exceptional comfort and stability. The "footprint" of the flexible honeycomb cells stimulates blood flow to enhance circulation. Adjustable velcro straps can accommodate for kyphosis or posterior pelvic tilt. Optional lumbar support. Available for rigid and folding frame wheelchairs from 14 to 24 inches wide. Both the honeycomb and removable, breathable cover are machine washable.



STIMULITE® LUMBAR SUPPORT

Whether you're sitting in a chair or car, this sleek and breathable Stimulite Lumbar Support will improve your posture and relieve pressure while keeping you cool and dry. Features a removable, adjustable strap. Washable.



STIMULITE® HONEYCOMB SHEETS

SBS-1: 1/2" thick, bare honeycomb sheet
 SBS-2: 1/2" thick, single-faced sheet
 SBS-3: 1/2" thick sheet with fabric on one side
 SBS-4: 11/16" thick, single-faced sheet
 SBS-5: 1" thick, single-faced sheet
 SBS-6: 3/4" thick & 1/4" thick, single-faced, double honeycomb panel
 SBS-7: 3/4" thick, double-faced, firmer honeycomb sheet; ideal for correcting pelvic obliquities

Sheets are 24" wide x 24" long



STIMULITE® ADJUSTABLE POSITIONING WEDGES

Covered in Cordura fabric with Velcro for easy attachment to cushion bottom. Available in 1-inch and 2-inch thickness.

- 14" wide x 10" long
- 15-3/4" wide x 10" long



STIMULITE® ON TOP CUSHION COVER

Don't yet have a Stimulite cushion? Just put Stimulite on top of your current cushion and like magic, your cushion will now breathe to control heat and moisture. This removable cover features a 3/4-inch thick layer of Stimulite honeycomb that fits over any cushion. Machine washable.



WATERPROOF COVER

This slick, stretch nylon waterproof cover will keep your cushion dry in the rain, at the pool or in the shower and is impervious to bodily fluids. Fits over any cushion. Machine washable and dryer safe.



OUTDOOR COVER

Made of a durable, black power mesh fabric, Supracor's breathable outdoor cover features side vents and a slip-resistant bottom to meet the rigorous demands of sports, camping and everyday use. Machine washable and dryer safe.