

Q500® M/F/R with SEDEO® PRO




ORDER FORM

Billing Address Or Customer-ID			Shipping Address		Identical To Billing Address
C.-ID:			C.-ID:		
Name:			Name:		
Street:			Street:		
Town/PC:			Town/PC:		
E-Mail:			E-Mail:		
Phone:			Phone:		
Order Information			User Weight		
Order	Calculation	Built-4-Me®	User Weight Capacity: 160 kg		
Purchase Order No.:			User Weight: kg		
Art. No.: Q500 M: 2250050M-000, Q500 F: 2250050F-000, Q500 R: 2250050R-000					
Order Date:	Day	Month	Year		

Client Measurements

a) Top Of Head To Seat	mm
b) Top Of Shoulder To Seat	mm
c) Underarm To Seat	mm
d) Elbow To Seat	mm
e) Top Of Pelvis To Seat	mm
f) Shoulder Width	mm
g) Trunk Width	mm
i) Hip Width	mm
j) Actual Sitting Depth	mm
m) Foot Length	mm
n) Lower Leg Length	mm
o) Actual Seat Width	mm



Seating Specifications

All Bases Prepared For DAHL Docking Station		Seat Type	SEDEO PRO
Typical WC Weight	Without Batteries 140 kg	Seat Size Range	Width 380 - 560 mm, Depth 400 - 560 mm
	With Batteries 165 kg	Seat To Floor Height Range	415 mm - 490 mm

MWD Base Specifications

User Weight Capacity	160 kg	Drive Wheels	13" Drive Wheels (Black)
Base Size	Width Min. 620 mm	Tires / Castors	7" Front / 7" Rear (Black)
Drive Control	VR2, Option R-Net	Charger	8 A Dual Mode Charger
Driving Range	Up To 40 km	Battery Size	60 And 80 Ah Batteries
Drive Motors	4-Pole 6 kph, Options 8, 10 And 12.5 kph	Colour Accents	7 Colour Accents
Ground Clearance	70 mm	Frame Colour	Black

FWD Base Specifications

User Weight Capacity	160 kg	Drive Wheels	13" / 14" Drive Wheels (Black)
Base Size	Width Min. 610 mm	Tires / Castors	9" Rear (Black), Options 10"
Drive Control	VR2, Option R-net	Charger	8 A Dual Mode Charger
Driving Range	Up To 40 km	Battery Size	60 And 80 Ah Batteries
Drive Motors	4-Pole 6 kph, Options 8, 10 And 12.5 kph	Colour Accents	7 Colour Accents
Ground Clearance	65 mm	Frame Colour	Black

RWD Base Specifications

User Weight Capacity	160 kg	Drive Wheels	13" / 14" Drive Wheels (Black)
Base Size	Width Min. 610 mm	Tires / Castors	9" Front (Black), Options 10"
Drive Control	VR2, Option R-net	Charger	8 A Dual Mode Charger
Driving Range	Up To 40 km	Battery Size	60 And 80 Ah Batteries
Drive Motors	4-Pole 6 kph, Options 8, 10 And 12.5 kph	Colour Accents	7 Colour Accents
Ground Clearance	65 mm	Frame Colour	Black

Base

Q500 Drive Bases

		£
QOP000023	Q500 M (MWD Base Chair)	5,400
QOP000024	Q500 F (FWD Base Chair)	5,600
QOP000029	Q500 R (RWD Base Chair)	5,600

Motor Packages

		£
QOP180009	4-Pole Motor, 6 kph / 4 mph	0
QOP180010	4-Pole Motor, 8 kph / 5 mph	162
QOP180007	4-Pole Motor, 10 kph / 6.5 mph (1) (2)	324
QOP180008	4-Pole Motor, 12.5 kph / 8 mph (1) (3)	405

Tracking

		£
QOP110009	Traction Control System For R-Net (GYRO Technology)	324
QOP110014	Omit Traction Control System For R-Net	0

Drive Wheels

MWD/FWD/RWD Base		
QOP070065	Drive Wheel, 13", Pneumatic, Knobbly, Black	0
QOP070069	Drive Wheel, 13", Solid, Knobbly, Black	0
FWD/RWD Base		
QOP070032	Drive Wheel, 14", Pneumatic, V-Type, Black	0
QOP070031	Drive Wheel, 14", Solid, V-Type, Black	0
QOP070034	Drive Wheel, 14", Pneumatic, Knobbly, Black	0
QOP070033	Drive Wheel, 14", Solid, Knobbly, Black	0

Castor Wheels

MWD Base		
QOP080213	Castor Wheel, 7", Solid, Black	0
FWD Base		
QOP080210	Castor Wheel, 9", Solid, Black	0
QOP080211	Castor Wheel, 9", Pneumatic, Black	0
RWD Base		
QOP080210	Castor Wheel, 9", Solid, Black	0
QOP080211	Castor Wheel, 9", Pneumatic, Black	0
QOP080212	Castor Wheel, 10", Pneumatic, Black	0

Mudguards (FWD/RWD Only)

		£
QOP010007	Mudguard, Drive Wheel, 14"	97
QOP010008	Mudguard, Castor Wheel, 9"	76

Kerb Climber (RWD Only)

		£
QOP090038	Kerb Climber	184

(1) Highly Recommended With 80 Ah Battery

(2) If FWD Selected: Highly Recommended With Traction Control System

(3) Not available with FWD

Base

Colour Accents

£

QOP140005	Apple Green (RAL 110 70 77)	0
QOP140004	Black (RAL 9005)	0
QOP140003	Blue (RAL 260 40 45)	0
QOP140006	Orange (RAL 060 60 80)	0
QOP140002	Red (RAL 040 40 67)	0
QOP140007	Silver (RAL 9006)	0
QOP140001	White (RAL 9010)	0

LED Lights & Indicators

£

QOP110201	LED Light And Indicator StVZO Approved (1)	319
QOP110209	Omit Light And Indicator	0

Seat Heights (2)

£

SH415	415 mm (3) (4)	0
SH425	425 mm (5)	0
SH44	440 mm (4)	0
SH45	450 mm (5)	0
SH465	465 mm	0
SH49	490 mm	0

Batteries & Chargers

£

Batteries

QOP160003	Battery, 50 Ah, Gel, Maintenance-Free	0
QOP160001	Battery, 60 Ah, Gel, Maintenance-Free	0
QOP160005	Battery, 80 Ah, Gel, Maintenance-Free (6)	130
QOP160009	Omit Batteries	0

Chargers

QOP160101	Charger, Automatic, 220 V / 8 A	0
QOP160100	Charger, Automatic, 220 V / 8 A, UK	0
QOP160103	Charger, Automatic, 220 V / 10 A, UK	184
QOP160104	Charger, Automatic, 220 V / 10 A	184
QOP160109	Omit Charger	0

(1) QOP110030 Lighting Module VR2 Or QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required

(2) Measured From Front Of Seat Pan To Floor. Obtain Quote To Verify Seat To Floor Heights

(3) Only With 60 Ah Batteries

(4) Not With QOP130009 Seat Elevator/Tilt Module, Powered

(5) Only For QOP130009 Seat Elevator (300 mm) / Tilt Module (Tilt 0° - 30°), Powered

(6) Highly Recommended With Battery Charger 220 V / 10 A

Powered Seating

Power Seating Options

		£
QOP130011	Fixed Tilt Module (Adjustable 0°, 3°, 6°, 9°, Factory Setting 6°)	0
QOP130004	Posterior Tilt, Powered, 0° - 30° (1)	891
QOP130012	Posterior Tilt, Powered, -5° - 25° (1) (2)	999
QOP130007	Posterior Tilt, Powered, 0° - 50° (1)	1107
QOP130005	Posterior Tilt, Powered, -5° - 45° (1) (2)	1215
QOP130009	Seat Elevator (300 mm) / Tilt Module (Tilt 0° - 30°), Powered (1)	1620

Backrest Frame Type

		£
QOP030083	Backrest, Tool Angle Adjustable, 90° - 120°	0
QOP030084	Backrest, Manual Recline, User Adjustable, 90° - 130°, Release Lever, Right (3)	248
QOP030085	Backrest, Manual Recline, User Adjustable, 90° - 130°, Release Lever, Left (3)	248
QOP030086	Backrest, Powered Recline, 90° - 130° (1)	454
QOP030087	Backrest, Powered Recline, 90° - 150° With 100 mm Mechanical Anti-Shear (1)	810

Backrest Frame Option

		£
QOP030098	Backrest, Flip Down, Quick Release	103

Back Cushion Type

		£
SEDEO PRO (See Section A)		
QOP030104	SEDEO PRO, Back Cushion	0
JAY J3 (See Section B)		
QOP030054	JAY J3, Back Cushion	0
Interfaces 		
QOP030106	SEDEO PRO, Interface Only	0
QOP030107	SEDEO PRO, Interface Only For Back Cushion, Flexible	0
QOP030055	JAY J3, Interface Only	0
QOP030209	Omit Back Cushion	0

A – SEDEO PRO Back Size

		£
Back Cushion Width		
BW38	380 mm, For Effective Back Width 380 mm	0
BW42	420 mm, For Effective Back Width 380 - 420 mm	0
BW46	460 mm, For Effective Back Width 420 - 460 mm	0
BW50	500 mm, For Effective Back Width 460 - 500 mm	0
BW54	540 mm, For Effective Back Width 500 - 540 mm	0
BW56	560 mm, For Effective Back Width 520 - 560 mm	0
Back Cushion Height		
BH44	440 mm, For Effective Back Height 440 - 500 mm (4)	0
BH50	500 mm, For Effective Back Height 500 - 560 mm	0
BH56	560 mm, For Effective Back Height 560 - 620 mm	0

A – SEDEO PRO Back Cushion

		£
QOP030109	SEDEO PRO, Back Cushion, Profiled	113
QOP030113	SEDEO PRO, Back Cushion, Profiled, Adjustable Straps	211
QOP030114	SEDEO PRO, Back Cushion, Deep Contour, Passive	270
QOP030115	SEDEO PRO, Back Cushion, Flexible, Deep Contour, Passive, Adjustable Straps	351
QOP030116	SEDEO PRO, Back Cushion, Flexible, Deep Contour, Active, Adjustable Straps	351
QOP030117	SEDEO PRO, Back Cushion, Flexible, Extra Deep Contour, Active, Adjustable Straps	378
	SEDEO PRO, Back Cushion, Adaptive (Choose B4Me Requirements Separately)	

(1) Not With QOP110036 Remote Joystick, VR2. In Case Of R-Net: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L

(2) Not With Seat Height 415 Or 440 mm

(3) Not On Same Joystick Side

(4) Only With QOP030116 SEDEO PRO, Back Cushion, Flexible, Deep Contour, Active, Adjustable Straps & QOP030117 SEDEO PRO, Back Cushion, Flexible, Extra Deep Contour, Active, Adjustable Straps

Powered Seating

A – SEDEO PRO Back Cushion Cover

£

QOP030059	Back Cushion Cover, Incontinence Protection	0
QOP030108	Back Cushion Cover, Comfort	0

A – SEDEO Lateral Thoracic Supports

£

QOP030607	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Right	292
QOP030608	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Left	292
QOP030609	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Short Pad, 90 x 120 mm, Right	292
QOP030610	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Short Pad, 90 x 120 mm, Left	292

B – JAY J3 Back Size

£

JAY030111	JAY J3 Back Size, 410 mm	0
JAY030112	JAY J3 Back Size, 460 mm	0
JAY030113	JAY J3 Back Size, 510 mm	0

B – JAY J3 Back, Shallow Contour

£

JAY030012	JAY J3 Back, Shallow Contour, Upper Thoracic, Short, 420 mm	572
JAY030019	JAY J3 Back, Shallow Contour, Upper Thoracic, Medium, 460 mm	572
JAY030025	JAY J3 Back, Shallow Contour, Upper Thoracic, Tall, 510 mm	572
JAY030015	JAY J3 Back, Shallow Contour, Shoulder Height, Short, 530 mm	572
JAY030022	JAY J3 Back, Shallow Contour, Shoulder Height, Medium, 570 mm	572
JAY030028	JAY J3 Back, Shallow Contour, Shoulder Height, Tall, 610 mm	572

B – JAY J3 Back, Mid Contour

£

JAY030013	JAY J3 Back, Mid Contour, Upper Thoracic, Short, 420 mm	572
JAY030020	JAY J3 Back, Mid Contour, Upper Thoracic, Medium, 460 mm	572
JAY030026	JAY J3 Back, Mid Contour, Upper Thoracic, Tall, 510 mm	572
JAY030016	JAY J3 Back, Mid Contour, Shoulder Height, Short, 530 mm	572
JAY030023	JAY J3 Back, Mid Contour, Shoulder Height, Medium, 570 mm	572
JAY030029	JAY J3 Back, Mid Contour, Shoulder Height, Tall, 610 mm	572

B – JAY J3 Back, Mid Deep Contour

£

JAY030045	JAY J3 Back, Mid Deep Contour, Upper Thoracic, Short, 420 mm	627
JAY030042	JAY J3 Back, Mid Deep Contour, Upper Thoracic, Medium, 460 mm	627
JAY030046	JAY J3 Back, Mid Deep Contour, Upper Thoracic, Tall, 500 mm	627
JAY030043	JAY J3 Back, Mid Deep Contour, Shoulder Height, Short, 530 mm	627
JAY030047	JAY J3 Back, Mid Deep Contour, Shoulder Height, Medium, 570 mm	627
JAY030048	JAY J3 Back, Mid Deep Contour, Shoulder Height, Tall, 610 mm	627
JAY030050	JAY J3 Back, Mid Deep Contour, Extra Deep, Upper Thoracic, Medium, 460 mm	627

B – JAY J3 Back, Deep Contour

£

JAY030014	JAY J3 Back, Deep Contour, Upper Thoracic, Short, 420 mm	627
JAY030021	JAY J3 Back, Deep Contour, Upper Thoracic, Medium, 460 mm	627
JAY030027	JAY J3 Back, Deep Contour, Upper Thoracic, Tall, 500 mm	627
JAY030017	JAY J3 Back, Deep Contour, Shoulder Height, Short, 530 mm	627
JAY030024	JAY J3 Back, Deep Contour, Shoulder Height, Medium, 570 mm	627
JAY030030	JAY J3 Back, Deep Contour, Shoulder Height, Tall, 610 mm	627
JAY030051	JAY J3 Back, Deep Contour, Extra Deep, Upper Thoracic, Medium, 460 mm	627

(1) Only Use With SEDEO Headrest Mount

(2) Only Use With Whitmyer Headrest Mount

(3) Choose A Whitmyer Headrest Mount And Headrest Receiver In Case You Would Like To Add A Head Control In The Special Control Section

(4) Only With QOP030109 SEDEO PRO, Back Cushion, Profiled & QOP030113 SEDEO PRO, Back Cushion, Profiled, Adjustable Straps & QOP030114 SEDEO PRO, Back Cushion, Deep Contour, Passive

(5) Not With Seat Width 380 & 400 mm

Powered Seating

B – JAY J3 Back, Spine Align Fit Kit

£

Spine Align Mini Fit Kit I

JAY090059	JAY J3 Spine Align Mini Fit Kit I, 360 mm	31
JAY090029	JAY J3 Spine Align Mini Fit Kit I, 410 mm	31
JAY090030	JAY J3 Spine Align Mini Fit Kit I, 460 mm	31
JAY090031	JAY J3 Spine Align Mini Fit Kit I, 510 mm	31

Spine Align Mini Fit Kit II

JAY090060	JAY J3 Spine Align Mini Fit Kit II, 360 mm	47
JAY090032	JAY J3 Spine Align Mini Fit Kit II, 410 mm	47
JAY090033	JAY J3 Spine Align Mini Fit Kit II, 460 mm	47
JAY090034	JAY J3 Spine Align Mini Fit Kit II, 510 mm	47

Positioning Boxes

JAY090035	JAY J3 Positioning Box For Back	151
JAY090036	JAY J3 Hardware Box For Back	474

B – JAY J3 Back, Chest Strap (1)

£

JAY090040	JAY J3 Chest Strap, Centre Open Harness, 4 Point Strap, Medium	157
JAY090041	JAY J3 Chest Strap, Centre Open Harness, 4 Point Strap, Large	157
JAY090042	JAY J3 Chest Strap, Centre Open Harness, 4 Point Strap, Extra Large	157

B – JAY J3 Lateral Thoracic Pads, Right

£

JAY090020	JAY J3 Lateral Thoracic Pad, (LxH) 100 x 100 mm, Right	46
JAY090021	JAY J3 Lateral Thoracic Pad, (LxH) 100 x 127 mm, Right	46
JAY090022	JAY J3 Lateral Thoracic Pad, (LxH) 100 x 152 mm, Right	46
JAY090023	JAY J3 Lateral Thoracic Pad, (LxH) 152 x 100 mm, Right	46
JAY090024	JAY J3 Lateral Thoracic Pad, (LxH) 152 x 127 mm, Right	46
JAY090025	JAY J3 Lateral Thoracic Pad, (LxH) 152 x 152 mm, Right	46

B – JAY J3 Lateral Thoracic Pads, Left

£

JAY090011	JAY J3 Lateral Thoracic Pad, (LxH) 100 x 100 mm, Left	46
JAY090012	JAY J3 Lateral Thoracic Pad, (LxH) 100 x 127 mm, Left	46
JAY090013	JAY J3 Lateral Thoracic Pad, (LxH) 100 x 152 mm, Left	46
JAY090014	JAY J3 Lateral Thoracic Pad, (LxH) 152 x 100 mm, Left	46
JAY090015	JAY J3 Lateral Thoracic Pad, (LxH) 152 x 127 mm, Left	46
JAY090016	JAY J3 Lateral Thoracic Pad, (LxH) 152 x 152 mm, Left	46

B – JAY J3 Lateral Thoracic Brackets

£

JAY090026	JAY J3 Lateral Thoracic, Fixed Bracket, Right	46
JAY090017	JAY J3 Lateral Thoracic, Fixed Bracket, Left	46
JAY090028	JAY J3 Lateral Thoracic, Bracket, Swing Away, Right	87
JAY090019	JAY J3 Lateral Thoracic, Bracket, Swing Away, Left	87

(1) Only With QOP020518 Seat Cushion Cover, Incontinence Protection, Only For NL

Powered Seating

Headrest Pad

£

QOP030430	SEDEO, Headrest Pad, Small, (LxH) 240 x 150 mm (1)	76
QOP030431	SEDEO, Headrest Pad, Small, With Cover, (LxH) 240 x 150 mm (1)	86
QOP030432	SEDEO, Headrest Pad, Flexible, (LxH) 400 x 110 mm (1)	130
WHT030401	Whitmyer, Headrest Pad, Plush, (LxH) 203 x 127 mm (2)	108
WHT030402	Whitmyer, Headrest Pad, Plush, (LxH) 254 x 127 mm (2)	119
WHT030403	Whitmyer, Headrest Pad, Plush, (LxH) 356 x 127 mm (2)	130
WHT030404	Whitmyer, Headrest Pad, Adjust-A-Plush, Small, (LxH) 355 x 76 mm (2)	173
WHT030406	Whitmyer, Headrest Pad, Adjust-A-Plush, Large, (LxH) 460 mm x 127 mm (2)	194

Headrest Mount (3)

£

QOP030433	SEDEO, Headrest Mount	108
QOP030434	SEDEO, Headrest Mount, With Quick Release	119
WHT030407	Whitmyer, Headrest Mount, LINX	130
WHT030408	Whitmyer, Headrest Mount, LINX²	130
WHT030409	Whitmyer, Headrest Mount, AXYS	108
WHT030410	Whitmyer, Headrest Mount, AXYS, Offset	119
WHT030411	Whitmyer, Headrest Mount, COBRA XTRA Flip Back	162

Headrest Receiver

£

QOP030119	Headrest Receiver, Back Cushion Mount (4)	0
QOP030120	Headrest Receiver, Backrest Frame Mount	0

Push Handles

£

QOP030204	Push Handles, Height-Adjustable	162
-----------	---------------------------------	-----

Seat Rail

£

QOP020085	Seat Rail, Short, For Seat Depth 400 - 520 mm	0
QOP020069	Seat Rail, Long, For Seat Depth 460 - 560 mm (5)	0

Powered Seating

Seat Cushion Type

£

SEDEO (See Section C)

QOP020027	SEDEO, Seat Cushion	0
-----------	---------------------	---

JAY (See Section D)

QOP020515	JAY, Seat Cushion	0
-----------	-------------------	---

Omit (Select Seat Width And Depth in Section C)

QOP020070	Omit Seat Cushion	0
-----------	-------------------	---

C – SEDEO Seat Size

£

Seat Cushion Width

SW38	380 mm, For Effective Seat Width 380 mm	0
SW42	420 mm, For Effective Seat Width 380 - 420 mm	0
SW46	460 mm, For Effective Seat Width 420 - 460 mm	0
SW50	500 mm, For Effective Seat Width 460 - 500 mm	0
SW54	540 mm, For Effective Seat Width 500 - 540 mm	0
SW56	560 mm, For Effective Seat Width 520 - 560 mm	0

Seat Cushion Depth

SD40	400 mm, Comes With Cushion Size 460 mm	0
SD42	420 mm, Comes With Cushion Size 460 mm	0
SD44	440 mm, Comes With Cushion Size 460 mm	0
SD46	460 mm, Comes With Cushion Size 500 mm	0
SD48	480 mm, Comes With Cushion Size 500 mm	0
SD50	500 mm, Comes With Cushion Size 500 mm	0
SD52	520 mm, Comes With Cushion Size 560 mm	0
SD54	540 mm, Comes With Cushion Size 560 mm	0
SD56	560 mm, Comes With Cushion Size 560 mm	0

C – SEDEO Seat Cushions

£

QOP020073	SEDEO, Seat Cushion, Profiled, Height 70 mm	119
QOP020075	SEDEO, Seat Cushion, Support, Height 90 mm	173
QOP020076	SEDEO, Seat Cushion, Support Memory, Height 90 mm	205
QOP020078	SEDEO, Seat Cushion, Stabilizer, Height 90 mm	227
	SEDEO, Seat Cushion, Adaptive (Choose B4Me Requirements Separately)	

C – SEDEO Seat Cushion Cover

£

QOP020518	Seat Cushion Cover, Incontinence Protection	0
QOP020080	Seat Cushion Cover, Comfort	0

C – SEDEO Seat Boards

£

QOP020082	SEDEO, Seat Board, Flat, Depth-Adjustable, Cushions Guide At Rear	108
QOP020083	SEDEO, Seat Pan, 50 mm, Depth-Adjustable, Cushions Guide At Rear	119

Powered Seating

D – JAY Seat Size

£

Seat Cushion Width

SW40	400 mm	0
SW42	420 mm	0
SW44	440 mm	0
SW46	460 mm	0
SW48	480 mm	0
SW50	500 mm, Comes With Cushion Size 500 mm	0
SW52	520 mm, Comes With Cushion Size 500 mm	0
SW54	540 mm, Comes With Cushion Size 560 mm	0
SW56	560 mm, Comes With Cushion Size 560 mm	0

Seat Cushion Depth

SD40	400 mm	0
SD42	420 mm	0
SD44	440 mm	0
SD46	460 mm	0
SD48	480 mm	0
SD50	500 mm	0
SD52	520 mm	0
SD54	540 mm (1)	0
SD56	560 mm (1)	0

D – JAY Seat Cushions

£

JAY020005	JAY Easy Fluid With Flat Base	308
JAY020003	JAY Soft Combi P	77
JAY020006	JAY Easy Visco With Flat Base	132
JAY020010	JAY Balance, Air, Single Valve, Microclimatic Cover, Until Seat Width 500 mm	503
JAY020012	JAY Balance, Air, Single Valve, Microclimatic Cover, From Seat Width 560 mm	503
JAY020014	JAY Balance, Air, Single Valve, Microclimatic Cover, With Positioning Kit, Until Seat Width 500 mm	568
JAY020016	JAY Balance, Air, Single Valve, Microclimatic Cover, With Positioning Kit, From Seat Width 560 mm	568
JAY020011	JAY Balance, Fluid, Microclimatic Cover, Until Seat Width 500 mm	468
JAY020013	JAY Balance, Fluid, Microclimatic Cover, From Seat Width 560 mm	468
JAY020015	JAY Balance, Fluid, Microclimatic Cover, With Positioning Kit, Until Seat Width 500 mm	533
JAY020017	JAY Balance, Fluid, Microclimatic Cover, With Positioning Kit, From Seat Width 560 mm	533
JAY020018	JAY Xtreme Active	506
JAY020023	JAY Leg Support Kit For The Xtreme Active (2)	39
JAY020019	JAY J2 Cushion, Until Seat Width 460 mm Or Seat Depth 480 mm	492
JAY020021	JAY J2 Cushion, From Seat Width 480 mm Or Seat Depth 500 mm	492
JAY020020	JAY J2 Cushion, Deep Contour, Until Seat Width 480 mm	525
JAY020022	JAY J2 Cushion, Deep Contour, From Seat Width 500 mm	525

Knee & Thigh Supports

£

QOP020342	Knee Support, Lateral, 120 x 120 mm, Swing Away, Right (3)	173
QOP020341	Knee Support, Lateral, 120 x 120 mm, Swing Away, Left (3)	173
QOP020344	Thigh Support, Lateral, Small, 120 x 120 mm, Swing Away, Right (4)	216
QOP020343	Thigh Support, Lateral, Small, 120 x 120 mm, Swing Away, Left (4)	216
QOP020347	Thigh Support, Lateral, Large, 180 x 120 mm, Swing Away, Right (4)	216
QOP020346	Thigh Support, Lateral, Large, 180 x 120 mm, Swing Away, Left (4)	216

(1) Only With Seat Width 560 mm

(2) Only For Cushions With Seat Depth 420 - 500 mm

(3) Not With QOP050015-17 SEDEO PRO, Legrest, Swing Away, QOP050080 SEDEO PRO, Hanger Receivers Only, QOP090122 Transfer Bar

(4) Not With Seat To Floor Height 415 mm, Not With QOP040035 Side Guard, Wedge Shaped, Padded

Armrest Types 		£
QOP040019	Armrest, Single Post, Flip Back (1)	0
QOP040021	Armrest, Single Post, Removable (1)	0
QOP040123	Armrest, Reclining, Flip Back	81

Arm Pads 		£
QOP040023	Arm Pad, Integral Foam, 300 mm, Right	0
QOP040024	Arm Pad, Integral Foam, 300 mm, Left	0
QOP040025	Arm Pad, Integral Foam, 400 mm, Right	0
QOP040026	Arm Pad, Integral Foam, 400 mm, Left	0
QOP040027	Arm Pad, Integral Foam, 450 mm, Right (2)	0
QOP040028	Arm Pad, Integral Foam, 450 mm, Left (2)	0
QOP040029	Omit Arm Pad, Right	0
QOP040030	Omit Arm Pad, Left	0

Comfort Arm Pads		
QOP040204	Arm Pad, Comfort, Padded, PU Coated, 400 mm, Right	124
QOP040203	Arm Pad, Comfort, Padded, PU Coated, 400 mm, Left	124

Hemiplegia Arm Pads		
QOP040031	Arm Pad, Hemiplegia, Integral Foam, Right (3)	221
QOP040032	Arm Pad, Hemiplegia, Integral Foam, Left (4)	221

Side Guard Pads For Armrest £

Side Guard Pads For Armrest, Single Post, Flip Back / Removable		
QOP040033	Side Guard (5)	0
QOP040034	Side Guard, Integral Foam (6)	0
QOP040035	Side Guard, Wedge Shaped, Padded (7)	216

Side Guard Pads For Armrest, Reclining, Flip Back		
QOP040036	Side Guard, Mounted On Seat Rail Bracket (8)	0

Centre Mount Legrest

QOP050086	Centre Mount Legrest, Tool Adjustable, With Flip Up And Angle-Adjustable Footboard (9) (12)	0
-----------	---	---

Pad Options		
QOP050087	Centre Calf Support, Integral Foam	30
QOP050092	Calf Pad, Integral Foam, Pair	0

Fixed Centre Mount Legrest

QOP050140	Fixed Centre Mount Legrest (10) (13)	0
QOP050072	Footplate, Dual, Height- And Angle-Adjustable, 260 x 150 mm	0
QOP050075	Calf Pad, Padded, Pair	0

Power Centre Mount Legrest

QOP050076	Centre Mount Legrest, Powered, 90° - 180°, 175 mm Extension (10) (11)	1188
QOP050071	Centre Mount Legrest, Powered, 90° - 180°, 175 mm Extension, Footplate To Floor (10) (11)	1404

Footplates		
QOP050072	Footplate, Dual, Height- And Angle-Adjustable, 260 x 150 mm	0
QOP050073	Footplate, Dual, Height- And Angle-Adjustable, Large, 310 x 170 mm	0

Pad Options		
QOP050075	Calf Pad, Padded, Pair	0

(1) Not With Seat To Floor Height 415 mm

(2) With QOP040019 Flip Back Armrest, Seat Width +20 mm

(3) Not With QOP110052 Joystick Mount, Right, QOP090040 Tray Table, Swing Away, QOP090046 Tray Table,

Half-Side, Swing Away, QOP090110 Rear View Mirror, Right

(4) Not With QOP110053 Joystick Mount, Left, QOP090040 Tray Table, Swing Away, QOP090046 Tray Table,

Half-Side, Swing Away, QOP090118 Rear View Mirror, Left

(5) Not With QOP040123 Armrest, Reclining, Flip Back

(6) Seat Width -20 mm Not With QOP040123 Armrest, Reclining, Flip Back

(7) Seat Width -80 mm, Only With QOP040021 Armrest, Single Post, Removable

(8) Only With QOP040123 Armrest, Reclining, Flip Back

(9) In Combination With Lift Or Tilt, Choose QOP050087 Centre Calf Support, Integral Foam Or QOP050092 Calf Pad,

Integral Foam, Pair

(10) Not With RWD

(11) Not With QOP110036 Remote Joystick, VR2; QOP110037 Remote Joystick, VR2, 2 Actuators & QOP110040

Remote Joystick, VR2, 2 Actuators And Lights. In Case Of R-Net: QOP110032 Actuator Control & Lights Module,

R-Net, ISM-6L

(12) Limited Knee Angle Possibilities With RWD/HWD

(13) In Combination With Lift Or Tilt, Choose QOP050075 Calf Pad, Padded, Pair

Powered Seating

SEDEO PRO Swing Away Legrest

£

QOP050015	SEDEO PRO, Legrest, 70°, Swing Away	0
QOP050016	SEDEO PRO, Legrest, Manual Elevating, Swing Away	292
QOP050017	SEDEO PRO, Legrest, Power Elevating, Swing Away	540
QOP050080	SEDEO PRO, Hanger Receivers Only	0
QOP050000	Omit Legrest	0

Footplates

QOP050088	Footplate, For Seat Width 385 - 400 mm	0
QOP050089	Footplate, For Seat Width 420 - 460 mm	0
QOP050090	Footplate, For Seat Width 480 - 520 mm	0
QOP050091	Footplate, For Seat Width 540 - 560 mm	0
QOP050081	Footplates, Coupled (1) (2)	0

Support Options

QOP050082	Calf Support, Integral Foam, Depth-Adjustable, Pair	86
QOP050083	Calf Support, Padded, Depth-Adjustable, Pair	103
QOP050084	Heel Loop	16
QOP050085	Calf Strap (2)	0
QOP050200	Knee Pad, Integral Foam	54

Amputee Support (3)

QOP050013	Amputee Support, Padded, Right	108
QOP050014	Amputee Support, Padded, Left	108

Positioning Supports & Accessories

£

QOP090128	Positioning Belt, Carseat Buckle, Short, 50 mm x 1700 mm	0
QOP090129	Positioning Belt, Carseat Buckle, Long, 50 mm x 2700 mm	22
QOP090110	Rear View Mirror, Right (4)	43
QOP090118	Rear View Mirror, Left (4)	43
QOP090009	Crutch Holder (5)	54
QOP090119	Brackets For Bag On Backrest Frame (5)	27
QOP090120	Bag For Armrest (6)	86
QOP090122	Transfer Bar (7)	205
QOP030605	Abduction Pommel (8)	205

Tray Table (Mounted Opposite Joystick)

QOP090040	Tray Table, Width-Adjustable, Swing Away	270
QOP090046	Tray Table, Half-Side, Swing Away (9)	216

Tray Table Accessories

QOP090125	Partial Padding For Tray Table	16
QOP090126	Raised Edges, Acrylic For Tray Table	22
QOP090127	Elbow Pads, Padded, Pair For Tray Table	32

(1) Choose Also A Footplate

(2) Only For QOP050015 SEDEO PRO, Legrest, 70°, Swing Away

(3) Extra On A Set Of Swing Away Legrests, Attached Lose

(4) Not In Combination With Same Entry Side QOP090040 Tray Table, Swing Away, QOP090046 Tray Table, Half-Side, Swing Away

(5) Not With QOP030086 Powered Recline, 90° - 130°, QOP030087 Powered Recline, 90° - 150° With 100 mm Mechanical Anti-Shear

(6) Not With QOP040123 Armrest, Reclining, Flip Back

(7) Not With QOP020342/QOP020341 Knee Support, Lateral, Swing Away, Right/Left &

QOP050015-17 SEDEO PRO, Legrest, Swing Away & QOP050080 SEDEO PRO, Hanger Receivers Only

(8) Not With (Power) Centre Mount Legrests

(9) Not With QOP090125-127 Tray Table Accessories

Joystick / Electronics**Remote Joystick – VR2**

£

QOP110036	Remote Joystick, VR2	0
QOP110037	Remote Joystick, VR2, 2 Actuators	0
QOP110040	Remote Joystick, VR2, 2 Actuators And Lights	0

Joystick Mounting

£

QOP110050	Joystick Mounting, Fixed	0
QOP110056	Joystick Mounting, Parallel, Swing Away, Capsulated	65

Joystick Mounting Location

£

QOP110052	Joystick Mounting Location, Right	0
QOP110053	Joystick Mounting Location, Left	0

Joystick Knobs (Optional)

£

QOP090505	Joystick Knob, Golf Ball, Ø 42 mm	32
QOP090506	Joystick Knob, T-Bar, 82 mm Wide	32
QOP090507	Joystick Knob, Soft Ball, Ø 70 mm	32
QOP090508	Joystick Knob, Ball, Ø 35 mm	32
QOP090045	Joystick Knob, Allround, Ball, Ø 48 mm	32
QOP090509	Joystick Knob, Chin Cup, Ø 37 mm	32
QOP090512	Joystick Knob, Stick, 80 x 15 mm	32
QOP090511	Joystick Knob, Mushroom, Ø 45 mm	32

Power Module – VR2

£

QOP110002	Power Module, VR2, 90 A	0
-----------	-------------------------	---

Lighting Module – VR2

£

QOP110030	Lighting Module, VR2 (1)	0
-----------	--------------------------	---

Attendant Controls

£

QOP110021	Attendant Control, Dual	346
QOP110020	Joystick Mounted As Attendant Control (2)	211

Attendant Mounting Location

QOP110025	Attendant Mounting Location, Right	0
QOP110026	Attendant Mounting Location, Left	0

(1) If LED Light And Indicator Is Chosen, Then This Is Required
 (2) Choose Also Remote Joystick

Joystick / Electronics**Remote Joystick – R-Net**

£

QOP110015	Remote Joystick, R-Net, LED, Actuators	378
QOP110016	Remote Joystick, R-Net, LED, Actuators And Lights	378
QOP110017	Remote Joystick, R-Net, Colour LCD	486
QOP110018	Remote Joystick, R-Net, Colour LCD, Heavy Duty	972
QOP110013	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared Control Included	594
QOP110023	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared And Bluetooth Control Included	756
QOP110019	Omit Remote Joystick, R-Net	0

Joystick Mounting

£

QOP110050	Joystick Mounting, Fixed	0
QOP110056	Joystick Mounting, Parallel, Swing Away, Capsulated	65

Joystick Mounting Location

£

QOP110052	Joystick Mounting Location, Right	0
QOP110053	Joystick Mounting Location, Left	0

Display Language Selection

£

QOP119001	English	0
-----------	---------	---

Joystick Knobs (Optional)

£

QOP090505	Joystick Knob, Golf Ball, Ø 42 mm	32
QOP090506	Joystick Knob, T-Bar, 82 mm Wide	32
QOP090507	Joystick Knob, Soft Ball, Ø 70 mm	32
QOP090508	Joystick Knob, Ball, Ø 35 mm	32
QOP090045	Joystick Knob, Allround, Ball Ø 48 mm	32
QOP090509	Joystick Knob, Chin Cup, Ø 37 mm	32
QOP090512	Joystick Knob, Stick, 80 x 15 mm	32
QOP090511	Joystick Knob, Mushroom, Ø 45 mm	32

Power Modules – R-Net

£

QOP110003	Power Module, R-Net, 90 A (1)	0
QOP110004	Power Module, R-Net, 120 A	108

Seating / Lighting Module Control – R-Net

£

QOP110032	Actuator Control & Lights Module, R-Net, ISM-6L	0
-----------	---	---

Attendant Controls

£

QOP110021	Attendant Control, Dual	346
QOP110020	Joystick Mounted As Attendant Control (1)	211

Attendant Mounting Location

QOP110025	Attendant Mounting Location, Right	0
QOP110026	Attendant Mounting Location, Left	0

(1) Choose Also Remote Joystick

Joysticks / Electronics

Table Top Controls

£

QOP110143	Table Top Control, LCD, Centre Tray Mounted	1728
-----------	---	------

Tray Table Mounting Location

QOP110080	Tray Table Mounting Location, Left	0
QOP110081	Tray Table Mounting Location, Right	0

#1 Ctrl+5 Switch Box Options

£

QOP110043	#1 Ctrl+5 Switch Box With Free Assignable Functions	702
-----------	---	-----

#1 Ctrl+5 Switch Box Mounting Location

QOP110180	#1 Ctrl+5 Switch Box Mounting Location, Behind Joystick (1)	0
QOP110181	#1 Ctrl+5 Switch Box Mounting Location, Opposite Joystick	0

#1 Ctrl+5 Switch Type

QOP110175	#1 Ctrl+5 Switch Type With Push Buttons	0
QOP110176	#1 Ctrl+5 Switch Type With Toggle Switches	0

Bluetooth Mouse Modules / Electronics – R-Net

£

QOP090068	Bluetooth Mouse Module, R-Net	232
QOP090141	2nd Bluetooth Mouse Module, R-Net (Android And Windows)	232
QOP090142	iOS Bluetooth Mouse Module, R-Net (Apple, iOS)	232

Special Controls

OMNI2 Input Control Module – R-Net (2)

£

QOP110116	OMNI2 Input Control Module, R-Net (3)	946
QOP110117	OMNI2 Input Control Module, R-Net, Scanning Mode (3)	946

OMNI2 Module Mount (4)

£

QOP110119	OMNI2 Module, Fixed Mount, Right	0
QOP110120	OMNI2 Module, Fixed Mount, Left	0
QOP110123	OMNI2 Module, Swing Away Mount, Right	65
QOP110124	OMNI2 Module, Swing Away Mount, Left	65

(1) Not Possible With QOP110056 Joystick Mounting, Parallel Swing Away, Capsulated (Here We Can't Swing Away)

(2) Choose OMNI Mount And Input Device in Sections Below And Choose Switches in Section E

(3) Display With Built In Infrared & Armrest Mount, Required For Specialty Input Devices

(4) If Used With Standard Joystick Selection, It Must Be Mounted Opposite Side

Hand Control ⁽¹⁾

£

Proportional Input Devices

QOP110301	Hand Control, Switch-It MicroGuide	2430
QOP110302	Hand Control, Switch-It MicroPilot	2700
QOP110303	Hand Control, Switch-It TouchDrive 2	3564
QOP110304	Hand Control, Switch-It VersaGuide	1026
QOP110305	Hand Control, Switch-It VersaGuide EZ ⁽²⁾	1026
QOP110306	Hand Control, Allround Joystick	864
QOP110307	Hand Control, Allround Joystick Light ⁽³⁾	1026
QOP110308	Hand Control, Allround Joystick Heavy Duty	1620
QOP110309	Hand Control, Allround Joystick, R-Net Bus Enabled ⁽⁴⁾	918
QOP110310	Hand Control, Allround Joystick Light, R-Net Bus Enabled ⁽⁴⁾	1080
QOP110311	Hand Control, Allround Joystick Heavy Duty, R-Net Bus Enabled ⁽⁴⁾	1728

Hand Control Mounting Platforms

£

QOP110342	Hand Control Mounting, Link-It To The Space Disk, Hand Tray ⁽⁵⁾	540
-----------	--	-----

**Hand Control Mounting **

QOP110350	Hand Control Mounting, Link-It To The Armrest, Fixed	162
QOP110353	Hand Control Mounting, Link-It To The Armrest, Swing Away	270

Hand Control Mounting Location

QOP110362	Hand Control Mounting Location, Right	0
QOP110361	Hand Control Mounting Location, Left	0

**Head Control ⁽⁶⁾ **

£

Proportional Input Devices

QOP110370	Head Control, Switch-It Dual Pro, 3-Way Proportional, Facial Pads	4644
QOP110371	Head Control, Switch-It Dual Pro, 3-Way Proportional, Spot Pads	4644
QOP110372	Head Control, Switch-It RIM, 3-Way Proportional	3456

Switched Input Devices

QOP110401	Head Control, 3-Switch Head Array, Wings, Short, Curved	3672
QOP110402	Head Control, 3-Switch Head Array, Wings, Short, Straight	3672
QOP110403	Head Control, 3-Switch Head Array, Wings, Long, Curved	3672
QOP110404	Head Control, 3-Switch Head Array, Wings, Long, Straight	3672
QOP110405	Head Control, 4-Switch Head Array, Wings, Long, Straight	4104
QOP110406	Head Control, 5-Switch Head Array, Wings, Long, Straight	4428
QOP110069	Head Control, 5-Switch Head Array, Wings, Long, Straight	4428

(1) For Switch Selection See Section E

(2) Not With QOP090512 Joystick Knob, Stick, 80 x 15 mm, QOP090511 Joystick Knob, Mushroom, Ø 45 mm, QOP090506 Joystick Knob, T-Bar, 82 mm Wide

(3) Only Use With Default Joystick Knob Or Joystick Knob, Allround, Ball Ø 48 mm

(4) OMNI Module Not Required

(5) Only For Micro Series Joysticks

(6) For Switch Selection See Section E. Choose A Whitmyer Headrest Mount On Previous Pages

Chin Control

£

Proportional Input Devices

QOP110470	Chin Control, Switch-It MicroGuide	2430
QOP110471	Chin Control, Switch-It MicroPilot	2700
QOP110472	Chin Control, Switch-It VersaGuide	864
QOP110473	Chin Control, Switch-It VersaGuide EZ	1026
QOP110474	Chin Control, Allround Joystick	864
QOP110475	Chin Control, Allround Joystick Light	1026
QOP110476	Chin Control, Allround Joystick, R-Net Bus Enabled (1)	918
QOP110477	Chin Control, Allround Joystick Light, R-Net Bus Enabled (1)	1080

Chin Control Mounting

QOP110480	Chin Control Mounting, Link-It To Chin Harness, Small (2)	486
QOP110481	Chin Control Mounting, Link-It To Chin Harness, Large (2)	540
QOP110483	Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Right	432
QOP110484	Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Left	432
QOP110485	Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Right	2160
QOP110486	Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Left	2160
QOP110487	Chin Control Mounting, Link-It To Headrest, Swing Away, Left	378
QOP110488	Chin Control Mounting, Link-It To Headrest, Swing Away, Right	378

Link-It Extension Kit

QOP110590	Link-It Extension Kit, Base Module, 100 mm Extension Tube, 5 Cable Clips	130
QOP110591	Link-It Extension Kit, Base Module, 150 mm Extension Tube, 5 Cable Clips	162

Sip 'N Puff (3)

£

QOP110510	Sip 'N Puff With Link-It Hardware, Swing Away Mount To The Headrest (4)	540
-----------	---	-----

(1) OMNI Module Not Required

(2) If Switches Required Only Allround Joystick With Satellite Switches Available

(3) For Switch Selection See Section E

(4) Headrest To Be Selected

Special Controls

E – Switches ⁽¹⁾

£

QOP110520	Buddy Button – Black	94
QOP110521	Buddy Button – Red	94
QOP110522	Buddy Button – Green	94
QOP110523	Egg Switch – Black	103
QOP110524	MicroLite Switch – Black	151
QOP110526	Piko Button – Black	94
QOP110527	Piko Button – Blue	94
QOP110528	Piko Button – Green	94
QOP110529	Piko Button – Red	94
QOP110530	Piko Button – Yellow	94
QOP110534	Twister Switch – Black	151
QOP110535	Twister Switch – Green	151
QOP110536	Twister Switch – Red	151

Switch Mounting

QOP110560	Switch Mounting, Link-It To The Hand Control ⁽²⁾	162
QOP110566	Switch Mounting, Link-It To The Armrest (Opposite Side As The Hand Control) ⁽³⁾	162
QOP110561	Switch Mounting, Link-It To The Headrest	216
QOP110563	Switch Mounting, Link-It To The Seat Rail ⁽⁴⁾	162
QOP110565	Switch Mounting, Link-It To Chin Boom ⁽⁵⁾	162
QOP110569	Omit Switch Mounting ⁽⁶⁾	0

Switch Mounting Location

QOP110572	Switch Mounting Location, Right	0
QOP110571	Switch Mounting Location, Left	0

Satellite Switches ⁽⁷⁾

QOP110537	Satellite Switch – Black	162
QOP110538	Satellite Switch – Green	162
QOP110539	Satellite Switch – Red	162

Multi-Function Switches ⁽⁸⁾

QOP110550	Mini 2 Button (Stereo Plug) ⁽⁹⁾	297
QOP110551	Mini 2 Button (Dual Mono Plug) ⁽⁹⁾	297
QOP110552	Mini 4 Button InLine (DB9 Connector) ⁽⁹⁾	499
QOP110554	4-Way Toggle Switched Joystick	499

Accessories

QOP110580	Mono Extension Cable (3.5 mm Jack Socket) (2 m)	108
QOP110581	Stereo Extension Cable (3.5 mm Jack Socket) (2 m)	140
QOP110582	Extension Cable DB9 (1 m)	130
QOP110583	Y-Cable, 2 Female Mono To 1 Male Stereo	162

Safety Devices

£

QOP090147	Remote Stop Switch (30 m Activation Range, Fit Via OBC Port) ⁽¹⁰⁾	756
-----------	--	-----

(1) Select 1 Or 2 Switches And Mounting Positions Accordingly. Single Function For Special Input Devices

(2) Not With HD Allround Joystick Or Space Disk Mount

(3) Armrest Rail Mount

(4) Knee, If 2 Switches Selected Only Left And Right Mount

(5) No Left / Right Switch Selection Required

(6) Switches Attached Lose

(7) Only For Allround Joystick, Maximum 2 Selectable

(8) Used In Combination With Assignable Button Programming

(9) Attached Lose, When Ordered. No Link-It Fixation Possible

(10) Not With VR2

Drive Wheels & Castor Wheels & Mudguards & Kerb Climber**◀ BACK**

// QOP070069



// QOP080213



// QOP010007



// QOP010008



// QOP90038

LED Lights & Chargers**◀ BACK**

// QOP110209



// QOP110209



// QOP160101



// QOP160104

SEDEO PRO Back Cushions & Interfaces**◀ BACK**

// QOP030109



// QOP030113



// QOP030114



// QOP030115



// QOP030116



// QOP030117

Headrest Pads & Headrest Mounts

◀ BACK




// QOP030430 & QOP030431



// WHT030401-WHT030403



// WHT030404-WHT030406



// QOP030433



// WHT030408



// WHT030409

SEDEO Seat Cushions

◀ BACK



// QOP020073



// QOP020075



// QOP020076



// QOP020078

JAY Seat Cushions[◀ BACK](#)

// JAY020005



// JAY020003



// JAY020006



// JAY020010



// JAY020011



// JAY020018



// JAY020019



// JAY020020

Armrest Types & Arm Pads[◀ BACK](#)

// QOP040019



// QOP040021



// QOP040123



// QPP040204



// QOP040027



// QOP040023

Power Centre Mount Legrest & Centre Mount Legrest[◀ BACK](#)


// QOP050076 & QOP050071



// QOP050086



// QOP050140

SEDEO PRO Swing Away Legrest & Support Options[◀ BACK](#)

// QOP050015



// QOP050016



// QOP050082

Remote Joystick– R-Net & Joystick Mounting[◀ BACK](#)

// QOP110017



// QOP110013



// QOP110018



// QOP110056

Remote Joystick – VR2[◀ BACK](#)

// QOP110036



// QOP110037

Joystick Knobs[◀ BACK](#)

// QOP090505



// QOP090045



// QOP090507



// QOP090508



// QOP090509



// QOP090512



// QOP090511



// QOP090506

Electronic Programmers – VR2[◀ BACK](#)

// QOP090063

Electronic Programmers – R-Net[◀ BACK](#)

// QOP090066



// QOP090064

Power Modules – R-Net & Seating / Lighting Module Control – R-Net[◀ BACK](#)

// QOP110004

Attendant Controls[◀ BACK](#)

// QOP110021

#1 Ctrl+5- Switch Box Options & OMNI Input Control Module – R-Net[◀ BACK](#)

// QOP110175



// QOP110176



// QOP090068



// QOP110101

Head Control

[◀ BACK](#)


// QOP110401 - QOP110406



// QOP110370 // QOP110371

Chin Control & Hand Control

[◀ BACK](#)


// QOP110470 // QOP110301



// QOP110471 // QOP110302

// QOP110472 // QOP110473 //
QOP110304 // QOP110305// QOP110306 - QOP110309
// QOP110474 - QOP110477// QOP110475 // QOP110477
// QOP110307 // QOP110310

Table Top Controls

[◀ BACK](#)


// QOP110143

Switches

◀ BACK



// QOP110521



// QOP110523



// QOP110524



// QOP110526 - QOP110530



// QOP110537 - QPP110539



// QOP110534 - QOP110536



// QOP110551



// QOP110552

Switch Mounting

◀ BACK



// QOP110526 // QOP110566

// QOP110566 // QOP110526
// QOP110529// QOP110560 // QOP110572
// QOP110013 // QOP110056
// QOP110529 // QOP110526

// QOP110563

Hand Control Mounting

◀ BACK



// QOP110353 // QOP110302
// QOP110362



// QOP110353 // QOP110302
// QOP110362 // QOP110342



// QOP110302



// QOP110353 // QOP110361
// QOP110524 // QOP110526
// QOP110560



// QOP110353 // QOP110306
// QOP110362 // QOP110538
// QOP110539

Chin Control Mounting

◀ BACK



// QOP110480 // QOP110481
// QOP110470 // QOP110471



// QOP110483 // QOP110484



// QOP110483 // QOP110484



// QOP110485 // QOP110486



// QOP110485 // QOP110474



// QOP110487 // QOP110471
// QOP110526 // QOP110529
// QOP110561 // QOP110572
// QOP110571



// QOP110487 // QOP110474

Comments For Customisations



Q500 M/F/R with SEDEO PRO

£

Base Price	
Total Of Additional Equipment	
Packaging And Transport	
Exclusive VAT	
VAT	
Total Price	

Date

Signature

Day Month Year

**Sunrise Medical**

Thorns Road, Brierley Hill, West Midland
DY5 2LD

Phone: +44 (0) 845 605 66 88

Fax: +44 (0) 845 605 66 89

E-Mail: enquiries@sunmed.co.uk

www.SunriseMedical.co.uk