

Deliver to:	
Name:	
Street:	
Town:	Postal code:
E-Mail:	
Tel:	Fax:



Invoice to:	
Name:	
Street:	
Town:	Postal code:
E-Mail:	
Tel:	Fax:

Order date:
Purchase order:
Number required:

Marked for:
Sunrise order no:

**Base Price      £4,000**

**Max. user weight: 34kg**



### TECHNICAL DATE OVERVIEW

#### Standard Features

##### Seating System - Base Configuration

- Stepless seat width, depth, back height and lower leg length
- Angle adjustable footrest
- 80° - 180° backrest recline (lie-flat)
- Canopy
- Air exchange spacer fabric with incontinence protection
- 2 point padded positioning belt
- Lateral thoracic supports (pair)
- Head supports (pair)
- Lateral thigh supports (pair)
- Medial thigh support (Pommel)

##### Frame - Base Configuration

- Foot operated brake
- 12" rear tires & 8" front castors
- Front castor suspension & swivel locks
- Quick release wheels
- Angle adjustable push handle
- Forward & rearward facing seat
- Under carriage storage
- Transit tie-down brackets (Transit Approved)

#### FRAME VERSION

199BS9	●	Frame, silver	(incl. Storage basket and canopy; Base configuration see above)
199KJ11	●	Transit Option	

#### SEATING SYSTEM

199KM01	●	Comfort Seating	(Base configuration see above)
199KM02	●	Recline Feature	

#### SEAT SIZE (Pre set)

<input type="radio"/>	Small
<input type="radio"/>	Medium
<input type="radio"/>	Large

cm	Seat Width	Seat Depth	Legrest Length	Back Height	Lateral Position
Small	24	20	15	41	13
Medium	27	27	20	49	13
Large	30	34	23	57	20

**CS note - please type 'small', 'medium' or 'large' into 'Add note' section on US configurator**

<b>Seat Depth:</b>	20 – 33cm (excl. Upholstery)
<b>Seat Width:</b>	24 – 30cm, with insert up to 15.5cm (excl. Upholstery)
<b>Back Height:</b>	41 – 61cm (excl. Upholstery)
<b>Seat Height:</b>	57.5 – 64cm (3 steps possible 57.5 – 60 – 66.4cm)
<b>Leg Length:</b>	15.5 – 31cm
<b>Angle Footplate:</b>	50°
<b>Overall Width</b>	61.5
<b>Overall Length</b>	110cm
<b>Overall Height</b>	81 - 122cm
<b>Weight Seat:</b>	9.8kg
<b>Weight Base:</b>	12.7kg

# ZIPPIE Voyage

● = Standard    ○ = No cost option    □ = Cost option

£

UPHOLSTERY (Select upholstery colour)			
● Seat upholstery with incontinence protection, removable, breathable, spacer fabric, black w/ dots			
199KN01	○ Blue trim on upholstery		
199KN02	○ Red trim on upholstery		
199KN03	○ Beige trim on upholstery		
199KN04	○ Grey trim on upholstery		
POSITIONING - LATERALS			
199KO01	○ Laterals, padded, spacer fabric, height and width-adjustable, small	(10.5cm x 15 cm; pair)	
199KO02	○ Laterals, padded, spacer fabric, height and width-adjustable, large	(11cm x 20 cm; pair)	
POSITIONING - HIP GUIDES			
199KP01	○ Hip Guides, padded, spacer fabric, width-adjustable	(pair)	
199KP03	□ Seat insert (select with small seat for seat width down to 15.5cm)		30
POSITIONING - HEAD SUPPORT			
199VJ01	○ Head support, padded, spacer fabric, height and width-adjustable, small	(10.5cm x 15 cm; pair)	
199VJ02	○ Head support, padded, spacer fabric, height and width-adjustable, large	(11cm x 20 cm; pair)	
POSITIONING - POMMEL			
199KQ02	○ Pommel, padded, spacer fabric, triangular		
199KQ01	○ Pommel, padded, spacer fabric, rectangular		
GRAB BAR			
199KR01	● Grab Bar, padded, removable (not to be use in conjunction with tray)		
TYRES			
<i>Standard with quick release rear and front wheels, and front wheel suspension with caster swivel locks</i>			
199WQ07	● 12" rear wheels (Rubber Solid rubber tyres with maintenance free inserts) and 8" front castors - solid rubber tyres		
FOOTREST			
199F329	● Footboard, platform, composite, angle-adjustable, flip-up, heel loop, black		
FOOTREST ACCESSORIES			
199A68	□ Shoe Holders, 6" long	(pair - see back page 2)	90
199A69	□ Shoe Holders, 8" long	(pair - see back page 2)	90
199A79	□ Foot Straps with Covers, 16" length - <b>only for use with above shoe holders</b>	(pair - see back page 2)	45
199A80	□ Foot straps with Covers, 18" length - <b>only for use with above shoe holders</b>	(pair - see back page 2)	45
JAY ANTERIOR TRUNK SUPPORTS			
<i>(see back page 2 for details)</i>			
Anterior Trunk Support - (please select style - opening type)			
199JY01	□ Chest strap, lateral opening, dynamic (ATS8100 Series)		120
199JY06	□ Chest strap, central opening, dynamic (ATS8800 Series)		120
Anterior Trunk Support - (please select size)			
199JZ05	○ Early Intervention (select for small seat size)	(see backpage 2)	
199JZ06	○ Xtra-Small (select for medium seat size)	(see backpage 2)	
199JZ01	○ Small (select for large seat size)	(see backpage 2)	
JAY ANTERIOR PELVIC BELTS			
<i>(see back page 2 for details)</i>			
Standard 2 Point Belt - Push Button, small (1")			
199ZZ10	○ no Pad (APB5110)		
199ZZ11	○ with Pad, small (4.5" long pads) (APB5115)		
199ZZ12	○ with Pad, medium (5.5" long pads) (APB5116)		
Standard 2 Point Belt - Push Button large (1 1/2")			
199ZZ14	○ no Pad (APB7210)		
199ZZ15	○ with Pad, large (7.5" long pads) (APB7212)		
V-Contour Pelvic 4 Point Belt			
199ZZ23	□ with Pad, small size, small Push Button (1" wide buckle) (APB3716)		120
199ZZ24	□ with Pad, medium size, large Push Button (1.5" wide buckle) (APB3627)		120

# ZIPPIE Voyage

● = Standard    ○ = No cost option    □ = Cost option

£

ACCESSORIES				£
199LF1	<input type="checkbox"/>	Hardware for LTV ventilators, dovetail slot, vertical frame mounted	(max equipment weight limit 7 Kg)	260
199LF11	<input type="checkbox"/>	Hardware for Trilogy ventilators, dovetail slot, vertical frame mounted	(max equipment weight limit 7 Kg)	260
199LF2	<input type="checkbox"/>	Battery and Vent/Utility underframe tray, central mounted, incl. 1 Utility Hook	(max equipment weight limit 23 Kg)	335
199LF3	<input type="checkbox"/>	Utility platform, front mounted	(max equipment weight limit 23 Kg)	205
199LF4	<input type="checkbox"/>	Rear add-on Accessory platform (Glider board)	(max equipment weight limit 20 Kg)	120
199A713	<input type="checkbox"/>	Bottle holder, Oxygen Tank	(max equipment weight limit 4.5 Kg)	120
199A714	<input type="checkbox"/>	IV/O2 Combination		240
199LF5	<input type="checkbox"/>	Cup holder	(max. load 0,5 kg)	20
199LF7	<input type="checkbox"/>	Utility Hook	(Attention: risk of tipping)	20
	<input type="checkbox"/>	Tray, push-fit, height-adjustable (not in combination with Grab bar)	(angle adjustable hardware included)	
199A91	<input type="checkbox"/>	Light Grey		235
199A92	<input type="checkbox"/>	Clear Poly		235

**Maximum weight limit (maximum user weight plus equipment) is 77kg**




**Table of maximum weight with different accessories - common combinations**

Maximum user weight	Storage Basket	Dovetail Vent	Central Underframe tray	Front utility mount	IV pole	Oxygen Holder	Total user & equipment weight (must not exceed 77kg)
34kg	7kg	7kg	23kg	23kg	1kg	4.5kg	
●	●						41 kg
●		●		●			64 kg
●			●				57 kg
●		●		●	●	●	69.5 kg
●			●		●	●	62.5 kg

## COMMENTS FOR SPECIAL BUILD:



<b>Zippie Voyage - Base Price</b>	£4,000
<b>Option price</b>	
<b>Total Price</b>	



**Measurements Zippie Voyage**

<b>Seat Depth:</b>	20 – 33cm (excl. Upholstery)
<b>Seat Width:</b>	24 – 30cm, with insert up to 15.5cm (excl. Upholstery)
<b>Back Height:</b>	41 – 61cm (excl. Upholstery)
<b>Seat Height:</b>	57.5 – 64cm (3 steps possible 57.5 – 60 – 66.4cm)
<b>Leg Length:</b>	15.5 – 31cm
<b>Angle Footplate:</b>	50°
<b>Overall Width</b>	61.5
<b>Overall Length</b>	110cm
<b>Overall Height</b>	81 - 122cm
<b>Weight Seat:</b>	9.8kg
<b>Weight Base:</b>	12.7kg

**Accessories for Posture, Function and Medical Equipment**



Dovetail Vertical Vent Mount (LTV/Trillogy)



Battery and Vent/Utility Underframe Tray



IV Support Pole



O2 Tank Holder



Upper Extremity Support Tray

**Lifestyle and Functional Accessories**



Front Mount Utility Platform



Rear Add On Accessory Platform



Cup Holder



Utility Hook



**Anterior Trunk Support - Style**

Lateral opening Dynamic (ATS8100)  
(Butterfly style)



	DYNAMIC	A	B
EI	ATS8107	6.5"	9.75"
XS	ATS8100	7.5"	12.5"
S	ATS8101	8.5"	15"

Centre Opening Dynamic (ATS8800 Series)



	DYNAMIC	A	B
EI	ATS8807	5"	9.75"
XS	ATS8800	5.75"	11"
S	ATS8801	6.5"	12.25"

**JAY Anterior Pelvic Belts**


Standard 2 Point Belt - Push Button 1"  
No Pad (APB5110)  
4.5" Pad (APB5115)  
5.5" Pad (APB5116)



1" WEBBING	NON-PADDED	4.5" PAD	5.5" PAD
Single Pull	APB5110	APB5115	APB5116

1.5" WEBBING	NON-PADDED	7.5" PAD
Single Pull	APB7210	APB7212

V-Contour Pelvic 4 Point Belt  
Push Button Varies



**Rectangular Pommel**



**Triangular Pommel**

