

Deliver to:	
Name:	
Street:	
Town:	Postal code:
E-Mail:	
Tel:	Fax:



Invoice to:	
Name:	
Street:	
Town:	Postal code:
E-Mail:	
Tel:	Fax:

Order date:
Purchase order:
Number required:

Marked for:
Sunrise order no:

Base Price £4,160

Max. user weight: 34kg



TECHNICAL DATE OVERVIEW

Standard Features

Seating System - Base Configuration

- Stepless seat width, depth, back height and lower leg length
- Angle adjustable footrest
- 80° - 180° backrest recline (lie-flat)
- Canopy
- Air exchange spacer fabric with incontinence protection
- 2 point padded positioning belt
- Lateral thoracic supports (pair)
- Head supports (pair)
- Lateral thigh supports (pair)
- Medial thigh support (Pommel)

Frame - Base Configuration

- Foot operated brake
- 12" rear tires & 8" front castors
- Front castor suspension & swivel locks
- Quick release wheels
- Angle adjustable push handle
- Forward & rearward facing seat
- Under carriage storage
- Transit tie-down brackets (Transit Approved)

FRAME VERSION

199BS9	●	Frame, silver	(incl. Storage basket and canopy; Base configuration see above)	
199KJ11	●	Transit Option		

SEATING SYSTEM

199KM01	●	Comfort Seating	(Base configuration see above)	
199KM02	●	Recline Feature		

SEAT SIZE (Pre set)

☐	Small	
☐	Medium	
☐	Large	

cm	Seat Width	Seat Depth	Legrest Length	Back Height	Lateral Position
Small	24	20	15	41	13
Medium	27	27	20	49	13
Large	30	34	23	57	20

CS note - please type 'small', 'medium' or 'large' into 'Add note' section on US configurator

Seat Depth:	20 – 33cm (excl. Upholstery)
Seat Width:	24 – 30cm, with insert up to 15.5cm (excl. Upholstery)
Back Height:	41 – 61cm (excl. Upholstery)
Seat Height:	57.5 – 64cm (3 steps possible 57.5 – 60 – 66.4cm)
Leg Length:	15.5 – 31cm
Angle Footplate:	50°
Overall Width	61.5
Overall Length	110cm
Overall Height	81 - 122cm
Weight Seat:	9.8kg
Weight Base:	12.7kg

= Standard = No cost option = Cost option			£
UPHOLSTERY (Select upholstery colour)			
Seat upholstery with incontinence protection, removable, breathable, spacer fabric, black w/ dots			
199KN01	Blue trim on upholstery		
199KN02	Red trim on upholstery		
199KN03	Beige trim on upholstery		
199KN04	Grey trim on upholstery		
POSITIONING - LATERALS			
199KO01	Laterals, padded, spacer fabric, height and width-adjustable, small (10.5cm x 15 cm; pair)		
199KO02	Laterals, padded, spacer fabric, height and width-adjustable, large (11cm x 20 cm; pair)		
POSITIONING - HIP GUIDES			
199KP01	Hip Guides, padded, spacer fabric, width-adjustable (pair)		
199KP03	Seat insert (select with small seat for seat width down to 15.5cm)	30	
POSITIONING - HEAD SUPPORT			
199VJ01	Head support, padded, spacer fabric, height and width-adjustable, small (10.5cm x 15 cm; pair)		
199VJ02	Head support, padded, spacer fabric, height and width-adjustable, large (11cm x 20 cm; pair)		
POSITIONING - POMMEL			
199KQ02	Pommel, padded, spacer fabric, triangular		
199KQ01	Pommel, padded, spacer fabric, rectangular		
GRAB BAR			
199KR01	Grab Bar, padded, removable (not to be use in conjunction with tray)		
TYRES Standard with quick release rear and front wheels, and front wheel suspension with caster swivel locks			
199WQ07	12" rear wheels (Rubber Solid rubber tyres with maintenance free inserts) and 8" front castors - solid rubber tyres		
FOOTREST			
199F329	Footboard, platform, composite, angle-adjustable, flip-up, heel loop, black		
FOOTREST ACCESSORIES			
199A68	Shoe Holders, 6" long (pair - see back page 2)	94	
199A69	Shoe Holders, 8" long (pair - see back page 2)	94	
199A79	Foot Straps with Covers, 16" length - only for use with above shoe holders (pair - see back page 2)	47	
199A80	Foot straps with Covers, 18" length - only for use with above shoe holders (pair - see back page 2)	47	
JAY ANTERIOR TRUNK SUPPORTS (see back page 2 for details)			
Anterior Trunk Support - (please select style - opening type)			
199JY01	Chest strap, lateral opening, dynamic (ATS8100 Series)	125	
199JY06	Chest strap, central opening, dynamic (ATS8800 Series)	125	
Anterior Trunk Support - (please select size)			
199JZ05	Early Intervention (select for small seat size) (see backpage 2)		
199JZ06	Xtra-Small (select for medium seat size) (see backpage 2)		
199JZ01	Small (select for large seat size) (see backpage 2)		
JAY ANTERIOR PELVIC BELTS (see back page 2 for details)			
Standard 2 Point Belt - Push Button, small (1")			
199ZZ10	no Pad (APB5110)		
199ZZ11	with Pad, small (4.5" long pads) (APB5115)		
199ZZ12	with Pad, medium (5.5" long pads) (APB5116)		
Standard 2 Point Belt - Push Button large (1 1/2")			
199ZZ14	no Pad (APB7210)		
199ZZ15	with Pad, large (7.5" long pads) (APB7212)		
V-Contour Pelvic 4 Point Belt			
199ZZ23	with Pad, small size, small Push Button (1" wide buckle) (APB3716)	125	
199ZZ24	with Pad, medium size, large Push Button (1.5" wide buckle) (APB3627)	125	

= Standard

= No cost option

= Cost option

£

ACCESSORIES						
199LF1	☐	Hardware for LTV ventilators, dovetail slot, vertical frame mounted			(max equipment weight limit 7 Kg)	271
199LF11	☐	Hardware for Trilogy ventilators, dovetail slot, vertical frame mounted			(max equipment weight limit 7 Kg)	271
199LF2	☐	Battery and Vent/Utility underframe tray, central mounted, incl. 1 Utility Hook			(max equipment weight limit 23 Kg)	349
199LF3	☐	Utility platform, front mounted			(max equipment weight limit 23 Kg)	214
199LF4	☐	Rear add-on Accessory platform (Glider board)			(max equipment weight limit 20 Kg)	125
199A713	☐	Bottle holder, Oxygen Tank			(max equipment weight limit 4.5 Kg)	125
199A714	☐	IV/O2 Combination				250
199LF5	☐	Cup holder			(max. load 0,5 kg)	21
199LF7	☐	Utility Hook			(Attention: risk of tipping)	21
	☐	Tray, push-fit, height-adjustable (not in combination with Grab bar)			(angle adjustable hardware included)	
199A91	☐	Light Grey				245
199A92	☐	Clear Poly				245

Maximum weight limit (maximum user weight plus equipment) is 77kg




Table of maximum weight with different accessories - common combinations

Maximum user weight	Storage Basket	Dovetail Vent	Central Underframe tray	Front utility mount	IV pole	Oxygen Holder	Total user & equipment weight (must not exceed 77kg)
34kg	7kg	7kg	23kg	23kg	1kg	4.5kg	
•	•						41 kg
•		•		•			64 kg
•			•				57 kg
•		•		•	•	•	69.5 kg
•			•		•	•	62.5 kg

COMMENTS FOR SPECIAL BUILD:

built4me

Zippie Voyage - Base Price	£4,160
Option price	
Total Price	



Measurements Zippie Voyage

Seat Depth:	20 – 33cm (excl. Upholstery)
Seat Width:	24 – 30cm, with insert up to 15.5cm (excl. Upholstery)
Back Height:	41 – 61cm (excl. Upholstery)
Seat Height:	57.5 – 64cm (3 steps possible 57.5 – 60 – 66.4cm)
Leg Length:	15.5 – 31cm
Angle Footplate:	50°
Overall Width	61.5
Overall Length	110cm
Overall Height	81 - 122cm
Weight Seat:	9.8kg
Weight Base:	12.7kg

Accessories for Posture, Function and Medical Equipment



Dovetail Vertical Vent Mount (LTV/Trilogy)



Battery and Vent/Utility Underframe Tray



IV Support Pole



O2 Tank Holder



Upper Extremity Support Tray

Lifestyle and Functional Accessories



Front Mount Utility Platform



Rear Add On Accessory Platform



Cup Holder



Utility Hook

Anterior Trunk Support - Style
Lateral opening Dynamic (ATS8100)
(Butterfly style)



	DYNAMIC	A	B
EI	ATS8107	6.5"	9.75"
XS	ATS8100	7.5"	12.5"
S	ATS8101	8.5"	15"

Centre Opening Dynamic (ATS8800 Series)



	DYNAMIC	A	B
EI	ATS8807	5"	9.75"
XS	ATS8800	5.75"	11"
S	ATS8801	6.5"	12.25"

JAY Anterior Pelvic Belts
Standard 2 Point Belt - Push Button 1"
No Pad (APB5110)
4.5" Pad (APB5115)
5.5" Pad (APB5116)



1" WEBBING	NON-PADDED	4.5" PAD	5.5" PAD
Single Pull	APB5110	APB5115	APB5116
1.5" WEBBING	NON-PADDED	7.5" PAD	
Single Pull	APB7210	APB7212	

V-Contour Pelvic 4 Point Belt
Push Button Varies

