



Supracor Stimulite Cushion Comparison Table

Model	Profile	Contour	Thickness	Weight	Best For	Key Features
Classic	Standard	Flat	2.75"	3 lbs	General use, moderate risk users	3-layer honeycomb, internal bolsters, washable, TPM(R)
Classic XS	Standard	Flat	2.75"	3 lbs	Sensitive skin, high skin protection	Softer top layer, flat for transfers, washable
Contoured	Standard	Deep	3.5"	3.25 lbs	Paraplegics, Quadraplegic	Supportive, comfort, adductors/abductor, TPM(R)
Contoured XS	Standard	Deep	3.5"	3.25 lbs	Sensitive skin, high positioning needs	Softer top layer, better positioning, washable
Slimline(TM)	Low	Mild	2.75"	3 lbs	Active users, lower risk for pressure ulcers	Lightweight, subtle contour, washable
Slimline(TM) XS	Low	Mild	2.75"	3 lbs	Active users with sensitive skin	Softer top layer, low profile, washable
Sport	Low	Flat	Varies	Light	Wheelchair athletes, active users	Breathable, mesh cover, fire-safe, washable
Corbee Paediatric	Low	Flat	2"	Light	Children, small-framed users	Small size, pediatric size, TPM(R), washable
Adjustable Contoured	Custom	Adjustable	3.5"	Varies	Users with complex positioning needs	Adjustable inserts, high customization, washable
Silver	Standard	Flat	2.75"	3 lbs	Users concerned with infection control	Anti-fungal, anti-bacterial silver-infused honeycomb, washable
Paediatric	Low	Flat	2"	Light	Children	Scaled-down cushion, washable