



Feature	Stimulite Corbee Paediatric	Stimulite Paediatric Contoured
Target User	Children who use wheelchairs or powerchairs	Children who use wheelchairs or powerchairs
Risk Level	Low risk of pressure injury	High risk of pressure injury
Cushion Type	Flat, flexible cushion for everyday use	Contoured design for complex postural and pressure care needs
Support Level	General support and comfort	Advanced postural support and positioning
Breathability	High – ventilated honeycomb design	High – ventilated honeycomb design
Pressure Relief	Gentle relief suitable for low-risk users	Deep pressure redistribution for high-risk users
Best For	Daily use at home, school, and light activities	Clinical environments, complex seating needs, and long sitting durations
Cover	Lightweight, machine washable	Durable, machine washable
Weight	Very lightweight and easy to handle	Slightly heavier due to contoured support structure
Use Environments	Home, school, travel	Home, school, hospitals, and specialist care